



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Power Retail Prescription Form

SEPTEMBER 2024
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX40 ULTRA

Standard Features

- Max user weight 136kg (21 Stone)
- Crash tested ISO 7176 -Part 19
- Ultra Low Maxx seating by Motion Concepts
- Seat widths of 405 mm - 505 mm / 16" - 20" & 480 mm - 555 mm / 19" - 22"
- Seat depth of 405 mm - 505 mm / 16" - 20" & 480 mm - 580 mm / 19" - 23"
- Backrest heights of 560mm - 660mm in 25mm increments
- Ultra Rail Mounted Flip back Cantilever Maxx Tilt Arm 230 - 320 mm (9-12.5")
- Fixed STF height and seat angle setting
- Fixed Angle Backrest with Rehab backrest
- Manual Centre Mount legrest
- LiNX control REM211 with standard bracket
- Choice of 5 Shroud colours, 5 Rim insert colours
- 9" (225 mm) x 2.80 (70 mm) Castor, pneumatic
- 14" (355 mm) x 3" (75 mm), Pneumatic drive wheels
- CTC Suspension Control - Traction - Comfort
- 6 km/h (4mph) 4 pole and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £9,095



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
JAA0001	AVIVA RX 40 Base Model Always with Performance Package: 120A Power Module + minimum 60 Ah Batteries	£9,095
JAX1300	Special Extended lead time: + 10 days	○
BASE WIDTH*		
JAL2810	CHASSIS Narrow (560 mm/22")	○
JAA2830	CHASSIS Wide (635 mm/25")	○
SEATING TYPE*		
JAA2020	Ultra Low Maxx Powered Positioning System	◆

SEAT ADJUSTMENTS*			
JAA2099		Fixed STF height and seat angle setting	(NC 2857, 2054, 2055, 2056, 1401) ◆
JAA2052		Ultra Low Maxx CG tilt (0 to 50°)	£1,322
JAA2053		Ultra Low Maxx CG tilt (0 to 45°) combined with 300 mm seat lifter increased STF Height by 4 mm	(NC 2313, 2314, 2315, 2263, 2264, 2261, 2262, 1494) £2,359
SEAT HEIGHT SETTING*			
JAA2856		Seat-to-Floor Height 420 mm	○
JAA2857		Seat-to-Floor Height 435 mm	(NC 2055, 2056, 2099) ○
JAA2858		Seat-to-Floor Height 460 mm	(NC 2055, 2056) ○
JAA2859		Seat-to-Floor Height 485 mm	(NC 2055, 2056) ○
ANTERIOR AND POSTERIOR TILT*			
JAA2057		Ultra Low Maxx CG Tilt Standard Mount No anterior or posterior tilt	◆
JAX2054		Ultra Low Maxx CG tilt 5° Fixed Anterior Assist mount Max tilt angle reduced by 5° Front STF height as defined above With centre-mounted legrests only Front STF height as defined above	(NC 2313, 2314, 2315, 2857, 2858, 2859, 2099) £349
JAX2055		Ultra Low Maxx CG tilt 10° Fixed Anterior Assist mount Max tilt angle reduced by 10° Front STF height as defined above With centre-mounted legrests only Front STF height as defined above	(NC 2312, 2313, 2314, 2315, 2857, 2858, 2859, 2099) £434
JAX2056		Ultra Low Maxx CG tilt 5° Fixed Posterior-tilt Assist mount Max tilt angle increased by 5°. Front STF height increased by 6 cm over chosen SH setting Front STF height increased by 6 cm over chosen seat height setting	(NC 2313, 2314, 2315, 0592, 2857, 2858, 2859, 2099) £349
BACKREST RECLINER*			
JAA2069		Fixed Angle Setting (NO POWERED RECLINE)	(NC 2096, 2097, 2098, 2263, 2264, 2245) ◆
JAA2060		Powered Recline 90 to 170° , Standard Mount	(NC 2510, 2511, 2512, 2513, 2514, 2515, 2516, 2540, 2541, 2542, 2543, 2261, 2262) £878
JAA2063		Powered Recline, 8° Precline mount, 82° to 162°	(NC 2510, 2511, 2512, 2513, 2514, 2515, 2516, 2540, 2541, 2542, 2543, 2261, 2262) £915
JAA2066		Powered Recline, 30° Precline mount, 60° to 140°	(NC 2510, 2511, 2512, 2513, 2514, 2515, 2516, 2540, 2541, 2542, 2543, 2261, 2262) £915
SEAT FRAME WIDTH RANGE*			
JAA2021		405 mm - 505 mm (increments of 25 mm) / 16" - 20"	(NC 2031) ○
JAA2003		480 mm - 555 mm (increments of 25 mm) / 19" - 22" Additional Delivery Time + 15 days	(NC 2026, 2027) ○
SEAT FRAME DEPTH RANGE*			
JAA2008		405 mm - 505 mm (increments of 25 mm) / 16" - 20"	(NC 2041, 2042) ○
JAA2024		480 mm - 580 mm (increments of 25 mm) / 19" - 23" Additional Delivery Time + 15 days	(NC 2036, 2037, 2810) ○
SEAT WIDTH*			
JAA2025		405 mm / 16" With centre-mount legrests only	(NC 2041, 2042, 0572, 0577, 0579, 0584, 2701, 2703, 2705) ○
JAA2026		431 mm / 17"	(NC 2041, 2042, 2701, 2703, 2705) ○
JAA2027		456 mm / 18"	(NC 2701, 2703, 2705) ○
JAA2028		480 mm / 19"	(NC 2704, 2706, 2550) ○
JAA2029		505 mm / 20"	(NC 2702, 2704, 2706) ○
JAA2030		530 mm / 21"	(NC 2035, 2036, 2145, 2702, 2704, 2706) ○
JAA2031		555 mm / 22"	(NC 2035, 2036, 2145, 2702, 2704, 2706) ○

SEAT DEPTH SETTING*				
JAA2035		405 mm / 16"	(NC 2030, 2031, 2701, 2702, 2703, 2704, 2705, 2706)	<input type="radio"/>
JAA2036		430 mm / 17"	(NC 2030, 2031)	<input type="radio"/>
JAA2037		455 mm / 18"		<input type="radio"/>
JAA2038		480 mm / 19"		<input type="radio"/>
JAA2039		505 mm / 20"		<input type="radio"/>
JAA2040		530 mm / 21"		<input type="radio"/>
JAA2041		555 mm / 22"	(NC 2025, 2026, 2145)	<input type="radio"/>
JAA2042		580 mm / 23"	(NC 2025, 2026, 2145, 2146)	<input type="radio"/>
SEAT CUSHION*				
JAA2144		Matrx Vi Cushion [Mesh fabric]		£360
JAA2145		Matrx Flo-Tech Lite Black smooth cover, 75 mm thick Up to 127 kg	(NC 2042, 2030, 2031)	£117
JAA2146		Matrx Flo-tech Contour Black smooth cover, 100 mm thick Up to 191 kg		£165
JAA2147		Matrx Libra Cushion 90 mm thick Standard: up to 136 kg Extended lead time: + 4 days Additional Delivery Time + 4 days		£388
JAA2140		Without seat cushion		<input type="radio"/>
BACK CANE HEIGHT - STRAIGHT OR ANGLED BACK CANES*				
JAA2510		Straight back cane height 455 mm / 18"	(NC 2542, 2543, 2050, 2051, 2320, 2551, 2552, 2553, 2554, 2555, 2556, 2557, 2177, 2558, 2559, 2178, 2174, 2175, 2560, 2561, 2562, 2563)	<input type="radio"/>
JAA2511		Straight back cane height 505 mm / 20"	(NC 2543, 2051, 2320, 2552, 2553, 2555, 2557, 2177, 2563, 2559, 2178, 2175, 2561, 2562)	<input type="radio"/>
JAA2512		Straight back cane height 555 mm / 22"	(NC 2320, 2553, 2555, 2177, 2563, 2178, 2175, 2561, 2562)	<input type="radio"/>
JAA2513		Straight back cane height 610 mm / 24"		<input type="radio"/>
JAA2514		Angled (5°) back cane height 505 mm / 20"	(NC 2543, 2051, 2320, 2552, 2553, 2555, 2557, 2177, 2563, 2559, 2178, 2175, 2561, 2562)	<input type="radio"/>
JAA2515		Angled (5°) back cane height 555 mm / 22"	(NC 2320, 2553, 2555, 2177, 2563, 2178, 2175, 2561, 2562)	<input type="radio"/>
JAA2516		Angled (5°) back cane height 610 mm / 24"		<input type="radio"/>
JAX3007		Please select If powered recline is required		<input type="radio"/>
ANGLE SETTING FOR STRAIGHT BACK CANES*				
JAA2517		75°	(NC 2515, 2516)	<input type="radio"/>
JAA2518		79°	(NC 2515, 2517)	<input type="radio"/>
JAA2519		85°	(NC 2515, 2518)	<input type="radio"/>
JAA2520		90°	(NC 2515, 2519)	<input type="radio"/>
JAA2521		95°	(NC 2515, 2520)	<input type="radio"/>
JAA2522		101°	(NC 2515, 2521)	<input type="radio"/>
JAA2523		106°	(NC 2515, 2522)	<input type="radio"/>
JAA2524		112°	(NC 2515, 2523)	<input type="radio"/>
JAA2525		116°	(NC 2515, 2524)	<input type="radio"/>
JAX3006		Please select if 5° Angle Back post or powered recline chosen		<input type="radio"/>

ANGLE SETTING FOR 5° ANGLED BACK CANES*				
JAA2526		81°	(NC 2511, 2512, 2513)	○
JAA2527		84°	(NC 2511, 2512, 2514)	○
JAA2528		91°	(NC 2511, 2512, 2515)	○
JAA2529		96°	(NC 2511, 2512, 2516)	○
JAA2530		101°	(NC 2511, 2512, 2517)	○
JAA2531		106°	(NC 2511, 2512, 2518)	○
JAA2532		111°	(NC 2511, 2512, 2519)	○
JAA2533		117°	(NC 2511, 2512, 2520)	○
JAA2534		121°	(NC 2511, 2512, 2521)	○
JAX3002		Straight Back post choosen or powered recline choosen		○
BACKREST SUPPORT*				
JAA2540		Standard Rehab Back Pan and Back Upholstery height 355 mm / 14"	(NC 2063, 2066)	○
JAA2541		Standard Rehab Back Pan and Back Upholstery height 405 mm / 16"	(NC 2063, 2066, 2510)	○
JAA2542		Standard Rehab Back Pan and Back Upholstery height 455 mm / 18"	(NC 2063, 2066, 2510, 2511, 2514)	○
JAA2543		Standard Rehab Back Pan and Back Upholstery height 505 mm / 20"	(NC 2063, 2066, 2510, 2511, 2512, 2514, 2515)	○
JAA2050		High Back, height 405 mm / 16" [Mesh fabric]		○
JAA2051		High Back, height 455 mm / 18" [Mesh fabric]	(NC 2511, 2514)	○
JAA2320		High Back, height 505 mm / 20" [Mesh fabric]	(NC 2511, 2512, 2514, 2515)	○
JAA2550		Matrx Elite E2 Back, height 355mm / 14" / 16 - 18" wide	(NC 2029)	£460
JAA2551		Matrx Elite E2 Back, height 405 mm / 16"		£460
JAA2552		Matrx Elite E2 Back, height 455 mm / 18"	(NC 2511, 2514)	£460
JAA2553		Matrx Elite E2 Back, height 505 mm / 20"	(NC 2511, 2512, 2514, 2515)	£460
JAA2554		Matrx Elite HD Back, height 405 mm /16" Extended lead time: + 4 days Additionnal Delivery Time + 4 days		£509
JAA2555		Matrx Elite HD Back, height 505 mm / 20" Extended lead time: + 4 days Additionnal Delivery Time + 4 days	(NC 2511, 2512, 2514, 2515)	£509
JAA2556		Matrx Elite TR Back, height 405 mm / 16" [Mesh fabric]		£496
JAA2557		Matrx Elite TR Back, height 455 mm / 18" [Mesh fabric]	(NC 2511, 2514)	£496
JAA2177		Matrx Elite TR Back, height 505 mm / 20" [Mesh fabric]	(NC 2511, 2512, 2514, 2515)	£496
JAA2563		Matrx Elite TR HD Back, height 505 mm / 20" [Mesh fabric] Extended lead time: + 4 days Additionnal Delivery Time + 4 days	(NC 2511, 2512, 2514, 2515)	£550
JAA2558		Matrx Elite Deep E2 Back, height 405 mm / 16" [Mesh fabric]		£550
JAA2559		Matrx Elite Deep E2 Back, height 455 mm / 18" [Mesh fabric]	(NC 2511, 2514)	£550
JAA2178		Matrx Elite Deep E2 Back, height 505 mm / 20" [Mesh fabric]	(NC 2511, 2512, 2514, 2515)	£550
JAA2174		Matrx PB Back, height 405 mm / 16" [Mesh fabric]		£496
JAA2175		Matrx PB Back, height 505 mm / 20" [Mesh fabric]	(NC 2511, 2512, 2514, 2515)	£496
JAA2560		Matrx PB Deep Back, height 405mm / 16"		£509
JAA2561		Matrx PB Deep Back, height 505 mm / 20"	(NC 2511, 2512, 2514, 2515)	£509
JAA2562		Matrx Elite Deep HD Back, height 505 mm / 20" Extended lead time: + 4 days Additionnal Delivery Time + 4 days	(NC 2511, 2512, 2514, 2515)	£613
JAA2570		Without Backrest Support, incl. Matrx interface	(NC 2322, 2323, 2324, 2325, 2326, 2327, 2230, 2247)	○

BACKREST HEIGHT (NON-RECLINE)*			
JAA2504		455 mm / 18" ONLY WITH: Back height 355mm/14", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2542, 2543, 2050, 2051, 2320, 2551, 2552, 2553, 2554, 2555, 2556, 2557, 2177, 2558, 2559, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) ○
JAA2505		485 mm / 19" ONLY WITH: Back height 355mm/14", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2542, 2543, 2050, 2051, 2320, 2551, 2552, 2553, 2554, 2555, 2556, 2557, 2177, 2558, 2559, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) ○
JAA2506		505 mm / 20" ONLY WITH: Back height 405mm/16", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2542, 2543, 2051, 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562, 2563, 2570) ○
JAA2501		535 mm / 21" ONLY WITH: Back height 405mm/16", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2542, 2543, 2051, 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562, 2563, 2570) ○
JAA2070		555 mm / 22" ONLY WITH Back height 455mm/18", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2541, 2543, 2050, 2320, 2550, 2551, 2553, 2554, 2555, 2556, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) ○
JAA2071		585 mm / 23" ONLY WITH Back height 455mm/18", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2541, 2543, 2050, 2320, 2550, 2551, 2553, 2554, 2555, 2556, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) ○
JAA2072		610 mm / 24" ONLY WITH Back height 505mm/20", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2541, 2542, 2050, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560, 2570) ○
JAA2073		635 mm / 25" ONLY WITH Back height 505mm/20", 100-125 mm/4-5" GAP IS RECOMMENDED FOR TILT-ONLY SYSTEMS (NON-RECLINE)	(NC 2541, 2542, 2050, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560, 2570) ○
JAX3020		Powered Recline option selected	○
BACKREST HEIGHT (ESR RECLINE SYSTEMS)*			
JAA2506		505 mm / 20" ONLY WITH: Elite Back, height 355mm/14"(2550), MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2051, 2320, 2551, 2552, 2553, 2554, 2555, 2556, 2557, 2177, 2558, 2559, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) ○
JAA2501		535 mm / 21" ONLY WITH: Elite Back, height 355mm/14"(2550), MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2051, 2320, 2551, 2552, 2553, 2554, 2555, 2556, 2557, 2177, 2558, 2559, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) ○
JAA2070		555 mm / 22" ONLY WITH: Back height 405mm/16", MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562, 2563, 2570) ○
JAA2071		585 mm / 23" ONLY WITH: Back height 405mm/16", MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2320, 2550, 2552, 2553, 2555, 2557, 2177, 2559, 2178, 2175, 2561, 2562, 2563, 2570) ○
JAA2072		610 mm / 24" ONLY WITH Back height 455mm/18", MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2320, 2550, 2551, 2553, 2554, 2555, 2556, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563) ○
JAA2073		635 mm / 25" ONLY WITH Back height 455mm/18", MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2320, 2550, 2551, 2553, 2554, 2555, 2556, 2177, 2558, 2178, 2174, 2175, 2560, 2561, 2562, 2563, 2570) ○
JAA2074		660 mm / 26" ONLY WITH Back height 505mm/20", MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2051, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560, 2570) ○
JAA2502		685 mm / 27" ONLY WITH Back height 505mm/20", MIN. 150-175 mm / 6-7" GAP IS RECOMMENDED FOR ESR RECLINE SYSTEMS	(NC 2051, 2550, 2551, 2552, 2554, 2556, 2557, 2558, 2559, 2174, 2560, 2570) ○
JAA3020		NON -RECLINE BACKREST SELECTED	(NC 2063, 2066, 2570) ○
ARMREST*			
JAA2307		Ultra Rail Mounted Flip back Cantilever Maxx Tilt Arm 230 - 320 mm (9-12.5") Not compatible with reclining backrest	(NC 2063, 2066, 2308, 2309) ◆
JAA2095		Cantilever Flip-Back Arm, width, height and angle adjustable, 240 - 330 mm (9.5" - 13")	(NC 2308, 2309) ○
JAA2096		Cantilever Flip-Back Arm, width, height and angle adjustable, 340 - 405 mm (13.5-16")	(NC 2308, 2309) ○
JAA2097		Two-post Recline Arm, Flip Back, adjustable height 230 - 330 mm (9-13")	○
JAA2098		Two-post Recline Arm, Flip Back, adjustable height 320 - 405 mm (12.5-16")	○

ARMREST ACCESSORIES			
JAA2308		1" Offset Armrest Spacer (pair), for two-post recline only (NC 2096, 2307)	£143
JAA2309		1/2" Offset Armrest Spacer (pair), for two-post recline only (NC 2096, 2307)	£143
ARMPAD LH*			
JAA2012		Desk arm pad with startex cover visco foam [255mm]; 40 mm thick	◆
JAA2013		Full length arm pad with startex cover visco foam [355mm], 40 mm thick	○
JAA2014		Ergonomic arm trough; 30 mm thick (NC 1405, 2311, 2310, 2347, 2348, 2349, 2350)	£117
JAA2015		without LH arm pad	○
ARMPAD RH*			
JAA2016		Desk arm pad with startex cover visco foam [255mm]; 40 mm thick	◆
JAA2017		Full length arm pad with startex cover visco foam [355mm]; 40 mm thick	○
JAA2018		Ergonomic arm trough, 30 mm thick (NC 1405, 2311, 2310, 2351, 2352, 2353, 2354)	£117
JAA2019		without RH arm pad	○
ARMPAD ACCESSORIES & MOUNTING OPTIONS LH			
JAA2080		Flat Hand Pad Left (NC 2013, 2015)	£117
JAA2081		Ergonomic Mounting Hardware to mount Otto Bock (EMO) hand pad options Left (NC 2013, 2015)	£52
ARMPAD ACCESSORIES & MOUNTING OPTIONS RH			
JAA2083		Flat Hand Pad Right (NC 2017, 2019)	£117
JAA2084		Ergonomic Mounting Hardware to mount Otto Bock (EMO) hand pad options Right (NC 2017, 2019)	£52
MULTI-AXIS UPPER EXTREMITY SUPPORTS LH			
JAA2082		Multi Axis Upper Extremity Support (MACES) Left, lowest armpad height will be increased by 43 mm	£295
MULTI-AXIS UPPER EXTREMITY SUPPORTS RH			
JAA2085		Multi Axis Upper Extremity Support (MACES) Right, lowest armpad height will be increased by 43 mm	£295

Height setting: ARMREST LH SEAT PAN TO TOP OF THE TUBE*

O001	<input type="checkbox"/>	Factory setting	<input type="radio"/>
O002	<input type="checkbox"/>	Customized	<input type="radio"/>
O003	<input type="checkbox"/>	230 mm	<input type="radio"/>
O004	<input type="checkbox"/>	240 mm	<input type="radio"/>
O005	<input type="checkbox"/>	250 mm	<input type="radio"/>
O006	<input type="checkbox"/>	260 mm	<input type="radio"/>
O007	<input type="checkbox"/>	270 mm	<input type="radio"/>
O008	<input type="checkbox"/>	280 mm	<input type="radio"/>
O009	<input type="checkbox"/>	290 mm	<input type="radio"/>
O010	<input type="checkbox"/>	300 mm	<input type="radio"/>
O011	<input type="checkbox"/>	310 mm	<input type="radio"/>
O012	<input type="checkbox"/>	320 mm	<input type="radio"/>
O013	<input type="checkbox"/>	330 mm	<input type="radio"/>
O014	<input type="checkbox"/>	340 mm	<input type="radio"/>
O015	<input type="checkbox"/>	350 mm	<input type="radio"/>
O016	<input type="checkbox"/>	360 mm	<input type="radio"/>
O017	<input type="checkbox"/>	370 mm	<input type="radio"/>
O018	<input type="checkbox"/>	380 mm	<input type="radio"/>
O019	<input type="checkbox"/>	390 mm	<input type="radio"/>
O020	<input type="checkbox"/>	405 mm	<input type="radio"/>

Height setting: ARMREST RH SEAT PAN TO TOP OF THE TUBE*

O021	<input type="checkbox"/>	Factory setting	<input type="radio"/>
O022	<input type="checkbox"/>	Customized	<input type="radio"/>
O023	<input type="checkbox"/>	230 mm	<input type="radio"/>
O024	<input type="checkbox"/>	240 mm	<input type="radio"/>
O025	<input type="checkbox"/>	250 mm	<input type="radio"/>
O026	<input type="checkbox"/>	260 mm	<input type="radio"/>
O027	<input type="checkbox"/>	270 mm	<input type="radio"/>
O028	<input type="checkbox"/>	280 mm	<input type="radio"/>
O029	<input type="checkbox"/>	290 mm	<input type="radio"/>
O030	<input type="checkbox"/>	300 mm	<input type="radio"/>
O031	<input type="checkbox"/>	310 mm	<input type="radio"/>
O032	<input type="checkbox"/>	320 mm	<input type="radio"/>
O033	<input type="checkbox"/>	330 mm	<input type="radio"/>
O034	<input type="checkbox"/>	340 mm	<input type="radio"/>
O035	<input type="checkbox"/>	350 mm	<input type="radio"/>
O036	<input type="checkbox"/>	360 mm	<input type="radio"/>
O037	<input type="checkbox"/>	370 mm	<input type="radio"/>
O038	<input type="checkbox"/>	380 mm	<input type="radio"/>
O039	<input type="checkbox"/>	390 mm	<input type="radio"/>
O040	<input type="checkbox"/>	405 mm	<input type="radio"/>

ELBOW BLOCK LH PADS			
JAA2347		Elbow Block XS (110 mm wide x 83 mm high)	£553
JAA2348		Elbow Block S (135 mm wide x 90 mm high)	£553
JAA2349		Elbow Block M (145 mm wide x 95 mm high)	£553
JAA2350		Elbow Block L (165 mm wide x 110 mm high)	£553
ELBOW BLOCK LH MOUNTING HARDWARE			
JAA2361		Vertical 3"/ Horizontal 5"	○
JAA2362		Vertical 3"/ Horizontal 7"	○
JAA2363		Vertical 5"/ Horizontal 5"	○
JAA2364		Vertical 5"/ Horizontal 7"	○
ELBOW BLOCKS RH PADS			
JAA2351		Elbow Block XS (110 mm wide x 83 mm high)	£553
JAA2352		Elbow Block S (135 mm wide x 90 mm high)	£553
JAA2353		Elbow Block M (145 mm wide x 95 mm high)	£553
JAA2354		Elbow Block L (165 mm wide x 110 mm high)	£553
ELBOW BLOCK RH MOUNTING HARDWARE			
JAA2371		Vertical 3"/ Horizontal 5"	○
JAA2372		Vertical 3"/ Horizontal 7"	○
JAA2373		Vertical 5"/ Horizontal 5"	○
JAA2374		Vertical 5"/ Horizontal 7"	○
HIP SUPPORTS LH			
JAA2184		Hip Support LH, Quick Release, pad size 140 x 100 mm (NC 2702, 2703, 2704, 2705, 2706)	£252
JAA2114		Hip Support LH, Quick Release, pad size 180 x 150 mm (NC 2702, 2703, 2704, 2705, 2706)	£252
JAA2186		Hip Support LH, Quick Release, pad size 230 x 100 mm (NC 2702, 2703, 2704, 2705, 2706)	£252
JAA2188		Hip Support LH, Quick Release, pad size 330 x 100 mm (NC 2702, 2703, 2704, 2705, 2706)	£252
HIP SUPPORTS RH			
JAA2185		Hip Support RH, Quick Release, pad size 140 x 100 mm (NC 2702, 2703, 2704, 2705, 2706)	£252
JAA2115		Hip Support RH, Quick Release, pad size 180 x 150 mm (NC 2702, 2703, 2704, 2705, 2706)	£252
JAA2187		Hip Support RH, Quick Release, pad size 230 x 100 mm (NC 2702, 2703, 2704, 2705, 2706)	£252
JAA2189		Hip Support RH, Quick Release, pad size 330 x 100 mm (NC 2702, 2703, 2704, 2705, 2706)	£252
MOUNTING HARDWARE FOR HIP SUPPORTS			
JAA2164		Mounting hardware for hip supports (6,1" / 155 mm high)	○
JAA2165		Mounting hardware for hip supports (8,1" / 206 mm high)	○

HEADREST*			
JAA2092		Onyx Auto Style headrest with standard mounting hardware (NC 2230, 2247, 2570)	£128
JAA2323		Elan Standard Pad 150 mm / 6" wide w/standard cover (ESP26), complete with Elan Headrest Hardware (MEHW) (NC 2230, 2247, 2570)	£294
JAA2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26), complete with Elan Headrest Hardware (MEHW) (NC 2230, 2247, 2570)	£294
JAA2322		Elan Standard Pad 355 mm / 14" wide w/standard cover (ESP26), complete with Elan Headrest Hardware (MEHW) (NC 2230, 2247, 2570)	£294
JAA2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23), complete with Elan Headrest Hardware (MEHW) (NC 2230, 2247, 2570)	£294
JAA2326		Elan Occipital Pad 305x200 mm / 12"x8" w/standard cover complete with Elan Headrest Hardware (MEHW) (NC 2230, 2247, 2570)	£294
JAA2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover complete with Elan Headrest Hardware (MEHW) (NC 2230, 2247, 2570)	£396
JAA1495		Head Control or equivalent Headrest, provided independently from this order According to the requirements of ISO 7176-19, a headrest is compulsory for using the wheelchair as vehicle seat	○
LATERAL SUPPORTS LH			
JAA2343		Lateral Support Maxx Style Quick Release, pad size [110 mm wide x 83 mm high] LH (NC 2559, 2178, 2560, 2561, 2562)	£427
JAA2345		Lateral Support Maxx Style Quick Release, pad size [135 mm wide x 90 mm high] LH (NC 2559, 2178, 2560, 2561, 2562)	£427
JAA2043		Lateral Support Maxx Style Quick Release, pad size [145 mm wide x 95 mm high] LH (NC 2559, 2178, 2560, 2561, 2562)	£427
JAA2045		Lateral Support Maxx Style Quick Release, pad size [165 mm wide x 110 mm high] LH (NC 2559, 2178, 2560, 2561, 2562)	£427
LATERAL SUPPORTS RH			
JAA2344		Lateral Support Maxx Style Quick Release, pad size [110 mm wide x 83 mm high] RH (NC 2559, 2178, 2560, 2561, 2562)	£427
JAA2346		Lateral Support Maxx Style Quick Release, pad size [135 mm wide x 90 mm high] RH (NC 2559, 2178, 2560, 2561, 2562)	£427
JAA2044		Lateral Support Maxx Style Quick Release, pad size [145 mm wide x 95 mm high] RH (NC 2559, 2178, 2560, 2561, 2562)	£427
JAA2046		Lateral Support Maxx Style Quick Release, pad size [165 mm wide x 110 mm high] RH (NC 2559, 2178, 2560, 2561, 2562)	£427
LATERAL SUPPORT ACCESSORIES			
JAA2049		Additional telescopic link for Maxx Style lateral support, 25 mm, Pair (NC 2559, 2178, 2560, 2561, 2562)	£169
CENTRE MOUNTED LEGREST*			
JAA0586		Centre mount legrest, angle settable Length: 310 - 410 mm (NC 0517, 1244, 0579, 0584, 0572, 0577)	◆
JAA0589		LNx Power centre mount legrest, length compensating; Length: 340 - 410 mm (NC 2262, 2263, 2264, 1244, 0579, 0584, 0572, 0577)	£1,640
JAX0592		LNx Power centre mount legrest with telescoping foot platform (lowers to the ground; Length: 380 - 419 mm) Selectable length depends on STF height. Custom-made, requires customized setting! Custom-made, requires customized setting! (NC 0516, 0517, 2056, 2261, 2262, 2263, 2264, 2365, 2366, 1244, 0579, 0584, 0572, 0577)	£2,912
JAX3003		Without Centre Mount Legrests	○

FOOTPLATES / FOOT PLATFORM*			
JAA2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5") (NC 3003)	£103
JAA0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5") (NC 3003)	○
JAA0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5") (NC 3003)	£135
JAA2366		Angle adjustable foot platform S (305 W x 215 mm D / 12"x8.5") (NC 0586, 3003)	£270
JAA0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75") (NC 0586, 3003)	£270
JAX3005		Please select if Centre mount footPlate/ platform is not required. (NC 0589)	○
Height setting: CENTRE MOUNT LEGREST, OPTIONS 0586 & 0589 . SEAT PAN TO FOOTPLATE*			
0041		Factory setting	○
0042		Customized	○
0043		310 mm	○
0044		325 mm	○
0045		335 mm	○
0046		345 mm	○
0047		355 mm	○
0048		364 mm	○
0049		375 mm	○
0050		385 mm	○
0051		410 mm	○
0077		Please select if CENTRE MOUNT LEGREST OPTIONS 0586 & 0589 is not required	○
Height setting: CENTRE MOUNT LEGREST with TELESCOPING FOOT PLATFORM . SEAT PAN TO FOOTPLATE*			
0052		Factory setting	○
0053		Customized	○
0054		380 mm	○
0055		393 mm	○
0056		406 mm	○
0057		419 mm	○
0076		Please select if CENTRE MOUNT LEGREST with TELESCOPING FOOT PLATFORM is not required	○
SIDE MOUNTED HANGERS LH*			
JAA0539		Hanger ADE 0° - 80° LH (NC 2055, 2261, 2262, 2263, 2264)	£299
JAA0572		Pivot Plus Manual Elevating Legrest LH	£655
JAA0577		Pivot Plus Power Elevating Legrest LH (NC 2261, 2262, 2263, 2264)	£1,038
JAA0535		Without Hanger LH	○
SIDE MOUNTED HANGERS RH*			
JAA0549		Hanger ADE 0° - 80° RH (NC 2261, 2262, 2263, 2264)	£299
JAA0579		Pivot Plus Manual Elevating Legrest RH	£655
JAA0584		Pivot Plus Power Elevating Legrest RH (NC 2261, 2262, 2263, 2264)	£1,038
JAA0545		Without Hanger RH	○
INDIVIDUAL FOOTPLATE LH*			
JAA0509		Angle Adjustable Footplate LH	£16
JAA0511		Angle and depth adjustable Footplate LH	£60
JAX3012		Please select if Side mounted Legrest are not required	○

INDIVIDUAL FOOTPLATE RH*

JAA0506	<input type="checkbox"/>	Angle Adjustable Footplate RH	£16
JAA0508	<input type="checkbox"/>	Angle and depth adjustable Footplate RH	£60
JAX3013	<input type="checkbox"/>	Please select if Side mounted Legrest are not required	<input type="radio"/>

Height setting: VARI LEGREST LH SEAT PAN TO FOOTPLATE*

0080	<input type="checkbox"/>	Factory Setting	<input type="radio"/>
0081	<input type="checkbox"/>	Customized	<input type="radio"/>
0082	<input type="checkbox"/>	270 mm	<input type="radio"/>
0083	<input type="checkbox"/>	280 mm	<input type="radio"/>
0084	<input type="checkbox"/>	290 mm	<input type="radio"/>
0085	<input type="checkbox"/>	300 mm	<input type="radio"/>
0086	<input type="checkbox"/>	310 mm	<input type="radio"/>
0087	<input type="checkbox"/>	320 mm	<input type="radio"/>
0088	<input type="checkbox"/>	330 mm	<input type="radio"/>
0089	<input type="checkbox"/>	340 mm	<input type="radio"/>
0090	<input type="checkbox"/>	350 mm	<input type="radio"/>
0091	<input type="checkbox"/>	360 mm	<input type="radio"/>
0092	<input type="checkbox"/>	370 mm	<input type="radio"/>
0093	<input type="checkbox"/>	380 mm	<input type="radio"/>
0094	<input type="checkbox"/>	390 mm	<input type="radio"/>
0095	<input type="checkbox"/>	400 mm	<input type="radio"/>
0096	<input type="checkbox"/>	410 mm	<input type="radio"/>
0097	<input type="checkbox"/>	420 mm	<input type="radio"/>
0098	<input type="checkbox"/>	430 mm	<input type="radio"/>

Height setting: VARI LEGREST RH SEAT PAN TO FOOTPLATE*

0099	<input type="checkbox"/>	Factory Setting	<input type="radio"/>
0100	<input type="checkbox"/>	Customized	<input type="radio"/>
0101	<input type="checkbox"/>	270 mm	<input type="radio"/>
0102	<input type="checkbox"/>	280 mm	<input type="radio"/>
0103	<input type="checkbox"/>	290 mm	<input type="radio"/>
0104	<input type="checkbox"/>	300 mm	<input type="radio"/>
0105	<input type="checkbox"/>	310 mm	<input type="radio"/>
0106	<input type="checkbox"/>	320 mm	<input type="radio"/>
0107	<input type="checkbox"/>	330 mm	<input type="radio"/>
0108	<input type="checkbox"/>	340 mm	<input type="radio"/>
0109	<input type="checkbox"/>	350 mm	<input type="radio"/>
0110	<input type="checkbox"/>	360 mm	<input type="radio"/>
0111	<input type="checkbox"/>	370 mm	<input type="radio"/>
0112	<input type="checkbox"/>	380 mm	<input type="radio"/>
0113	<input type="checkbox"/>	390 mm	<input type="radio"/>
0114	<input type="checkbox"/>	400 mm	<input type="radio"/>
0115	<input type="checkbox"/>	410 mm	<input type="radio"/>
0116	<input type="checkbox"/>	420 mm	<input type="radio"/>
0117	<input type="checkbox"/>	430 mm	<input type="radio"/>

Height setting: PIVOT PLUS LEGREST LH SEAT PAN TO FOOTPLATE*			
0058		Factory setting	○
0059		Customized	○
0060		402 mm	○
0061		414 mm	○
0062		426 mm	○
0063		438 mm	○
0064		450 mm	○
0065		462 mm	○
0066		474 mm	○
0080		Please select if PIVOT PLUS LEGREST is not required	○
Height setting: PIVOT PLUS LEGREST RH SEAT PAN TO FOOTPLATE*			
0067		Factory setting	○
0068		Customized	○
0069		402 mm	○
0070		414 mm	○
0071		426 mm	○
0072		438 mm	○
0073		450 mm	○
0074		462 mm	○
0075		474 mm	○
0079		Please select if PIVOT PLUS LEGREST is not required	○
CONTROLS*			
JAA2203		Remote LiNX REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators	(NC 1437, 2231, 2232, 2229, 2233, 2230, 2247, 2257, 2240, 2243, 2271, 2255, 2256, 2270, 2272, 1908, 1909, 2209, 2210) ◆
JAA2204		Remote LiNX REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights	(NC 2232, 2229, 2233, 2230, 2247, 2257, 2240, 2243, 2271, 2255, 2256, 2270, 2272, 2209, 2210) £151
JAA2205		Rehab Remote LiNX REM400, touch-screen user interface Configurable jack sockets	(NC 2210) £341
JAA2206		Rehab Remote LiNX REM500 without joystick for alternative access users Requires a user input/alternative joystick: 2231, 2232, 2229, 2230, 2247, 2234	(NC 1207, 1208, 1707, 1708, 1709, 1710) £776
REM400 TOGGLE OPTION			
JAA2270		Assignable Toggle Switch Module f. REM400 (DLX-SW400TGL)	(NC 2204, 2206) £316

CONTROL BRACKETRY*				
JAA1207		Standard Bracket RH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404, 2311)	○
JAA1208		Standard Bracket LH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405, 2310)	○
JAA1707		Swing Away Bracket RH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404, 2311)	£133
JAA1708		Swing Away Bracket LH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405, 2310)	£133
JAA1709		Matrx Resolve Swing in/out bracket RH	(NC 2283, 2287, 2207, 2265, 2275, 2235, 2267, 2241, 1404)	£219
JAA1710		Matrx Resolve Swing in/out bracket LH	(NC 2284, 2288, 2208, 2266, 2276, 2236, 2268, 2242, 1405)	£219
JAA2209		Rem500 Swing away RH, multi-adjustable	(NC 2204, 2205, 2285, 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404, 2311)	○
JAA2210		Rem500 Swing away LH, multi-adjustable	(NC 2204, 2205, 2286, 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1405, 2310)	○
JAA2075		Swing Away + height / width adjust. Mount on RH for user, Quad link design		£244
JAA2076		Swing Away + height / width adjust. Mount on LH for user, Quad link design		£244
JAA2213		Centre Mount, swing-away to RH, Nucleus	(NC 2283, 2287, 2207, 2265, 2275, 2217, 1404, 1405, 2310, 2311, 1415)	£498
JAA2214		Centre Mount, swing-away to LH, Nucleus	(NC 2284, 2288, 2208, 2266, 2276, 2218, 1404, 1405, 2310, 2311, 1339)	£498
JAA2227		Bracket RH, mounted for Attendant	(NC 1221, 2206, 1287, 2219)	£64
JAA2228		Bracket LH, mounted for Attendant	(NC 1222, 2206, 1286, 2219)	£64
JAX2239		Individual mounting option per request (Special) Requires additional information/specification		POA
DUAL CONTROL				
JAA2219		Compact Dual Control LiNX (DLX-ACU200)	(NC 1222, 2227, 2228)	£395
INVACARE LiNX IDC DUAL CONTROL				
JAA2245		Invacare LiNX IDC Dual control	(NC 2260, 2069)	£750
DUAL CONTROL BRACKETRY				
JAA2259		Bracket RH, Attendant	(NC 2245)	○
JAA2260		Bracket LH, Attendant	(NC 2245)	○
INVACARE LiNX IDC DUAL CONTROL BRACKETRY				
JAA1221		Bracket RH, detachable	(NC 2219, 2227)	○
JAA1222		Bracket LH, detachable	(NC 2219, 2228)	○
LiNX COMPACT REMOTE				
JAA2231		Proportional Compact Remote DLX-CR401	(NC 2204)	£719
JAA2232		Proportional Compact Remote Low Force DLX-CR-400LF	(NC 2204)	£960
LiNX COMPACT REMOTE BRACKET				
JAA2207		RH Mount	(NC 1707, 2209, 2213, 2217, 1404, 2311)	○
JAA2208		LH Mount	(NC 1708, 2210, 2214, 2218, 1405, 2310)	○
JAA2265		Centre Mount, swing-away to RH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!	(NC 1707, 2209, 2213, 2217, 1404, 1405, 2310, 2311, 1415)	£498
JAA2266		Centre Mount, swing-away to LH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!	(NC 1708, 2210, 2214, 2218, 1404, 1405, 2310, 2311, 1339)	£498
JAX2215		Individual mounting option per request (Special) Requires additional information/specification		POA
CHIN CONTROL				
JAA2229		Proportional Chin Joystick (ASL138) Selected headrest always with STD mounting hardware	(NC 2203, 2204, 2230, 2247, 2233, 2257, 2240)	£1,734
UNI9250		Mo-vis Micro joystick		POA
UNI9251		Mo-vis Multi joystick		POA

CHIN CONTROL BRACKETRY			
JAA1289		Swing Away Bracket RH, side manual (ASL628) Requires Elan Headrest	○
JAA1288		Swing Away Bracket LH, side manual (ASL628) Requires Elan Headrest	○
JAX1287		Swing Away Bracket RH, powered	(NC 2259, 1221, 1415) £2,183
JAX1286		Swing Away Bracket LH, powered	(NC 2260, 1222, 1339) £2,183
LINX HEAD CONTROL			
JAA2230		Head Control with ultrasonic proximity sensors (ATOM, ASL104) with standard wings M4	(NC 2204, 2229, 2233, 2092, 2323, 2324, 2322, 2325, 2326, 2327, 2257, 2240, 2570) £3,719
JAA2247		Head Control with ultrasonic proximity sensors (PROTON, ASL104P) with PROTON wings	(NC 2204, 2229, 2233, 2092, 2323, 2324, 2322, 2325, 2326, 2327, 2257, 2240, 2570) £4,356
LINX FOUR -SWITCH PROXIMITY ARRAY			
JAA2257		Four-Switch Proximity Array (ASL106)	(NC 2204, 2229, 2230, 2247, 2233, 2240) £5,760
JAA2275		Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1707, 2209, 2213, 2217, 1404, 1405, 2310, 2311, 1415) £498
JAA2276		Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1708, 2210, 2214, 2218, 1404, 1405, 2310, 2311, 1339) £498
JAX2269		Individual mounting option per request (Special) Requires additional information/specification	POA
LINX WIRELESS ACCESSORIES FOR ATOM			
JAA2249		Wireless Triple Switch Receiver (ASL557-3)	£341
JAA2250		Wireless Mouse Emulator (ASL558)	£817
LINX MICRO EXTREMITY CONTROL			
JAA2233		MEC Proportional Joystick (Micro Extremity Control ASL130)	(NC 2204, 2229, 2230, 2247, 2257) £3,612
LINX HEAVY DUTY REMOTE			
JAA2240		Proportional Heavy Duty Remote Joystick w mode switch (ASL133)	(NC 2204, 2229, 2230, 2247, 2257) £2,013
LINX MEC /HEAVY DUTY REMOTE BRACKET			
JAA2285		Tray Mount RH (100x150 mm)	(NC 1707, 2209, 2213, 2217, 1415, 1404, 1405, 2310, 2311) ○
JAA2286		Tray Mount LH (100x150 mm)	(NC 1708, 2210, 2214, 2218, 1339, 1404, 1405, 2310, 2311) ○
JAA2283		Centre Mount, swing-away to RH, Nucleus	(NC 1707, 2209, 2213, 2217, 1415, 1404, 1405, 2310, 2311) £498
JAA2284		Centre Mount, swing-away to LH, Nucleus	(NC 1708, 2210, 2214, 2218, 1339, 1404, 1405, 2310, 2311) £498
JAA2287		Centre Mount, swing-away to RH, Nucleus Incl. Tray ASL603, 100x150 mm	(NC 1707, 2209, 2213, 2217, 1415, 1404, 1405, 2310, 2311) £498
JAA2288		Centre Mount, swing-away to LH, Nucleus Incl. Tray ASL603, 100x150 mm	(NC 1708, 2210, 2214, 2218, 1339, 1404, 1405, 2310, 2311) £498
JAX2289		Individual mounting option per request (Special) Requires additional information/specification	POA
LINX PEDIATRIC COMPACT			
JAA2243		Proportional Pediatric Compact (ASL132), incl. Input Module	(NC 2203, 2204, 2255, 2256) £2,586
LINX PEDIATRIC COMPACT BRACKETRY			
JAA2235		Tray Mount RH (100x150 mm)	(NC 1707, 2209, 2213, 2217, 1415, 1404, 1405, 2310, 2311) ○
JAA2236		Tray Mount LH (100x150 mm)	(NC 1708, 2210, 2214, 2218, 1339, 1404, 1405, 2310, 2311) ○
JAA2267		Centre Mount, swing-away to RH, Nucleus	(NC 1707, 2209, 2213, 2217, 1415, 1404, 1405, 2310, 2311) £498
JAA2268		Centre Mount, swing-away to LH, Nucleus	(NC 1708, 2210, 2214, 2218, 1339, 1404, 1405, 2310, 2311) £498
JAA2241		Centre Mount, swing-away to RH, Nucleus Incl. Tray ASL603, 100x150 mm	(NC 1707, 2209, 2213, 2217, 1415, 1404, 1405, 2310, 2311) £498
JAA2242		Centre Mount, swing-away to LH, Nucleus Incl. Tray ASL603, 100x150 mm	(NC 1708, 2210, 2214, 2218, 1339, 1404, 1405, 2310, 2311) £498
JAX2277		Individual mounting option per request (Special) Requires additional information/specification	POA

LINX SIP N PUFF HEAD CONTROL			
JAA2255		Sip'n Puff Head Control (ASL109), includes Input Module and Lip Switch (ASL 301 LIP) Combines simple sip and puff controls with head movements	(NC 2203, 2204, 2092, 2323, 2324, 2322, 2325, 2326, 2327, 2243) £5,562
LINX SIP N PUFF CONTROL & INPUT MODULE			
JAA2256		Sip'n Puff Control includes Sip'n Puff kit, Lip Switch ASL301LIP and Input module	(NC 2203, 2204, 2243) £930
LINX INTERFACE OPTIONS			
JAA2271		Input Module DLX-IN500-A (sub D, Sip'n Puff)	(NC 2255, 2256, 2203, 2204) £786
JAA2273		IN200 Module (Jack port for external power button + control inputs via DB9)	£241
JAA2272		LiNX Output Module (DLX-OUT500)	(NC 2204) £947
JAA1681		5:1 Switch Adapter (9-pin sub D to 5 jack ports)	£121
SPECIAL ALTERNATIVE JOYSTICKS & INPUTS			
JAX2216		Special alternative Joystick/Input and bracket as per request	POA
LINX 10-WAY SWITCH			
JAA2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access	£505
LINX 10-WAY SWITCH BRACKET			
JAA2217		Bracket RH, mounted for user	(NC 2207, 2287, 2283, 2275, 2213, 2209, 2265, 1404, 2311) ○
JAA2218		Bracket LH, mounted for user	(NC 2208, 2288, 2284, 2286, 2276, 2214, 2210, 2266, 1405, 2310) ○
JAA2226		For attendant, mounted on back shroud	£64
TOGGLES / BUTTONS			
JAA2261		Stereo Toggle - direct control for tilt-only	(NC 0584, 0589, 0592, 2053, 2060, 2063, 2066) £291
JAA2262		Stereo Button - direct control for tilt-only	(NC 0584, 0589, 0592, 2053, 2060, 2063, 2066) £278
JAA2263		Four-Way Toggle - direct control for tilt & recline	(NC 2069, 0577, 0584, 0589, 0592) £362
JAA2264		Four-Way Button - direct control for tilt & recline	(NC 2069, 0577, 0584, 0589, 0592) £387
TOGGLES / BUTTONS MOUNT			
JAA2279		RH Mount (of Toggle/Button)	○
JAA2280		LH Mount (of Toggle/Button)	○
JAA2278		Without mounting	○
LINX USB CHARGER PORT			
JAA2253		USB Charger Port DLX-USB02-A	(NC 2702, 2703, 2704, 2705, 2706) £165
LINX USB CHARGER PORT BRACKET			
JAA2281		RH Mount (of Charger Port)	○
JAA2282		LH Mount (of Charger Port)	○
REMOTE STOP SWITCH			
JAA2258		Remote Stop Switch (ASL504)	£1,050
FURTHER LINX COMPONENTS / ACCESSORIES			
JAA2300		To be enclosed only	○
JAX2301		To be connected and programmed (Special)	POA
BUDDY BUTTONS			
JAA1685		Red Buddy Button	£88
JAA1686		Yellow Buddy Button	£88
JAA1687		Green Buddy Button	£88
JAA1690		Blue Buddy Button	£88
EGG SWITCH			
JAA2254		Egg Switch, black	£88

PIKO BUTTONS			
JAA1711		Red Piko Button	£158
JAA1712		Yellow Piko Button	£158
JAA1713		Green Piko Button	£158
JAA1714		Orange Piko Button	£158
JAA1715		Blue Piko Button	£158
JAA1716		Black Piko Button	£158
BODYPOINT JOYSTICK KNOBS			
JAA2721		3 inch U-Shaped, Joystick Handle, Bodypoint PC101A	£65
JAA2722		4 inch U-Shaped, Joystick Handle, Bodypoint PC102A	£65
JAA2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£65
JAA2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£65
SHROUD COLOUR*			
JAA2871		Mid range - Black/Anthracite Shroud	○
JAA2872		Mid range - Silver Shroud	○
JAA2873		Mid range - Bengal Red Shroud	○
JAA2874		Mid range - Mercury Blue Shroud	○
JAA2875		Mid range - Midnight Purple Shroud	○
JAX0299		Premium - Individual Shroud Colour RAL number required Extended lead lime: + 15 days	(NC 2883, 2884, 2885, 2886, 2888, 2889, 2890, 2891) £214
JAA2865		Premium - Racing Green Shroud	(NC 2883, 2884, 2885, 2886, 2887, 2889, 2890, 2891) £104
JAA2866		Premium - Deep Ocean Blue Satin Shroud	(NC 2883, 2884, 2885, 2886, 2887, 2888, 2890, 2891) £104
JAA2867		Premium - Frosty Blue Shroud	(NC 2883, 2884, 2885, 2886, 2887, 2888, 2889, 2891) £104
JAA2868		Premium - Gold Shroud RAL1036	(NC 2883, 2884, 2885, 2886, 2887, 2888, 2889, 2890) £104
RIM COLOUR*			
JAA2882		Mid range - Black	○
JAA2883		Mid range - Silver	○
JAA2884		Mid range - Bengal Red	○
JAA2885		Mid range - Mercury Blue	○
JAA2886		Mid range - Midnight Purple	○
JAX2887		Premium - Individual Rim Insert Colour RAL number required Extended lead lime: + 15 days Price Included with Premium Shroud Kit	(NC 2872, 2873, 2874, 2875, 2865, 2866, 2867, 2868) ○
JAA2888		Premium - Racing Green Rim Inserts Price Included with Premium Shroud Kit	(NC 2872, 2873, 2874, 2875, 0299, 2866, 2867, 2868) ○
JAA2889		Premium - Deep Ocean Blue Satin Rim Inserts Price Included with Premium Shroud Kit	(NC 2872, 2873, 2874, 2875, 0299, 2865, 2867, 2868) ○
JAA2890		Premium - Frosty Blue Rim Inserts Price Included with Premium Shroud Kit	(NC 2872, 2873, 2874, 2875, 0299, 2865, 2866, 2868) ○
JAA2891		Premium - Gold Rim Inserts RAL1035 Price Included with Premium Shroud Kit	(NC 2872, 2873, 2874, 2875, 0299, 2865, 2866, 2867) ○
FORKS*			
JAA2631		Single sided fork	£352
JAA2632		Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆

FRONT WHEEL*			
JAA0835		9" (225 mm) x 2.80 (70 mm) Castor, pneumatic	◆
JAA0836		9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	○
JAA0838		9" (225 mm) x 2.80 (70 mm) Castor, puncture protected Gel-sealed	○
JAA0828		8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre (NC 1909)	○
CASTOR LOCKS			
JAA1490		Castor Lock	£197
DRIVE WHEEL*			
JAA0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	◆
JAA0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	○
JAA0706		DW 14" (355 mm) x 3" (75 mm), Puncture protected Gel-sealed	○
C.T.C SUSPENSION*			
JAA2870		CTC Suspension Control - Traction - Comfort	◆
SUSPENSION SETTING*			
JAA1910		Soft (suggested for up to 70 kg User Weight)	○
JAA1912		Medium (suggested for up to 90 kg User Weight)	○
JAA1911		Firm (suggested for up to 110 kg User Weight)	○
JAA1915		Very firm (suggested for up to 136 kg User Weight)	○
MOTOR TYPE*			
JAA1294		Indoor setting, reduced speed approx. 3 km/h	○
JAA1282		6 km/h	◆
JAA1908		10 km/h Package UK only Always combined with LED 1437 (NC 1437, 1254, 2837)	£543
JAA1909		12 km/h Ultra High Speed Package UK only Always combined with LED 1437 (NC 1437, 1254, 2837, 0828)	£653
BATTERY TYPE*			
JAA1259		60 Ah Gel	◆
JAA1260		73,5 Ah Gel (NC 1450)	£92
G TRAC			
JAA1215		G-Trac Module	£324
CHARGER*			
JAA1343		OptiCharge Charger 8A (IP44)	◆
JAA1313		OptiCharge Charger 10A (IP44)	£151
ANTI-TIPPERS*			
JAA1205		Anti-tipper	◆
LIGHTS			
JAA1437		LED Light system with Indicators 30 Lux (NC 1908, 1909)	£479
JAA2837		LED Light system extra bright with indicators 60 Lux (NC 1908, 1909)	£143
JAA1254		Without Light system (NC 1909)	○
KERB RIDER			
JAA1244		Kerb Rider With park position feature (NC 0586, 0589, 0592, 1294)	£212

TRANSPORTATION KIT*			
JAA1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory) Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering	£586
JAA1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) (Adequate headrest is compulsory)	◆
JAA1401		6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19)) required for variants with combined Lift & Tilt module (Adequate headrest is compulsory)	(NC 2052) ○
EASY ACCESS FRONT TIE DOWNS			
JAA2894		Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down For transportworthy versions (1494/1401) only	£49
LAP BELT*			
JAA2113		Lap belt with Auto Lock Flexible mount in side rails	◆
POSITIONING BELT			
JAA2701		BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm	(NC 2185, 2186, 2187, 2188, 2189, 2253, 2114, 2115, 2025, 2026, 2027) £77
JAA2702		BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm	(NC 2185, 2186, 2187, 2188, 2189, 2253, 2114, 2115, 2028, 2029, 2030, 2031) £77
JAA2703		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm	(NC 2185, 2186, 2187, 2188, 2189, 2253, 2114, 2115, 2025, 2026, 2027) £79
JAA2704		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm	(NC 2185, 2186, 2187, 2188, 2189, 2253, 2114, 2115, 2028, 2029, 2030, 2031) £79
JAA2705		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm	(NC 2185, 2186, 2187, 2188, 2189, 2253, 2114, 2115, 2025, 2026, 2027) £79
JAA2706		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm	(NC 2185, 2186, 2187, 2188, 2189, 2253, 2114, 2115, 2028, 2029, 2030, 2031) £79
REAR VIEW MIRROR			
JAA1415		Rearview mirror RH	(NC 2265, 2267, 2275) £44
JAA1339		Rearview mirror LH	(NC 2266, 2268, 2276) £44
OTHER OPTIONS			
JAA1438		Luggage rack	£329
JAA1416		Fender for Castors	£88
JAA1417		Fender for Drive wheels	£169
JAA1404		Tray, Swing Away RH	(NC 1707, 2207, 2209, 2210, 2214, 2218, 2235, 2236, 2267, 2268, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2241, 2242) £167
JAA1405		Tray, Swing Away LH	(NC 1708, 2208, 2209, 2210, 2214, 2218, 2235, 2236, 2267, 2268, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2241, 2242) £167
JAA2310		Tray Flip and fold away 1/2 lap tray LH (14" W x 10" D)	(NC 1707, 2207, 2209, 2210, 2214, 2218, 2235, 2236, 2267, 2268, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2241, 2242) £738
JAA2311		Tray Flip and fold away 1/2 lap tray RH (14" W x 10" D)	(NC 1708, 2208, 2209, 2210, 2214, 2218, 2235, 2236, 2267, 2268, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2241, 2242) £738
JAX2312		Vent Box, Low Profile Recline, Universal	(NC 2054, 2055, 2056) POA
JAX2314		Vent Box, Low Profile Recline, Trilogy	(NC 2054, 2055, 2056) POA
JAX2315		Vent Box, Low Profile Recline, LTV – Easy Fit	(NC 2054, 2055, 2056) POA
JAX2313		Oxygen Holder, Low Profile Recline	(NC 2054, 2055, 2056) POA
TOTAL PRICE			

Non Conforming Options:

JAL2810	CHASSIS Narrow (560 mm/22")
JAA2099	Fixed STF height and seat angle setting
JAA2052	Ultra Low Maxx CG tilt (0 to 50°)
JAA2053	Ultra Low Maxx CG tilt (0 to 45°) combined with 300 mm seat lifter increased STF Height by 4 mm
JAA2857	Seat-to-Floor Height 435 mm
JAA2858	Seat-to-Floor Height 460 mm
JAA2859	Seat-to-Floor Height 485 mm
JAX2054	Ultra Low Maxx CG tilt 5° Fixed Anterior Assist mount Max tilt angle reduced by 5° Front STF height as defined above With centre-mounted legrests only Front STF height as defined above
JAX2055	Ultra Low Maxx CG tilt 10° Fixed Anterior Assist mount Max tilt angle reduced by 10° Front STF height as defined above With centre-mounted legrests only Front STF height as defined above
JAX2056	Ultra Low Maxx CG tilt 5° Fixed Posterior-tilt Assist mount Max tilt angle increased by 5°. Front STF height increased by 6 cm over chosen SH setting Front STF height increased by 6 cm over chosen seat height setting
JAA2069	Fixed Angle Setting (NO POWERED RECLINE)
JAA2060	Powered Recline 90 to 170° , Standard Mount
JAA2063	Powered Recline, 8° Precline mount, 82° to 162°
JAA2066	Powered Recline, 30° Precline mount, 60° to 140°
JAA2025	405 mm / 16" With centre-mount legrests only
JAA2026	431 mm / 17"
JAA2027	456 mm / 18"
JAA2028	480 mm / 19"
JAA2029	505 mm / 20"
JAA2030	530 mm / 21"
JAA2031	555 mm / 22"
JAA2035	405 mm / 16"
JAA2036	430 mm / 17"
JAA2037	455 mm / 18"
JAA2041	555 mm / 22"
JAA2042	580 mm / 23"
JAA2145	Matrx Flo-Tech Lite Black smooth cover, 75 mm thick Up to 127 kg
JAA2146	Matrx Flo-tech Contour Black smooth cover, 100 mm thick Up to 191 kg
JAA2510	Straight back cane height 455 mm / 18"
JAA2511	Straight back cane height 505 mm / 20"
JAA2512	Straight back cane height 555 mm / 22"
JAA2513	Straight back cane height 610 mm / 24"
JAA2514	Angled (5°) back cane height 505 mm / 20"
JAA2515	Angled (5°) back cane height 555 mm / 22"
JAA2516	Angled (5°) back cane height 610 mm / 24"
JAA2517	75°
JAA2518	79°
JAA2519	85°
JAA2520	90°
JAA2521	95°

JAA2522	101°
JAA2523	106°
JAA2524	112°
JAA2540	Standard Rehab Back Pan and Back Upholstery height 355 mm / 14"
JAA2541	Standard Rehab Back Pan and Back Upholstery height 405 mm / 16"
JAA2542	Standard Rehab Back Pan and Back Upholstery height 455 mm / 18"
JAA2543	Standard Rehab Back Pan and Back Upholstery height 505 mm / 20"
JAA2050	High Back, height 405 mm / 16" [Mesh fabric]
JAA2051	High Back, height 455 mm / 18" [Mesh fabric]
JAA2320	High Back, height 505 mm / 20" [Mesh fabric]
JAA2550	Matrx Elite E2 Back, height 355mm / 14" / 16 - 18" wide
JAA2551	Matrx Elite E2 Back, height 405 mm / 16"
JAA2552	Matrx Elite E2 Back, height 455 mm / 18"
JAA2553	Matrx Elite E2 Back, height 505 mm / 20"
JAA2554	Matrx Elite HD Back, height 405 mm /16" Extended lead time: + 4 days Additionnal Delivery Time + 4 days
JAA2555	Matrx Elite HD Back, height 505 mm / 20" Extended lead time: + 4 days Additionnal Delivery Time + 4 days
JAA2556	Matrx Elite TR Back, height 405 mm / 16" [Mesh fabric]
JAA2557	Matrx Elite TR Back, height 455 mm / 18" [Mesh fabric]
JAA2177	Matrx Elite TR Back, height 505 mm / 20" [Mesh fabric]
JAA2563	Matrx Elite TR HD Back, height 505 mm / 20" [Mesh fabric] Extended lead time: + 4 days Additionnal Delivery Time + 4 days
JAA2558	Matrx Elite Deep E2 Back, height 405 mm / 16" [Mesh fabric]
JAA2559	Matrx Elite Deep E2 Back, height 455 mm / 18" [Mesh fabric]
JAA2178	Matrx Elite Deep E2 Back, height 505 mm / 20" [Mesh fabric]
JAA2174	Matrx PB Back, height 405 mm / 16" [Mesh fabric]
JAA2175	Matrx PB Back, height 505 mm / 20" [Mesh fabric]
JAA2560	Matrx PB Deep Back, height 405mm / 16"
JAA2561	Matrx PB Deep Back, height 505 mm / 20"
JAA2562	Matrx Elite Deep HD Back, height 505 mm / 20" Extended lead time: + 4 days Additionnal Delivery Time + 4 days
JAA2570	Without Backrest Support, incl. Matrx interface
JAA2307	Ultra Rail Mounted Flip back Cantilever Maxx Tilt Arm 230 - 320 mm (9-12.5") Not compatible with reclining backrest
JAA2096	Cantilever Flip-Back Arm, width, height and angle adjustable, 340 - 405 mm (13.5-16")
JAA2097	Two-post Recline Arm, Flip Back, adjustable height 230 - 330 mm (9-13")
JAA2098	Two-post Recline Arm, Flip Back, adjustable height 320 - 405 mm (12.5-16")
JAA2308	1" Offset Armrest Spacer (pair), for two-post recline only
JAA2309	1/2" Offset Armrest Spacer (pair), for two-post recline only
JAA2013	Full length arm pad with startex cover visco foam [355mm], 40 mm thick
JAA2015	without LH arm pad
JAA2017	Full length arm pad with startex cover visco foam [355mm]; 40 mm thick
JAA2019	without RH arm pad
JAA2347	Elbow Block XS (110 mm wide x 83 mm high)
JAA2348	Elbow Block S (135 mm wide x 90 mm high)
JAA2349	Elbow Block M (145 mm wide x 95 mm high)

JAA2350	Elbow Block L (165 mm wide x 110 mm high)
JAA2351	Elbow Block XS (110 mm wide x 83 mm high)
JAA2352	Elbow Block S (135 mm wide x 90 mm high)
JAA2353	Elbow Block M (145 mm wide x 95 mm high)
JAA2354	Elbow Block L (165 mm wide x 110 mm high)
JAA2114	Hip Support LH, Quick Release, pad size 180 x 150 mm
JAA2186	Hip Support LH, Quick Release, pad size 230 x 100 mm
JAA2188	Hip Support LH, Quick Release, pad size 330 x 100 mm
JAA2185	Hip Support RH, Quick Release, pad size 140 x 100 mm
JAA2115	Hip Support RH, Quick Release, pad size 180 x 150 mm
JAA2187	Hip Support RH, Quick Release, pad size 230 x 100 mm
JAA2189	Hip Support RH, Quick Release, pad size 330 x 100 mm
JAA2092	Onyx Auto Style headrest with standard mounting hardware
JAA2323	Elan Standard Pad 150 mm / 6" wide w/standard cover (ESP26), complete with Elan Headrest Hardware (MEHW)
JAA2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26), complete with Elan Headrest Hardware (MEHW)
JAA2322	Elan Standard Pad 355 mm / 14" wide w/standard cover (ESP26), complete with Elan Headrest Hardware (MEHW)
JAA2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23), complete with Elan Headrest Hardware (MEHW)
JAA2326	Elan Occipital Pad 305x200 mm / 12"x8" w/standard cover complete with Elan Headrest Hardware (MEHW)
JAA2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover complete with Elan Headrest Hardware (MEHW)
JAA0586	Centre mount legrest, angle settable Length: 310 - 410 mm
JAA0589	LNX Power centre mount legrest, length compensating; Length: 340 - 410 mm
JAX0592	LNX Power centre mount legrest with telescoping foot platform (lowers to the ground; Length: 380 - 419 mm) Selectable length depends on STF height. Custom-made, requires customized setting! Custom-made, requires customized setting!
JAX3003	Without Centre Mount Legrests
JAA2365	Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")
JAA0516	Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5)
JAA2366	Angle adjustable foot platform S (305 W x 215 mm D / 12"x8.5")
JAA0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
JAA0572	Pivot Plus Manual Elevating Legrest LH
JAA0577	Pivot Plus Power Elevating Legrest LH
JAA0579	Pivot Plus Manual Elevating Legrest RH
JAA0584	Pivot Plus Power Elevating Legrest RH
JAA2203	Remote LiNX REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators
JAA2204	Remote LiNX REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights
JAA2205	Rehab Remote LiNX REM400, touch-screen user interface Configurable jack sockets
JAA2206	Rehab Remote LiNX REM500 without joystick for alternative access users Requires a user input/alternative joystick: 2231, 2232, 2229, 2230, 2247, 2234
JAA2270	Assignable Toggle Switch Module f. REM400 (DLX-SW400TGL)
JAA1207	Standard Bracket RH
JAA1208	Standard Bracket LH
JAA1707	Swing Away Bracket RH
JAA1708	Swing Away Bracket LH

JAA1709	Matrx Resolve Swing in/out bracket RH
JAA1710	Matrx Resolve Swing in/out bracket LH
JAA2209	Rem500 Swing away RH, multi-adjustable
JAA2210	Rem500 Swing away LH, multi-adjustable
JAA2213	Centre Mount, swing-away to RH, Nucleus
JAA2214	Centre Mount, swing-away to LH, Nucleus
JAA2227	Bracket RH, mounted for Attendant
JAA2228	Bracket LH, mounted for Attendant
JAA2219	Compact Dual Control LiNX (DLX-ACU200)
JAA2245	Invacare LiNX IDC Dual control
JAA2259	Bracket RH, Attendant
JAA2260	Bracket LH, Attendant
JAA1221	Bracket RH, detachable
JAA1222	Bracket LH, detachable
JAA2231	Proportional Compact Remote DLX-CR401
JAA2232	Proportional Compact Remote Low Force DLX-CR-400LF
JAA2207	RH Mount
JAA2208	LH Mount
JAA2265	Centre Mount, swing-away to RH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!
JAA2266	Centre Mount, swing-away to LH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!
JAA2229	Proportional Chin Joystick (ASL138) Selected headrest always with STD mounting hardware
JAX1287	Swing Away Bracket RH, powered
JAX1286	Swing Away Bracket LH, powered
JAA2230	Head Control with ultrasonic proximity sensors (ATOM, ASL104) with standard wings M4
JAA2247	Head Control with ultrasonic proximity sensors (PROTON, ASL104P) with PROTON wings
JAA2257	Four-Switch Proximity Array (ASL106)
JAA2275	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JAA2276	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JAA2233	MEC Proportional Joystick (Micro Extremity Control ASL130)
JAA2240	Proportional Heavy Duty Remote Joystick w mode switch (ASL133)
JAA2285	Tray Mount RH (100x150 mm)
JAA2286	Tray Mount LH (100x150 mm)
JAA2283	Centre Mount, swing-away to RH, Nucleus
JAA2284	Centre Mount, swing-away to LH, Nucleus
JAA2287	Centre Mount, swing-away to RH, Nucleus Incl. Tray ASL603, 100x150 mm
JAA2288	Centre Mount, swing-away to LH, Nucleus Incl. Tray ASL603, 100x150 mm
JAA2243	Proportional Pediatric Compact (ASL132), incl. Input Module
JAA2235	Tray Mount RH (100x150 mm)
JAA2236	Tray Mount LH (100x150 mm)
JAA2267	Centre Mount, swing-away to RH, Nucleus
JAA2268	Centre Mount, swing-away to LH, Nucleus
JAA2241	Centre Mount, swing-away to RH, Nucleus Incl. Tray ASL603, 100x150 mm

JAA2242	Centre Mount, swing-away to LH, Nucleus Incl. Tray ASL603, 100x150 mm
JAA2255	Sip'n Puff Head Control (ASL109), includes Input Module and Lip Switch (ASL 301 LIP Combines simple sip and puff controls with head movements)
JAA2256	Sip'n Puff Control includes Sip'n Puff kit, Lip Switch ASL301LIP and Input module
JAA2271	Input Module DLX-IN500-A (sub D, Sip'n Puff)
JAA2272	LiNX Output Module (DLX-OUT500)
JAA2217	Bracket RH, mounted for user
JAA2218	Bracket LH, mounted for user
JAA2261	Stereo Toggle - direct control for tilt-only
JAA2262	Stereo Button - direct control for tilt-only
JAA2263	Four-Way Toggle - direct control for tilt & recline
JAA2264	Four-Way Button - direct control for tilt & recline
JAA2253	USB Charger Port DLX-USB02-A
JAA2872	Mid range - Silver Shroud
JAA2873	Mid range - Bengal Red Shroud
JAA2874	Mid range - Mercury Blue Shroud
JAA2875	Mid range - Midnight Purple Shroud
JAX0299	Premium - Individual Shroud Colour RAL number required Extended lead lime: + 15 days
JAA2865	Premium - Racing Green Shroud
JAA2866	Premium - Deep Ocean Blue Satin Shroud
JAA2867	Premium - Frosty Blue Shroud
JAA2868	Premium - Gold Shroud RAL1036
JAA2883	Mid range - Silver
JAA2884	Mid range - Bengal Red
JAA2885	Mid range - Mercury Blue
JAA2886	Mid range - Midnight Purple
JAX2887	Premium - Individual Rim Insert Colour RAL number required Extended lead lime: + 15 days Price Included with Premium Shroud Kit
JAA2888	Premium - Racing Green Rim Inserts Price Included with Premium Shroud Kit
JAA2889	Premium - Deep Ocean Blue Satin Rim Inserts Price Included with Premium Shroud Kit
JAA2890	Premium - Frosty Blue Rim Inserts Price Included with Premium Shroud Kit
JAA2891	Premium - Gold Rim Inserts RAL1035 Price Included with Premium Shroud Kit
JAA0828	8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre
JAA1294	Indoor setting, reduced speed approx. 3 km/h
JAA1908	10 km/h Package UK only Always combined with LED 1437
JAA1909	12 km/h Ultra High Speed Package UK only Always combined with LED 1437
JAA1437	LED Light system with Indicators 30 Lux
JAA2837	LED Light system extra bright with indicators 60 Lux
JAA1254	Without Light system
JAA1244	Kerb Rider With park position feature
JAA1450	Preparation for DAHL-Docking station (Adequate headrest is compulsory) Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering
JAA1494	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) (Adequate headrest is compulsory)

JAA1401	6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19)) required for variants with combined Lift & Tilt module (Adequate headrest is compulsory)
JAA2701	BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm
JAA2702	BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm
JAA2703	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm
JAA2704	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm
JAA2705	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm
JAA2706	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm
JAA1415	Rearview mirror RH
JAA1339	Rearview mirror LH
JAA1404	Tray, Swing Away RH
JAA1405	Tray, Swing Away LH
JAA2310	Tray Flip and fold away 1/2 lap tray LH (14" W x 10" D)
JAA2311	Tray Flip and fold away 1/2 lap tray RH (14" W x 10" D)
JAX2312	Vent Box, Low Profile Recline, Universal
JAX2314	Vent Box, Low Profile Recline, Trilogy
JAX2315	Vent Box, Low Profile Recline, LTV – Easy Fit
JAX2313	Oxygen Holder, Low Profile Recline



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):