

## **Certificate in relation to air transportation by passengers with walking disabilities**

We hereby certify that the our product

### **e-motion M25**

including the lithium-ion batteries fulfil the requirements of

### **IATA chapter 2.3.2.4**

and therefore may be carried as checked baggage without restriction.

The built-in lithium-ion batteries of the e-motion M25, which are not intended for removal by the user, fulfil the requirements of the tests listed in the UN Manual of Tests and Criteria" Part III subsection 38.3.

The e-motion M25 can be set into a dedicated "flight mode". In this mode the system is completely deactivated. While in this mode the e-motion M25 can only be switched on by a coded procedure. Unintended switching on by luggage is not possible. Furthermore, the Bluetooth device of the system is completely switched off.

**We strongly recommend contacting the airline sufficiently in advance to request approval for transport.**

This certificate does not apply to damaged or defective batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31<sup>st</sup> December 2022**.

Albstadt, 11<sup>th</sup> of November 2021

Alber GmbH



Hansjörg Reiner  
Managing Director