

Support Belts

SystemRoMedic™



Support for every step

Sitting transfers, sit-to-stand, standing and walking



FlexiBelt

Versatile support belt for sitting transfers and active assistance – designed for comfort and a secure grip in every situation.

FlexiBelt supports sit-to-stand transfers and sitting transfers from wheelchair to bed, toilet or car. It can also be used when standing and walking. The belt is worn around the waist or hips and features multiple sturdy handles – both vertical and horizontal – giving caregivers a firm, secure grip from any position.

The classic FlexiBelt features a soft polyester velvet inner surface for all-day comfort. Durable nylon outer material ensures long product life. Anti-slip treated inner surface keeps the belt firmly in place. Reinforced seams in exposed areas and a sturdy finished edge increase longevity.

> Ergonomic handles

Vertical and horizontal handles provide a firm grip from all positions.

> Quick-connect buckle

Strong, safe and easy to open, close and tighten.

> Tightening strap

Sturdy and easy to adjust for a secure fit.



> Soft inner surface

Polyester velvet lining for user comfort during extended wear.

> Reinforced seams

Extra stitching in exposed areas for maximum durability.

> Anti-slip treated

Inner surface prevents the belt from moving during transfer.

EasyBelt

Robust support belt for sitting transfers and active assistance – with moisture-repellent material for use in hygiene-sensitive environments.

EasyBelt supports sit-to-stand transfers and sitting transfers from wheelchair to bed, toilet or car. It can also be used when standing and walking. The belt is worn around the waist or hips and features multiple sturdy handles – both vertical and horizontal – giving caregivers a firm, secure grip in any situation.

EasyBelt's inner surface is made of moisture-repellent, PU-coated nylon – ideal for use where moisture, dampness or special hygiene requirements are present. Available in sizes XS to L to suit a wide range of users.



> Ergonomic handles

Vertical and horizontal handles provide a firm grip from all positions.

> Moisture-repellent inner surface

PU-coated nylon – suitable for hygiene-sensitive environments.

> Available XS–L

Four sizes to suit waist measurements from 45 to 160 cm.

> Tightening strap

Sturdy and easy to adjust for a secure fit.

> Quick-connect buckle

Strong, safe and easy to open, close and tighten.

FlexiMove Mini

Support for sit-to-stand and positioning in bed but also for sitting transfers.

FlexiMove has multiple areas of application. For example as support when helping someone get up from a bed or chair, when moving higher up in bed or when transferring to and from a wheelchair – in that case, often on a sliding board. FlexiMove can also be used for positioning in combination with EasySlide (50 x 60) and WendyLean.

FlexiMove is available in two variants with different inside material. One variant has an inside made of soft and comfortable polyester velvet. The other variant has an inside made of nylon which is moisture repellent and resists soiling. The second variant can easily be cleaned using disinfectant wipes.

> Two different surface materials available

Soiling resistant and moisture-repellent nylon or soft and comfortable polyester velvet.

> Gripping options

Both short and long handles enable more gripping options.



> Anti-slip function

Both versions provide anti-slip function to prevent the device from sliding against the body.

ReadyBelt

For support during mobility training and transfers.

ReadyBelt is a belt that is designed to provide the user with support and help during mobility training and different transfer situations. With ReadyBelt, the caregiver can easily and in an ergonomically correct manner assist the user during rising and during transfers, for example, between bed and wheelchair and between wheelchair and toilet or car.

Several ergonomically designed handles provide the caregiver with a secure grip, allowing the ideal working stance in every situation.

ReadyBelt is not to be used as an assistive device for lifting.

> Sizing & fit

Available in different sizes, suitable for most users.

> Application

Placed around the user's waist or hips. Secured with a clasp that is easy to open, close and tighten.

> Comfort & safety

Soft and comfortable to wear. Non-slip inner surface prevents the belt from sliding up.



Technical data ▼

	Size ▼	Measurement ▼	Waist measurement ▼	Handles ▼
FlexiBelt	S	105 cm	60–100 cm	7 pcs
	M	120 cm	70–120 cm	9 pcs
	L	135 cm	100–160 cm	11 pcs
	XL	150 cm	115–175 cm	13 pcs
EasyBelt	XS	75 cm	45–70 cm	5 pcs
	S	105 cm	60–100 cm	7 pcs
	M	120 cm	70–120 cm	9 pcs
	L	135 cm	100–160 cm	11 pcs
FlexiMove		20 x 60 cm		
ReadyBelt	S	105 x 10 cm		7 pcs
	M	120 x 10 cm		9 pcs
	L	140 x 10 cm		11 pcs

Invacare International GmbH
 Neuhofweg 51, 4147 Aesch BL, Switzerland
globalcomms@dhcare.com
www.dhcare.com

Invacare International GmbH is a member of the DHCare group of companies. Invacare International GmbH reserves the right to modify existing product models at any time. Product specifications may vary across regions. ©2026 Invacare International GmbH. All rights reserved. Support Belts_EU_14072026

