

Support Belts

SystemRoMedic™



Support for every step

Sitting transfers, sit-to-stand, standing and walking



FlexiBelt

Versatile support belt for sitting transfers and active assistance — designed for comfort and a secure grip in every situation.

FlexiBelt supports sit-to-stand transfers and sitting transfers from wheelchair to bed, toilet or car. It can also be used when standing and walking. The belt is worn around the waist or hips and features multiple sturdy handles — both vertical and horizontal — giving caregivers a firm, secure grip from any position.

The classic FlexiBelt features a soft polyester velvet inner surface for all-day comfort. Durable nylon outer material ensures long product life. Anti-slip treated inner surface keeps the belt firmly in place. Reinforced seams in exposed areas and a sturdy finished edge increase longevity.

> Ergonomic handles

Vertical and horizontal handles provide a firm grip from all positions.

> Quick-connect buckle

Strong, safe and easy to open, close and tighten.

> Tightening strap

Sturdy and easy to adjust for a secure fit.



> Soft inner surface

Polyester velvet lining for user comfort during extended wear.

> Reinforced seams

Extra stitching in exposed areas for maximum durability.

> Anti-slip treated

Inner surface prevents the belt from moving during transfer.

EasyBelt

Robust support belt for sitting transfers and active assistance — with moisture-repellent material for use in hygiene-sensitive environments.

EasyBelt supports sit-to-stand transfers and sitting transfers from wheelchair to bed, toilet or car. It can also be used when standing and walking. The belt is worn around the waist or hips and features multiple sturdy handles — both vertical and horizontal — giving caregivers a firm, secure grip in any situation.

EasyBelt's inner surface is made of moisture-repellent, PU-coated nylon — ideal for use where moisture, dampness or special hygiene requirements are present. Available in sizes XS to L to suit a wide range of users.



> Ergonomic handles

Vertical and horizontal handles provide a firm grip from all positions.

> Moisture-repellent inner surface

PU-coated nylon — suitable for hygiene-sensitive environments.

> Available XS-L

Four sizes to suit waist measurements from 45 to 160 cm.

> Tightening strap

Sturdy and easy to adjust for a secure fit.

> Quick-connect buckle

Strong, safe and easy to open, close and tighten.

Technical data ▼

	Size ▼	Measurement ▼	Waist measurement ▼	Handles ▼
EasyBelt	XS	75 cm	45–70 cm	5 pcs
	S	105 cm	60–100 cm	7 pcs
	M	120 cm	70–120 cm	9 pcs
	L	135 cm	100–160 cm	11 pcs
FlexiBelt	S	105 cm	60–100 cm	7 pcs
	M	120 cm	70–120 cm	9 pcs
	L	135 cm	100–160 cm	11 pcs
	XL	150 cm	115–175 cm	13 pcs