

 **TOP END**®



Invacare® Top End® Force™ -3 Handcycle

The **Invacare® Top End® Force™-3 Handcycle** is a great handcycle at a great price. This adjustable handcycle arrives fully assembled and is now a 30 speed handcycle. The Force-3 handcycle is an optimal recreation handcycle with function, comfort and performance—made right here in the United States.

Features

No Sizing Required

Stealth black powder coat paint

60" Wheelbase

Adjustable back angle, seat and footrests

30 speed drive train

170 mm cranks

Cable driven upper derailleur

Adjustable crank height: diagonal and vertical

Alex 26-559 wheels

1° camber

Clincher tires

Seat and back cushions

Chest and footrest safety straps

Drafting bumper

Threaded axles

Cable springs

Frictionless ovalized ergonomic hand pedals

Safety flag

Chain guard

Reflectors

Locking parking brake

250 lb. weight capacity

Model #: TE10039 (Non-disc brakes)



1. Wraparound frame with adjustable seat upholstery
2. Cable driven shifting for upper chainrings
3. Front disk brake and welded footrest clamps



Model #: TE10040 (Disc brakes)

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End

www.topendwheelchair.com
Sports and Recreation
Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(800) 532-8677
Fax: (727) 522-1007

©2018

Invacare Corporation.

All rights reserved.

Trademarks are identified

by the symbols ™ and ®.

All trademarks are owned

by or licensed to Invacare

Corporation unless

otherwise noted. Shimano,

Rapid Fire and Tektro

are registered trademarks

of Shimano Inc. Avid is a

registered trademark of

SRAM Corporation.

BionX is a registered

trademark of BionX

International Corporation.

Specifications are subject

to change without

notification.

Form No. 18-181

180619

Options

Alignment gauge

Water bottle

Cordless computer

Helmet

Safety lights

V-Crank upgrade

Bike rack

Indoor training roller

Knobby tires

 **TOP END**®