

# Invacare® | *Cadiz shower stool and chair*




- A high quality and cost effective solution for showering and dressing.
- Easy to handle seating height adjustment from 360 to 510 mm - in 13 steps
- Ergonomic seat surface
- Stable and secure footprint, due to angled rubber feet
- High load capacity: 143 kg (22.5 stone)

Invacare *Cadiz shower stool*

## Technical data

	<i>Invacare Cadiz stool</i>
Shower stool w/o backrest (pk 2)	1544791
<b>Dimensions</b>	
Seat width (A)	510 mm
Seat depth	310 mm
Seat height (B)	360 - 510 mm
Total width (C)	440 mm
<b>More info</b>	
Load capacity	143 kg (22.5 stone)
Weight	1,8 kg
Colour	White



Warranty 2 years






- Easy to handle height adjustment from 360 to 510 mm - in 13 steps
- Ergonomic shaped backrest
- Stable and secure footprint, due to angled rubber feet
- Easy to clean

Invacare **Cadiz shower chair**

**Technical data**

	<i>Invacare</i> <b>Cadiz chair</b>
Shower stool with backrest (pk 2)	1544829

<b>Dimensions</b>	
Seat width (A)	510 mm
Seat depth	310 mm
Seat height	360 - 510 mm
Backrest width	390 mm
Backrest height	230 mm
Total width (C)	440 mm
Total depth	415 - 460 mm
Total height (B)	790 - 940 mm



<b>More info</b>	
Load capacity	143 kg (22.5 stone)
Weight	2,6 kg
Colour	White
Warranty	2 Years