

Invacare Colibri®

Feel the Freedom



Stylish, lightweight and transportable

Designed with Invacare Litelock™ technology



Yes, you can.®

Feel the difference

Our scooters have been designed together with healthcare experts, to provide unrivalled standards of quality, safety and comfort.

Through continuous dedicated research and an unrivalled experience in designing for mobility, at Invacare we pride ourselves in scooters that improve quality of life. Being able to regain confidence, trust and independence from your scooter is important.

From folding transportable designs, cushioned seat and adjustable seat risers, every detail of our scooters is considered and refined to make getting about feel great.

Quality



Market-leading build quality with exceptional attention to detail.

Safety



Compliant with the latest European Medical Device Standards.

Comfort



Adaptive and adjustable design to provide lasting comfort.



Feel in control with Invacare Colibri®

Introducing the Invacare Colibri®, our lightest and most transportable scooter. Designed to complement your style with a choice of colours and improve your independence.

Compact and easy to drive, this scooter has been designed with a patented Litelock™ system making it easy to transport into the back of your car. The swivel seat makes getting on and off your scooter simple. The Colibri® is

incredibly simple to start and steer. It has been designed to give you a comfortable and safe drive with the ability to manoeuvre in tight spaces.

up to
8km/h
max speed

up to
16km
max distance
on a single
charge

136 kg
max user
weight



Easy transportation with the Invacare Litelock™ system



Easy to use and transport

- Travel easily in style
- Ergonomic steering for added control
- Simple to steer and drive



Safety in mind

- Reflectors
- Additional optional extras available



Built around comfort

- High quality cushioned seat
- Adjustable seat height
- Flip up armrests and swivel seat makes getting on and off easy

Our **lightest** and
most **transportable**
scooter

+ **-** **Easy
battery
charge**

in situ or remove
from scooter



Adjustable steering
for greater comfort



Swivel seat for easy transfer



Soft, padded seat for enhanced comfort



Adjustable steering for improved seating position



Add to your scooter

Tailor the Invacare Colibri® to your taste and needs. Turn over to explore the optional extras available.

Add to your scooter

We have a wide range of optional extras available, so you can customize your Invacare Colibri® scooter. Ask your retailer to discover the full range.



Additional battery pack



Large mirror



Crutch holder



Optional size wheelsets



'Easy fit' suspension kit



Interchangeable shroud panels

Technical data

For more information about this product, including the user manual, please visit our website www.invacare.eu.com



Backrest width ▼

465 mm



Seat depth ▼

400 mm



Seat to floor height ▼

370 - 420 mm



Backrest height ▼

385 mm



Driving unit width ▼

505 mm



Total length ▼

1010 mm



Empty weight ▼

44 kg (12 Ah)
49 kg (18 Ah)



Max. user weight ▼

136 kg



Battery capacity ▼

2 x 12 Ah
2 x 18 Ah



Turning diameter ▼

2200 mm



Max. climbable obstacle height ▼

45 mm



Max. safe slope ▼

6° / 10.5%



Motor capacity ▼

200 W



Max. driving distance ▼

Up to
16 km



Max. driving speed ▼

8 km/h

Colours*



Sapphire Blue



Jasper Red



Jade Green



Black Vinyl

* Scooter colours may vary slightly. Please check with your retailer.

Upholstery

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland



Yes, you can.®