



reddot award 2014  
winner



### The colourful route to an effortless, safe drive

The 4 wheel *Invacare Colibri* is a stylish, simple and colourful micro lightweight scooter, designed for those who enjoy a totally independent lifestyle.

Thanks to the unique patented *Invacare LiteLock™* system the scooter can easily be taken apart without the need for tools and fits neatly into any car boot. Plus, with its small footprint the *Invacare Colibri* offers greater manoeuvrability in tight, or limited spaces.

Trusted Solutions, Passionate People



Yes, you can.®





**Features and Options**



**LiteLock system - easy to take apart and reassemble**

The ingenious and unique patented *Invacare LiteLock™* system automatically releases once the battery box has been removed. It takes less than a minute to take the scooter apart - without tools. Its front section can be simply lifted and placed in the boot of a car, along with the other sections. It's just as easy to reassemble and the scooter is secure.



**Effortless Dismantling**

It takes less than a minute to take the *Invacare Colibri* apart ready for transportation.



**Option of a bigger wheel set and battery pack**

Perfect for travelling longer distances, on bumpier terrain, such as cobbled streets. Choose where to charge your battery: directly on the *Colibri* or detached from it.



**Increased comfort**

The soft, padded seat provides comfort during driving. For increased comfort there is also the optional 'easy fit' seat suspension kit.

**Technical data**

465 mm	400 mm	370 - 420 mm	385 mm	505 mm
1010 mm	44 kg (12 Ah) 49 kg (18 Ah)	136 kg	2 x 12 Ah 2 x 18 Ah	2200 mm
45 mm	6° / 10.5%	200 W	Up to 16 km	8 km/h



**LiteLock™ system**

Invacare Limited  
Pencoed Technology Park -  
Pencoed - Bridgend CF35 5AQ -  
United Kingdom  
Tel: +44 1 656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS -  
County Dublin - Ireland  
Tel. +353 1 8107084  
Fax +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

© 2022 Invacare International Sàrl All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Colibri - UK - 04/2022 - 1576204

**Shroud colours**




Sapphire Blue



Jasper Red



Jade Green



Marble White



Diamond Silver

<sup>1/</sup> Compliant with the directive on the Restriction of the use of certain Hazardous Substances in electrical and electronic equipment - RoHS (2011/65/EC)



Follow us on Twitter:  
[@Think\\_Mobility](https://twitter.com/Think_Mobility)



Like us on Facebook:  
Invacare UK



Follow us on YouTube:  
Invacare Ltd - UK



**Yes, you can.®**