

# Comet<sup>®</sup> Pro

Feel the Freedom



A premium scooter for advanced mobility

giving exceptional power and performance



Yes, you can.™

# Feel the difference

Our scooters have been designed together with healthcare experts, to provide unrivalled standards of quality, safety and comfort.

Through continuous dedicated research and an unrivalled experience in designing for mobility, at Invacare we pride ourselves in scooters that improve quality of life. Being able to regain confidence, trust and independence from your scooter is important.

From ergonomic steering controls to the adaptable cushioned seats and adjustable seat risers, every detail of our scooters is considered and refined to make getting about feel great.

## Quality



Market-leading build quality with exceptional attention to detail.

## Safety



Compliant with the latest European Medical Device Standards.

## Comfort



Adjustable design to provide lasting comfort.



# Feel in control with Comet® Pro

The Comet Pro is our most advanced outdoor scooter. Designed with exceptional features and options to take you wherever you want to go in comfort, safety and style. It's a design to provide drive performance and total control without compromise.


Comet Pro features exceptional driving features and modern aesthetics to deliver premium performance and control. The 2 x 12 V / 75 Ah batteries allow users to enjoy challenging journeys beyond the city's limits.

It also has large front and rear wheels, combined with a superior suspension system that allows safe negotiation of kerbs and obstacles smoothly and comfortably.

Up to  
15 km/h  
max speed

Up to  
55 km  
max distance  
on a single  
charge

160 kg  
max user  
weight



## Improved performance

- Easy drive direction to keep you in a straight line at high speeds
- Advanced suspension ensures a smooth and comfortable drive over uneven terrain
- Ergonomic steering for added control
- Designed for all-terrain performance
- Anti-splash guards



## Safety in mind

- Automatic speed reduction helps reduce speed around corners
- Two-step disengaging lever to prevent free-wheeling
- Brake lights and indicators
- Headlights for improved visibility



## Built around comfort

- Adjustable steering column for personalised positioning
- High quality cushioned seats
- Adjustable seat height
- Flip up armrests and swivel seat makes getting on and off easy



## Auto-speed reduction

reduces speed around corners for added safety



## Optional connectivity

Designed with Zeta™ connectivity, an innovative technology allowing access to your mobility scooter's battery and journey data via an easy-to-use mobile app.



## Adjustable steering column

for responsive control and improved comfort



Improved suspension



Reinforced seat column



Height and width adjustable armrest



## Add to your scooter

Tailor the Comet® Pro to your taste and needs. Turn over to explore the optional extras available.

## Add to your scooter

We have a wide range of optional extras available, so you can customise your Comet Pro scooter. Ask your retailer to discover the full range.



Rear basket



Side & rear view mirrors



Front & rear bumper



Mobile phone holder



Crutch holder



Seatbelt



Cupholder



Interchangeable shroud panels

## Technical data

For more information about this product, including the user manual, please visit our website [www.invacare.eu.com](http://www.invacare.eu.com)



Backrest width ▼

510 mm



Seat depth ▼

Deluxe: 470 mm  
Premium: 460 mm



Seat to floor height ▼

Min. 440 mm  
Max. 510 mm



Backrest height ▼

Deluxe: 490 mm  
Premium: 630 mm



Driving unit width ▼

665 - 685 mm



Total length ▼

1450 mm



Empty weight ▼

136 kg



The heaviest part ▼

63 kg



Max. user weight ▼

160 kg



Battery capacity ▼

2 x 12 V/75 Ah (C20)  
AGM



Turning diameter ▼

2750 mm



Max. climbable obstacle height ▼

100 mm



Ground clearance ▼

100 mm



Max. safe slope ▼

10°



Motor capacity ▼

12 km/h: 550 W / 1500 W<sup>1</sup>  
15 km/h: 550 W / 1600 W<sup>1</sup>



Max. driving distance ▼

12 km/h: 55 km  
15 km/h: 55 km



Max. driving speed ▼

12 / 15 km/h

<sup>1</sup> Peak performance

## Colours\*



Marble White



Alium Blue



Diamond Silver



Garnet Red



Black Vinyl

\*Scooter colours may vary slightly. Please check with your retailer.

## Upholstery

**Tel: +41 61 487 70 70**  
hqeurope@invacare.com  
[www.invacare.eu.com](http://www.invacare.eu.com)

Invacare International GmbH  
Neuhofweg 51  
4147 Aesch BL  
Switzerland



Yes, you can.®