

# Comet® Ultra

Feel the Freedom



**A premium scooter with extra weight capacity**  
for exceptional power and performance



**Yes, you can.**



# Feel the difference

Our scooters have been designed together with healthcare experts, to provide unrivalled standards of quality, safety and comfort.

Through continuous dedicated research and an unrivalled experience in designing for mobility, at Invacare we pride ourselves in scooters that improve quality of life. Being able to regain confidence, trust and independence from your scooter is important.

From ergonomic steering controls to the adaptable cushioned seats and adjustable seat risers, every detail of our scooters is considered and refined to make getting about feel great.

## Quality



Market-leading build quality with exceptional attention to detail.

## Safety



Compliant with the latest European Medical Device Standards.

## Comfort



Adjustable design to provide lasting comfort.





# Feel in control with Comet® Ultra

Ideal for those with larger body shapes, the Comet® Ultra offers a secure, supportive and powerful driving experience without compromising on style or functionality.

Packed with all of the safety features and driving characteristics seen throughout the Comet® family, the Comet® Ultra will get you to

your destination quickly and efficiently, in complete comfort and control.

Up to  
10 km/h  
max speed

Up to  
49 km  
max distance  
on a single  
charge

220 kg  
max user  
weight



## Improved performance

- Easy drive direction to keep you in a straight line at high speeds
- Advanced suspension ensures a smooth and comfortable drive over uneven terrain
- Ergonomic steering for added control
- Designed for all-terrain performance
- Anti-splash guards



## Safety in mind

- Automatic speed reduction for safe cornering
- Two-step disengaging lever to prevent free-wheeling
- Brake lights and indicators
- Headlights for improved visibility



## Built around comfort

- High quality cushioned seat
- Adjustable steering column
- Adjustable seat height
- Flip-up armrests and swivel seat makes getting on and off easy



## Auto-speed reduction

reduces speed  
around corners for  
added safety



## Optional connectivity

Designed with Zeta™ connectivity, an innovative technology allowing access to your mobility scooter's battery and journey data via an easy-to-use mobile app.



## Adjustable steering column for responsive control and improved comfort



Improved suspension



Reinforced seat column



Height and width adjustable reinforced armrest for extra support



## + Add to your scooter

Tailor the Comet® Ultra to your taste and needs. Turn over to explore the optional extras available.



## Add to your scooter

We have a wide range of optional extras available, so you can customise your Comet® Ultra scooter. Ask your retailer to discover the full range.



Rear basket



Side & rear view mirrors



Front & rear bumper



Mobile phone holder



Crutch holder



Seatbelt



Cupholder



Interchangeable shroud panels

## Technical data

For more information about this product, including the user manual, please visit our website [www.invacare.eu.com](http://www.invacare.eu.com).



Backrest width ▼

660 mm



Seat depth ▼

535 mm



Seat to floor height ▼

Min. 440 mm  
Max. 510 mm



Backrest height ▼

515 mm



Driving unit width ▼

665 - 685 mm



Total length ▼

1450 mm



Empty weight ▼

148 kg



The heaviest part ▼

approx. 63 kg



Max. user weight ▼

220 kg



Battery capacity ▼

2 x 12 V/75 Ah  
(C20) AGM



Turning diameter ▼

2750 mm



Max. climbable obstacle height ▼

100 mm



Ground clearance ▼

100 mm



Max. safe slope ▼

1 °



Motor capacity ▼

650 W / 1800 W <sup>1</sup>



Max. driving distance ▼

49 km



Max. driving speed ▼

10 km/h

<sup>1</sup> Peak performance

## Colours\*



Marble White



Alum Blue



Black Vinyl

\*Scooter colours may vary slightly. Please check with your retailer.

## Upholstery

Tel: +41 61 487 70 70  
[hqeurope@invacare.com](mailto:hqeurope@invacare.com)  
[www.invacare.eu.com](http://www.invacare.eu.com)

Invacare International GmbH  
Benkenstrasse 260  
4108 Witterswil  
Switzerland



Yes, you can.®