



Product Service

# CERTIFICATE

No. Z1 09 11 52355 016

**Holder of Certificate:** **Küschall AG**  
 Benkenstrasse 260  
 4108 Witterswil  
 SWITZERLAND

**Production Facility(ies):** 52355

**Certification Mark:**



**Product:** **Active wheelchairs**

**Model(s):** **Ultra-Light / Compact**

<b>Parameters:</b>		Ultra-Light	Compact
Total length incl. footrests		730-1100 mm	740-1160 mm
Total length without footrests		----	570-850 mm
Wheelchair width		490-750 mm	490-750 mm

- Further parameters see enclosure

**Tested according to:** DIN EN 12182:1999  
 DIN EN 12183:1999

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition the certification holder must not transfer the certificate to third parties. See also notes overleaf.

**Test report no.:** 71358709

**Date,** 2009-11-25 (Torsten Zimmer)  
 Page 1 of 2





Enclosure to Certificate No.: Z1 09 11 52355 016

**Model designation for the active wheelchairs Ultra-Light / Compact**

	<b>Ultra-Light</b>	<b>Compact</b>
Total length incl. foot rests	730 - 1100 mm	740 - 1160 mm
Total length without foot rests	----	570 - 850 mm
Wheelchair width	490 - 750 mm	
Wheelchair height	690 - 1010 mm	670 - 1010 mm
Seat width	280, 300, 320, 340, 360, 380, 400, 420 mm (Single diagonal strut) 440, 460, 480, 500 mm (Double diagonal strut)	
Seat depth	320 - 500 mm	
Seat height at the front edge	430 - 530 mm	370 - 530 mm
Seat angle	-2° - 16°	
Width of the backrest	280 - 500 mm	
Height of the backrest	300 - 510 mm	
Angle of the backrest	82° - 102°	
Max. user weight	100 kg (Single diagonal strut) 130 kg (Double diagonal strut)	
Driving wheels		
- Size	22, 24, 25 and 26 x 1" - 1 3/4"	
- Tyres	<input checked="" type="checkbox"/> air-filled	<input checked="" type="checkbox"/> accident proof
- Adjustability	<input checked="" type="checkbox"/> yes	<input type="checkbox"/> no
- Kind of fixation	<input checked="" type="checkbox"/> quick release	<input type="checkbox"/> rigid
- Adjustment of the camber	0° / 3°	
Front wheels		
- Size	3, 4, 5, 6 and 7 x 3/4" - 1 1/4"	
- Tyres	<input checked="" type="checkbox"/> air-filled	<input checked="" type="checkbox"/> accident proof
- Adjustability	<input checked="" type="checkbox"/> yes	<input type="checkbox"/> no
- Kind of fixation	<input type="checkbox"/> quick release	<input checked="" type="checkbox"/> rigid
Brake system	Knee lever brakes	
Frame material	Aluminium	
Arm rests	Height adjustable and removable	
Foot rests	Height adjustable, foot plate(s) can be swivelled upwards	Height adjustable, can be turned aside and removable, foot plate(s) can be swivelled upwards
Anti tip device	existing	

Department: MHS 5 / me  
Date: 2009-11-25