

Certificate in terms of transportation

We hereby confirm that the removable battery packs used in our products

e-fix (E25, E26, Esprit Action 4 NG)^{P(small),M(extended)}
e-fix Esprit Action (E34)^{K,T}
e-fix eco (E35)^T
viamobil eco (V14)^P
viaplus (V12)^P
scalamobil (S30, S31, S34, S35, S36, S38, S39)^P
Adventure (A10)^{P,M}

are containing non-spillable lead-acid batteries (VRLA) –with referenced to:

UN Nr. 2800 „Batteries, gel, non-spillable, electric storage“

IATA: Non-spillable batteries are not subject to DG regulations since they meet the requirements of Packing Instructions 872 of Special Provision A67.

IMDG: According IMDG Provision 238 the batteries are not subject to the provisions of the IMDG code.

ADR/RID: New and spent batteries are exempt from all ADR/RID (Special Provision 598)

US DOT: Non-spillable batteries are not subject to DG regulations , since they meet the requirements of 49 CFR 173.159(d).

Albstadt, 2021-06-21
Alber GmbH



Hansjörg Reiner
Managing Director

P – Panasonic LC Series, M – Multipower MP Series (Jewo), K – Kung Long VRLA, T – Tianneng TNE