

Propad Leg Trough



Shaped leg support to facilitate
comfort and recovery



The *Invacare Propad Leg Trough* has been designed in consultation with clinicians to provide maximum support for the whole leg.

Offers comfort and support to facilitate recovery

Provides optimum support and protection for the whole leg, reducing pressure and facilitating recovery.

Relieves pressure on vulnerable heels

Elevates the heels and allows them to be suspended over the edge of the trough, providing pressure relief for the vulnerable heel area.

Features and Options



Water resistant cover

Two-way stretch, vapour permeable and anti-static polyurethane (PU) water resistant cover.



Available in three sizes

Wide range of patients can be accommodated.



Cover includes a zip, and is removable

Facilitates laundering of the cover.



Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Length	Width	Height	Total product weight	Max. cleaning temperature
<i>Propad Leg Trough</i>	510, 610, 660 mm	110 mm	115 mm	0.8 kg	80°

Manufactured to European Standards:
 Interior Foam: BS 5852: 1982 Crib 5, BSEN597-1,2
 Cover: BS7175: 1989 Crib 5

Colour

Blue



QR Code

Invacare Ltd
 Pencoed Technology Park - Pencoed
 Bridgend CF35 5AQ - United Kingdom
 Tel: +44 1 656 776 222
 E-mail: uk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
 Unit 5 - Seatown Business Campus Seatown Road - SWORDS
 County Dublin - Ireland
 Tel. +353 1 8107084
 Email: Ireland@invacare.com
www.invacare.ie



Yes, you can.®