

Softform Heelpad



Pressure reducing Heelpad



The *Invacare Softform Heelpad* is designed to reduce pressure on the vulnerable heel area. The controlled volume of gel in the sacs allows pressure displacement and even weight distribution. Pressure on a patient's heels are kept to a minimum when the heels are in contact with the gel sacs. The *Heelpad* was designed to be used in conjunction with the *Invacare Odstock Wedge*, but can also be used independently as an effective pressure reducing device.



Reduces pressure on the heels

The *Invacare Softform Heelpad* enables vulnerable heels to be immersed into independent gel sacs, offering significant pressure reduction and comfort to users.

Reduces the requirement for patient handling

When used in conjunction with the *Odstock Wedge*, the *Heelpad* will provide improved support for the lower limbs.

Features and Options



Single and Double Heelpads available

Single and Double formats available to meet specific patient requirements.



Heelpads supplied with three platilon covers

Platilon protective covers are disposable, and facilitate infection control.








Designed to be used in conjunction with the *Invacare Odstock Wedge*

Used in conjunction, these two products combine to significantly reduce pressure on the vulnerable heel area.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local *Invacare* website.

Width	Length	Height	Total product weight	Max. cleaning temperature
				
620 mm	480 mm	50 mm, 32 mm	3.58 kg	80°

Softform Heelpad

Manufactured to European Standards:
Interior Foam: BS 5852: 1982 Crib 5, BSEN 597-1,2
Cover: BS 7175: 1989 Crib 5

Colour

Fuchsia



QR Code

Invacare Ltd
Pencoed Technology Park - Pencoed
Bridgend CF35 5AQ - United Kingdom
Tel: +44 1 656 776 222
E-mail: uk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS
County Dublin - Ireland
Tel: +353 1 8107084
Email: Ireland@invacare.com
www.invacare.ie

