

Softform Flexipad



Range of highly versatile support pads



Invacare **Flexipads** are highly versatile pads filled with silicone gel, and can be used to position and support limbs, joints, the torso or head.

Facilitates positioning

Invacare **Flexipads** can be used to position a patient's limbs for treatment.

Facilitates support

Can provide supplementary short-term support for patients experiencing difficulty in maintaining posture.

Features and Options



Single Flexipad

Fitted with a strap for fastening to chairs, the **Invacare Flexipad** increases the surface area under elbows, arm or back.



Double Flexipad

Two gel sacs with a foam core – highly versatile in a number of applications, including supporting large bony prominences.








Flexipad Ultra

A lighter, more flexible format, available with or without straps.



Technical data


For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Width	Length	Height	Total product weight	Max. cleaning temperature
 220 mm	 330 mm	 50 mm, 32 mm	 0.8 kg	 80°

Softform Flexipad

Manufactured to European Standards:
Cover: BS 7175: 1989 (Crib 5)

Colours



QR Code

Invacare Ltd
Pencoed Technology Park - Pencoed
Bridgend CF35 5AQ - United Kingdom
Tel +44 1 656 776 222
E-mail: uk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus - Seatown Road -
SWORDS, County Dublin - Ireland
Tel. +353 1 8107084
Email: Ireland@invacare.com
www.invacare.ie

