



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

# Manual Retail Prescription Form

APRIL 2024  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Rea Dahlia 30

### Standard Features

- Max User Weight 135KG/21 Stone
- Crash Tested ISO 7176-Part 19
- 16" Transit Propulsion with Attendant Hub Brake and Solid Black Tyres
- Or Self Propel 20", 22" or 24" Self Propel Wheels
- Choice of Seat Widths 390mm, 440mm or 490mm (+50mm or +100mm Width Adj)
- Choice of 8 Seat Heights 330mm - 500mm (13" - 20")
- Depth Adjustable Seat 430mm - 500mm (17" - 20")
- Without Seat Cushion, Supplied with Seat Board
- Flex 3 Backrest Plate, Width and Height Adjustable, with Neckrest
- Passad 2 Backrest Cushion with Multi-stretch Black PU coated fabric Cover and Push Handles Connected
- Assistant Operated Gas Piston Seat Tilt +30° and Backrest Recline 30° with Tilt Lock
- Fixed 80° Legrests Swing In Swing Out. Height, Depth and Angle Adjustable Footplates
- Height Adjustable Ducktail Pad Armrests with Autolock
- Black Frame Colour, Solid Black Front Castors, User Manoeuvred Tyre Brakes
- Pelvic Belt, Adj both sides, Anti Tippers, Front Transportation Brackets

**Standard Base price: £2,550**



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description		
<b>BASE PRODUCT*</b>			
BVA0010	DAHLIA 30 °		£2,550
BVA8000	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE		POA
<b>PROPULSION*</b>			
BVA0020	Transit 16" wheels with hub brakes	(NC 0140, 0150, 0160, 0190, 1555, 1565, 1575, 1584, 1588)	<input type="radio"/>
BVA0030	Self Propel 24"/22"/20" wheels with user brakes	(NC 0090, 0180, 1545, 1588)	<input type="radio"/>
BVA0042	One Arm Drive Right with Quick Release Rear Wheels	(NC 0130, 0140, 0150, 0180, 0190, 1545, 1555, 1584, 1586)	£431
BVA0044	One Arm Drive Left with Quick Release Rear Wheels	(NC 0130, 0140, 0150, 0180, 0190, 1545, 1555, 1584, 1586)	£431

SEAT WIDTH*				
BVA0050		390 mm/15" (Small chassis)	(NC 0463, 0464, 0465, 0466)	○
BVA0060		440 mm/17" (Medium chassis)	(NC 0461, 0465, 0466)	○
BVA0070		490 mm/19" (Large chassis)	(NC 0461, 0462, 0463)	○
SEAT WIDTH EXTENSION*				
BVA0080		+ 50 mm pre-adjustment on armrests and legrests		○
BVA0090		+ 100 mm pre-adjustment on armrests and legrests		○
BVA0100		Without seat width extension, set at 0mm adjustment but can be adjusted wider if required	(NC 0460, 0470)	○
SEAT DEPTH*				
BVA0110		430-500 mm / 17" - 20"		◆
SEAT HEIGHT*				
BVA0130		350 mm/ 14"	(NC 0225, 0637, 0670, 1435, 1455, 1475, 1495, 1515, 1575)	○
BVA0140		380 mm/ 15"	(NC 0200, 0220, 0227, 0230, 0637, 0670, 1435, 1440, 1495, 1515, 1545, 1565, 1575)	○
BVA0150		400 mm/ 16"	(NC 0200, 0230, 0670, 1435, 1440, 1495, 1515, 1545, 1575)	○
BVA0160		430 mm/ 17"	(NC 0200, 0230, 0670, 1435, 1440, 1495, 1515, 1545, 1555)	○
BVA0170		450 mm/ 18"	(NC 0200, 0210, 0225, 0720, 0770, 1435, 1440, 1555, 1565)	○
BVA0180		480 mm/ 19"	(NC 0200, 0210, 0225, 0227, 0720, 0770, 1435, 1440, 1455, 1475, 1555, 1565, 1575, 1725, 1727)	○
BVA0190		500 mm/ 20"	(NC 0200, 0210, 0225, 0227, 0720, 0770, 1435, 1440, 1455, 1475, 1545, 1565, 1575)	○
SEAT TILT ANGLE*				
BVA0200		Reduced Tilt: Positive Seat Angle 10°		○
BVA0210		Reduced Tilt: Positive Seat Angle 20°		○
BVA0220		Full Tilt: Positive Seat Angle 30°		○
BACKREST RECLINE ANGLE*				
BVA0225		Reduced Backrest Recline: Angle 10°		○
BVA0227		Reduced Backrest Recline: Angle 20°		○
BVA0230		Full Backrest Recline: Angle 30°		○
NEGATIVE SEAT ANGLE				
BVA0240		Negative Seat Angle -1 to aid standing transfers		○
SEAT CUSHION*				
BVA0270		Cushion Flosshape Black Multi-stretch PU coated fabric		£57
SWO0271		Cushion Flosshape Black Multi-stretch PU coated fabric seat upholstery Shorter depth 450mm deep. Only use with seat depth in shortest position		£57
BVA0274		Cushion Flosshape Black Spacer		£57
BVA0300		Cushion Vicair Adjuster O2 (10cm)		£547
BVA0320		Without Seat Cushion, Seat Plate only		◆
BACKREST TUBES*				
BVA0360		Straight Back Tubes (Flex 3 and Matrix)		◆
BACKREST RECLINE PIVOT POINT*				
BVA0370		Plus 50mm width extension, outrigs backrest pivot points by 5cm overall to allow a wider backrest to be fitted to a smaller chassis	(NC 0390, 0400, 0470)	£84
BVA1111		Without Backrest Pivot Point Width Extension		◆

<b>BACKREST*</b>			
BVA0390		Matrx Elite E2 (includes hardware mounted to frame, backrest supplied in box)	(NC 0460, 0470, 0490, 0500, 0510, 0520, 0570, 1140, 1150, 1160, 1170, 1240, 1250, 1260, 1280) £494
BVA0400		Matrx Elite E2 Deep (includes hardware mounted to frame, backrest supplied in box)	(NC 0460, 0470, 0490, 0500, 0510, 0520, 0570, 1140, 1150, 1160, 1170, 1240, 1250, 1260, 1280) £546
BVA0430		Flex 3 Backrest, Fixed Plate (Width Adjustable)	(NC 0455, 0456, 0457, 0461, 0462, 0463, 0464, 0465, 0466, 0490) ◆
BVA0440		Flex 3 Backrest Tension (Width Adjustable)	(NC 0455, 0456, 0457, 0461, 0462, 0463, 0464, 0465, 0466, 0470) £102
BVA0450		Backrest Tubes Without Backrest	(NC 0455, 0456, 0457, 0460, 0461, 0462, 0463, 0464, 0465, 0466, 0470, 0490, 0500, 0510, 0520, 1080, 1100, 1110, 1140, 1150, 1160, 1170, 1240, 1250, 1260, 1280) ○
<b>MATRX BACKREST HEIGHT</b>			
BVA0455		Matrx Backrest Height 410 mm / 16"	○
BVA0456		Matrx Backrest Height 460 mm / 18"	○
BVA0457		Matrx Backrest Height 510 mm / 20"	○
<b>MATRX BACKREST WIDTH</b>			
BVA0461		Matrx Backrest Width 380mm / 15" (Chest width 290mm / 11 1/2")	○
BVA0462		Matrx Backrest Width 410mm / 16" (Chest width 315mm / 12 1/2")	○
BVA0463		Matrx Backrest Width 430mm / 17" (Chest width 340mm / 13")	○
BVA0464		Matrx Backrest Width 460mm / 18" (Chest width 370mm / 14 1/2")	○
BVA0465		Matrx Backrest Width 480mm / 19" (Chest width 390mm / 15")	○
BVA0466		Matrx Backrest Width 510mm / 20" (Chest width 420mm / 16 1/2")	○
<b>BACKREST WIDTH EXTENSION*</b>			
BVA0460		+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required	○
BVA0470		+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required	(NC 1725, 1727) ○
BVA0480		Without Backrest width pre-adjustment, set at 0mm but can be adjusted if required	○
<b>BACKREST CUSHION*</b>			
BVA0490		Cover for tension adjustable upholstery	○
BVA0500		Cushion, Laguna (Lateral Support)	(NC 2222) ○
BVA0510		Cushion, Mistral 2 (Waist Support)	(NC 0575, 2222) £144
BVA0520		Cushion, Passad 2 (Shoulder Support)	(NC 2222) ◆
BVA0550		Without Backrest Cover/ Cushion	(NC 0570, 0575) ○
<b>BACKREST UPHOLSTERY*</b>			
BVA0570		Black Multi-stretch Black PU coated fabric	◆
BVA0575		Black Spacer TR36	○
BVA2222		No upholstery, select with Matrx and Without backrest options	○
<b>PUSH HANDLES*</b>			
BVA0580		Push Bar, Angle Adjustable	£102
BVA0600		Push Handles connected	◆
<b>SEAT TILT*</b>			
BVA0620		Assistant operated gas piston seat tilt	(NC 3333) ◆
BVA0632		User operated seat tilt, with height adjustable armrest. Limited to 10 degrees of tilt	(NC 0915, 0920, 0930, 0950, 0975, 0980, 0990, 1010, 1750, 1755, 1756, 1757, 1760, 1765, 1770, 1775, 1780, 1790, 1850, 3333) £84
BVA0637		Electrical operated seat and backrest tilt	(NC 0640, 0670) £616

<b>BACKREST REGULATION*</b>			
BVA0640		Assistant Operated gas piston for backrest recline	◆
BVA3333		Only select if option 0637 Electrical operated is required	○
<b>CENTRAL LEGREST</b>			
BVA0670		Central Legrest, includes calfpad and footrest. Also select options 0740, 0790, 4444, 5555	(NC 0681, 0690, 0692, 0700, 0710, 0720, 0730, 0750, 0760, 0770, 0780, 0805, 0855, 0865, 1040, 1310, 1320, 1330, 1340, 1350, 1351, 1352, 1353) £303
<b>LEGREST HANGER ATTACHMENT*</b>			
BVA0680		Standard position for legrest mounting	◆
SWO0681		Standard position with 0-30° Legrest Abduction function	(NC 0670, 0805) £100
BVA0690		High position 50mm higher for leg length compensation	○
BVA0692		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest	(NC 0710, 0720, 0730, 0760, 0770, 0780) ○
<b>LEGREST HANGER LH*</b>			
BVA0700		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH	(NC 0770, 0855, 1330, 8888) ○
BVA0710		Legrest Hanger fixed 80 deg, swing in/swing out LH	(NC 1302, 1304, 1310) ◆
BVA0720		Legrest Hanger fixed 90 deg, swing in/swing out LH	(NC 0750, 0865, 1302, 1304, 1310, 1495, 1515) ○
BVA0730		Legrest Amputee LH	(NC 0805, 0825, 0855, 0865, 1302, 1304, 1310, 1330, 1340, 1350) £137
BVA0740		Without Legrest, LH	(NC 0805, 0825, 0855, 0865, 1302, 1304, 1310, 1330, 1340, 1350) ○
<b>LEGREST HANGER RH*</b>			
BVA0750		Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support with locking clip RH	(NC 0855, 1330, 8888) ○
BVA0760		Legrest Hanger fixed 80 deg, swing in/swing out RH	(NC 1312, 1314, 1320) ◆
BVA0770		Legrest Hanger fixed 90 deg, swing in/swing out RH	(NC 0865, 1312, 1314, 1320, 1495, 1515) ○
BVA0780		Legrest Amputee RH	(NC 0805, 0845, 0855, 0865, 1312, 1314, 1320, 1330, 1340, 1350) £137
BVA0790		Without Legrest, RH	(NC 0805, 0845, 0855, 0865, 1312, 1314, 1320, 1330, 1340, 1350) ○
<b>FOOTREST LH*</b>			
BVA0805		Footrest one piece, depth and angle adjustable short tube	(NC 0845, 1310, 1320, 1351, 1352, 1353) £118
BVA0825		Footrest lower part, depth and angle adjustable with heel strap, LH	(NC 1330) ○
BVA4444		Without footrest, select with options Central legrest, Amputee legrest or Without legrest	○
<b>FOOTREST RH*</b>			
BVA0845		Footrest lower part, depth and angle adjustable with heel strap, RH	(NC 1330) ○
BVA5555		Without footrest, select with options Central legrest, Amputee legrest, Without legrest or One piece footrest	○
<b>FOOTREST TUBE LENGTH*</b>			
BVA0855		Short Tube 24CM	○
BVA0865		Long Tube 36CM	○
BVA6666		No Tube length. Select with Central Legrest, Amputee legrest or Without legrest	○

ARMREST LH*			
BVA0900		Armrest Black PVC Moulded Narrow Pad LH, Height Adjustable from 200 mm - 320 mm with Autolock. Adjustment taken from seat plate to top of armpad	(NC 1756) ○
BVA0910		Armrest Ducktail Pad LH, Height Adjustable from 200 mm - 320mm with Autolock	(NC 1756) ◆
BVA0915		Armrest Extended Height, Ducktail pad LH, Adjusts from 255 mm - 355 mm with Autolock	(NC 1756) ○
BVA0920		Armrest Wide Pad LH, Height Adjustable from 200 mm - 320mm with autolock.	(NC 1756) ○
BVA0930		Armrest Hemiplegic LH, Height Adjustable from 200 mm - 320 mm with autolock	(NC 1750, 1755, 1756, 1760, 1780, 1790) £179
BVA0950		Without Armrest LH Select when choosing BVA1756 Half Tray LH, Narrow Armrest is included	(NC 1020, 1750, 1755, 1760, 1780, 1790) ○
ARMREST RH*			
BVA0960		Armrest Black PVC Moulded Narrow Pad RH, Height Adjustable from 200 mm - 320 mm with Autolock. Adjustment taken from seat plate to top of armpad	(NC 1757) ○
BVA0970		Armrest Ducktail Pad RH, Height Adjustable from 200 mm - 320mm with Autolock	(NC 1757) ◆
BVA0975		Armrest Extended Height, Ducktail pad RH, Adjusts from 255 mm - 355 mm with Autolock	(NC 1757) ○
BVA0980		Armrest Wide Pad RH, Height Adjustable from 200 mm - 320mm with autolock.	(NC 1757) ○
BVA0990		Armrest Hemiplegic RH, Height Adjustable from 200 mm - 320 mm with autolock	(NC 1020, 1750, 1755, 1760, 1770, 1775, 1780, 1790) £179
BVA1010		Without Armrest RH Select when choosing BVA1757 Half Tray RH, Narrow Armrest is included	(NC 1020, 1750, 1755, 1760, 1770, 1775, 1780, 1790) ○
ARMREST OPTIONS			
BVA1020		Pad for sideguard, removable with hook and loop fastening	£51
SEAT OPTIONS			
BVA1030		Incontinence Cover	£40
BVA1040		Pommel Abduction Cushion	(NC 1060) £118
BVA1050		Seat Plate Extension Plates. Width, Forward, Rear or Assymetric (Supplied as a kit of 2 metal plates unmounted)	(NC 1060) £42
BVA1060		2 cm Increased seat height plate. Allows an increased 20mm of seat height by adding a foam plate to the seat board with hook and loop	£64
BVA1070		Tilt Scale for seat tilt	£44
HEADREST OPTIONS*			
BVA1080		Headrest, Sideway Adjustable	(NC 1087, 1088, 1089, 1110) ○
BVA1100		Headrest with cheek support, sideway adjustable	(NC 1087, 1088, 1089, 1090, 1110) £107
BVA1081		Matrx Elan with Standard pad (254mm / 10")	(NC 1090, 1110, 9999) £124
BVA1082		Matrx Elan with Standard pad (355mm / 14")	(NC 1090, 1110, 9999) £124
BVA1083		Matrx Elan with Occipital pad (229mm / 9")	(NC 1090, 1110, 9999) £124
BVA1084		Matrx Elan with Occipital pad (305mm / 12")	(NC 1090, 1110, 9999) £124
BVA1085		Matrx Elan with 4-Point pad (229mm W x 177mm H / 9" W x 7" H)	(NC 1090, 1110, 9999) £232
BVA1086		Matrx Elan with 4-Point pad (280mm W x 254mm H / 11" W x 10" H)	(NC 1090, 1110, 9999) £232
BVA9991		Without headrest	(NC 1087, 1088, 1089) ○

<b>MATRIX HEADREST HARDWARE*</b>				
BVA1087	<input type="checkbox"/>	Matrix Elan hardware	(NC 1090, 1110, 1130)	£184
BVA1088	<input type="checkbox"/>	Matrix Elan extended hardware, (Extra link)	(NC 1090, 1110, 1130)	£284
BVA1089	<input type="checkbox"/>	Matrix Elan LOXX Hardware	(NC 1090, 1110, 1130)	£166
BVA9999	<input type="checkbox"/>	Without Elan headrest hardware (Only select with option 1080, 1100 Headrest and without headrest)		○
<b>OTHER HEADREST OPTIONS</b>				
BVA1090	<input type="checkbox"/>	Cover for headrest	(NC 1110)	£31
BVA1110	<input type="checkbox"/>	Neckrest		◆
BVA1130	<input type="checkbox"/>	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)		£48
BVA9992	<input type="checkbox"/>	Without neckrest		○
<b>TRUNK SUPPORTS LH</b>				
BVA1140	<input type="checkbox"/>	Trunk Support Multi Adjustable with standard arm LH	(NC 1260)	○
BVA1160	<input type="checkbox"/>	Trunk Support Multi Adjustable with extended height arm LH	(NC 1260)	○
BVA1240	<input type="checkbox"/>	Trunk Support Multi Adjustable Swingaway LH	(NC 1260)	£113
<b>TRUNK SUPPORTS RH</b>				
BVA1150	<input type="checkbox"/>	Trunk Support Multi Adjustable with standard arm RH	(NC 1260)	○
BVA1170	<input type="checkbox"/>	Trunk Support Multi Adjustable with extended height arm RH	(NC 1260)	○
BVA1250	<input type="checkbox"/>	Trunk Support Multi Adjustable Swingaway RH	(NC 1260)	£113
<b>BACKREST OPTIONS</b>				
BVA1260	<input type="checkbox"/>	Trunk Support Attachment. If trunk supports are chosen, the attachment bracket is already included		○
BVA1280	<input type="checkbox"/>	Wedges (One Pair)		£11
BVA1300	<input type="checkbox"/>	Tilt Scale for backrest recline		£44
<b>ELEVATING LEGREST OPTIONS*</b>				
BVA1302	<input type="checkbox"/>	Calf pad cover, for LH elevating legrest		○
BVA1304	<input type="checkbox"/>	Knee pad for LH elevating legrest		◆
BVA1312	<input type="checkbox"/>	Calf pad cover, for RH elevating legrest		○
BVA1314	<input type="checkbox"/>	Knee pad for RH elevating legrest		◆
BVA8888	<input type="checkbox"/>	Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)		○
<b>LEGREST OPTIONS</b>				
BVA1310	<input type="checkbox"/>	Footplate Extender LH	(NC 1340)	£34
BVA1320	<input type="checkbox"/>	Footplate Extender RH	(NC 1340)	£34
BVA1330	<input type="checkbox"/>	Calf Strap (for One piece footrest)		○
BVA1340	<input type="checkbox"/>	Cover Footplate, Black Plush Fabric		£54
BVA1350	<input type="checkbox"/>	Calf Strap Padded		£73
<b>BODYPOINT ANKLE HUGGERS</b>				
BVA1351	<input type="checkbox"/>	Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)		£114
BVA1352	<input type="checkbox"/>	Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)		£114
BVA1353	<input type="checkbox"/>	Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)		£114
<b>CHASSIS COLOUR*</b>				
BVA1405	<input type="checkbox"/>	Black Frame		◆
<b>CASTOR ATTACHMENT*</b>				
BVA1410	<input type="checkbox"/>	Fixed		◆

<b>FRONT CASTORS*</b>			
BVA1435		100x27 mm (4" x 1") Softroll Black Tyre For Seat Height 330mm	○
BVA1440		125x27 mm (5" x 1") Medium Hard Black Tyre For Seat Height 350 mm	○
BVA1455		140x45 mm (5 1/2" x 1 3/4") Soft Black Tyre For Seat Height 380-450 mm	£21
BVA1475		150X30 mm (6" x 1 1/4") Medium Hard Black Tyre For Seat Height 380-450 mm	○
BVA1495		200x32 mm (8" x 1 1/4") Black Solid Tyre For Seat Height 450 - 500 mm	○
BVA1515		200x45 mm (8" x 1 3/4") Black Solid Tyre For Seat Height 450-500 mm	◆
<b>REAR WHEEL POSITION*</b>			
BVA1534		Standard	◆
<b>REAR WHEEL SIZE*</b>			
BVA1545		16" / 405 mm Transit rear wheel with fixed attachment	(NC 1584, 1588, 1605, 1640, 1650, 1660, 1710, 1725, 1727, 7777) ○
BVA1555		20" / 515 mm Self Propel rear wheel with QR attachment	(NC 1586, 1588, 1625) ○
BVA1565		22" / 560 mm Self Propel rear wheel with QR attachment	(NC 7777) ○
BVA1575		24" / 600 mm Self Propel rear wheel with QR attachment	(NC 7777) ○
<b>REAR WHEEL TYPE*</b>			
BVA1584		Standard self propel with colourless hub, clear anodized rim and Stainless Steel spokes	(NC 1728) ○
BVA1586		Hub Brake Self Propel or Transit 16" Transit with Black Plastic Hub and Rim 22" and 24" / 560 mm and 600 mm Self Propel with colourless hub, clear anodized rim and Stainless Steel spokes	○
BVA1588		One Arm Drive Black hub and Zinc Plated Spokes with OAD Handrim System	(NC 1605, 1650, 1660, 1710, 1728) ○
<b>REAR WHEEL TYRE*</b>			
BVA1605		Pneumatic Profile Black 1" 3/8 tyres	○
BVA1625		Solid Black Tyres 1"3/8	◆
BVA1640		Pneumatic low profile Black 1" tyres	○
<b>HANDRIM</b>			
BVA1650		Aluminium handrim	○
BVA1660		Plastic coated handrim	£48
<b>SPOKEGUARD</b>			
BVA1710		Frost spokeguard	£166
<b>BRAKE*</b>			
BVA1720		User manoeuvred tyre brakes, includes extended handles	◆
BVA1725		One Arm Brake LH	£126
BVA1727		One Arm Brake RH	£126
<b>ATTENDANT HUB BRAKE</b>			
BVA1728		Attendant Hub Brake	£140

<b>PELVIC BELT*</b>			
BVA1730		Pelvic belt	◆
BVA1731		Bodypoint 2-Point Centre-pull Push Button Medium (Mounted to chair)	£83
BVA1732		Bodypoint 2-Point Centre-pull Push Button Large (Mounted to chair)	£83
BVA1733		Bodypoint 4-Point Centre-pull Push Button Medium (Mounted to chair)	£98
BVA1734		Bodypoint 4-Point Centre-pull Push Button Large (Mounted to chair)	£98
<b>BODYPOINT HARNESS</b>			
BVA1736		Bodypoint Pivotfit Harness Medium (Supplied with fixing kit, unmounted)	£109
BVA1737		Bodypoint Pivotfit Harness Large (Supplied with fixing kit, unmounted)	£109
<b>OTHER OPTIONS</b>			
BVA1750		Table Tray without elbow pads includes mounting kit (NC 1765, 1790)	£138
BVA1755		Table Tray with elbow Pads includes mounting kit (NC 1765, 1790)	£138
BVA1756		Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included) (NC 1760, 1790)	£103
BVA1757		Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included) (NC 1760, 1790)	£103
BVA1770		Swing away table without elbow pads (NC 1765, 1780, 1790)	£152
BVA1775		Swing away table with elbow pads (NC 1765, 1780, 1790)	£152
BVA1760		Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray) (NC 1790)	£51
BVA1765		Half table tray padded cushion (NC 1790)	£46
BVA1780		Locking device for table tray (NC 1790)	£15
BVA1790		Mounting kit for table (if ordered separately)	£9
BVA1820		Pump	£22
BVA1830		Cane Holder	£49
BVA1840		Clinic Kit	£55
BVA1850		Harness Attachment Bracket Seat (Pair) on seat	£23
BVA1860		Harness Attachment Bracket Upper	£23
BVA1870		Kerb Lever (NC 1900)	£42
BVA1900		Accessories Tray to suit Ventilators etc. Positioned at the rear of the frame between the rear wheels	£269
BVA1930		Restrictor Clip kit, Plastic clips which attach to the gas struts to limit the seat tilt or backrest recline	£8
<b>TRANSPORTATION BRACKETS*</b>			
BVA1890		Transportation Brackets	◆
BVA7777		Without Transportation Brackets, only select with 20" Self Propel wheels	○

<b>TOTAL PRICE</b>	
--------------------	--



**Non Conforming Options:**

BVA0090	+ 100 mm pre-adjustment on armrests and legrests
BVA0130	350 mm/ 14"
BVA0140	380 mm/ 15"
BVA0150	400 mm/ 16"
BVA0160	430 mm/ 17"
BVA0180	480 mm/ 19"
BVA0190	500 mm/ 20"
BVA0200	Reduced Tilt: Positive Seat Angle 10°
BVA0210	Reduced Tilt: Positive Seat Angle 20°
BVA0220	Full Tilt: Positive Seat Angle 30°
BVA0225	Reduced Backrest Recline: Angle 10°
BVA0227	Reduced Backrest Recline: Angle 20°
BVA0230	Full Backrest Recline: Angle 30°
BVA0390	Matrx Elite E2 (includes hardware mounted to frame, backrest supplied in box)
BVA0400	Matrx Elite E2 Deep (includes hardware mounted to frame, backrest supplied in box)
BVA0455	Matrx Backrest Height 410 mm / 16"
BVA0456	Matrx Backrest Height 460 mm / 18"
BVA0457	Matrx Backrest Height 510 mm / 20"
BVA0461	Matrx Backrest Width 380mm / 15" (Chest width 290mm / 11 1/2")
BVA0462	Matrx Backrest Width 410mm / 16" (Chest width 315mm / 12 1/2")
BVA0463	Matrx Backrest Width 430mm / 17" (Chest width 340mm / 13")
BVA0464	Matrx Backrest Width 460mm / 18" (Chest width 370mm / 14 1/2")
BVA0465	Matrx Backrest Width 480mm / 19" (Chest width 390mm / 15")
BVA0466	Matrx Backrest Width 510mm / 20" (Chest width 420mm / 16 1/2")
BVA0460	+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required
BVA0470	+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required
BVA0490	Cover for tension adjustable upholstery
BVA0500	Cushion, Laguna (Lateral Support)
BVA0510	Cushion, Mistral 2 (Waist Support)
BVA0520	Cushion, Passad 2 (Shoulder Support)
BVA0570	Black Multi-stretch Black PU coated fabric
BVA0575	Black Spacer TR36
BVA2222	No upholstery, select with Matrx and Without backrest options
BVA0637	Electrical operated seat and backrest tilt
BVA0640	Assistant Operated gas piston for backrest recline
BVA3333	Only select if option 0637 Electrical operated is required
BVA0670	Central Legrest, includes calfpad and footrest. Also select options 0740, 0790, 4444, 5555
SWO0681	Standard position with 0-30° Legrest Abduction function
BVA0690	High position 50mm higher for leg length compensation
BVA0692	90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest
BVA0700	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH
BVA0710	Legrest Hanger fixed 80 deg, swing in/swing out LH
BVA0720	Legrest Hanger fixed 90 deg, swing in/swing out LH

BVA0730	Legrest Amputee LH
BVA0750	Elevating Legrest, angle adjustable 0 to 75°swing In/swing out with knee pad and calf support with locking clip RH
BVA0760	Legrest Hanger fixed 80 deg, swing in/swing out RH
BVA0770	Legrest Hanger fixed 90 deg, swing in/swing out RH
BVA0780	Legrest Amputee RH
BVA0805	Footrest one piece, depth and angle adjustable short tube
BVA0825	Footrest lower part, depth and angle adjustable with heel strap, LH
BVA0845	Footrest lower part, depth and angle adjustable with heel strap, RH
BVA0855	Short Tube 24CM
BVA0865	Long Tube 36CM
BVA0915	Armrest Extended Height, Ducktail pad LH, Adjusts from 255 mm - 355 mm with Autolock
BVA0920	Armrest Wide Pad LH, Height Adjustable from 200 mm - 320mm with autolock.
BVA0930	Armrest Hemiplegic LH, Height Adjustable from 200 mm - 320 mm with autolock
BVA0950	Without Armrest LH Select when choosing BVA1756 Half Tray LH, Narrow Armrest is included
BVA0975	Armrest Extended Height, Ducktail pad RH, Adjusts from 255 mm - 355 mm with Autolock
BVA0980	Armrest Wide Pad RH, Height Adjustable from 200 mm - 320mm with autolock.
BVA0990	Armrest Hemiplegic RH, Height Adjustable from 200 mm - 320 mm with autolock
BVA1010	Without Armrest RH Select when choosing BVA1757 Half Tray RH, Narrow Armrest is included
BVA1020	Pad for sideguard, removable with hook and loop fastening
BVA1040	Pommel Abduction Cushion
BVA1060	2 cm Increased seat height plate. Allows an increased 20mm of seat height by adding a foam plate to the seat board with hook and loop
BVA1080	Headrest, Sideway Adjustable
BVA1100	Headrest with cheek support, sideway adjustable
BVA1087	Matrx Elan hardware
BVA1088	Matrx Elan extended hardware, (Extra link)
BVA1089	Matrx Elan LOXX Hardware
BVA9999	Without Elan headrest hardware (Only select with option 1080, 1100 Headrest and without headrest)
BVA1090	Cover for headrest
BVA1110	Neckrest
BVA1130	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)
BVA1140	Trunk Support Multi Adjustable with standard arm LH
BVA1160	Trunk Support Multi Adjustable with extended height arm LH
BVA1240	Trunk Support Multi Adjustable Swingaway LH
BVA1150	Trunk Support Multi Adjustable with standard arm RH
BVA1170	Trunk Support Multi Adjustable with extended height arm RH
BVA1250	Trunk Support Multi Adjustable Swingaway RH
BVA1260	Trunk Support Attachment. If trunk supports are chosen, the attachment bracket is already included
BVA1280	Wedges (One Pair)
BVA1302	Calf pad cover, for LH elevating legrest
BVA1304	Knee pad for LH elevating legrest
BVA1312	Calf pad cover, for RH elevating legrest
BVA1314	Knee pad for RH elevating legrest
BVA8888	Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)

BVA1310	Footplate Extender LH
BVA1320	Footplate Extender RH
BVA1330	Calf Strap (for One piece footrest)
BVA1340	Cover Footplate, Black Plush Fabric
BVA1350	Calf Strap Padded
BVA1351	Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)
BVA1352	Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)
BVA1353	Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)
BVA1435	100x27 mm (4" x 1") Softroll Black Tyre For Seat Height 330mm
BVA1440	125x27 mm (5" x 1") Medium Hard Black Tyre For Seat Height 350 mm
BVA1455	140x45 mm (5 1/2" x 1 3/4") Soft Black Tyre For Seat Height 380-450 mm
BVA1475	150X30 mm (6" x 1 1/4") Medium Hard Black Tyre For Seat Height 380-450 mm
BVA1495	200x32 mm (8" x 1 1/4") Black Solid Tyre For Seat Height 450 - 500 mm
BVA1515	200x45 mm (8" x 1 3/4") Black Solid Tyre For Seat Height 450-500 mm
BVA1545	16" / 405 mm Transit rear wheel with fixed attachment
BVA1555	20" / 515 mm Self Propel rear wheel with QR attachment
BVA1565	22" / 560 mm Self Propel rear wheel with QR attachment
BVA1575	24" / 600 mm Self Propel rear wheel with QR attachment
BVA1584	Standard self propel with colourless hub, clear anodized rim and Stainless Steel spokes
BVA1586	Hub Brake Self Propel or Transit 16" Transit with Black Plastic Hub and Rim 22" and 24" / 560 mm and 600 mm Self Propel with colourless hub, clear anodized rim and Stainless Steel spokes
BVA1588	One Arm Drive Black hub and Zinc Plated Spokes with OAD Handrim System
BVA1605	Pneumatic Profile Black 1" 3/8 tyres
BVA1625	Solid Black Tyres 1"3/8
BVA1640	Pneumatic low profile Black 1" tyres
BVA1650	Aluminium handrim
BVA1660	Plastic coated handrim
BVA1710	Frost spokeguard
BVA1725	One Arm Brake LH
BVA1727	One Arm Brake RH
BVA1728	Attendant Hub Brake
BVA1750	Table Tray without elbow pads includes mounting kit
BVA1755	Table Tray with elbow Pads includes mounting kit
BVA1756	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)
BVA1757	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)
BVA1770	Swing away table without elbow pads
BVA1775	Swing away table with elbow pads
BVA1760	Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)
BVA1765	Half table tray padded cushion
BVA1780	Locking device for table tray
BVA1790	Mounting kit for table (if ordered separately)
BVA1850	Harness Attachment Bracket Seat (Pair) on seat
BVA1900	Accessories Tray to suit Ventilators etc. Positioned at the rear of the frame between the rear wheels
BVA7777	Without Transportation Brackets, only select with 20" Self Propel wheels





For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):