

# Dyna-Form<sup>®</sup> Maximus

Bariatric foam mattress



**Comfort and pressure redistribution**

for high-risk patients up to 413 kg



# High-risk bariatric pressure care



## Bariatric pressure care

- > Finest quality Combustion Modified Ether (CME) foam construction supports the weight of patients and prevents bottoming out, while delivering pressure redistribution.
- > Specifically designed for bariatric patients considered at High Risk of pressure ulcer development.
- > Weight limit of 413 kg



## Cover and hygiene

- > Multi-stretch, vapour permeable cover reduces shear and friction forces and helps control infection.
- > Cover fabric is impregnated with an anti-fungal agent and can be laundered up to 80 °C.
- > Hypochlorite compatible at 1,000 ppm (0.1%), to be used in line with local infection control guidelines.



## Durable and adaptable

- > Toughened, non-slip polyurethane base gives a robust surface against hydrolysis and strikethrough, and prevents slippage.
- > Supplied at 880 mm wide as standard, and can be made to the dimensions of any bariatric bed frame.
- > Supplied with a 4-year warranty.

## Technical data ▼



Height ▼

150 mm



Width ▼

880 mm\*



Length ▼

1980 mm\*



Max. User Weight ▼

413 kg



Total weight ▼

17 kg



Max. Cleaning Temperature ▼

80 °C

\* Specific sizes available on request.

Manufactured to: BS EN 597-1 Cig, BS EN 597-2 Match and BS 7177 (ignition Source 5)

Essential nursing care is pivotal in pressure injury management. This mattress may contribute to the outcome of a pressure injury prevention plan that includes education, clinical judgment and action planning based on vulnerability, which are fundamental factors in the prevention of pressure injuries.

Contact us:  
Invacare International GmbH  
Neuhofweg 51  
4147 Aesch BL, Switzerland  
globalcomms@dhcare.com  
[www.dhcare.com](http://www.dhcare.com)



Legal Manufacturer:  
DHG Ltd.  
Withey Court  
Western Industrial Estate  
Caerphilly, United Kingdom  
CF83 1BF



Invacare International GmbH is a member of the DHCare group of companies. Invacare International GmbH reserves the right to modify existing product models at any time. Product specifications may vary across regions.  
©2026 Invacare International GmbH. All rights reserved. Dyna-Form\_Maximus\_EU\_25092026