

EasyGlide

Transfer boards



Safe moving & handling

Effortless sitting transfers for greater independence



Feel safe and comfortable with EasyGlide

EasyGlide is a sliding board designed for safe, effortless sitting transfers to and from wheelchairs, beds, toilets and car seats.

Available in two models and multiple sizes, EasyGlide is made of strong polyethylene. The longer versions are ideal for wheelchair-to-car transfers. The flexible three-part model features foldable sides to protect clothing from wheelchair wheels.

Rounded corners and a smooth, low-friction surface make positioning and removal effortless. Double anti-slip tape on the underside keeps the board securely in place during transfer.



Easy to transport

Lightweight and compact – fits in a bag or wheelchair pocket.



Very convenient

Quick-guide on the underside for immediate use in any situation.



Easy to clean

All surfaces can be cleaned and disinfected with 70% ethanol or 45% isopropanol.



- > **EasyGlide 5020**
(600 × 330 mm)
Three-part, foldable
Foldable long sides protect clothing from wheelchair wheels during transfer.



- > **EasyGlide 5050**
(750 × 330 mm)
Extended length
Longer format ideal for wheelchair-to-car transfers – same foldable design.



- > **EasyGlide 5060**
(600 × 250 mm)
Slim & versatile
Narrower one-piece board for standard transfers between bed, chair and toilet.



- > **EasyGlide 5070**
(450 × 200 mm)
Compact format
The smallest model – easy to carry and ideal for short-distance transfers.

Technical data ▼



Length ▼



Width ▼



Height ▼

Unit ▼

EasyGlide 5020	600 mm	330 mm	5 mm	1
EasyGlide 5050	750 mm	330 mm	6 mm	1
EasyGlide 5060	600 mm	250 mm	5 mm	1
EasyGlide 5070	450 mm	200 mm	5 mm	1