

EasySlide

SystemRoMedic™



Safe and easy

Reduces friction for gentle, effortless transfers



Trusted since 1984

EasySlide is the original SystemRoMedic transfer aid – used in hospitals and nursing homes worldwide for over 40 years.

Its unique Ultra Low Friction (ULF) material on the inner surface enables caregivers to transfer users on virtually any surface with minimal effort and no lifting. Available in several models, materials and sizes, EasySlide suits

a wide range of transfer situations. Classic models are made of microfibre with polyester padding and ULF lining. Most sizes can be fitted with a replaceable or disposable protective cover for improved hygiene.

EasySlide 190 x 60

Designed for supine transfers to and from beds, x-ray tables, operating tables, shower trolleys and stretchers — and equally effective for turning and repositioning in bed.

> ULF sliding material

Ultra Low Friction inner layer reduces resistance on any surface.

> Seat belt included

Safely helps to stabilise the user during transfer.

> Soft padded surface

Gentle on skin — comfortable for users and caregivers alike.

> Protective cover available

Replaceable polyester cover with carbon fibre reduces static electricity.

> Disposable covers available

Single-use plastic covers for maximum hygiene and less laundry.



Ultra Low Friction

Unique ULF material reduces friction on any surface — enabling gentle transfers with minimal effort and no lifting.



Multiple models & sizes

Available in microfibre and disposable versions to suit any transfer situation or hygiene requirement.

EasySlide disposable

EasySlide disposable is a thin, pliable and versatile sliding mat on a roll — up to 60 individual mats per pack. Suitable for turning, repositioning, supine and seated transfers, leg lifting, dressing assistance and sling application. Each roll comes in a hygienic wall-mounted dispenser for easy access.

> Up to 60 mats per roll

Always fresh, clean and ready — supplied in practical two-packs.

> Tubular or open

Use as a tube or separate into a flat sheet.

> Perforated every 50 cm

Simply unroll and tear off the required length.

> Eco-friendly material

Tough, durable and not water-soluble — ideal for showering situations.

> Instructions on every mat

User guide printed directly on each sheet.

> Wall dispenser available

Hygienic storage within reach — also available as a steel tabletop dispenser.



EasySlide protective cover

Fits most EasySlide models. Made of polyester with carbon fibre to reduce static electricity. Easily washable and reusable.




> **Carbon fibre lining**

Reduces static electricity for smoother, safer transfers.

EasySlide disposable cover

Disposable plastic covers for EasySlide reduce laundry and improve infection control. Made of bio-degradable polyethylene.



> **Biodegradable polyethylene**

Hygienic, eco-conscious and low-maintenance.

DHCare Ltd.
4D/E, The Ark, 201 Talgarth Road,
Hammersmith, London, W6 8BJ, UK
www.dhcare.com

©2026 DHCare Ltd. All rights reserved. DHCare Ltd reserves the right to modify existing models at any time. Specifications may vary across regions. EasySlide Transfer Aids_EU_29052026

