



**Yes, you can:**

Tel: 01656 776222  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

# Power Retail Prescription Form

JANUARY 2025  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Esprit Action

### Standard Features

- Max User Weight : 125kg (19.7 Stone)
- Crash tested ISO 7176 -Part 19\*
- Double crossbar
- Seat widths 380mm to 605mm (15" - 24")
- Seat depths 400mm, 450mm or 500mm (16", 18" or 20")
- Seat heights 460mm or 485mm (18" or 19")
- Nylon seat upholstery
- Fixed 7° 430mm backrest
- Lumbar Support (+1 STRAP) BI Material Padded backrest
- 80° swing in / swing out hangers
- Flip up footplates including heel strap
- Removable full length armrests
- Choice of 4 standard colours (with 9 special colour options)
- Anti tippers
- Buckle belt
- 
- \* Battery MUST be removed when used as a seat in a vehicle



**Standard Base price: £3,706**

Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
SPE0060	<input type="checkbox"/> ESPRIT ACTION ADULT (DOUBLE CROSSBAR)	£3,706
<b>CHAIR TYPE*</b>		
SPE0100	<input type="checkbox"/> SPECIAL	POA
SPE0100	<input type="checkbox"/> BASE MODEL ESPRIT ACTION	◆
<b>CHASSIS*</b>		
SPE0210	<input type="checkbox"/> DOUBLE CROSSBAR (SW 380 MM TO 605 MM / SD FROM 400 MM TO 500 MM) USER MAX WEIGHT = 125Kg	◆

BRACKET TYPE*			
SPE0300	<input type="checkbox"/>	STANDARD SWING AWAY BRACKET	◆
SPE0310	<input type="checkbox"/>	ALBER E35 SWIVEL ARM	£360
USER REMOTE*			
SPE0400	<input type="checkbox"/>	USER REMOTE CONTROL RH	○
SPE0410	<input type="checkbox"/>	USER REMOTE CONTROL LH	○
DUAL REMOTE			
SPE0440	<input type="checkbox"/>	INTUITIVE DUAL REMOTE CONTROL RH	£853
SPE0450	<input type="checkbox"/>	INTUITIVE DUAL REMOTE CONTROL LH	£853
SPE0460	<input type="checkbox"/>	INTUITIVE SINGLE DUAL REMOTE CONTROL (ONE HANDLE) (ADJUSTABLE RH TO LH SIDE)	£608
SPEED*			
SPE0500	<input type="checkbox"/>	6 KM/H	◆
BATTERY*			
SPE0600	<input type="checkbox"/>	BATTERY 7,65AH/36V (NON SPILLABLE LEAD ACID BATTERY)	◆
POWER OPTIONS			
SPE0610	<input type="checkbox"/>	ADDITIONAL BATTERY 7,65AH/36V	£271
SPE0630	<input type="checkbox"/>	REMOTE PROTECTOR	£62
SEAT WIDTH*			
SPE1040	<input type="checkbox"/>	380 MM SEAT WIDTH (NC 1250)	○
SPE1050	<input type="checkbox"/>	405 MM SEAT WIDTH (NC 1250)	○
SPE1060	<input type="checkbox"/>	430 MM SEAT WIDTH	○
SPE1070	<input type="checkbox"/>	455 MM SEAT WIDTH	○
SPE1080	<input type="checkbox"/>	480 MM SEAT WIDTH	○
SPE1090	<input type="checkbox"/>	505 MM SEAT WIDTH	○
SPE1100	<input type="checkbox"/>	555 MM SEAT WIDTH (REINFORCED SEAT UPHOLSTERY) (NC 2000, 2010, 2030, 2200, 2210, 2220, 2230, 5020, 5040, 5080, 8200)	£103
SPE1110	<input type="checkbox"/>	605 MM SEAT WIDTH (REINFORCED SEAT UPHOLSTERY) (NC 2000, 2010, 2030, 2200, 2210, 2220, 2230, 5020, 5040, 5060, 5070, 5080, 8200)	£103
SEAT DEPTH*			
SPE1210	<input type="checkbox"/>	400 MM SEAT DEPTH (NC 5090)	○
SPE1220	<input type="checkbox"/>	400 TO 450 MM ADJUSTABLE SEAT DEPTH SET AT 430 MM	○
SPE1230	<input type="checkbox"/>	450 MM SEAT DEPTH (NC 5080, 5090)	○
SPE1240	<input type="checkbox"/>	450 TO 500 MM ADJUSTABLE SEAT DEPTH (NC 5080)	£30
SPE1250	<input type="checkbox"/>	500 MM SEAT DEPTH (NC 5080)	○
SEAT HEIGHT*			
SPE1300	<input type="checkbox"/>	460 MM SEAT HEIGHT	○
SPE1310	<input type="checkbox"/>	485 MM SEAT HEIGHT	○
SEAT UPHOLSTERY TYPE*			
SPE1410	<input type="checkbox"/>	NYLON WITH LOOP BANDS STRIPS SEAT UPHOLSTERY (DOUBLE CROSSBAR ONLY)	◆

BACKREST TYPE*			
SPE2000		FIXED 7° BACKREST 2 HEIGHTS 400/430 MM (DELIVERED 430 MM)	◆
SPE2010		FIXED 7° BACKREST 510 MM HEIGHT	○
SPE2020		FOLDING BACKREST 430 MM HEIGHT	£32
SPE2030		FOLDING BACKREST 510 MM HEIGHT	£32
SPE2040		ANGLE ADJUSTABLE AND FOLD DOWN (-15° TO +15°) FIXED 7° HEIGHT 425-450 MM / 17"-18" SET AT 425 MM / 17" (2 POSITIONS; 25 MM / 1" STEPS THROUGH UPHOLSTERY)	£84
SPE2050		RECLINING/FOLDING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	(NC 2210, 2220, 2230) £147
SPE2060		RECLINING (RATCHET) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	£124
SPE2070		RECLINING (RATCHET) 0/30° BACKREST 3 HEIGHTS 470/490/510 MM (DELIVERED 510 MM)	£124
SPE2080		RECLINING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	£51
SPE2090		RECLINING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 470/490/510 MM (DELIVERED 510 MM)	£64
BACKREST UPHOLSTERY TYPE*			
SPE2200		MATRIX ELITE E2 BACK WITH HARDWARE 405 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS	(NC 2310, 2320, 5140) £557
SPE2210		MATRIX ELITE E2 BACK WITH HARDWARE 455 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS	(NC 2310, 2320, 5140) £557
SPE2220		MATRIX ELITE DEEP E2 BACK WITH HARDWARE 405 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS	(NC 2310, 2320, 5140) £647
SPE2230		MATRIX ELITE DEEP E2 BACK WITH HARDWARE 455 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS	(NC 2310, 2320, 5140) £647
SPE2300		LUMBAR SUPPORT (+1 STRAP) BI MATERIAL PADDED COVER	◆
SPE2310		CONTOUR ADJUSTABLE (4+1 STRAPS) BI MATERIAL COVER	○
SPE2320		CONTOUR ADJUSTABLE (FULL STRAPS) BI MATERIAL COVER	○
HANGERS RH*			
SPE3010		80° LONG, UL FROM 340 TO 470MM	○
SPE3020		ANGLE ADJUSTABLE (0° TO 75°) RH SWING IN/SWING OUT HANGER, UL FROM 340 TO 500MM (INCLUDING CALFPAD)	(NC 5150) £88
SPE3030		AMPUTEE LEGREST RH	(NC 3210, 3220) £129
HANGERS LH*			
SPE3110		80° LONG, UL FROM 340 TO 470MM	○
SPE3120		ANGLE ADJUSTABLE (0° TO 75°) LH SWING IN/SWING OUT HANGER, UL FROM 340 TO 500MM (INCLUDING CALFPAD)	(NC 5150) £88
SPE3130		AMPUTEE LEGREST LH	(NC 3300, 3310) £129
FOOTPLATES RH*			
SPE3200		FLIP UP RH FOOTPLATE INCLUDING HEEL STRAP	○
SPE3210		FLIP UP ANGLE AND DEPTH ADJUSTABLE RH FOOTPLATE, INCLUDING HEEL STRAP	£42
SPE3220		FLIP UP ONE PIECE FOOTPLATE, DEPTH & ANGLE ADJUSTABLE	(NC 3310) £76
SPZ7000		Amputee legrest selected	(NC 3010, 3020) ○

FOOTPLATES LH*			
SPE3300		FLIP UP LH FOOTPLATE, INCLUDING HEEL STRAP	○
SPE3310		FLIP UP ANGLE AND DEPTH ADJUSTABLE LH FOOTPLATE, INCLUDING HEEL STRAP	£42
SPZ7001		Amputee legrest or footboard selected	○
ARMRESTS RH*			
SPE4000		REMOVABLE RH DESK, ADJ IN HEIGHT (235-280 for FSTF 460 & 210-280 for FSTF 485MM) DEPTH (+35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	○
SPE4010		REMOVABLE RH FULL LENGTH, ADJ IN HEIGHT (210 TO 280 MM) DEPTH (+35 MM BACKWARD) & WIDTH (+25MM OUTWARDS)	◆
SPE4020		REMOVABLE RH ARMPAD HEMI, BLACK PU TR26, MULTIPOSITION AUTOLOCK SAFETY DEVICE INCLUDED, ADJ IN HEIGHT (210 TO 345 MM) DEPTH (35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	£32
ARMRESTS LH*			
SPE4100		REMOVABLE LH DESK, ADJ IN HEIGHT (235-280 for FSTF 460 & 210-280 for FSTF 485MM) DEPTH (+35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	○
SPE4110		REMOVABLE LH FULL LENGTH, ADJ IN HEIGHT (210 TO 280 MM) DEPTH (35 MM BACKWARD) & WIDTH (+25MM OUTWARDS)	◆
SPE4120		REMOVABLE LH ARMPAD HEMI, BLACK PU TR26, MULTIPOSITION AUTOLOCK SAFETY DEVICE INCLUDED, ADJ IN HEIGHT (210 TO 345 MM) DEPTH (35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	£32
SEAT CUSHION*			
SPE5010		BASIC 50 MM FOAM CUSHION - 50 MM HEIGHT, NYLON & ANTI-SLIP COVER	○
SPE5020		MATRX FLO-TECH LITE CUSHION - 78 MM HEIGHT	£117
SPE5040		MATRX CONTOUR VISCO NG CUSHION - 100 MM HEIGHT MAXIMUM USER WEIGHT 100 KG	£164
SPE5060		MATRX VI CUSHION - 82 MM HEIGHT, ADDITIONAL DELIVERY TIME + 5 DAYS	£359
SPE5070		MATRX LIBRA CUSHION - 70 MM HEIGHT, ADDITIONAL DELIVERY TIME + 5 DAYS	£387
SEAT CUSHION DEPTH*			
SPE5080		FOR 400 MM SEAT DEPTH	○
SPE5085		FOR 450 MM SEAT DEPTH	○
SPE5090		FOR 500 MM SEAT DEPTH	○
SEATING OPTIONS			
SPE5100		NECKREST BLACK PU TR26	£127
SPE5110		HEADREST BLACK PU TR26	£153
SPE5120		HEADREST WITH CHEEK SUPPORT BLACK PU TR26 (INCLUDING BACKBRACE EXCEPT FOR LAGUNA)	£177
SPE5140		LATERAL POSITIONING PADS (PAIR)	£59
SPE5150		ADJUSTABLE CALF STRAP	○
BACK BRACE			
SPE5170		BACK BRACE FIXED	○

<b>COLOURS*</b>			
SPE6000	<input type="checkbox"/>	ANTHRACITE BLACK	○
SPE6015	<input type="checkbox"/>	DARK PURPLE	○
SPE6035	<input type="checkbox"/>	NAVY BLUE	○
SPE6040	<input type="checkbox"/>	SLATE GREY	○
SWO6925	<input type="checkbox"/>	WATER BLUE - 3 week lead time	○
SWO0806	<input type="checkbox"/>	Electric Green - 3 week lead time	£164
SWO6033	<input type="checkbox"/>	Red Matt - 3 week lead time	£164
SWO0607	<input type="checkbox"/>	Dark Red - 3 week lead time	£164
SWO0515	<input type="checkbox"/>	Pearl Grey - 3 week lead time	£164
SWO0647	<input type="checkbox"/>	Intense Yellow - 3 week lead time	£164
SWO0648	<input type="checkbox"/>	Vibrant Orange - 3 week lead time	£164
SWO0651	<input type="checkbox"/>	Polar Blue - 3 week lead time	£164
SWO0652	<input type="checkbox"/>	Candy Pink - 3 week lead time	£164
SPE6050	<input type="checkbox"/>	SPECIAL COLOUR (ON DEMAND ONLY)	POA
<b>ACCENT COLOUR*</b>			
SPE6100	<input type="checkbox"/>	BLACK	○
SPE6120	<input type="checkbox"/>	LIME	○
<b>FRONT WHEELS*</b>			
SPE7000	<input type="checkbox"/>	200 X 50 MM (8" X 2") BLACK SOLID TYRE	£47
SPE7010	<input type="checkbox"/>	200 X 50 MM (8" X 2") BLACK PNEUMATIC TYRE	◆
<b>REAR WHEELS*</b>			
SPE7100	<input type="checkbox"/>	355 X 70 MM (14" X 27") BLACK SOLID WITH PUNCTURE PROOF TYRE	£82
SPE7110	<input type="checkbox"/>	355 X 70 MM (14" X 27") BLACK PNEUMATIC TYRE	◆
<b>ANTI TIPPER*</b>			
SPE8000	<input type="checkbox"/>	REAR ANTI-TIPPERS, FIXED, PAIR	◆
SPE8010	<input type="checkbox"/>	REAR ANTI-TIPPERS REMOVABLE WITH JACK UP, PAIR	£276
<b>LAP BELT*</b>			
SPE8020	<input type="checkbox"/>	POSTURE BELT	◆
<b>TRANSPORTATION BRACKET*</b>			
SPE8030	<input type="checkbox"/>	TRANSPORTATION BRACKET	◆
<b>OTHER OPTIONS</b>			
SPE8100	<input type="checkbox"/>	CANE HOLDER	£20
SPE8200	<input type="checkbox"/>	ACRYLIC TRAY FLIP UP	£308

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

SPE1250	500 MM SEAT DEPTH
SPE2000	FIXED 7° BACKREST 2 HEIGHTS 400/430 MM (DELIVERED 430 MM)
SPE2010	FIXED 7° BACKREST 510 MM HEIGHT
SPE2030	FOLDING BACKREST 510 MM HEIGHT
SPE2200	MATRX ELITE E2 BACK WITH HARDWARE 405 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS
SPE2210	MATRX ELITE E2 BACK WITH HARDWARE 455 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS
SPE2220	MATRX ELITE DEEP E2 BACK WITH HARDWARE 405 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS
SPE2230	MATRX ELITE DEEP E2 BACK WITH HARDWARE 455 MM BACK HEIGHT ADDITIONAL DELIVERY TIME + 5 DAYS
SPE2310	CONTOUR ADJUSTABLE (4+1 STRAPS) BI MATERIAL COVER
SPE2320	CONTOUR ADJUSTABLE (FULL STRAPS) BI MATERIAL COVER
SPE3010	80° LONG, UL FROM 340 TO 470MM
SPE3020	ANGLE ADJUSTABLE (0° TO 75°) RH SWING IN/SWING OUT HANGER, UL FROM 340 TO 500MM (INCLUDING CALFPAD)
SPE3210	FLIP UP ANGLE AND DEPTH ADJUSTABLE RH FOOTPLATE, INCLUDING HEEL STRAP
SPE3220	FLIP UP ONE PIECE FOOTPLATE, DEPTH & ANGLE ADJUSTABLE
SPE3300	FLIP UP LH FOOTPLATE, INCLUDING HEEL STRAP
SPE3310	FLIP UP ANGLE AND DEPTH ADJUSTABLE LH FOOTPLATE, INCLUDING HEEL STRAP
SPE5020	MATRX FLO-TECH LITE CUSHION - 78 MM HEIGHT
SPE5040	MATRX CONTOUR VISCO NG CUSHION - 100 MM HEIGHT MAXIMUM USER WEIGHT 100 KG
SPE5060	MATRX VI CUSHION - 82 MM HEIGHT, ADDITIONAL DELIVERY TIME + 5 DAYS
SPE5070	MATRX LIBRA CUSHION - 70 MM HEIGHT, ADDITIONAL DELIVERY TIME + 5 DAYS
SPE5080	FOR 400 MM SEAT DEPTH
SPE5090	FOR 500 MM SEAT DEPTH
SPE5140	LATERAL POSITIONING PADS (PAIR)
SPE5150	ADJUSTABLE CALF STRAP
SPE8200	ACRYLIC TRAY FLIP UP



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

**Special Instructions**

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_

Description of special request (please include relevant option code or part number):  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_