



Yes, you can:

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Power Invacare Ireland Prescription Form

07 JANUARY 2026
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Esprit Action Junior

Standard Features

- Max User Weight 80kg (12.6 Stone)
- Crash tested ISO7176 - 19*
- Single crossbar (Width Adjustable)
- Seat widths 305mm to 380mm 12" - 15")
- Seat depth 325mm (13") to 375mm (15") adjustable (set at 350mm)
- Seat height 460mm (18")
- Nylon seat upholstery
- FIXED 7° BACKREST 430mm backrest
- Standard nylon sling type backrest upholstery
- 80° swing in / swing out legrest hangers
- Flip up angle and depth adjustable footplates
- Height adjustable, removable arm armrests
- Choice of 3 standard colours (with 10 Special colour options)
-
- * Battery MUST be removed when used as a seat in a vehicle

Standard Base price: €4,959



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
SPD0050	ESPRIT ACTION JUNIOR (SINGLE CROSSBAR) (NC 2300)	€4,959
CHAIR TYPE*		
SPDSPEC	SPECIAL	POA
SPD0100	BASE MODEL ESPRIT ACTION	◆
CHASSIS*		
SPD0200	SINGLE CROSSBAR (SW 305 MM TO 380 MM / SD 325 MM TO 375 MM) USER MAX WEIGHT = 80Kg	◆

BRACKET TYPE*			
SPD0300	<input type="checkbox"/>	STANDARD SWING AWAY BRACKET	○
SPD0310	<input type="checkbox"/>	ALBER E35 SWIVEL ARM	€473
USER REMOTE*			
SPD0400	<input type="checkbox"/>	USER REMOTE CONTROL RH	○
SPD0410	<input type="checkbox"/>	USER REMOTE CONTROL LH	○
DUAL REMOTE			
SPD0440	<input type="checkbox"/>	INTUITIVE DUAL REMOTE CONTROL RH	€1,119
SPD0450	<input type="checkbox"/>	INTUITIVE DUAL REMOTE CONTROL LH	€1,119
SPD0460	<input type="checkbox"/>	INTUITIVE SINGLE DUAL REMOTE CONTROL (ONE HANDLE) (ADJUSTABLE RH TO LH SIDE)	€798
SPEED*			
SPD0500	<input type="checkbox"/>	6 KM/H	◆
BATTERY*			
SPD0600	<input type="checkbox"/>	BATTERY 7,65AH/36V (NON SPILLABLE LEAD ACID BATTERY)	◆
POWER OPTIONS			
SPD0610	<input type="checkbox"/>	ADDITIONAL BATTERY 7,65AH/36V	€356
SPD0630	<input type="checkbox"/>	REMOTE PROTECTOR	€82
SEAT WIDTH*			
SPD1000	<input type="checkbox"/>	305 MM SEAT WIDTH	○
SPD1010	<input type="checkbox"/>	330 MM SEAT WIDTH	○
SPD1020	<input type="checkbox"/>	355 MM SEAT WIDTH	○
SPD1030	<input type="checkbox"/>	380 MM SEAT WIDTH	○
SEAT DEPTH*			
SPD1200	<input type="checkbox"/>	325 TO 375 MM ADJUSTABLE SEAT DEPTH SET AT 350 MM	○
SEAT HEIGHT*			
SPD1300	<input type="checkbox"/>	460 MM SEAT HEIGHT	○
SPD1310	<input type="checkbox"/>	485 MM SEAT HEIGHT	○
SEAT UPHOLSTERY TYPE*			
SPD1400	<input type="checkbox"/>	NYLON WITH LOOP BANDS STRIPS SEAT UPHOLSTERY, WIDTH & DEPTH ADJUSTABLE (SINGLE CROSSBAR ONLY)	◆
BACKREST TYPE*			
SPD2000	<input type="checkbox"/>	FIXED 7° BACKREST 2 HEIGHTS 400/430 MM (DELIVERED 430 MM)	◆
SPD2020	<input type="checkbox"/>	FOLDING BACKREST 430 MM HEIGHT	€42
SPD2040	<input type="checkbox"/>	ANGLE ADJUSTABLE AND FOLD DOWN (-15° TO + 15°) FIXED 7° HEIGHT 425-450 MM / 17"-18" SET AT 425 MM / 17" (2 POSITIONS; 25 MM / 1" STEPS THROUGH UPHOLSTERY)	€127
SPD2050	<input type="checkbox"/>	RECLINING/FOLDING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	€84
SPD2060	<input type="checkbox"/>	RECLINING (RATCHET) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	€163
SPD2080	<input type="checkbox"/>	RECLINING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	€67
BACKREST UPHOLSTERY TYPE*			
SPD2300	<input type="checkbox"/>	LUMBAR SUPPORT (+1 STRAP) BI MATERIAL PADDED COVER	◆
SPD2320	<input type="checkbox"/>	CONTOUR ADJUSTABLE (FULL STRAPS) BI MATERIAL COVER	○

HANGERS RH*			
SPD3000	<input type="checkbox"/>	80° SHORT, UL FROM 200 TO 330MM	◆
HANGERS LH*			
SPD3100	<input type="checkbox"/>	80° SHORT, UL FROM 200 TO 330MM	◆
FOOTPLATES RH*			
SPD3200	<input type="checkbox"/>	FLIP UP RH FOOTPLATE INCLUDING HEEL STRAP	○
SPD3210	<input type="checkbox"/>	FLIP UP ANGLE AND DEPTH ADJUSTABLE RH FOOTPLATE, INCLUDING HEEL STRAP	€56
SPD3220	<input type="checkbox"/>	FLIP UP ONE PIECE FOOTPLATE, DEPTH & ANGLE ADJUSTABLE (NC 3310)	€100
FOOTPLATES LH*			
SPD3300	<input type="checkbox"/>	FLIP UP LH FOOTPLATE, INCLUDING HEEL STRAP	○
SPD3310	<input type="checkbox"/>	FLIP UP ANGLE AND DEPTH ADJUSTABLE LH FOOTPLATE, INCLUDING HEEL STRAP	€56
SPD9925	<input type="checkbox"/>	No LH footplate required due to FLIP UP ONE PIECE BEING SELECTED (NC 3210)	○
ARMRESTS RH*			
SPD4000	<input type="checkbox"/>	REMOVABLE RH DESK, ADJ IN HEIGHT (235-280 for FSTF 460 & 210-280 for FSTF 485MM) DEPTH (+35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	○
SPD4010	<input type="checkbox"/>	REMOVABLE RH FULL LENGTH, ADJ IN HEIGHT (210 TO 280 MM) DEPTH (+35 MM BACKWARD) & WIDTH (+25MM OUTWARDS)	◆
ARMRESTS LH*			
SPD4100	<input type="checkbox"/>	REMOVABLE LH DESK, ADJ IN HEIGHT (235-280 for FSTF 460 & 210-280 for FSTF 485MM) DEPTH (+35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	○
SPD4110	<input type="checkbox"/>	REMOVABLE LH FULL LENGTH, ADJ IN HEIGHT (210 TO 280 MM) DEPTH (35 MM BACKWARD) & WIDTH (+25MM OUTWARDS)	◆
SEAT CUSHION*			
SPD5000	<input type="checkbox"/>	ESSENTIAL CUSHION 375 MM DEPTH	○
SEATING OPTIONS			
SPD5110	<input type="checkbox"/>	HEADREST BLACK PU TR26	€202
SPD5120	<input type="checkbox"/>	HEADREST WITH CHEEK SUPPORT BLACK PU TR26 (INCLUDING BACKBRACE EXCEPT FOR LAGUNA)	€233
SPD5140	<input type="checkbox"/>	LATERAL POSITIONING PADS (PAIR)	€78
SPD5150	<input type="checkbox"/>	ADJUSTABLE CALF STRAP	○
BACK BRACE			
SPD5160	<input type="checkbox"/>	BACK BRACE WIDTH ADJUSTABLE	○

COLOURS*			
SPD6000	<input type="checkbox"/>	ANTHRACITE BLACK	○
SPD6015	<input type="checkbox"/>	DARK PURPLE	○
SPD6020	<input type="checkbox"/>	WATER BLUE	○
SWO6529	<input type="checkbox"/>	POP BLUE - 3 week lead time	€215
SWO6259	<input type="checkbox"/>	SLATE GREY - 3 week lead time	€215
SWO0806	<input type="checkbox"/>	Electric Green - 3 week lead time	€215
SWO6033	<input type="checkbox"/>	Red Matt - 3 week lead time	€215
SWO0607	<input type="checkbox"/>	Dark Red - 3 week lead time	€215
SWO0515	<input type="checkbox"/>	Pearl Grey - 3 week lead time	€215
SWO0647	<input type="checkbox"/>	Intense Yellow - 3 week lead time	€215
SWO0648	<input type="checkbox"/>	Vibrant Orange - 3 week lead time	€215
SWO0651	<input type="checkbox"/>	Polar Blue - 3 week lead time	€215
SWO0652	<input type="checkbox"/>	Candy Pink - 3 week lead time	€215
SPD6050	<input type="checkbox"/>	SPECIAL COLOUR (ON DEMAND ONLY)	POA
ACCENT COLOUR*			
SPD6100	<input type="checkbox"/>	BLACK	○
SPD6120	<input type="checkbox"/>	LIME	○
FRONT WHEELS*			
SPD7000	<input type="checkbox"/>	200 X 50 MM (8" X 2") BLACK SOLID TYRE	€62
SPD7010	<input type="checkbox"/>	200 X 50 MM (8" X 2") BLACK PNEUMATIC TYRE	◆
REAR WHEELS*			
SPD7100	<input type="checkbox"/>	355 X 70 MM (14" X 27") BLACK SOLID WITH PUNCTURE PROOF TYRE	€108
SPD7110	<input type="checkbox"/>	355 X 70 MM (14" X 27") BLACK PNEUMATIC TYRE	◆
ANTI TIPPER*			
SPD8000	<input type="checkbox"/>	REAR ANTI-TIPPERS, FIXED, PAIR	◆
SPD8010	<input type="checkbox"/>	REAR ANTI-TIPPERS REMOVABLE WITH CRUTCHES (E25), PAIR	€362
LAP BELT*			
SPD8020	<input type="checkbox"/>	POSTURE BELT	◆
TRANSPORTATION BRACKET*			
SPD8030	<input type="checkbox"/>	TRANSPORTATION BRACKET	◆
OTHER OPTIONS			
SPD8100	<input type="checkbox"/>	CANE HOLDER	€26
SPD8200	<input type="checkbox"/>	ACRYLIC TRAY FLIP UP	€403

TOTAL PRICE	
--------------------	--

Non Conforming Options:

SPD2300 LUMBAR SUPPORT (+1 STRAP) BI MATERIAL PADDED COVER

SPD3210 FLIP UP ANGLE AND DEPTH ADJUSTABLE RH FOOTPLATE, INCLUDING HEEL STRAP

SPD3310 FLIP UP ANGLE AND DEPTH ADJUSTABLE LH FOOTPLATE, INCLUDING HEEL STRAP



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):