


Essential Care

Quality castellated foam mattress

A comfortable pressure-redistributing support surface for people at risk of developing pressure injuries.




A cost-effective solution

For comfort and pressure redistribution



Yes, you can.®

Affordable quality mattress



Pressure Redistribution

- › Castellations help redistribute weight over a greater area.
- › Durable castellated foam for immersion and envelopment to provide optimal pressure redistribution.
- › Excellent recovery properties ensuring good pressure redistribution.



Shear & Friction Reduction

- › Horizontal castellations help reducing shear forces when moving up and down the bed.
- › Multi-stretch, water resistant, polyurethane fabric cover helps further reducing shear and friction forces.
- › Smooth cover supports in the reduction of frictional forces.



Comfort

- › Castellations designed to improve comfort and support.
- › Three-sided zipped cover with white inner liner for easy strikethrough inspection and removal for laundering.
- › Suitable for multiple care environments.

Technical data ▼

For more information about this product, including the user manual, please visit our website: www.invacare.co.uk



Sizes available ▼



Max. User Weight ▼



Max. Cleaning Temperature ▼

| | Sizes available ▼ | Max. User Weight ▼ | Max. Cleaning Temperature ▼ |
|----------------------|---|--------------------|-----------------------------|
| 120 mm height | 1780 × 780 mm 1950 × 880 mm | 110 kg | 95° |
| 150 mm height | 1950 × 880 mm 2000 × 880 mm 2050 × 880 mm | 140 kg | 95° |

Manufactured to: BS EN 597-1 Cig, BS EN 597-2 Match and BS 7177 (ignition Source 5)

Essential nursing care is pivotal in pressure injury management. This mattress may positively contribute to the outcome of a pressure injury prevention plan that includes education, clinical judgment and action planning based on vulnerability, which are fundamental factors in the prevention of pressure injuries. A range of assessment scales can be used as a formal method of assessing risk from pressure injury development, and should be used in conjunction with an informal assessment (informed nursing judgement) as part of a comprehensive and individualised plan.

T: +44 1656 776222
uk@invacare.com
www.invacare.co.uk
www.invacare.ie

Invacare UK & Ireland
Unit 4 Pencoed Technology
Park
Bridgend CF35 5AQ



Yes, you can.®



© 2024 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. Essential Care - UK - 11/2024