



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

MARCH 2021
Customer Ref:

Customer Ref\*, Contact No\*, Delivery Address\*, Account No\*, Order No\*, Order Qty\*, Order date\*

Fox Modulite

Standard Features

- Modulite highly adjustable seating system
Telescopic adjustable seat widths from 330mm - 530mm (13" - 21")
Seat depths 410mm - 510mm stepless adjustment (16" - 20")
Tension adjustable backrest
Standard T-armrest, detachable with short armpads
Manual settable seat Tilt (from 0° to 8°)
Manual settable backrest Recline (0° - 30°, in 3.75° steps)
Standard 80° legrest hangers with fixed flip-up footplates
LiNX control REM110 with standard bracket
Choice of 5 colour options
Solid Grey castors 200mm x 50mm (8" x 2"), Silver Rim
Grey drive wheels, puncture proof solid inlay, (12 1/2" x 2 1/4"), Silver Rim (aluminium)
6 km/h (4mph) motors with 40Ah Gel battery
8A (200 - 250V) charger
4-point tie down system
Lapbelt with autolock (adjustable both sides)
Wheel lock / parking brake
Max User Weight 127Kg (20 stone)



Standard Base price: £3,135

Please make your selection with a cross in the box [ ] Free of charge [O] standard feature [◆] POA special request

Customer weight = KG

Table with columns: Code, Description, Price. Rows include BASE PRODUCT\* (AXB2100 FOX MODULITE £3,135), CHAIR TYPE\* (AXB0000 Fox Modulite ◆, AXZ0000 Special Fox Modulite POA), SEATING TYPE\* (AXB2100 Modulite seating unit ◆)

SEAT FRAME AND SEAT PLATE*			
AXB2111	<input type="checkbox"/>	Telescopic seat frame	○
SEAT WIDTH*			
AXZ1001	<input type="checkbox"/>	330mm (13")	○
AXB2121	<input type="checkbox"/>	380mm (15")	(NC 2701, 2703, 2705) ○
AXB2122	<input type="checkbox"/>	430mm (17")	(NC 2701, 2703, 2705) ○
AXB2123	<input type="checkbox"/>	480mm (19")	(NC 2702, 2704, 2706) ○
AXB2124	<input type="checkbox"/>	530mm (21")	(NC 2174, 2175, 2176, 2177, 2178, 2702, 2704, 2706) ○
SEAT DEPTH*			
AXB2132	<input type="checkbox"/>	410mm (16")	(NC 2196, 2197, 2701, 2702, 2703, 2704, 2705, 2706) ○
AXB2134	<input type="checkbox"/>	460mm (18")	○
AXB2136	<input type="checkbox"/>	510mm (20")	○
SEAT OPTION			
AXZSWO	<input type="checkbox"/>	H Frame interface kit	£444
SEAT CUSHION*			
AXB2142	<input type="checkbox"/>	Matrx Flo-tech Lite	£112
AXB2147	<input type="checkbox"/>	Matrx Libra	(NC 2155) £389
AXB2140	<input type="checkbox"/>	Without seat cushion	◆
BACKREST SUPPORT*			
AXB2163	<input type="checkbox"/>	Short 480mm	(NC 2172, 2173) ○
AXB2162	<input type="checkbox"/>	Tall 540mm	◆
PUSH HANDLES FOR BACKPOSTS*			
AXB2107	<input type="checkbox"/>	WITH push handles	◆
AXB2108	<input type="checkbox"/>	WITHOUT push handles	(NC 1495, 2219, 2220, 2227, 2228) £36
BACKREST UPHOLSTERY / CUSHION*			
AXB2171	<input type="checkbox"/>	Tension adjustable backrest	(NC 2152, 1496, 1445, 1446, 2155) ◆
AXB2179	<input type="checkbox"/>	Flex3 Tension adjustable	(NC 2151, 2199, 1446, 2155) £138
AXB2172	<input type="checkbox"/>	Flex3 backrest without cushion	(NC 2163, 2128) £223
AXB2173	<input type="checkbox"/>	Flex3 backrest with Laguna Backrest	(NC 2163, 2128) £362
AXB2174	<input type="checkbox"/>	Matrx PB 400mm	(NC 1445, 1446, 1496, 2124) £493
AXB2175	<input type="checkbox"/>	Matrx PB 500mm	(NC 1445, 1446, 1496, 2124) £493
AXB2176	<input type="checkbox"/>	Matrx Elite	(NC 1445, 1446, 1496, 2124) £493
AXB2177	<input type="checkbox"/>	Matrx Elite TR	(NC 1445, 1446, 1496, 2124) £493
AXB2178	<input type="checkbox"/>	Matrx Elite Deep	(NC 1445, 1446, 1496, 2124) £548
AXB2170	<input type="checkbox"/>	Without backrest upholstery	(NC 1445, 1446, 1496, 2322, 2323, 2324, 2325, 2326, 2327) ○
UPHOLSTERY TYPE*			
AXB2151	<input type="checkbox"/>	Black fabric	(NC 2179, 2170) ◆
AXB2152	<input type="checkbox"/>	Black multi stretch PU coated fabric	(NC 2171, 2170) ○
AXB2155	<input type="checkbox"/>	Black Spacer fabric	(NC 2147, 2171, 2179, 2170) ○
AXZ1002	<input type="checkbox"/>	Without	(NC 2171, 2179, 2172, 2173, 2174, 2175, 2176, 2177, 2178) ○
ARMREST*			
AXB2191	<input type="checkbox"/>	T-Armrest standard	◆
AXB2181	<input type="checkbox"/>	T-Armrest extended height	£36
AXB2183	<input type="checkbox"/>	In-rigged armrests	£25
ARMPADS*			
AXB2195	<input type="checkbox"/>	Short 325mm	◆
AXB2196	<input type="checkbox"/>	Long 375mm	(NC 2132) ○
AXB2197	<input type="checkbox"/>	Long ergonomic 397mm	(NC 2132) ○
AXB2199	<input type="checkbox"/>	Soft	(NC 2179) £86

<b>HEADREST / NECKREST - Adequate headrest option must be selected*</b>			
AXB1350	<input type="checkbox"/>	Neckrest	£72
AXB1432	<input type="checkbox"/>	Headrest	£203
AXB1427	<input type="checkbox"/>	Headrest With cheek support	£203
AXB2323	<input type="checkbox"/>	Elan standard pad 150mm (6") (NC 2170)	£293
AXB2324	<input type="checkbox"/>	Elan standard pad 255mm (10") (NC 2170)	£293
AXB2322	<input type="checkbox"/>	Elan standard pad 355mm (14" ) (NC 2170)	£293
AXB2325	<input type="checkbox"/>	Elan Occipital pad 230mm x 180mm (9" X 7") (NC 2170)	£293
AXB2326	<input type="checkbox"/>	Elan Occipital pad 305mm x 200mm (12" X 8") (NC 2170)	£293
AXB2327	<input type="checkbox"/>	Elan 4 point headrest pad 280mm x 255mm (11" X 10") (NC 2170)	£394
AXB1496	<input type="checkbox"/>	All bracketry as per option 1432 (NC 2170, 2171, 2174, 2175, 2176, 2177, 2178)	○
AXB1495	<input type="checkbox"/>	Push bar 22mm (NC 2108)	○
<b>LATERAL SUPPORT</b>			
AXB1446	<input type="checkbox"/>	Rigid pad (NC 2170, 2171, 2174, 2175, 2176, 2177, 2179)	£149
AXB1445	<input type="checkbox"/>	Multi-adjustable pad (NC 2170, 2171, 2174, 2175, 2176, 2177)	£149
<b>SEAT HEIGHT SETTING*</b>			
AXB0312	<input type="checkbox"/>	Low Seat to Floor height 480mm	○
AXB0344	<input type="checkbox"/>	High Seat to Floor height 510mm	○
<b>SEAT TILT*</b>			
AXB1002	<input type="checkbox"/>	Manually settable	◆
<b>BACKREST RECLINER*</b>			
AXB0975	<input type="checkbox"/>	Manual settable	◆
AXB0977	<input type="checkbox"/>	Powered	£421
<b>FOOTREST HANGER TYPE*</b>			
AXB0526	<input type="checkbox"/>	Standard 80° hangers	◆
AXB0529	<input type="checkbox"/>	Manual elevating legrest	£189
AXB0501	<input type="checkbox"/>	Without hanger (NC 0512)	○
<b>AMPUTEES SUPPORT</b>			
AXB0554	<input type="checkbox"/>	Amputee Support RH	£173
AXB0555	<input type="checkbox"/>	Amputee Support LH	£173
<b>FOOTPLATE TYPE*</b>			
AXB0512	<input type="checkbox"/>	Fixed Footplate (NC 0501)	◆
AXB0513	<input type="checkbox"/>	Angle and depth adjustable Footplate (NC 0501)	£79
ARZ1003	<input type="checkbox"/>	Without Footplates (NC 0526, 0529)	○
<b>CONTROLS*</b>			
AXB2201	<input type="checkbox"/>	REM110 (NC 0977)	◆
AXB2203	<input type="checkbox"/>	REM211	○
<b>CONTROL BRACKETS*</b>			
AXB1207	<input type="checkbox"/>	Standard Bracket RH (NC 1404)	○
AXB1208	<input type="checkbox"/>	Standard Bracket LH (NC 1405)	○
AXB1707	<input type="checkbox"/>	Swing away RH (NC 1404)	£118
AXB1708	<input type="checkbox"/>	Swing away LH (NC 1405)	£118
AXB2227	<input type="checkbox"/>	RH Bracket mounted for Attendant (NC 2108, 2219, 2220, 2128)	£56
AXB2228	<input type="checkbox"/>	LH Bracket mounted for Attendant (NC 2108, 2219, 2220, 2128)	£56
<b>DUAL CONTROL</b>			
AXB2220	<input type="checkbox"/>	Dual control REM050 (NC 2108, 2227, 2228)	£353
AXB2219	<input type="checkbox"/>	Compact (ACU200) (NC 2108, 2227, 2228)	£353

<b>DUAL CONTROL BRACKET</b>				
AXB1207	<input type="checkbox"/>	Standard Bracket RH	(NC 2227, 2228)	<input type="radio"/>
AXB1208	<input type="checkbox"/>	Standard Bracket LH	(NC 2227, 2228)	<input type="radio"/>
<b>BODYPOINT JOYSTICK KNOBS</b>				
AXB2721	<input type="checkbox"/>	3" U-shaped, joystick handle		£65
AXB2722	<input type="checkbox"/>	4" U-shaped, joystick handle		£65
AXB2723	<input type="checkbox"/>	Domed rubber, joystick handle		£65
AXB2724	<input type="checkbox"/>	Mushroom grooved, joystick handle		£65
<b>WHEEL BASE SETTING*</b>				
AXB1226	<input type="checkbox"/>	Recommended Factory Setting		◆
<b>SHROUD COLOUR*</b>				
AXB0201	<input type="checkbox"/>	Pearl Night Blue		<input type="radio"/>
AXB0233	<input type="checkbox"/>	Metallic Silver		<input type="radio"/>
AXB0228	<input type="checkbox"/>	Silky Matt Black		<input type="radio"/>
AXB0204	<input type="checkbox"/>	Transparent Cherry Red		<input type="radio"/>
AXB0222	<input type="checkbox"/>	White		<input type="radio"/>
<b>CASTORS*</b>				
AXB0802	<input type="checkbox"/>	200mm x 50mm (8" x 2") PU tyre Grey		<input type="radio"/>
AXB0803	<input type="checkbox"/>	200mm x 50mm (8" x 2") Solid inlay Grey		◆
<b>DRIVE WHEELS*</b>				
AXB0709	<input type="checkbox"/>	Puncture Proof Solid inlay (12 1/2 x 2 1/4") Grey		<input type="radio"/>
AXB0710	<input type="checkbox"/>	Puncture Proof Solid inlay (12 1/2 x 2 1/4") Grey		◆
<b>MOTOR TYPE*</b>				
AXB1294	<input type="checkbox"/>	3 km/h (2mph)		<input type="radio"/>
AXB1282	<input type="checkbox"/>	6 km/h (4mph)		◆
<b>BATTERY*</b>				
AXB1255	<input type="checkbox"/>	40 Ah Gel		◆
AXB1248	<input type="checkbox"/>	50 Ah AGM		£33
<b>CHARGER*</b>				
AXB1343	<input type="checkbox"/>	OptiCharge 8A (200 - 250V)		◆
AXB1313	<input type="checkbox"/>	OptiCharge 10A (200 - 250V)		£136
<b>WHEEL LOCK / PARKING BREAK*</b>				
AXB1206	<input type="checkbox"/>	Manual wheel locks		◆
<b>KERB RIDER</b>				
AXB1244	<input type="checkbox"/>	Kerbrider		£189
<b>TRANSPORTATION KIT*</b>				
AXB1494	<input type="checkbox"/>	4-Point tie-down		◆
<b>FIXATION BELT*</b>				
AXB1329	<input type="checkbox"/>	Lap Belt with auto-lock	(NC 2701, 2702, 2703, 2704, 2705, 2706)	◆
AXB2701	<input type="checkbox"/>	Evoflex pelvic stabilizer, size L	(NC 2121, 2122, 2132, 1329)	£84
AXB2702	<input type="checkbox"/>	Evoflex pelvic stabilizer, size M	(NC 2123, 2124, 2132, 1329)	£84
AXB2703	<input type="checkbox"/>	2-Point hip belt centre pull, size L	(NC 2121, 2122, 2132, 1329)	£77
AXB2704	<input type="checkbox"/>	2-Point hip belt centre pull, size M	(NC 2123, 2124, 2132, 1329)	£77
AXB2705	<input type="checkbox"/>	2-Point hip belt dual pull, size L	(NC 2121, 2122, 2132, 1329)	£77
AXB2706	<input type="checkbox"/>	2-Point Hip Belt Dual pull, size M	(NC 2123, 2124, 2132, 1329)	£77

**OTHER OPTIONS**

AXB2128		Bag hooks	(NC 2172, 2173, 2227, 2228)	£69
AXB1312		RH		£38
AXB1363		LH		£38
AXB1404		Tray, Swing Away RH		£196
AXB1405		Tray, Swing Away LH		£196
AXB1490		Castor lock levers		£29

**TOTAL PRICE**

**Non Conforming Options:**

AXB2121	380mm (15")
AXB2122	430mm (17")
AXB2123	480mm (19")
AXB2124	530mm (21")
AXB2132	410mm (16")
AXB2147	Matrx Libra
AXB2163	Short 480mm
AXB2108	WITHOUT push handles
AXB2171	Tension adjustable backrest
AXB2179	Flex3 Tension adjustable
AXB2172	Flex3 backrest without cushion
AXB2173	Flex3 backrest with Laguna Backrest
AXB2174	Matrx PB 400mm
AXB2175	Matrx PB 500mm
AXB2176	Matrx Elite
AXB2177	Matrx Elite TR
AXB2178	Matrx Elite Deep
AXB2170	Without backrest upholstery
AXB2151	Black fabric
AXB2152	Black multi stretch PU coated fabric
AXB2155	Black Spacer fabric
AXB2196	Long 375mm
AXB2197	Long ergonomic 397mm
AXB2199	Soft
AXB2323	Elan standard pad 150mm (6")
AXB2324	Elan standard pad 255mm (10")
AXB2322	Elan standard pad 355mm (14" )
AXB2325	Elan Occipital pad 230mm x 180mm (9" X 7")
AXB2326	Elan Occipital pad 305mm x 200mm (12" X 8")
AXB2327	Elan 4 point headrest pad 280mm x 255mm (11" X 10")
AXB1496	All bracketry as per option 1432
AXB1495	Push bar 22mm
AXB1446	Rigid pad
AXB1445	Multi-adjustable pad
AXB0977	Powered
AXB0526	Standard 80° hangers
AXB0529	Manual elevating legrest
AXB0501	Without hanger
AXB0512	Fixed Footplate
AXB2227	RH Bracket mounted for Attendant
AXB2228	LH Bracket mounted for Attendant
AXB2220	Dual control REM050
AXB2219	Compact (ACU200)

AXB1329	Lap Belt with auto-lock
AXB2701	Evoflex pelvic stabilizer, size L
AXB2702	Evoflex pelvic stabilizer, size M
AXB2703	2-Point hip belt centre pull, size L
AXB2704	2-Point hip belt centre pull, size M
AXB2705	2-Point hip belt dual pull, size L
AXB2706	2-Point Hip Belt Dual pull, size M
AXB2128	Bag hooks
AXB1404	Tray, Swing Away RH
AXB1405	Tray, Swing Away LH



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

**Special Instructions**

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):