

**matrix**<sup>®</sup>  
SEATING SERIES

# **Flo-tech™ Contour Visco**



Contoured Visco memory foam that  
conforms to the shape of the body



The *Invacare Matrix Flo-tech Contour Visco* cushion is moulded from Visco-elastic foam. Its contoured surface area conforms to the body, improving posture, increasing comfort and maximising the surface area in contact with the individual, thereby increasing pressure redistribution.

### Deep moulded low memory foam cushion

The visco-elastic foam conforms to the body. By maximising the contact area it redistributed pressure and delivers exceptional levels of comfort. The deep foam makes this cushion particularly suitable for heavier and taller people.

### Contoured seating surface with pommel and side support

Offers increased lateral support and pelvic stability to enhance posture.

## Features and Options



### Ingress resistant cover

Cover is two-way stretch, water-resistant and vapour-permeable. Cover features an anti-slip base with carry handle for easy transportation.



### Available in 58 sizes

The *Invacare Matrix Flo-tech Contour Visco* is available in a range of 58 sizes, from 14-20" widths and 13-21" depths. (See size chart for details).









### Four foam densities available

Available in four densities to accommodate a range of client weights - Soft (up to 51 kg), Medium (up to 76 kg), Firm (up to 108 kg) and Very Firm (up to 127 kg).

## Technical data


For more comprehensive pre-sales information about this product, including the product's user manual, please visit [www.invacare.co.uk](http://www.invacare.co.uk)

Various sizes	Height	Total product weight	Max. user weight <sup>1</sup>	Max. cleaning temperature	Warranty
					
Up to 58 sizes	120 mm	1.6 kg	soft: 51kg (8 stone) medium: 76kg (12 stone) firm: 108kg (17 stone) very firm: 127kg (20 stone)	80°	3 years

1. Guide maximum user weight limits are set assuming that the cushion selected is of appropriate width and depth for the intended user/s.

Manufactured to European Standards: Interior Foam: BS EN 1021-1, BS EN 1021-2, BS 5852: 1982 Crib 5

## Cover options



Invacare Limited  
Pencoed Technology Park - Pencoed  
Bridgend CF35 5AQ - United Kingdom  
Tel: +44 1656 776 222  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS  
County Dublin - Ireland  
Tel. +353 1 8107084  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

