



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Powerchairs Retail Prescription Form

FEBRUARY 2021
 Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

Fox in a Box

Standard Features

- Max user weight 127Kg (20 stone)
- Crash tested ISO7176-Part 19
- Seat width 400mm (16") or 450mm (18")
- Seat depth Set at 460mm (18") adjustable to 370mm (14.5")
- Seating Black Nylon sling seat
- Backrest Black Nylon sling upholstery
- Backrest height 490mm (19")
- Detachable T-armrest with width adjustment
- Long Ergonomic armpad (397mm)
- Low seat to floor height 480mm adjustable to 510mm
- Manually settable seat Tilt (0° / 2° / 4° / 6° / 8°) factory setting: 4°
- Backrest Recliner, settable in steps (85° / 95° / 105° / 115°)
- Legrest hangers Standard hangers 80°
- Footrests Fixed angle, flip-up footplates
- Control system LiNX drive only Simply SMART
- Shroud colours Pearly Night Blue or Cherry Red
- Castor Grey 200mm x 50mm (8" x 2") solid inlay, Silver Rim
- Drive wheel puncture proof solid Inlay, (12 1/2" x 2 1/4"), Silver Rim
- Speed setting 6km/h (4mph)
- 34Ah battery (range 27km)
- Lap belt with auto lock (adjustable both sides)
- Manual wheel locks
- Charger 8A (200 - 250V)
- 4-Point tie-down



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

| Code | Description | |
|-------------------------------|---|--------|
| FOX 400mm x 400mm RED | | |
| 1640339 | <input type="checkbox"/> Fox 16" x 16" Red | £2,514 |
| FOX 450mm x 430mm RED | | |
| 1640345 | <input type="checkbox"/> Fox 18" x 17" Red | £2,514 |
| FOX 400mm x 400mm BLUE | | |
| 1640346 | <input type="checkbox"/> Fox 16" x 16" Blue | £2,514 |
| FOX 450mmx 430mm BLUE | | |
| 1640347 | <input type="checkbox"/> Fox 18" x 17" Blue | £2,514 |

| | |
|--------------------|--|
| TOTAL PRICE | |
|--------------------|--|