

Tel: 01 810 7084 Fax: 01 810 7085  
Email: ordersireland@invacare.com  
Online Spares available at: www.invacare.ie

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :

**Invacare® FOX (Modulite™ Seating System)**  
Ireland  
Max User Weight 127Kg (20 stone)

**Standard Features**



Invacare® **Modulite™**  
An evolution in seating designed around you.

- **Modulite** highly adjustable seating system
- Telescopic adjustable Seat Widths from 330mm - 530mm (13" - 21")
- Seat Depths 410mm - 510mm Stepless Adjustment (16" - 20")
- Tension adjustable backrest
- Standard T-armrest, detachable with Short armpads
- Manual settable seat Tilt (from 0° to 8°)
- Manual settable backrest Recline (0° - 30° in 3.75° steps)
- Standard 80° Legrest hangers with fixed flip-up footplates
- **LiNX** Control REM110 with standard bracket
- Choice of 5 colour options
- Solid Grey castors 200mm x 50mm (8" x 2"), Silver Rim
- Grey Drive Wheels, Puncture Proof Solid Inlay, (12 1/2" x 2 1/4"), Silver Rim (aluminium)
- 6km/h (4mph) Motors with 40Ah Gel Battery
- 8A (200 - 250V) Charger
- 4-Point tie-down system
- Lapbelt with autolock (adjustable both sides)
- Wheel lock / parking brake









**Standard Base Price : €3,837**

\* Configurations may vary from the image shown


	<input checked="" type="checkbox"/>	Please make your selection with a cross in the box	◆ Standard feature	○ Free of Charge		
Customer weight =		Kg				
Code		Description		Ireland		
AXB0000		Fox Basic Model		◆		
AXZ0000		Fox Special Basic Model		POA		
<b>SEATING TYPE : MODULITE MULTI ADJUSTABLE SEATING SYSTEM</b>						
AXB2100	<input checked="" type="checkbox"/>	Modulite Seating Unit		€3,837		
<b>SEAT FRAME AND SEAT PLATE</b>						
AXB2111	<input checked="" type="checkbox"/>	Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)		◆		
<b>MODULITE SEAT SIZE SETTING CHART:</b>						
		<b>SEAT WIDTHS</b>				
		330mm* (13")	AXB2121 380mm (15")	AXB2122 430mm (17")	AXB2123 480mm (19")	AXB2124 530mm (21")
<b>SEAT DEPTHS</b>	AXB2132 (16") 410mm	○	○	○	○	○
	AXB2134 (18") 460mm	○	○	○	○	○
	AXB2136 (20") 510mm	○	○	○	○	○



\*330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)  
 Seat depth 410mm: cannot combine T-Armrests (2191) with flexible lapbelt mount (2113)  
 Modulite Teen seating available, seat width and depth from 305mm (12")

<b>SEAT OPTION</b>			
AXBSWO	"H" Frame Interface Kit		€460
<b>SEAT CUSHION (sizes according to seat frame settings)</b>			
AXB2142	Matrx Flo-tech Lite Cushion		€116
AXB2147	Matrx Libra Cushion, 90mm thick	Not available with seat width 530mm with 410mm seat depth, or spacer fabric upholstery	€483
AXB2140	Without seat cushion		◆
<b>BACKREST SUPPORT</b>			
AXB2163	Backposts 480mm (Short)	Not available with flex3 adjustable backrests	○
AXB2162	Backposts 540mm (Tall)	Must be selected for flex3 adjustable backrest plates	◆
<b>PUSH HANDLES FOR BACKPOSTS</b>			
AXB2107	Backposts WITH push handles		◆
AXB2108	Backposts WITHOUT push handles Push handles will be cut down on the vertical section.	Not available with attendant controls, dual controls, matrx backrest's, headrest 1495 or tension adjustable backrest in combination with a headrest.	€44
<b>BACKREST UPHOLSTERY / CUSHION</b>			
AXB2171	Basic Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear, Black fabric only	Not available with laterals, headrest 1496	◆
AXB2179	Flex3 Tension Adjustable (540mm - 480mm) Black Dartex only	Not available with laterals, soft armpad, chin remote	€165
AXB2172	Flex3 Adjustable Backrest Plates (without Backrest Cushion)		€229
AXB2173	Flex3 Adjustable Backrest Plates (with Laguna Backrest Cushion)		€372
AXB2174	Matrx Posture Back, height 400mm (Mesh fabric)	Not available with laterals, 530mm seat width, headrest 1496	€519
AXB2175	Matrx Posture Back, height 500mm (Mesh fabric)	Not available with laterals, 530mm seat width, headrest 1496	€519
AXB2176	Matrx Elite Back, height 500mm (Mesh fabric)	Not available with laterals, 530mm seat width, headrest 1496, seat width of 380mm 2121 in combination with flip up, following armrests	€519
AXB2177	Matrx Elite TR Back, height 500mm (Mesh fabric)	Not available with laterals, 530mm seat width, headrest 1496	€519
AXB2178	Matrx Elite Deep Back, height 500mm (Mesh fabric)	Not available with laterals, headrest 1496	€672
AXB2170	Without Backrest Upholstery	Not available with laterals, headrest 1496	○
<b>UPHOLSTERY TYPE</b>			
AXB2151	Black fabric		◆
AXB2152	Black Dartex		○
AXB2155	Black Spacer fabric		○
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;"> <p><b>T-Armrest 2191</b></p>  </div> <div style="text-align: center;"> <p><b>In-rigged armrest 2183</b></p>  </div> <div style="text-align: center;"> <p><b>Short armpad</b></p>  </div> <div style="text-align: center;"> <p><b>Long armpad</b></p>  </div> <div style="text-align: center;"> <p><b>Long ergonomic armpad</b></p>  </div> <div style="text-align: center;"> <p><b>Soft armpad</b></p>  </div> </div>			
<b>ARMREST (Pair)</b>			
AXB2191	T-Armrest Standard Range of adjustment from seat plate 245mm to 315mm		◆
AXB2181	T-Armrest Extended height Range of adjustment from seat plate 295mm to 365mm		€44
AXB2183	In-rigged Armrests (reduces seat width by 25mm each side)		€26
<b>ARMPADS (Pair)</b>			
AXB2195	Short Armpad (325mm)		◆
AXB2196	Long Armpad (375mm)	Not available with seat depth 410mm	○
AXB2197	Long Ergonomic Armpad (397mm)	Not available with seat depth 410mm	○
AXB2199	Soft Armpad (REA Azalea Wide Pad)		€93

**SEAT ACCESSORIES**



<b>HEADREST / NECKREST - Adequate headrest option must be selected</b>			
AXB1350		Neckrest: height and depth adjustable	€77
AXB1432		Headrest: height, depth and sideways adjustable	€215
AXB1427		Headrest with cheek support: height and depth adjustable	€215
AXB2323		Elan Standard Pad 150mm (6") wide with standard cover complete with Elan Headrest Hardware	Not available with head controls €369
AXB2324		Elan Standard Pad 255mm (10") wide with standard cover complete with Elan Headrest Hardware	Not available with head controls €369
AXB2322		Elan Standard Pad 355mm (14") wide with standard cover complete with Elan Headrest Hardware	Not available with head controls €369
AXB2325		Elan Occipital Pad 230mm x 180mm (9" x 7") with standard cover complete with Elan Headrest Hardware	Not available with head controls €369
AXB2326		Elan Occipital Pad 305mm x 200mm (12" x 8") with standard cover complete with Elan Headrest Hardware	Not available with head controls €369
AXB2327		Elan 4 Point Headrest Pad 280mm x 255mm (11" x 10") with standard cover complete with Elan Headrest Hardware	Not available with head controls €497
AXB1496		All bracketry as option 1432 without Headrest cushion	○
AXB1495		Push bar dia. 22mm in conjunction with telescopic seat frame / Backposts, No Headrest - no mounting hardware.	Not available without push handles 2108 ○



**Rigid pad  
1446**



**Multi-adjustable pad  
1445**

<b>LATERAL SUPPORT (Pair)</b>			
AXB1446		Rigid pad (Azalea type, Dartex)	€154
AXB1445		Multi-adjustable pad (Azalea type, Dartex)	€154

**SEAT ADJUSTMENTS**

<b>SEAT HEIGHT SETTING</b>			
AXB0312		Low Seat To Floor height 480mm at 4° factory setting (-30mm / + 30mm depending on Tilt settings)	○
AXB0344		High Seat To Floor height 510mm at 4° factory setting (- 30mm / + 30mm depending on Tilt settings)	○

<b>SEAT TILT</b>			
AXB1002	X	Manually settable Seat Tilt (0° / 2° / 4° / 6° / 8°), Factory setting: 4°	◆

<b>BACKREST RECLINER</b>			
AXB0975		Manual Settable Recliner (0° - 30° in 3.75° steps) Foldable	◆
AXB0977		Powered Recliner (0° - 30°) Foldable	€307

**LEGREST / FOOTREST**


<b>FOOTREST HANGER Type, Swing away and detachable (Pair)</b>			
AXB0526		Standard 80° Hangers	◆
AXB0529		Manual Elevating Legrests	€195
AXB0501		Without Hanger	○

<b>AMPUTEE SUPPORTS (Replaces legrest hanger)</b>			
AXB0554		Amputee Support RH	€183
AXB0555		Amputee support LH	€183

<b>FOOTPLATE TYPE (Pair)</b>			
AXB0512		Fixed Angle, Flip-up Footplates	◆
AXB0513		Angle and Depth Adjustable Footplates (Adapt)	€83



**LINX PRIMARY CONTROLS**



<b>LINX PRIMARY REMOTE &amp; BRACKET</b>			
AXB2201		REM110 Drive only	◆
AXB2203		REM211 Drive & Actuators	○
AXB1207		Standard Bracket RH	○
AXB1208		Standard Bracket LH	○
AXB1707		Swing away RH	€77
AXB1708		Swing away LH	€77
AXB2227		Bracket RH mounted for Attendant	€17
AXB2228		Bracket LH mounted for Attendant	€17

<b>LINX DUAL CONTROL &amp; BRACKET: CAN ONLY BE SELECTED IN ADDITION TO LINX USER CONTROL</b>			
AXB2220		Dual Control REM050 Drive only	€424
AXB2219		Compact Dual Control (ACU200)	€415
AXB1207		Standard Bracket RH	○
AXB1208		Standard Bracket LH	○
<div style="display: flex; justify-content: space-around; align-items: flex-start;"> <div style="text-align: center;"> <p><b>NEW 3" U-Shape</b></p> </div> <div style="text-align: center;"> <p><b>NEW 4" U-Shape</b></p> </div> <div style="text-align: center;"> <p><b>NEW Rubber dome</b></p> </div> <div style="text-align: center;"> <p><b>NEW Grooved Mushroom</b></p> </div> </div>			
<b>BODYPOINT JOYSTICK KNOBS</b>			
AXB2721		3" wide U-Shaped, Joystick Handle	€68
AXB2722		4" wide U-Shaped, Joystick Handle	€68
AXB2723		Wide Rubber Dome, Joystick Handle	€68
AXB2724		Grooved Mushroom Joystick Handle	€68
<b>DRIVE UNIT OPTIONS</b>			
<b>WHEELBASE SETTING</b>			
AXB1226	X	Recommended factory setting, based on seat depth and seat to floor height Wheel setting positions can be changed as desired after performing the appropriate risk assessment	◆
<b>SHROUD COLOUR</b>			
AXB0201		Pearl Night Blue	○
AXB0233		Metallic Silver	○
AXB0228		Silky Matt Black	○
AXB0204		Transparent Cherry Red	○
AXB0222		White	○
<b>CASTORS</b>			
AXB0802		Castor Grey 200mm x 50mm (8" x 2") PU Tyre, Black Rim	○
AXB0803		Castor Grey 200mm x 50mm (8" x 2") Solid Inlay, Silver Rim	◆
<b>DRIVE WHEELS</b>			
AXB0709		Drive Wheel Grey, Puncture Proof Solid Inlay, (12 1/2" x 2 1/4"), Black Rim (plastic)	○
AXB0710		Drive Wheel Grey, Puncture Proof Solid Inlay, (12 1/2" x 2 1/4"), Silver Rim (aluminium)	◆
<b>PERFORMANCE OPTIONS</b>			
<b>MOTOR TYPE</b>			
AXB1294		3 km/h (2mph)	○
AXB1282		6 km/h (4mph)	◆
<b>BATTERY</b>			
AXB1255		40Ah Gel	◆
AXB1248		50Ah AGM	€33
<b>PERFORMANCE ACCESSORIES</b>			
<b>CHARGER</b>			
AXB1343		OptiCharge Charger 8A (200 - 250V)	◆
AXB1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	€140
<b>SAFETY ACCESSORIES</b>			
<b>WHEEL LOCK / PARKING BRAKE</b>			
AXB1206	X	Manual Wheel Locks	◆
<b>KERB RIDER</b>			
AXB1244		Kerbrider	€195
<b>TRANSPORTATION KIT</b>			
AXB1494	X	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) <b>(Adequate headrest is recommended)</b>	◆
<b>LAP BELT</b>			
AXB1329	X	Lap Belt with Auto-Lock	◆
<b>BODYPOINT PELVIC BELTS (Includes mounting kit)</b>			
AXB2701		Evoflex Pelvic Stabilizer, Rehab Latch size L (padded length 620mm)	Not available with seat width 380mm, 430mm, seat depth 410mm, hip supports €107
AXB2702		Evoflex Pelvic Stabilizer, Rehab Latch size M (padded length 460mm)	Not available seat width 480mm, 530mm, seat depth 410mm, hip supports €107
AXB2703		2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)	Not available with seat width 380mm, 430mm, seat depth 410mm, hip supports €98
AXB2704		2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)	Not available seat width 480mm, 530mm, seat depth 410mm, hip supports €98
AXB2705		2-Point Hip Belt, Dual pull, Side release button, Cinch mount size L (padded length 620mm)	Not available with seat width 380mm, 430mm, seat depth 410mm, hip supports €98
AXB2706		2-Point Hip Belt, Dual pull, Side release button, Cinch mount size M (padded length 460mm)	Not available seat width 480mm, 530mm, seat depth 410mm, hip supports €98

**ACCESSORIES****BAG HOOKS**

AXB2128		Bag hooks (Max load up to 6 kg = impact on max user weight)	Not available with flex3 backrest 2172 / 2173, attendant brackets, dual controls	€86
<b>CANE HOLDER</b>				
AXB1312		Cane holder RH		€24
AXB1363		Cane holder LH		€24
<b>TRAY - mounted opposite side to controller</b>				
AXB1404		Tray, Swing Away RH		€203
AXB1405		Tray, Swing Away LH		€203
<b>CASTOR LOCKS</b>				
AXB1490		Castor Lock Levers		€30

**Total Price:****Special Instructions - Please supply the following information:****Name:****Email:****Phone:****Fax:****Description of special request (please include relevant option code or part number):**