

Tel: 01810 7084 Fax: 01 810 7085
Email: ordersireland@invacare.com
Online Spares available at: www.invacare.ie

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :

Invacare® FOX Modulite Teen
Ireland
Max User Weight 127Kg (20 stone)

Invacare® **Modulite™**
An evolution in seating designed around you.

Standard Features



- Telescopic adjustable Seat Widths from 305mm - 380mm (12" - 15")
- Seat Depths 305mm - 510mm Stepless Adjustment (12" - 20")
- Choice of backrest (Tension adjustable, Flex 3 Laguna, Flex 3 Mistral and Matrix)
- Detachable T-armrest, Height, Width and Depth adjustable with standard armpads
- Manually settable seat Tilt (from 0° to 8°) and backrest recline (0° - 30°)
- Standard 80° Hangers with fixed flip-up footplates
- **LiNX** Control System with standard bracket
- Choice of 5 colour options
- Solid Grey castors 200mm x 50mm (8" x 2"), Silver Rim
- Grey Drive Wheels, Puncture Proof Solid Inlay, (12 1/2" x 2 1/4"), Silver Rim (aluminium)
- 6 km/h (4mph) Motors with 40 Ah Gel Battery
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)
- Wheel lock / parking brake
- 8A (200 - 250V) Charger



* Configurations may vary from the image shown

Standard Base Price : €3,256

	X	Please make your selection with a cross in the box	◆ Standard feature	○ Free of Charge
		Customer weight =	Kg	
Code		Description	Ireland	
AXZ0000		Fox Modulite Teen	€3,256	
AXZ0000		Fox Modulite Teen Special Basic Model	POA	
SWO0025	X	Invacare ONE to ONE fitting service	◆	
SEATING TYPE : MODULITE MULTI ADJUSTABLE SEATING SYSTEM				
AXZ2100	X	Modulite Seating Unit	◆	
SEAT FRAME AND SEAT PLATE				
AXB2111	X	Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)	◆	

MODULITE SEAT SIZE SETTING CHART:

		SEAT WIDTHS		
		**305mm (12")	*330mm (13")	380mm (15")
SEAT DEPTH	(12") 305mm	○	○	○
	(13") 330mm	○	○	○
	(14") 355mm	○	○	○
	(15") 380mm	○	○	○
	(16") 410mm	○	○	○
	(18") 460mm	○	○	○
	(20") 510mm	○	○	○



**305mm (12") seat width achieved on a 380mm seat frame (2121) with inrigged armrest (2183) and additional i
*330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)

SEAT OPTION

AXBSWO	"H" Frame Interface Kit	Only Available in a 15" seat width	€ 460
--------	-------------------------	------------------------------------	-------

SEAT CUSHION (sizes according to seat frame settings)

AXB2142	Matrx Flo-tech Lite Cushion	Only Available in a 15" seat width	€116
AXB2147	Matrx Libra Cushion, 90mm thick	Only Available in a 15" seat width, Spacer fabric upholstery	€483

MATRIX MINI SERIES CUSHIONS

1626906	PS cushion with Kidabra print 310mm x 310mm, 12" x 12" MAX User Weight 68kg	€194
1626907	PS cushion with Kidabra print 310mm x 360mm, 12" x 14" MAX User Weight 68kg	€194
1626908	PS cushion with Kidabra print 310mm x 410mm, 12" x 16" MAX User Weight 68kg	€194
1627456	PS cushion with Kidabra print 330mm x 310mm, 13" x 12" MAX User Weight 68kg	€194
1627457	PS cushion with Kidabra print 330mm x 360mm, 13" x 14" MAX User Weight 68kg	€194
1627458	PS cushion with Kidabra print 330mm x 410mm, 13" x 16" MAX User Weight 68kg	€194
1626909	PS cushion with Kidabra print 360mm x 310mm, 14" x 12" MAX User Weight 68kg	€194
1626910	PS cushion with Kidabra print 360mm x 360mm, 14" x 14" MAX User Weight 68kg	€194
1626911	PS cushion with Kidabra print 360mm x 410mm, 14" x 16" MAX User Weight 68kg	€194
AXB2140	Without seat cushion	◆

BACKREST SUPPORT

AXB2162	Backposts 540mm (Tall)	Must be selected for Flex3 adjustable backrest plates	◆
AXB2163	Backposts 480mm (Short)	Not available with Flex3 adjustable backrests	○

PUSH HANDLES FOR BACKPOSTS

AXB2107	Backposts WITH push handles	◆	
AXB2108	Backposts WITHOUT push handles Push handles will be cut down on the vertical section.	Not available with Attendant controls, Dual controls, Matrix backrest's, Headrest 1495 or Tension adjustable backrest in combination with a Headrest.	€44














BACKREST UPHOLSTERY / CUSHION

AXB2171	Basic Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear, Black fabric only	Not available with laterals, headrest 1496 only available with 15" seat width	◆
AXB2172	Flex3 Adjustable Backrest Plates (without Backrest Cushion)	Only available in a 15" seat width	€229
AXB2173	Flex3 Adjustable Backrest Plates (with Laguna Backrest Cushion)	Only available in a 15" seat width	€372
SP1562143	Mistral Flex 3 Backrest shell - shorter & narrower than the Laguna Flex 3	Only available in a 15" seat width	€131.81
SP1431668-TR26	Mistral Backrest Cushion, Black Dartex	Only available in a 15" seat width	€139.34
SP1431668-TR35	Mistral Backrest Cushion, Black Plush	Only available in a 15" seat width	€151.94

MATRIX MINI SERIES BACKREST: Extended Mounting Pins Required with all Mini Backs

EMP	Extended Mounting Pins	
MINI1207	Mini Back with Kidabra print 310mm x 180mm, 12" x 7" MAX User Weight 68kg	€586
MINI1210	Mini Back with Kidabra print 310mm x 260mm, 12" x 10" MAX User Weight 68kg	€586
MINI1212	Mini Back with Kidabra print 310mm x 310mm, 12" x 12" MAX User Weight 68kg	€586
MINI1214	Mini Back with Kidabra print 310mm x 360mm, 12" x 14" MAX User Weight 68kg	€586
MINI1307	Mini Back with Kidabra print 330mm x 180mm, 13" x 7" MAX User Weight 68kg	€586
MINI1310	Mini Back with Kidabra print 330mm x 260mm, 13" x 10" MAX User Weight 68kg	€586
MINI1312	Mini Back with Kidabra print 330mm x 310mm, 13" x 12" MAX User Weight 68kg	€586
MINI1314	Mini Back with Kidabra print 330mm x 360mm, 13" x 14" MAX User Weight 68kg	€586
MINI1407	Mini Back with Kidabra print 360mm x 180mm, 14" x 7" MAX User Weight 68kg	€586
MINI1410	Mini Back with Kidabra print 360mm x 260mm, 14" x 10" MAX User Weight 68kg	€586
MINI1412	Mini Back with Kidabra print 360mm x 310mm, 14" x 12" MAX User Weight 68kg	€586
MINI1414	Mini Back with Kidabra print 360mm x 360mm, 14" x 14" MAX User Weight 68kg	€586

AXB2170		Without Backrest Upholstery	Not available with Laterals, Headrest 1496	○
---------	--	-----------------------------	--	---

BACKREST OPTIONS			
SP1571542 RH SP1571553 LH	Angle adjustable backrest support brackets	Not compatible with Powered Backrest recliner	€45.34
UPHOLSTERY TYPE			
AXB2151	Black fabric		◆
AXB2152	Black Dartex		○
AXB2155	Black Spacer fabric		○
ARMREST (Pair)			
	T-Armrest 2191		In-rigged armrest 2183
			Hip supports
			Short armpad
			Long armpad
AXB2191	T-Armrest Standard Range of adjustment from seat plate 245mm to 315mm		◆
AXB2181	T-Armrest Extended height Range of adjustment from seat plate 295mm to 365mm		€44
AXB2183	In-rigged Armrests (reduces seat width by 25mm each side)	Must Be selected with a 12" & 13" Seat Width	○
SP1582159	Inrigged Armrest Standard inrigged arm option MUST BE SELECTED with additional brackets.	Must Be selected with a 12" & 13" Seat Width, Not Available with hip support	€36.35
ARMPADS (Pair)			
AXB2195	Short Armpad (325mm)		◆
AXB2196	Long Armpad (375mm)	Not available with seat depth 410mm	○
SEAT ACCESSORIES			
HEADREST / NECKREST - Adequate headrest option must be selected			
	Headrest 1432		Matrix Elan Adjustment
			Elan Standard Pad, 2 sizes
			Elan Occipital Pad
			Elan 4 Point Pad
			Matrix Elan Vs Matrix mini Elan
ELAN SERIES HEADREST'S			
1603255	Elan Standard Pad - Small (160mm x 90mm) Standard cover		€162
1603254	Elan Standard Pad - Medium (260mm x 130mm) Standard cover		€162
1603291	Elan Occipital Pad (230mm x 180mm) Standard cover		€162
1629822	Elan 4 Point Headrest Pad (Small) Standard cover		€300
Please select Mounting hardware type for Elan headrest's			
1603250	Matrix Elan 6-point Hardware Infinite adjustments		€238
1603251	Matrix Elan 6-point Mini (paediatric) 100mm shorter post		€238
1603249	Matrix standard multi-axis hardware Limited adjustability		€238
1603252	Elan adapter for Flex 3 back		€20
STANDARD HEADREST OPTIONS			
AXB1432	Headrest: height, depth and sideways adjustable		€226
AXB1496	All bracketry as option 1432 without Headrest cushion		○
AXB1495	Push bar dia. 22mm in conjunction with telescopic seat frame / Backposts, No Headrest - no mounting hardware.		○
LATERAL SUPPORT (Pair) For Flex3 Backs			
	Lateral support 1446		Lateral multi-adjustable pad 1445
AXB1446	Rigid pad (Azalea type, Dartex)	Not available with Tension Backrest or Matrix Backrest's	€154
AXB1445	Multi-adjustable pad (Azalea type, Dartex)	Not available with Tension Backrest or Matrix Backrest's	€154
MATRIX MINI LATERAL SUPORTS			
1627466	Swing-away lateral Hardware (each)		€149
EPAD-XS	Extra small lateral pad (each)		€78
1627465	Chest harness interface kit		€56
SEAT ADJUSTMENTS			
SEAT HEIGHT SETTING			
AXB0312	Low Seat To Floor height 480mm at 4° factory setting (-30mm / + 30mm depending on Tilt settings)		○
AXB0344	High Seat To Floor height 510mm at 4° factory setting (- 30mm / + 30mm depending on Tilt settings)		○

SEAT TILT			
AXB1002	X	Manually settable Seat Tilt (0° / 2° / 4° / 6° / 8°), Factory setting: 4°	◆
BACKREST RECLINER			
AXB0975		Manual Settable Recliner (0° - 30° in 3.75° steps) Foldable	◆
AXB0977		Powered Recliner (0° - 30°) Foldable	€307
LEGREST / FOOTREST			
FOOTREST HANGER Type, Swing away and detachable (Pair)			
AXB0526		Standard Hangers 80°	◆
1606173-0601		Black short 90° Legrest Hanger x 2 - Action style.	€88.49
1606172-0601		Black short 80° Legrest Hanger x 2 - Action style.	€86.79
AXB0501		Without Hanger	○
NB if either of the above legrest options is selected , the following items will be required:			
SP1531195		Legrest support carrier LH	€34.04
SP1531196		Legrest support carrier RH	€34.04
AMPUTEE SUPPORTS (Replaces legrest hanger)			
AXB0554		Amputee Support RH	€183
AXB0555		Amputee support LH	€183
FOOTPLATE TYPE (Pair), Standard legrest options.			
AXB0512		Fixed Angle, Flip-up Footplates	◆
AXB0513		Angle And Depth Adjustable Footplates (Adapt)	€83
FOOTPLATE OPTIONS: for use with Short 80° & 90° Legrest			
1428660 LH 1428654 RH		Angle and depth adjustable Footplate pair. Including stem	€32.41
1530620 LH 1530614 RH		Fixed Footplate pair. Including stem	€17.72
1492150 LH 1492149 RH		Upmounted Footplate	€108.79



LiNX PRIMARY CONTROLS

Drive Only				Drive & Actuators				Dual Control				Compact Dual Control				
REM110				REM211				REM050				ACU200				
LiNX PRIMARY REMOTE & BRACKET																
AXB2201		REM110 Drive only														◆
AXB2203		REM211 Drive & Actuators														○
AXB1207		Standard Bracket RH														○
AXB1208		Standard Bracket LH														○
AXB1707		Swing Away RH														€77
AXB1708		Swing Away LH														€77
AXB2227		Bracket RH mounted for Attendant														€17
AXB2228		Bracket LH mounted for Attendant														€17
LiNX DUAL CONTROL & BRACKET: CAN ONLY BE SELECTED IN ADDITION TO LiNX USER CONTROL																
AXB2220		Dual Control REM050 Drive only														€424
AXB2219		Compact Dual Control (ACU200)														€415
AXB1207		Standard Bracket RH														○
AXB1208		Standard Bracket LH														○
DRIVE UNIT OPTIONS																
WHEELBASE SETTING																
AXB1226	X	Recommended factory setting, based on seat depth and seat to floor height Wheel setting positions can be changed as desired after performing the appropriate risk assessment														◆
SHROUD COLOUR																
AXB0201		Pearl Night Blue														○
AXB0233		Metallic Silver														○
AXB0228		Silky Matt Black														○
AXB0204		Transparent Cherry Red														○
AXB0222		White														○

CASTORS			
AXB0802		Castor Grey 200mm x 50mm (8" x 2") PU Tyre, Black Rim	○
AXB0803		Castor Grey 200mm x 50mm (8" x 2") Solid Inlay, Silver Rim	◆
DRIVE WHEELS			
AXB0709		Drive Wheel Grey, Puncture Proof Solid Inlay, (12 1/2" * 2 1/4"), Black Rim (plastic)	○
AXB0710		Drive Wheel Grey, Puncture Proof Solid Inlay, (12 1/2" * 2 1/4"), Silver Rim (aluminium)	◆
PERFORMANCE OPTIONS			
MOTOR TYPE			
AXB1294		3 km/h (2mph)	○
AXB1282		6 km/h (4mph)	◆
BATTERY			
AXB1248		50 Ah AGM	€ 33
AXB1255		40 Ah Gel	◆
PERFORMANCE ACCESSORIES			
CHARGER			
AXB1343		OptiCharge Charger 8A (200 - 250V)	◆
AXB1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	€ 140
SAFETY ACCESSORIES			
WHEEL LOCK / PARKING BRAKE			
AXB1206	X	Manual Wheel Locks	◆
KERB RIDER			
AXB1244		Kerbrider	€195
TRANSPORTATION KIT			
AXB1494	X	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is recommended)	◆
LAP BELT			
AXB1329	X	Lap Belt with Auto-Lock	◆
BODYPOINT PELVIC BELTS (Includes mounting kit)			
HB216-S38-B2		HB216 2-point belt Small x 380mm Pads, cinch mount	€84
HB216-M46-B2		HB216 2-point belt Medium x 460mm Pads, cinch mount	€84
HB219-S38-B2		HB219 2-point belt Small x 380mm Pads, cinch mount	€96
HB219-M46-B2		HB219 2-point belt Medium x 460mm Pads, cinch mount	€96
HB235-S32-B2		HB235 2-point belt Extra small x, 320mm Pads, cinch mount	€100
HB235-S38-B2		HB235 2-point belt Small x 380mm Pads, cinch mount	€100
EB235S		Evoflex Small pelvic stabilizer, push button	€120
EB225M		Evoflex Medium pelvic stabilizer, Latch	€120
EB225L		Evoflex Large pelvic stabilizer, Latch	€120
BODYPOINT HARNESSSES			
SH250XXS-B2		PivotFit XX-Small, front pull, cinch mount, dynamic trunk support	€143
SH250XS-B2		PivotFit X-Small, front pull, cinch mount, dynamic trunk support	€143
SH250S-B2		PivotFit Small, front pull, cinch mount, dynamic trunk support	€143
SH250M-B2		PivotFit Medium, front pull, cinch mount, dynamic trunk support	€143
SH250XXS-C2-B2		PivotFit XX-Small, shoulder harness, Plain top strap, Cinch mount	€152
SH250XS-C2-B2		PivotFit Small, shoulder harness, Plain top strap, Cinch mount	€152
SH250M-C2-B2		PivotFit Medium, shoulder harness, Plain top strap, Cinch mount	€152
SH260XXS-B2		PivotFit XX-small, rear pull, cinch mount, standard	€152
SH260XS-B2		PivotFit X-small, rear pull, cinch mount, standard	€152
SH260S-B2		PivotFit Small, rear pull, cinch mount, standard	€152
SH260M-B2		PivotFit Medium, rear pull, cinch mount, standard	€152
SH270XXS-B2		PivotFit XX-Small, Front-Pull cinch mount, dynamic trunk support	€152
SH270XS-B2		PivotFit X-Small, Front-Pull cinch mount, dynamic trunk support	€152
SH270S-B2		PivotFit Small, Front-Pull, cinch mount, dynamic trunk support	€152
SH270M-B2		PivotFit Medium, Front-Pull cinch mount, dynamic trunk support	€152
SH270XXS-C2-B2		PivotFit XXS, Front pull, Plain top strap, cinch mount	€158
SH270XS-C2-B2		PivotFit XS, Front pull, Plain top strap, cinch mount	€158
SH270S-C2-B2		PivotFit Small, Front pull, Plain top strap, cinch mount	€158
SH270M-C2-B2		PivotFit Medium, Front pull, Plain top strap, cinch mount	€158
SH280XXS-B2		PivotFit XX-Small, rear pull, cinch mount, dynamic	€158
SH280XS-B2		PivotFit X-Small, rear pull, cinch mount, dynamic	€158
SH280S-B2		PivotFit Small, rear pull, cinch mount, dynamic	€158
SH280M-B2		PivotFit Medium, rear pull, cinch mount, dynamic	€158
1627465		Chest harness interface kit	€48

ACCESSORIES

BAG HOOKS *NEW*



AXB2128		Bag hooks (Max load up to 6 kg = impact on max user weight)	Not available with Flex3 Backrest 2172/2173, Attendant Brackets, Dual Control	€72
CANE HOLDER				
AXB1312		Cane holder RH		€24
AXB1363		Cane holder LH		€24
TRAY - mounted opposite side to controller				
AXB1404		Tray, Swing Away RH		€203
AXB1405		Tray, Swing Away LH		€203
CASTOR LOCKS				
AXB1490		Castor Lock Levers		€30

Total Price:



Please refer to our Invacare Unique special order form, for more complex custom requirements. Available at www.invacare.co.uk

Special Instructions - Please supply the following information:

Name: _____ **Email:** _____

Phone: _____ **Fax:** _____

Description of special request (please include relevant option code or part number):