

Tel: 01 810 7084 Fax: 01 810 7085
Email: ordersireland@invacare.com

Online spares available at www.invacare.ie

Customer Ref :	Account No :
Contact No:	Order No :
Delivery Address :	Order Qty :
	Order Date :

Standard Features:

Passed crash test ISO 7176 pt 19
Folding Cross Brace and Chassis
Attendant Propulsion
Choice of 11 Seat widths from 255mm (10") to 510mm (20")
Choice of 9 Seat depths from 305mm (12") to 510mm (20")
Seat heights from 270mm (13") to 470mm (18.5")
Seat tilt from 0° to 45° - lock tilt feature
Solid Castors and Rear Wheel Tyres
Swing Away and Removable Legrest Hangers
Flip Up Footplates (height adjustable)
16 Frame Colours (including Black)
Occupant Controlled Push to Lock Brakes
Anti Tip Device complete with Tipping Lever

Standard Base Price : €3,071



p Please make your selection with a cross in the box : u Standard feature m Free of Charge (NC : 9999) Non Conformance with 9999

Code	Description	NC's (Non Conformance)	Retail
NXT-IVC-30304	X NXT Folding chair with seat tilt		€3,071
NXT-99999-1	X Transport package		◆
NXT Special Basic Model - PLEASE COMPLETE DETAILS AT THE END OF THE FORM			
LOWER EXTREMITY ACCOMODATION FRAME			
NXT-30324	LEAF Seat Frame Pair		€45
NXT-30325-L	LEAF Seat Frame Left		€23
NXT-30325-R	LEAF Seat Frame Right		€23
TILT ACTIVATION			
NXT-IVC-HandL	Single lever hand activated seat tilt		€133
NXT-IVC-FootP	Single pedal foot activated seat tilt		○
NXT-30274	Hand activated seat tilt and attendant drum brake (NC: NXT-IVC-HandL, NXT-IVC-FootP) (Not available with 12" Poly Mag or spoked wheels)		€175
SEAT DIMENSIONS			
SEAT WIDTH			
NXT-30250	255mm (10") seat width - plus 50mm (2") width adjustment	(NC : NXT-36020, NXT-36021)	○
NXT-30210	280mm (11") seat width - plus 50mm (2") width adjustment	(NC : NXT-36020, NXT-36021)	○
NXT-30220	305mm (12") seat width - plus 50mm (2") width adjustment	(NC : NXT-36020, NXT-36021)	○
NXT-30230	330mm (13") seat width - plus 50mm (2") width adjustment	(NC : NXT-36020, NXT-36021)	○
NXT-30240	355mm (14") seat width - plus 50mm (2") width adjustment	(NC : NXT-36021)	○
NXT-30290	380mm (15") seat width - plus 50mm (2") width adjustment		○
NXT-30260	405mm (16") seat width - plus 50mm (2") width adjustment		○
NXT-30261	430mm (17") seat width - plus 25mm (1") width adjustment		○
NXT-30262	460mm (18") seat width - no width adjustment		○
NXT-30263	485mm (19") seat width - no width adjustment		○
NXT-30264	510mm (20") seat width - no width adjustment		○
SEAT DEPTH			
NXT-30504	305mm (12") seat depth - adjustable from 12" to 16" (Small frame)		○
NXT-30505	330mm (13") seat depth - adjustable from 12" to 16" (Small frame)		○
NXT-30506	355mm (14") seat depth - adjustable from 14" to 18" (Medium frame)		○
NXT-30521	380mm (15") seat depth - adjustable from 14" to 18" (Medium frame)		○
NXT-30522	405mm (16") seat depth - adjustable from 16" to 20" (Long frame)		○

NXT-30523		430mm (17") seat depth - adjustable from 16" to 20" (Long frame)		○
NXT-30524		460mm (18") seat depth - adjustable from 16" to 20" (Long frame)		○
NXT-30527		485mm (19") seat depth - adjustable from 16" to 20" (Long frame)		○
NXT-30528		510mm (20") seat depth - adjustable from 16" to 20" (Long frame)		○
NXT-IVC-30525		2" extension (pair)		€46
CHASSIS				
FRAME COLOURS				
NXT-B		Black		○
NXT-BB		Blizzard Blue		○
NXT-BLUE BLAST		Blue Blast		○
NXT-BS		Blue Stardust		○
NXT-CP		Cosmic Purple		○
NXT-EG		Extreme Green		○
NXT-MT		Mystic Teal		○
NXT-OO		Octane Orange		○
NXT-PB		Pacific Blue		○
NXT-PP		Passion Purple		○
NXT-PINK		Pink Parfait		○
NXT-RR		Really Red		○
NXT-RS		Red Sunset		○
NXT-SG		Sparkle Green		○
NXT-T		Twilight (Sparkly Black)		○
NXT-Y		Yellow		○
PUSH HANDLES				
NXT-30518	X	Traditional bend push handles		◆
PUSH HANDLE HEIGHT				
NXT-30515		Short - 485mm (19")		○
NXT-30514		Medium - 585mm (23")		○
NXT-30516		Tall - 660mm (26")		○
PUSH HANDLE EXTENSION OPTIONS				
NXT - 30513		Width adjustable and removable		○
NXT - 30526		Stroller handle extension - width adjustable and removable		€114
NXT - 30517		Individual stroller handle extensions - (pair) removable		€4
LEGRESTS / FOOTRESTS				
LEGREST HANGER TYPES				
NXT-IVC-30280		Pop up and swing away 75° hanger - knee to heel length 150mm - 230mm (6" - 9") (seat rail to footplate)	(NC : NXT-36407)	○
NXT-IVC-30281		Pop up and swing away 90° hanger - knee to heel length 150mm - 230mm (6" - 9") (seat rail to footplate)	(NC : NXT-36403)	€192
NXT-IVC-30282		Locking swing away 75° hanger - knee to heel length 150mm - 230mm (6" - 9") (seat rail to footplate)	(NC : NXT-36407)	€192
NXT-IVC-30283		Locking swing away 90° hanger - knee to heel length 150mm - 230mm (6" - 9") (seat rail to footplate)	(NC : NXT-36403)	€192
NXT-IVC-30279		Swing away elevating hanger - knee to heel length 150mm - 230mm (6" - 9") (seat rail to footplate)		€17
NXT-IVC-33000		No legrest option		○
HANGER EXTENSIONS				
NXT-36402		110mm (4.5") Extension - knee to heel length up to 343mm (13.5") (seat rail to footplate)		€49
NXT-36403		165mm (6.5") Extension - knee to heel length up to 393mm (15.5") (seat rail to footplate)	Not available with 90° hangers	€61
NXT-36407		High mounted clamp - knee to heel length 50mm - 125mm (2" - 5") (seat rail to footplate)	Not available with 75° hangers	€4
FOOTREST TYPES				
NXT-IVC-17070		Composite flip up footplate		○
NXT-36020		Small aluminium angle adjustable flip up footplate	For chairs with a 14" seat width and above	€139
NXT-36022		Heel loop for footplate above (NXT-36020)		€46
NXT-36021		Large aluminium angle adjustable flip up footplate complete with heel loop	For chairs with a 15" seat width and above	€139
NXT-36023		Omit heel loop for footplate above (NXT-36021)		○
NXT-36029		Aluminium, angle adjustable and removable footboard		€139
NXT-IVC-36408-L		MAPS(MULTI AXIS POSITION)LEFT (Must choose aluminum footplates) NC with composite footrest plates or footboard		€40
NXT-IVC-36408-R		MAPS(MULTI AXIS POSITION) RIGHT (Must choose aluminum footplates) NC with composite footrest plates or footboard		€40
NXT-36404		No footplate option		○
REAR WHEELS, TYRES AND ACCESSORIES				
MAG REAR WHEELS SIZE AND TYRE TYPE (Complete with double locking safety axle receiver)				
NXT-IVC-31525		Ø300mm (12") Pneumatic mag - available with 3" - 7" front castors, 13" - 18" seat to ground heights		○
NXT-31526		Ø300mm (12") Solid tyre - available with 3" - 7" front castors, 13" - 16.5" seat to ground heights		€33
NXT-31523		Ø300mm (12") Semi Pneumatic mag - available with 3" - 7" front castors, 13" - 18" seat to ground heights		€152
NXT-31530		Ø410mm (16") Pneumatic mag - available with 3" - 5" front castors, 12.5" - 16.5" seat to ground heights		○
NXT-31533		Ø460mm (18") Pneumatic mag - available with 3" - 6" front castors, 13.5" - 17.5" seat to ground heights		○
NXT-IVC-31527		Ø510mm (20") Pneumatic mag complete with aluminium handrim - available with 5" or 6" front castors, 13" - 17" seat to ground heights		€51
NXT-IVC-31528		Ø560mm (22") Pneumatic mag complete with aluminium handrim - available with 7" front castors, 15.5" - 18" seat to ground heights		€51

NXT-IVC-31529	Ø610mm (24") Pneumatic mag complete with aluminium handrim - available with 7" front castors, 17" - 18.5" seat to ground heights Only on seat depths 17" and above	€51
SPOKED REAR WHEELS SIZE AND TYRE TYPE (Complete with double locking safety axle receiver)		
NXT-31540	Ø510mm (20") Pneumatic spoke (NC : NXT-31005)	○
NXT-31541	Ø560mm (22") Pneumatic spoke (NC : NXT-31005)	○
NXT-31542	Ø610mm (24") Pneumatic spoke (NC : NXT-31005)	○
REAR WHEEL ACCESSORIES		
NXT-31005	Plastic coated handrims Only available with mag rear wheels	€108
NXT-31000	Airless inserts (to create solid tyres)	○
FRONT CASTORS		
NXT-34888	75mm x 25mm (3" x 1") Solid castor	○
NXT-40000	100mm x 25mm (4" x 1") Solid castor	○
NXT-43000	100mm x 35mm (4" x 1.4") Poly aluminium mag	€2
NXT-IVC-34555	125mm x 25mm (5" x 1" Poly) Solid castor	€107
NXT-34222	125mm x 37.5mm (5" x 1.5") Solid castor	○
NXT-42000	125mm x 35mm (5" x 1.4") Poly aluminium mag	€2
NXT-IVC-34666	150mm x 25mm (6" x 1") Solid castor	€137
NXT-34999	150mm x 37.5mm (6" x 1.5") Solid castor	○
NXT-44000	150mm x 35mm (6" x 1.4") Poly aluminium mag	€2
NXT-34444	150mm x 50mm (6" x 2") Poly with tread	€2
NXT-34777	180mm x 25mm (7" x 1") Solid castor	○
NXT-34000	180mm x 37.5mm (7" x 1.5") Solid castor	€107
NXT-34200	180mm x 50mm (7" x 2") Poly	€107
SAFETY		
BRAKE TYPE		
NXT-31501	Push to lock brakes	◆
NXT-IVC-31506	Pull to lock brakes	€192
NXT-IVC-31801	Brake handle extension	€192
	(For hand activated seat tilt and attendant drum brake see option NXT-30274)	
ARMRESTS		
BACK POST MOUNTED ARMREST TYPES		
NXT-IVC-38204	Flip back standard 330mm (13" long) (NC : NXT-IVC-38214)	○
NXT-IVC-38202	Flip back long 405mm (16" long)	€146
NXT-38208	Flip back telescopic	€111
NXT-38201	Locking flip back standard 330mm (13" long) (NC : NXT-IVC-38214)	€34
NXT-38200	Locking flip back long 405mm (16" long)	€34
NXT-38207	Locking, flip back telescopic	€34
NXT-38209	No flip back armrest option	○
SEAT FRAME MOUNTED ARMREST TYPES		
NXT-38203	Single post, height adjustable and removable	€17
ARMREST PAD		
NXT-IVC-38213	Desk length	○
NXT-IVC-38214	Full length Only available on 16" armrests	€20
BELTS & HARNESSSES		
IVC-NXT-6B	Strap Risers Pair	€33
SFDI3 -S	Four Point Belt - Small	€6
SFDI3 -M	Four Point Belt - Medium	€6
SFDI3 -L	Four Point Belt - Large	€6
SB-3331-S	Centre pull/pls Push Button - Small	€1
SB-3331-M	Centre pull/pls Push Button - Medium	€1
SB-3331-L	Centre pull/pls Push Button - Large	€1
SFDI-6NNZ-S	Anti Trunk Support Narrow W/O Zip - Small	€126
SFDI-6NNZ-M	Anti Trunk Support Narrow W/O Zip - Medium	€126
SFDI-6NNZ-L	Anti Trunk Support Narrow W/O Zip - Large	€126
SFDI-6SNZ-S	Anti Trunk Support W/O ZIP - Small	€121
SFDI-6SNZ-M	Anti Trunk Support W/O ZIP - Medium	€121
SFDI-6SNZ-L	Anti Trunk Support W/O ZIP - Large	€121
SFDI-6SZ-S	Anti Trunk Support W/ZIP - Small	€126
SFDI-6SZ-M	Anti Trunk Support W/ZIP - Medium	€126
SFDI-6SZ-L	Anti Trunk Support W/ZIP - Large	€126
SFDI-9BB-S	Hi-Tops/Buckle -Small	€13
SFDI-9BB-M	Hi-Tops/Buckle -Medium	€13
SFDI-9BB-L	Hi-Tops/Buckle -Large	€13
SFDI-9BV-S	Hi-Tops/Velcro-Small	€13
SFDI-9BV-M	Hi-Tops/Velcro-Medium	€13
SFDI-9BV-L	Hi-Tops/Velcro-Large	€13
SAFETY OPTIONS		
NXT-31600	Rear anti tipper device, complete with stepping lever	○

ACCESSORIES		
NXT-TFBP	Back pack (black only)	€49
SEAT PAN OPTION		
NXT-IVC-10909	SEAT PAN - REMOVABLE W/ EZLITE	€140
TOTAL PRICE:		

Non Conforming Options:	
NXT-36403	165mm (6.5") Extension - knee to heel length up to 393mm (15.5") (seat rail to footplate)
NXT-36407	High mounted clamp - knee to heel length 50mm - 125mm (2" - 5") (seat rail to footplate)
NXT-36020	Small aluminium angle adjustable flip up footplate
NXT-36021	Large aluminium angle adjustable flip up footplate complete with heel loop
NXT-31005	Plastic coated handrims
NXT-IVC-38214	Full length

Special Instructions - Extended lead times will apply

Name:

Email:

Phone:



Description of special request (please include relevant option code):