

# ***Flo-tech<sup>®</sup> Lite Visco***




Low profile visco elastic low  
memory foam cushion offering  
exceptional comfort



Yes, you can.

The *Invacare Flo-tech Lite Visco* cushion is made of high-density visco elastic foam, allowing it to conform to the user's shape. This affords exceptional levels of comfort and increases pressure redistribution. This low profile cushion format is ideal for both chair and wheelchair use, facilitating side transfers.



## Features and Options



### Low memory foam designed to maximise body contact area

The visco-elastic foam conforms to the user's body, maximising the contact area, delivering exceptional levels of comfort and optimising pressure redistribution.

### Ingress Resistant Cover

Cover is two-way stretch, water-resistant and vapour-permeable. Cover features an anti-slip base with carry handle.

### Available in 61 sizes

The *Invacare Flo-tech Lite Visco* is available in a range of 61 sizes, from 250 x 250 to 480 x 530 mm. (See size chart for details).

### Three foam densities available






Available in three densities to accommodate a range of user weights – Soft (up to 32 kg), Medium (up to 51 kg) and Firm (up to 100 kg).

### Contoured seating surface with pommel and side support

Offers mild lateral support and pelvic stability to enhance the user's posture.

## Technical data


For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

Various sizes	Height	Total product weight	Max. user weight <sup>1</sup>	Max. cleaning temperature
				
A range of 61 sizes	78 mm	1.26 kg	100 kg	80°

Manufactured to European Standards:  
Interior Foam: BS EN 1021-1, BS EN 1021-2,  
BS 5852: 1982 Crib 5

1. Guide maximum user weight limits are set assuming that the cushion selected is of appropriate width and depth for the intended user/s.

## Colours



T: +44 1656 776222  
uk@invacare.com  
www.invacare.co.uk  
www.invacare.ie

Invacare UK & Ireland  
Unit 4 Pencoed  
Technology Park  
Bridgend CF35 5AQ



**Yes, you can.®**