

# *Proportional Chin Control*



Ingenious design for a perfect fit



Thanks to the new Proportional Joystick (ASL 138), the Invacare chin control has been designed to be more compact and refined. The reduced size of the joystick has allowed this control to fit perfectly on Modulite seating and **Ultra Low Maxx** power chairs. The sleek minimal design of the *Invacare* chin control is available with a manual or electric swing away function, to the left or to the right to suit individual constraints or needs.

#### Electric swing away chin control

- Horizontal electric swing away
- Allows the individual to easily adjust the height of the joystick, enabling them to adapt to their seating position throughout the day
- Retractable even when the chair is off
- Precise and discreet control

#### Manual swing away chin control

- Horizontal manual swing away
- Stable and robust
- The arm can be easily retracted by a carer



#### Egg switch and headrest included

The Egg switch allows access to the different functions of the powerchair. These are pre-programmed but can be changed. The choice of Matrix Elan headrest pads provides exceptional support and comfort.



#### Ideal proportional joystick

Exterior coating resistant to saliva. Deflection force: 100g.

Invacare® **Modulite™**

**ULTRALow maxx**  
By Motion Concepts

#### Easy to order

Available as a standard option on our range of power wheelchairs equipped with Ultra Low Maxx and Modulite seats.

For more information about this product, please visit [www.invacarelinx.co.uk](http://www.invacarelinx.co.uk)