



# Invacare Jasmine

## Product Data Sheet

---

### Product Range Available

Code: CJ1541963 - Description: Jasmine mobile hoist with 2 point 45cm spreader bar

Code: CJ1541965 – Description: Jasmine mobile hoist with 4 point 45cm spreader bar

All products are delivered with battery, charger lead and 2 point spreader bar as standard.  
Control boxes have Linak Jumbo Care™ functionality.

### Key Features

- Aluminium construction
- Ergonomic push bar
- 200kg (32 stone) SWL
- Wide base
- High lift
- Covered leg actuators
- Externally tested and verified to ISO 10535
- Electrical leg operation
- Warranty 5 years\*

### Technical Specification

- 200g (32 stone) SWL
- Overall product weight - 43kg
- Castors 120mm (2 braked castors at rear)
- Maximum leg opening – 1050 mm (internal)
- Closed leg dimension – 565 mm (internal)
- Closed leg dimension – 690 mm (external)
- Lift height range – 2120 mm
- Base length – 1340 mm
- Turning radius – 1530 mm
- Number of cycles per full charge – 40
- Battery – 2.9ah
- Emergency stop
- Emergency lower –Manual and electrical



- **Accessories:**

- 2 point and 4 point spreader bars in 550mm widths
- Homecare weigh scale with adaptor
- Jumbo Care service interface diagnostic tool
- Wide range of slings
- Compatible with other manufacturers slings - for sling compatibility see bulletin PU661

**Supporting Documentation**

- Jasmine range brochure
- Jasmine user manual
- Weigh scale log book
- Jumbo Care installation manual
- Jumbo Care service data tool
- Sling and Hoist Assessment Guidelines
- Sling Poster and A5 filo-fax guide

All above documentation can be found on [www.invacare.co.uk](http://www.invacare.co.uk)

Invacare's Business Development Managers can provide in field support and training on hoists and slings

\*5 years on control box, 3 years on frame and 1 year on battery and handset

---

**Manufacturers Contact Details**

Invacare Ltd, Unit 4 Pencoed Technology Park, Bridgend, CF35 5AQ

**Tel:** 01656 776222 **E-mail:** [UK@invacare.com](mailto:UK@invacare.com)

**Web:** [www.invacare.co.uk](http://www.invacare.co.uk) / [www.thinkindependence.co.uk](http://www.thinkindependence.co.uk)



safepatienthandling