



Yes, you can.®

## > *küschall*® *K-Junior*®

**Achieve your dreams in lightness and design**

*K-Junior*® is the new rigid wheelchair for children and teenies designed by *küschall*®. Its lightness and design makes the *K-Junior*® an attractive wheelchair for the trendy and active youth of today.



  
**JUNIOR**



High-mounted and retracted footplate



"Buggy" pushhandle




Seat length increase by 2.5 cm in one easy handling operation



Smiley and Teddy footplates

### Technical data




24 - 36 cm  
(in steps of 2 cm)



25 - 40 cm  
(in steps of 2,5 cm)



front: 38 - 45 cm  
rear: 33 - 47 cm  
(in steps of 1 cm)



27 - 48 cm  
(in steps of 1,5 cm)




15 - 29 cm  
(high mounted)  
30 - 41 cm  
(in steps of 1 cm)



variable,  
standard 4 cm




-12° / -8° / -4° /  
0° / 4° / 8° / 12°



SW +21cm (at 3°)  
(SW +17cm at 0°)



approx. 76 cm



approx. 9,2 kg



max. 60 kg




approx. 6 kg  
(without rear wheels)


### Frame colours




Glossy black



Royal red




Sapphire blue



Xanthic yellow




Cyber green




Silver grey


### CPS \* (Colour Power Set)




Techno anthracite




Apricot



Violet



B-blue



Lemon

### Upholstery



Nylon Black

### Enjoy the freedom thanks to best performance and stability

The **küschall® K-Junior's®** main purpose is to accompany the children of today through childhood and provide them with the right tools to express themselves in their very active daily lives.



### The perfect combination of lightness and quality

The lightness of the **K-Junior's®** (9,2kg in standard configuration) enables best handling and driving performance.

Keeping in mind that today's kids want to look and be cool, the trendy and dynamic design was created by keeping visible parts to a minimum.

\* CPS includes Starec hubs, 24" tyre rims, back angle adjustment plate, forks as well as castor housing and seat module caps.

All information quoted is believed to be correct at time of print. Invacare® reserves the right to alter product specifications without prior consultation.



INVACARE® Deutschland GmbH  
Kleiststraße 49 · 32457 Porta Westfalica  
Phone + 49 (0) 57 31 / 7 54-0 · Fax + 49 (0) 57 31 / 7 54-5 21 91  
E-mail info@invacare.de · [www.invacare.de](http://www.invacare.de)



K-Junior - EU - 06/06  
© küschall® AG, 2006. All Rights Reserved.