



Küschall®
ELEMENT

Ahead of
the curve



Küschall.com

Limited edition



IN COLLABORATION
WITH BREITFUSS

ELEMENT

An exclusive design, the new Küschall Element limited edition is a refined active wheelchair, designed to showcase the passion of Küschall.

Every curve and contour highlights our commitment to achieving excellence. Each handcrafted chair is a statement of performance, precision and personal expression.

In collaboration with Breitfuss - innovators in wood-based design, the Küschall Element is born from combined expertise and the vision to create something truly exceptional.

BOLD AND
DISTINCTIVE

NATURALLY
EXCEPTIONAL





Built for strength

Crafted from high quality natural European Ash wood – a trusted material for modern performance, the Küschall Element is precisely engineered for resilience and durability.



— Distinctly original

Bold and beautiful, the Küschall Element delivers active mobility with natural elegance. A chair that elevates confidence, celebrates individuality and makes every outing an occasion.





Precise geometry

An optimised wheelbase and seat position gives comfort, effortless manoeuvrability and control. The Küschall Element has a thoughtful design with a focus on performance.



A closer look

Easier handling

Two integrated handles on back and seat.

Large choice of handrims

Choose from ergonomic and lightweight handrims.

Performance wheel

24" Spinergy CLX as standard or optional Corima carbon with Schwalbe Marathon Plus tyres.

Breathable back support

Padded back support in black mesh fabric.

Secured positioning

A range of Bodypoint belts to suit individual needs.

Lightweight active brakes

Easy to operate with a clean, minimalistic design.

High strength core

Natural Ash frame, built for resilience and tested to endure.

Carbon footrest

Durable and lightweight.

[CLICK HERE TO WATCH THE VIDEO >>](#)



Technical specification

SEAT WIDTH

340 – 420 mm



SEAT DEPTH

410 – 460 mm



KNEE-TO-HEEL LENGTH

390 – 460 mm



BACKREST HEIGHT

300 – 360 mm



BACKREST ANGLE

85°



SEAT HEIGHT

Front: 480 – 500 mm



SEAT ANGLE

8° / 12°



FRAME ANGLE

90°



TOTAL WIDTH

655 mm for SW 420 mm



TOTAL LENGTH

865 mm for SD 460 mm



CENTER OF GRAVITY

95mm
(seat adjustable by
-15 mm / +15 mm)



STATIC STABILITY UPHILL

4.8° (with 12° seat angle)
9.1° (with 8° seat angle)



TOTAL WEIGHT UNLOADED

Approx. 12.9 kg
(with standard back
& seat cushion)



TRANSPORT WEIGHT

Approx. 8.2 kg
(without rear wheel,
back & seat cushion)



MAX. USER WEIGHT

100 kg



All measurements and weights indicated are based on seat width 400 mm and seat depth 400 mm in the lightest and most compact configuration.
Weight and measurements can change with the different configurations. For further information please check the user manual on our website.

Tel: +41 61 487 70 70
info@kuschall.com
www.kuschall.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland



©2025 Invacare International GmbH. All rights reserved. Kuschall® Element™ designed in collaboration with Breitfuss. EU_17122025

Kuschall®
UNLIMIT YOUR WORLD