

Tel: 01 810 7084 Fax: 01 810 7085
Email: ordersireland@invacare.com
Online Spares available at www.invacare.ie

JANUARY 2020

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :



Kite® (Modulite® Seating System)

Max User Weight 160Kg (25 stone)

Ireland

Standard Features



Invacare® Modulite®
An evolution in seating designed around you.



- Telescopic Seat Frame and Plate with Side Rails width 330mm to 530mm, (13" - 21") with an additional 50mm, 2" by arms
- Seat Depths 410mm - 510mm (16" - 20")
- Tension Adjustable Backrest cushion
- Standard T-Armrests, detachable with Short armpad
- Manual settable seat Tilt (0° - 12°)
- Manual settable backrest Recline (90° to 120°)
- Standard 80° hangers with fixed footplate
- LiNX Control REM110 with standard bracket
- Black or Arctic White Frame colour
- 7 Body Shroud colours, 7 side shroud colours and 3 Rim colours
- 250mm (10") Puncture proof PU-foam sealed, Black castors
- 350mm (14") Puncture proof solid inlay Black drive wheels
- Dual Swing Technology
- 6 km/h (4mph) 4 Pole Motors and 60Ah Gel Batteries
- 8A (200 - 250V) Charger
- Anti-tippers
- 4-Point tie-down system
- Lapbelt with autolock (adjustable both sides)



*Configurations may vary from the image shown

Standard Base Price €7,166

	X	Please make your selection with a cross in the box	◆ Standard feature	○ Free of Charge
Customer weight =		Kg		
Code	Description	Ireland		
Weight Range - Must Be selected to process order correctly				
AVB1920	User weight up to 125kg			○
AVB1921	User weight 126kg to 136kg			○
AVB1922	User weight 137kg to 150kg	Not available with Neckrest (1350), headrest sideways adjustable (1432), Headrest with cheek support (1427) with 30° COG tilt (1093) or seat lifter (1094)		
AVB1923	User weight 151kg to 160kg	Not available with Neckrest (1350), headrest sideways adjustable (1432), Headrest with cheek support (1427)with 30° COG tilt (1093)or seat lifter (1094)		
MODELS				
AVB0010	Kite Basic Model			◆
AVZ0010	Kite Special Basic Model PLEASE COMPLETE DETAILS AT THE END OF THE FORM			POA
SEATING TYPE				
AVB2100	Modulite Seating Unit			€7,166
AVZ2101	Modulite HD Seating Unit - 530 to 550mm seat widths	Must be selected for Seat Width 550mm, Only Available with Tension Adjustable Backrest and Tall Push handles. Black Nylon Fabric Only.		
				€7,585

SEAT FRAME AND SEAT PLATE					
AVB2111		Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)	◆		
AVZ2102		Telescopic seat frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest) 530mm to 550mm	○		
MODULITE SEAT SIZE SETTING CHART:					
					
		SEAT WIDTHS		SEAT WIDTHS	
		330mm* (13")	AVB2121 380mm (15")	AVB2122 430mm (17")	AVB2123 480mm (19")
SEAT DEPTHS	Seat Rail 410mm - 510mm	AVB2132 410mm (16")	○	○	○
		AVB2134 460mm (18")	○	○	○
		AVB2136 510mm (20")	○	○	○
		HD SEATING		SEAT WIDTHS	
SEAT DEPTHS	AVB2136 510mm (20")		○	○	
	AVZ2137 580mm (23")		○	○	
<p>*330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)</p> <ol style="list-style-type: none"> 1. Seat depth 410mm: cannot combine T-Armrests (2191) with flexible lapbelt mount 2. Seat width 530mm (2124) not possible in combination with following armrests (2192/3) 3. Modulite Teen seating available, seat width and depth from 305mm (12") 					
SEAT OPTION					
AVBSWO		"H" Frame Interface Kit	€460		
SEAT CUSHION (sizes according to seat frame settings)					
AVB2142		Matrix Flo-tech Lite	Not available with HD seating		
AVB2143		Matrix Flo-tech Contour Visco	Not available with 510mm seat depth, 380mm or 530mm seat widths, black spacer fabric, with HD seating		
AVB2147		Matrix Libra 90mm thick	Not available with combination 530mm seat with and 410mm seat depth, seat depth 560mm only available with 530mm seat width, black spacer fabric, with HD seating		
AVB2148		Seat Cushion (90mm high) (required for HD seating Only)	○		
AVB2140		Without seat cushion	Not available with HD seating		
BACKREST SUPPORT					
AVB2163		Backposts 480mm (Short)	Not available with flex3 adjustable backrest		
AVB2162		Backposts 540mm (Tall) (required for HD seating)	Must be selected for flex3 adjustable backrest		
PUSH HANDLES FOR BACKPOSTS					
AVB2107		Backposts WITH push handles (required for HD seating)	◆		
AVB2108		Backposts WITHOUT push handles, Push handles will be cut down on the vertical section.	Not available with attendant controls, dual controls, matrix backrest or tension adjustable backrest in combination with a headrest.		
BACKREST UPHOLSTERY / CUSHION					
AVB2171		Basic Tension Adjustable, (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear, Black fabric only	Not available with laterals, chin remote, headrest 1496, (required for HD seating) FOC with HD.		
AVB2179		Flex3 Tension Adjustable (540mm - 480mm) Black Dartex only	Not available with flip up or following armrests, chin control, soft armpads		
AVB2172		Flex3 Adjustable Backrest Plates (without Backrest Cushion)	○		
AVB2173		Flex3 Adjustable Backrest Plates with Laguna Backrest cushion	◆		
AVB2174		Matrix Posture Back height, 400mm (Mesh fabric)	Not available with laterals or 530mm (21") seat width, headrest option 1496		
AVB2175		Matrix Posture Back, height 500mm (Mesh fabric)	Not available with laterals or 530mm (21") seat width, headrest option 1496		
AVB2176		Matrix Elite Back, height 500mm (Mesh fabric)	Not available with laterals or 530mm (21") seat width, headrest option 1496, seat width of 380mm (15") 2121 in combination with flip up or following armrests		
AVB2177		Matrix Elite TR Back, height 500mm (Mesh fabric)	Not available with laterals or 530mm (21") seat width, headrest option 1496		
AVB2178		Matrix Elite Deep Back, height 500mm (Mesh fabric)	Not available with laterals or headrest option 1496		
AVB2170		Without Backrest Upholstery (Straight 1" diameter Back Posts/Back Frame only suitable for mounting specialist backrests)	Not available with laterals or headrest option 1496		
UPHOLSTERY TYPE					
AVB2151		Black fabric	Not available with HD seating		
AVB2152		Black Dartex	Not available with HD seating		
AVB2155		Black Spacer fabric	Not available with HD seating		
AVB2156		Black Nylon (required for HD seating Only)	○		
					

ARMREST (Pair)			
AVB2191	T-Armrest Standard Range of adjustment from seat plate 245mm to 315mm		◆
AVB2181	T-Armrest Extended height Range of adjustment from seat plate 295mm to 365mm		€46
AVB2183	In-rigged Armrests (reduces seat width by 25mm each side)	Not available with hip support, with HD seating	○
AVB2193	Flip-up Armrest (height and depth adjustable) (height adjustable from 230mm - 360mm)	Not available with laterals, in combination with 380mm seat width and matrix elite back, with HD seating	€51
AVB2192	Following Armrest (height adjustable from 230mm - 300mm)	Not available with laterals, in combination with 530mm seat plate and flex3 adjustable backrest plates, with HD seating	€149
AVB2194	Following Armrest (height adjustable from 300mm - 360mm)	Not available with laterals, in combination with 530mm seat plate and flex3 adjustable backrest plates, with HD seating	€149


ARMPAD (Pair)			
AVB2195	Short Armpad (325mm)		◆
AVB2196	Long Armpad (375mm)		○
AVB2197	Long Ergonomic Armpad (397mm)		○
AVB2199	Soft Armpad (REA Azalea Wide Pad)	Not available with HD seating	€98

HIP SUPPORTS, Multi adjustable LH (for flip up and following armrests only)			
AVB2184	SMALL 140mm x 100mm	Not available with HD seating	€259
AVB2186	MEDIUM 230mm x 100mm	Not available with HD seating	€259
AVB2188	LARGE 330mm x 100mm	Not available with HD seating	€259

HIP SUPPORTS, Multi adjustable RH (for flip up and following armrests only)			
AVB2185	SMALL 140mm x 100mm	Not available with HD seating	€259
AVB2187	MEDIUM 230mm x 100mm	Not available with HD seating	€259
AVB2189	LARGE 330mm x 100mm	Not available with HD seating	€259

EXTENDED HEIGHT (+2") MOUNTING HARDWARE FOR HIP SUPPORTS			
AVB2165	Extended height mounting hardware for hip supports (8.1" / 206mm high)		○

SEAT ACCESSORIES



HEADREST / NECKREST (Dartex) - Adequate headrest option must be selected

AVB1350	Neckrest: height and depth adjustable		€81
AVB1432	Headrest: height, depth and sideways adjustable		€198
AVB1427	Headrest with cheek support: height and depth adjustable		€226
AVB2323	Elan Standard Pad 150mm (6") wide with standard cover complete with Elan Headrest Hardware	Not available with tension adjustable backrests, head controls	€369
AVB2324	Elan Standard Pad 255mm (10") wide with standard cover complete with Elan Headrest Hardware	Not available with tension adjustable backrests, head controls	€369
AVB2322	Elan Standard Pad 355mm (14") wide with standard cover complete with Elan Headrest Hardware	Not available with tension adjustable backrests, head controls	€369
AVB2325	Elan Occipital Pad 230mm x 180mm (9" x 7") with standard cover complete with Elan Headrest Hardware	Not available with tension adjustable backrests, head controls	€369
AVB2326	Elan Occipital Pad 305mm x 200mm (12" x 8") with standard cover complete with Elan Headrest Hardware	Not available with tension adjustable backrests, head controls	€369
AVB2327	Elan 4 Point Headrest Pad 280mm x 255mm (11" x 10") with standard cover, complete with Elan Headrest Hardware	Not available with tension adjustable backrests, head controls	€497
AVB1496	All bracketry as option 1432 without Headrest cushion		€92
AVB1495	Push bar dia. 22mm in conjunction with telescopic seat frame / backposts, No Headrest - no mounting hardware.		○





LATERAL SUPPORT (PAIR)			
AVB1446	Rigid pad (Azalea type, Dartex)	€2,048	Not available with flex3 tension backrest, with HD seating
AVB1445	Multi-adjustable pad (Azalea type, Dartex)		Not available with HD seating






SEAT ADJUSTMENTS

SEAT LIFTER / SEAT TILT						
Modulite	Tilt Range	Seat Height				
		AVB0337 440mm	AVB0331 460mm	AVB0339 480mm	AVB0381 500mm	
AVB1091	Manually settable	0° - 12°				◆
AVB1092	Powered basic	0° - 20°				€452
AVB1093	Powered centre of gravity shift	0° - 30°				€478
AVB1094	Seat Lifter with electric tilt (COG Shift)	0° - 30°				€2,048

HD	Tilt Range	Seat Height				
		AVB0337 440mm	AVB0331 460mm	AVB0339 480mm	AVB0381 500mm	
AVB1091	Manually settable Seat Tilt	0° - 12°				◆
AVB1095	Powered Tilt	0° - 20°				€452
AVB1093	Powered Tilt with (COG shift)	0° - 30°				€478
AVB1094	Seat Lifter with electric Tilt (COG Shift)	0° - 30°				€2,048

BACKREST RECLINER				
		Recline Range		Selection
	AVB0975	Manually settable recline	90° - 120°	◆
	AVB0977	Powered Recliner with quick release for easy folding of backrest	90° - 120°	€318

LEGREST

	80° standard 0526		REA Style angle adjustable° 0529		VARI F 0503		VARI A 0533		ADE 0531
---	----------------------	---	-------------------------------------	---	----------------	--	----------------	---	-------------

FOOTREST HANGER TYPE / PAIRS







AVB0526	Standard 80° Hangers	Only with fixed footplate, adapt footplate, with HD seating	◆
AVB0529	Hanger, manually elevating	Only with fixed footplate, adapt footplate, with HD seating	€168
AVB0503	VARI F Angle adjustable with tools 0° - 70°, length 290mm - 460mm		€44
AVB0533	VARI A : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	Not available with fixed footplate	€118
AVB0531	ADE Powered : Angle adjustable powered by joystick, 0° - 80° length, 290mm - 460mm	Not available with fixed footplate or footboard	€653
AVB0501	Without Hanger		○

FOOTPLATE TYPE (Pair)

AVB0512	Fixed Footplate pair		◆
AVB0504	Angle adjustable Footplate pair		€29
AVB0513	Angle and depth adjustable Footplate pair (Adapt)		€118
AVB0522	Angle adjustable Footboard		€152



LINX PRIMARY CONTROLS (SECTION 1)













	Drive Only		Drive & Actuators		Drive, Actuators & Lights		Touch-screen & Joystick		REM400 with Toggle Switch Module		Touch-screen
REM110		REM211		REM216		REM400				REM500	




LINX PRIMARY REMOTE & BRACKET

AVB2201	REM110 Drive only	Not available with lights	◆
AVB2203	REM211 Drive & Actuators	Not available with lights	○
AVB2204	REM216 Drive, Actuators & Lights		€147
AVB2205	REM400, touch-screen interface (Built in Bluetooth HID)		€449
AVB2270	Assignable Toggle Switch Module for REM400 (as default set to left ON / OFF, right speed dial)		€263
AVB1207	Standard Bracket RH		○
AVB1208	Standard Bracket LH		○
AVB1707	Swing Away RH		€144
AVB1708	Swing Away LH		€144
AVB2213	Centre mount, swing away to RH, Nucleus NEW	Not available with a tray	€670
AVB2214	Centre mount, swing away to LH, Nucleus NEW	Not available with a tray	€670
AVB2227	Bracket RH mounted for Attendant		€65
AVB2228	Bracket LH mounted for Attendant		€65
AVZ2239	Individual mounting option per request (Special) - Requires additional information / specification		POA




LINX REM500 REMOTE & BRACKET (Secondary User input required)

AVB2206	REM500, without joystick for alternative access (Built in Bluetooth HID)		€825
AVB2209	Swing away RH, multi-adjustable		○
AVB2210	Swing away LH, multi-adjustable		○
AVB2213	Centre mount, swing away to RH, Nucleus	Not available with a tray	€670
AVB2214	Centre mount, swing away to LH, Nucleus	Not available with a tray	€670
AVZ2239	Individual mounting option per request (Special) - Requires additional information / specification		POA

	Dual Control		Compact Dual Control		LiNX IDC Dual Control
REM050		ACU200			
LiNX DUAL CONTROL & BRACKET: CAN ONLY BE SELECTED IN ADDITION TO LiNX USER CONTROL					
AVB2220	Dual Control REM050 as secondary remote				€444
AVB2219	Compact Dual Control (ACU200)				€423
AVB1207	Standard Bracket RH				○
AVB1208	Standard Bracket LH				○
INVACARE LiNX IDC DUAL CONTROL & BRACKET					
AVB2245	Invacare IDC Dual control				€678
AVB1221	IDC Dual control Detachable Bracket RH				○
AVB1222	IDC Dual control Detachable Bracket LH				○
If no specialist controls are required please move to the frame colour section on page 7					
LiNX SPECIALIST CONTROLS (SECTION 2)					
	Compact Remote 2231		Compact Remote Low Force 2232		Chin Joystick 2229
CR-400		CR-400LF			Head Control (ATOM)
					Head Control (PROTON)
LiNX - COMPACT REMOTES & BRACKET (Requires Primary Remote REM400 or REM500)					
AVB2231	Proportional Compact Remote CR400				€767
AVB2232	Proportional Compact Remote Low Force CR-400LF				€1,088
AVB2207	RH Mount				○
AVB2208	LH Mount				○
AVB2213	Centre Mount swing away to RH (Nucleus)			Not available with a tray	€670
AVB2214	Centre Mount swing away to LH (Nucleus)			Not available with a tray	€670
AVZ2215	Individual mounting option per request (Special) - Requires additional information / specification				POA
LiNX CHIN CONTROL & BRACKET (Requires Primary Remote REM400 or REM500)					
AVB2229	Proportional Chin Joystick (ASL138)			Not available with HD seating	€2,817
AVB1289	Swing Away Bracket RH side manual (ASL628) (Requires ELAN Headrest)				○
AVB1288	Swing Away Bracket LH side manual (ASL628) (Requires ELAN Headrest)				○
AVB1287	Swing Away Bracket RH side powered,			Has to be mounted on opposite side of attendant bracket, cane holder, mirror or tray	€2,268
AVB1286	Swing Away Bracket LH side powered,			Has to be mounted on opposite side of attendant bracket, cane holder, mirror or tray	€2,268
HEAD CONTROL (with ATOM technology) (Requires Primary Remote REM400 or REM500)					
AVB2230	Head Control with ultrasonic proximity sensors (ATOM, ASL104) with standard wings M3				€4,079
AVB2247	Head Control with ultrasonic proximity sensors (PROTON, ASL104P) with PROTON wings				€4,760
FOUR-SWITCH PROXIMITY ARRAY (with ATOM technology)					
AVZ2257	Four-Switch Proximity Array (ASL106)				€7,440
AVZ2275	Centre Mount, swing-away to RH, Nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)				€815
AVZ2276	Centre Mount, swing-away to LH, Nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)				€815
AVZ2215	Individual mounting option per request (Special) - Requires additional information / specification				POA
WIRELESS ACCESSORIES FOR ATOM					
AVB2249	Wireless Triple Switch Receiver (ASL557-3) (Sub D plug)				€383
AVB2250	Wireless Mouse Emulator (ASL558)				€919
	Micro Extremity Control		Heavy Duty Remote Joystick		10-Way Switch
					LiNX USB Charger Port

MICRO EXTREMITY CONTROL (Requires Primary Remote REM400 or REM500)		
AVB2233	MEC Proportional Joystick (Micro Extremity Control)	€3,837
AVB2235	Tray Mount RH (ASL603, 100mm x 150mm)	○
AVB2236	Tray Mount LH (ASL603, 100mm x 150mm)	○
AVB2213	Centre Mount, swing away to RH, Nucleus	€670
AVB2214	Centre Mount, swing away to LH, Nucleus	€670
AVB2241	Centre Mount, swing away to RH, Nucleus (incl. Tray ASL603, 100mm x 150mm)	€670
AVB2242	Centre Mount, swing away to LH, Nucleus (incl. Tray ASL603, 100mm x 150mm)	€670
AVZ2215	Individual mounting option per request (Special) - REQUIRES ADDITIONAL INFORMATION / SPECIFICATION	POA
HEAVY DUTY REMOTE JOYSTICK (Requires Primary Remote REM400 or REM500)		
AVB2240	Proportional Heavy Duty Remote Joystick with mode switch (ASL133)	€2,433
AVB2235	Tray Mount RH (ASL603, 100mm x 150mm)	○
AVB2236	Tray Mount LH (ASL603, 100mm x 150mm)	○
AVB2213	Centre Mount, swing away to RH, Nucleus	€670
AVB2214	Centre Mount, swing away to LH, Nucleus	€670
AVB2241	Centre Mount, swing away to RH, Nucleus (incl. Tray ASL603, 100mm x 150mm)	€670
AVB2242	Centre Mount, swing away to LH, Nucleus (incl. Tray ASL603, 100mm x 150mm)	€670
AVZ2215	Individual mounting option per request (Special) - REQUIRES ADDITIONAL INFORMATION / SPECIFICATION	POA
LINX 10-WAY SWITCH & BRACKET		
AVB2251	10-Way Switch DLX-SW10-A for direct Powered Seating Access (includes loom to connect + additional jack for mode switch)	€489
AVB2217	Bracket RH, mounted for user	○
AVB2218	Bracket LH, mounted for user	○
AVB2227	Bracket RH, mounted for attendant	○
AVB2228	Bracket LH, mounted for attendant	○
LINX ACCESSORIES		
AVB2252	Loom providing jack sockets for Drive and Mode (includes connection for 10-way Switch)	€100
LINX USB CHARGER PORT & BRACKET		
AVB2253	USB Charger Port DLX-USB02-A	€147
AVB2207	RH Mount	○
AVB2208	LH Mount	○
FURTHER LINX COMPONENTS / ACCESSORIES (SECTION 3)		
AVB2300	To be enclosed only - not mounted to the chair	○
AVZ2301	To be connected and programmed (Special)	POA
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">  <p>Buddy Button</p> </div> <div style="text-align: center;">  <p>Egg Switch</p> </div> <div style="text-align: center;">  <p>Piko Button</p> </div> </div>		
BUDDY BUTTONS (Accessory loom 2252 required if no REM400 or REM500 is selected)		
AVB1685	Red Buddy Button	€82
AVB1686	Yellow Buddy Button	€82
AVB1687	Green Buddy Button	€82
AVB1690	Blue Buddy Button	€82
EGG SWITCH (Accessory loom 2252 required if no REM400 or REM500 is selected)		
AVB2254	Egg Switch, black	€82
PIKO BUTTONS (Accessory loom 2252 required if no REM400 or REM500 is selected)		
AVB1711	Red Piko Button	€146
AVB1712	Yellow Piko Button	€146
AVB1713	Green Piko Button	€146
AVB1714	Orange Piko Button	€146
AVB1715	Blue Piko Button	€146
AVB1716	Black Piko Button	€146
TOGGLES / BUTTONS		
AVZ2261	Stereo Toggle - direct control for Tilt only	€248
AVZ2262	Stereo Button - direct control for Tilt only	€237
AVZ2263	Four-Way Toggle - direct control for Tilt & Recline	€310
AVZ2264	Four-Way Button - direct control for Tilt & Recline	€331
AVZ2207	RH Mount	○
AVZ2208	LH Mount	○
AVZ2225	Without Mounting	○
REMOTE STOP SWITCH		
AVZ2258	Remote Stop Switch (ASL504)	€970

NEW 3" U-Shape		NEW 4" U-Shape		NEW Rubber dome		NEW Grooved Mushroom	
BODYPOINT JOYSTICK KNOBS							
AVB2721		3" wide U-Shaped, Joystick Handle					€68
AVB2722		4" wide U-Shaped, Joystick Handle					€68
AVB2723		Wide Rubber Dome, Joystick Handle					€68
AVB2724		Grooved Mushroom Joystick Handle					€68
DRIVE UNIT OPTIONS							
FRAME COLOUR							
AVB0228		Frame colour Black					<input type="radio"/>
AVB0252		Frame colour Arctic White					<input type="radio"/>
BODY SHROUD COLOURS							
AVB0242		Body shroud colour Pearl Night Blue					<input type="radio"/>
AVB0243		Body shroud colour Metallic Silver					<input type="radio"/>
AVB0224		Body shroud colour Black					<input checked="" type="radio"/>
AVB0204		Body shroud colour Transparent Cherry Red					<input type="radio"/>
AVB0222		Body shroud colour Arctic White					<input type="radio"/>
AVB0248		Body shroud colour Grasshopper Green					<input type="radio"/>
AVB0279		Body shroud Sunshine Yellow					<input type="radio"/>
SIDE SHROUD COLOURS							
AVB0284		Side shroud colour Pearl Night Blue					<input type="radio"/>
AVB0285		Side shroud colour Metallic Silver					<input type="radio"/>
AVB0265		Side shroud colour Black					<input checked="" type="radio"/>
AVB0261		Side shroud colour Transparent Cherry Red					<input type="radio"/>
AVB0253		Side shroud colour Arctic White					<input type="radio"/>
AVB0250		Side shroud colour Grasshopper Green					<input type="radio"/>
AVB0254		Side shroud colour Sunshine Yellow					<input type="radio"/>
RIM COLOUR							
AVB0290		Rim colour White	Not available with grey front wheels				<input type="radio"/>
AVB0292		Rim colour Silver					<input checked="" type="radio"/>
AVB0291		Rim colour Black	Not available with grey front wheels				<input type="radio"/>
FRONT WHEEL							
AVB0804		254mm (10") Castor pneumatic, Black (non-marking)					<input type="radio"/>
AVB0805		254mm (10") Castor puncture proof solid inlay, Black (non-marking)					<input checked="" type="radio"/>
AVB0806		254mm (10") Castor puncture proof (sealed), Black (non-marking)					<input type="radio"/>
DRIVE WHEELS							
AVB0704		350mm (14") DW pneumatic, Black (non-marking)	Not available with grey front wheels				<input type="radio"/>
AVB0705		350mm (14") DW puncture proof solid inlay, Black (non-marking)	Not available with grey front wheels				<input type="radio"/>
AVB0706		350mm (14") DW puncture proof (sealed), Black (non-marking)	Not available with grey front wheels				<input type="radio"/>
PERFORMANCE OPTIONS							
DUAL SWING TECHNOLOGY							
AVB1200	<input checked="" type="checkbox"/>	Dual-swing technology suspension system					<input checked="" type="radio"/>
REAR SUSPENSION							
AVB1217		Rear suspension with soft shock absorbers recommended for user weight up to 90 kg	Not available with seat lifter with electrical tilt, with HD seating				<input type="radio"/>
AVB1218		Rear suspension with medium shock absorbers recommended for user weight 90 - 136 kg	Not available with seat lifter with electrical tilt, with HD seating				<input type="radio"/>
AVB1219		Rear suspension with hard shock absorbers recommended for user weight 136 - 160 kg; mandatory for lifter version option AVB1094					<input type="radio"/>
MOTOR TYPE							
AVB1282		6km/h (4mph)					<input checked="" type="radio"/>
AVB1908		10km/h (6mph) includes LED light system with indicators					€431
AVB1909		12km/h (8mph) includes LED light system with indicators					€431
BATTERY							
AVB1259		60Ah (Gel)					<input checked="" type="radio"/>
AVB1364		73.5Ah (Gel)					€93
PERFORMANCE ACCESSORIES							
G-TRAC							
AVB1215			G-trac ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain.				€824
CHARGER							
AVB1343		OptiCharge Charger 8A (200 - 250V)					<input checked="" type="radio"/>
AVB1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster					€158

ANTI-TIPPERS			
AVB1205	X	Standard anti-tipper	◆
LIGHT			
AVB1437		LED light system with Indicators	€353
AVB1254		Without Light system	◆
SAFETY ACCESSORIES			
KERB RIDER			
AVB1244		Kerb Rider	Not available with a footboard €227
TRANSPORTATION KIT			
AVB1494	X	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory)	◆
LAP BELT			
AVB1329		Lap belt with auto lock (standard mount on backrest)	◆
AVB2113		Lap belt with auto lock (flexible mount in side rails)	○
BODYPOINT PELVIC BELTS (Includes mounting kit)			
AVB2701		Evoflex Pelvic Stabilizer, Rehab Latch size L (padded length 620mm)	Not available with seat width 380mm, 430mm, seat depth 410mm, hip supports, with HD seating €107
AVB2702		Evoflex Pelvic Stabilizer, Rehab Latch size M (padded length 460mm)	Not available seat width 480mm, 530mm, seat depth 410mm, hip supports, with HD seating €107
AVB2703		2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)	Not available with seat width 380mm, 430mm, seat depth 410mm, hip supports, with HD seating €98
AVB2704		2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)	Not available seat width 480mm, 530mm, seat depth 410mm, hip supports, with HD seating €98
AVB2705		2-Point Hip Belt, Dual pull, Side release button, Cinch mount size L (padded length 620mm)	Not available with seat width 380mm, 430mm, seat depth 410mm, hip supports, with HD seating €98
AVB2706		2-Point Hip Belt, Dual pull, Side release button, Cinch mount size M (padded length 460mm)	Not available seat width 480mm, 530mm, seat depth 410mm, hip supports, with HD seating €98
REAR VIEW MIRROR - Must be selected with 10km/h & 12km/h motor			
AVB1415		Rear View Mirror RH	€42
AVB1339		Rear View Mirror LH	€42
ACCESSORIES			
BAG HOOKS			
			
AVB2128		Bag hooks (Max load up to 6 kg = impact on max user weight)	Not available with flex3 backrest 2172 / 2173, attendant brackets, dual control €86
LUGGAGE RACK			
AVB1438		Luggage rack	€339
CANE HOLDER (in combination with attendant or dual control the cane holder has to be mounted on the opposite side of the remote)			
AVB1312		Cane Holder RH	€42
AVB1363		Cane Holder LH	€42
FENDER (Black)			
AVB1416		Fender for castors	€104
AVB1417		Fender for rear wheels	€104
TRAY (PLEASE CHOOSE OPPOSITE HAND TO USER CONTROL)			
AVB1404		Swing Away RH	€167
AVB1405		Swing Away LH	€167
			Total Price :
			
<p>Need a different option not on the form? Take a look at our Powerchair Alternative Options form, available as standard Powerchair Alternative Option required:</p>			
			
Special Instructions - Please supply the following information:			
Name:		Email:	
Phone:		Fax:	
Description of special request (please include relevant option code or part number:			