



Yes, you can:

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Manual Invacare Ireland Prescription Form

07 JANUARY 2026
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Kuschall Champion 2.0

Standard Features

- Max User Weight 120kg (19 Stone)
- Lightweight Folding Fully Hydroformed Aluminium Wheelchair
- Multiple adjustments, Front, Rear Seat Heights and Centre of Gravity
- Aluminium Front Frame with 75° or 90° Front Frame Angles
- 24" Self Propel Q/R Wheels with Grey Aluminium Handrims
- and Pneumatic Tyres
- Choice of 7 Seat Widths 360mm to 480mm (14" - 19")
- Choice of 8 Seat Depths 340mm to 480mm (13" - 19")
- Choice of 5 Standard Frame Colours
- Range of Front Seat Heights 450mm to 540mm (17" - 21")
- Range of Rear Seat Heights 390mm to 480mm (15" - 19")
- Folding Backrest, Angle Adjustable and Soft Lock Feature
- Black Nylon Light Upholstery
- Angle and Height Adjustable Rigid Footrest, with Calf Strap
- Without Push Handles, Flip up Carbon Clothesguards
- Unfolding Help Cord. Black Comfort Front Castor
- Standard Brake Lever Bent
- Tool Kit



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
DHA0031	<input type="checkbox"/> Kuschall Champion 2.0	○
UNIQUE SPECIALS (For special modifications only) Additional lead time of 4 weeks		
DHA0000	<input type="checkbox"/> UNIQUE Special Request - For use when requesting any special options or builds	POA
UNIQUE CUSTOMISATION AND PERSONALISATION OPTIONS Additional 4 weeks lead time for Unique options		
UHA0001	<input type="checkbox"/> UNIQUE Personalisation - For use when selecting UHA customisation options	○

FRAME*			
DHA1150		Frame Angle 75° Aluminium	(NC 1350) €3,622
DHA1160		Frame Angle 90° Aluminium	(NC 1350) €3,622
DHA1170		Frame Angle 75° Aluminium - Tapered/Adducted Front Frame	(NC 1320, 1325, 7010) €3,622
DHA1180		Frame Angle 90° Aluminium - Tapered/Adducted Front Frame	(NC 1320, 1325, 7010) €3,622
CHAMPION SL HYBRID RIGID FRAME			
UHA1260		Complete frame with fully welded front seat height & centre of gravity rear wheel supports	(NC 6180) €578
SEAT WIDTH*			
DHA1320		360 mm	(NC 2710, 2748, 7000) ○
DHA1325		380 mm	(NC 7000) ○
DHA1330		400 mm	○
DHA1335		420 mm	○
DHA1340		440 mm	○
DHA1345		460 mm	○
DHA1350		480 mm	○
SEAT DEPTH*			
DHA1405		340 mm	(NC 2710, 2720, 2748, 2750, 2860, 7810, 7820, 7830) ○
DHA1415		360 mm	(NC 2710, 2720, 2748, 2860, 7810, 7820, 7830) ○
DHA1425		380 mm	(NC 2720) ○
DHA1435		400 mm	○
DHA1445		420 mm	○
DHA1460		440 mm	○
DHA1470		460 mm	○
DHA1480		480 mm	(NC 2710, 2748) ○
BACK HEIGHT*			
DHA1615		300 mm	(NC 2200, 2240, 2400) ○
DHA1620		315 mm	(NC 2200, 2240, 2400) ○
DHA1625		330 mm	(NC 2240) ○
DHA1630		345 mm	(NC 2240) ○
DHA1635		360 mm	○
DHA1640		375 mm	○
DHA1645		390 mm	○
DHA1650		405 mm	○
DHA1655		420 mm	○
DHA1660		435 mm	○
DHA1665		450 mm	○
DHA1670		465 mm	○

FRONT SEAT HEIGHT*				
DHA1740		450 mm	(NC 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960)	<input type="radio"/>
DHA1745		460 mm	(NC 1948, 1950, 1952, 1954, 1956, 1958, 1960)	<input type="radio"/>
DHA1750		470 mm	(NC 1950, 1952, 1954, 1956, 1958, 1960)	<input type="radio"/>
DHA1755		480 mm	(NC 1952, 1954, 1956, 1958, 1960)	<input type="radio"/>
DHA1760		490 mm	(NC 1954, 1956, 1958, 1960)	<input type="radio"/>
DHA1765		500 mm	(NC 1956, 1958, 1960)	<input type="radio"/>
DHA1770		510 mm	(NC 1958, 1960)	<input type="radio"/>
DHA1775		520 mm	(NC 1960)	<input type="radio"/>
DHA1780		530 mm	(NC 5000)	<input type="radio"/>
DHA1785		541 mm Only possible with 5" castors	(NC 5000, 5050, 5100, 5150, 5200, 5250)	<input type="radio"/>

REAR SEAT HEIGHT. On 75° frame, Rear Seat Height can be 0mm - 70mm lower then front seat height and 90° can be 30mm to 100mm lower*

DHA1815		390 mm Only possible with 22" & 24" rear wheels	(NC 6420, 6430)	<input type="radio"/>
DHA1820		400 mm Only possible with 22" & 24" rear wheels	(NC 6420, 6430)	<input type="radio"/>
DHA1825		410 mm	(NC 6430)	<input type="radio"/>
DHA1830		420 mm	(NC 6430)	<input type="radio"/>
DHA1835		430 mm		<input type="radio"/>
DHA1840		440 mm		<input type="radio"/>
DHA1845		450 mm		<input type="radio"/>
DHA1850		460 mm		<input type="radio"/>
DHA1855		470 mm Only possible with 22" & 24" rear wheels	(NC 6410)	<input type="radio"/>
DHA1860		480 mm Only possible with 22" & 24" rear wheels	(NC 6410)	<input type="radio"/>

LOWER LEG LENGTH. Cushion height must be deducted*

DHA1924		320 mm		<input type="radio"/>
DHA1926		330 mm		<input type="radio"/>
DHA1928		340 mm		<input type="radio"/>
DHA1930		350 mm		<input type="radio"/>
DHA1932		360 mm		<input type="radio"/>
DHA1934		370 mm		<input type="radio"/>
DHA1936		380 mm		<input type="radio"/>
DHA1938		390 mm		<input type="radio"/>
DHA1940		400 mm		<input type="radio"/>
DHA1942		410 mm		<input type="radio"/>
DHA1944		420 mm		<input type="radio"/>
DHA1946		430 mm		<input type="radio"/>
DHA1948		440 mm		<input type="radio"/>
DHA1950		450 mm		<input type="radio"/>
DHA1952		460 mm		<input type="radio"/>
DHA1954		470 mm		<input type="radio"/>
DHA1956		480 mm		<input type="radio"/>
DHA1958		490 mm		<input type="radio"/>
DHA1960		500 mm		<input type="radio"/>

BACKREST*

DHA2140		Angle Adjustable, Forward Folding with Soft Lock. Backrest remains locked for transportation		<input type="checkbox"/>
---------	--	---	--	--------------------------

BACKREST UPHOLSTERY*			
DHA2200		Black Standard Light Upholstery, without Tension Adjustment (NC 2241, 2242, 2243, 2244, 2245, 0002, 0003, 0004, 2430)	◆
DHA2220		Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener (NC 2241, 2242, 2243, 2244, 2245, 0002, 0003, 0004)	€160
DHA2230		Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener	€253
DHA2240		Without Backrest Upholstery Select if fitting a Matrx/Fixed backrest, no upholstery required (NC 2430, 2241, 2242, 2243, 2244, 2245, 0002, 0003, 0004)	○
BACK UPHOLSTERY STITCHING COLOUR			
UHA2241		Neon Yellow	€96
UHA2242		Neon Pink	€96
UHA2243		Blue	€96
UHA2244		White	€96
UHA2245		Red	€96
UHA0002		Your personal colour stitching. Add RAL colour of your choice in special comments box, for a personalised colour stitching on back upholstery.	€96
BACK UPHOLSTERY EMBROIDERY			
UHA0003		Your personal design embroidery. Supply your design with your order, subject to copyright. Design size up to 250mm x 250mm (NC 0004)	€256
BACK UPHOLSTERY FLOCKING			
UHA0004		Your personal flock design. Supply your design with your order, subject to copyright. Design size up to 280mm x 400mm	€127
PUSH HANDLES*			
DHA2400		Without Push Handles	◆
DHA2410		Push Handles Fixed Short/Mini Grip	€149
DHA2420		Push Handles Fixed Long Grip	○
DHA2430		Push Handles Height Adjustable, Rear Set (NC 3480)	€294
DHA2440		Push Handles Height Adjustable, Integrated in backrest	€225
DHA2450		Push Handles Foldable	€287
BACKREST ANGLE*			
DHA2510		76° angle from the seat	○
DHA2530		80.5° angle from the seat	○
DHA2550		85° angle from the seat	◆
DHA2570		89.5° angle from the seat	○
DHA2590		94° angle from the seat	○
SEAT UPHOLSTERY*			
DHA2600		Seat Upholstery with 1 central loop	○
DHA2601		Standard Seat Upholstery with 1 central loop and 2 side hook	◆
DHA2610		Tension Adjustable Seat Upholstery with one 5cm central band (Loop)	€54

SEAT CUSHION			
DHA2700		Standard Seat Cushion, 60mm high Double Foam with Visco Top Layer - Comfort/Low Risk	(NC 2810, 2820, 2830, 2840, 2850) €74
DHA2710		Flo-Tech Lite Cushion 78mm high- Comfort/Low Risk	(NC 2810, 2820, 2830, 2840, 2850) €168
DHA2720		Matrx Libra Cushion 83mm high, Black Smooth Cover - Very High Risk	€645
DHA2748		Contour Visco NG Cushion 100mm high - High Risk	(NC 2810, 2820, 2830, 2840, 2850) €236
DCA2750		Vicair Active O2 lightweight cushion 60mm high - High risk	(NC 2810, 2820, 2830, 2840, 2850, 2860) €731
SEAT BAG			
DHA2800		Seat Bag, Click in and removable, includes the fixation clips	(NC 7010) €85
MATRX LIBRA CUSHION ACCESSORIES			
DHA2810		Fluid supplement sac (Pair)	€123
DHA2820		Tapered wedge (Pair)	€28
DHA2830		Abductor	€58
DHA2840		Lateral side wedge	€46
DHA2850		Anterior front wedge	€43
SEAT RIGIDIZER			
DHA2860		Seat rigidizer	€79
FOOTREST*			
DHA3240		Rigid, Plastic Footplate, Set to front, Height and Angle Adjustable. Note: This does not flip up for transfer, but folds easy for transport and storage	◆
DHA3250		Rigid, Plastic Footplate, Set to rear, Height and Angle Adjustable. Note: This does not flip up for transfer, but folds easy for transport and storage	○
CALF STRAP*			
DHA3320		Calf Strap Mounted on Supporter	◆
DHA3330		Calf Strap Mounted on Frame	○
CLOTHESGUARDS/MUDGUARDS*			
DHA3405		Composite Clothesguards, Flip Back - (only possible with 24" and 25" rear wheels)	(NC 6430) ○
DHA3420		Carbon Clothesguards, Flip Up	◆
DHA3425		Composite Mudguards, Fixed (only possible with 24" and 25" rear wheels)	(NC 3490, 6430) €91
DHA3430		Carbon Mudguards, Fixed	(NC 3490) €192
DHA3440		Carbon Mudguards, Removable	(NC 3490) €327
DHA3450		Composite Mudguards, Removable (only possible with 24" and 25" rear wheels)	(NC 3490, 6430) €91
ARMREST			
DHA3480		Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included	€281
DHA3490		Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included	€109

FRAME COLOUR*				
DHA4000		Space Black (RAL 9005)	(NC 4066)	○
DHA4010		Anthracite Grey (RAL 7016)	(NC 4066)	○
DHA4020		Pure White (RAL 9016)	(NC 4066, 4220)	○
DHA4030		Carmine Red (RAL 3002)	(NC 4066)	○
DHA4040		Azure Blue (RAL 5017)	(NC 4066)	○
DHA4060		Flashy Yellow (RAL 1516) (Additional 15 day lead time)	(NC 4066)	€39
DHA4100		Black Glossy (RAL 9005) (Additional 15 day lead time)	(NC 4066)	€39
DHA4130		Silver Grey (RAL 9006) (Additional 15 day lead time)	(NC 4066)	€39
DHA4145		Pearl Night Blue (RAL 5026) (Additional 15 day lead time)	(NC 4066)	€39
DHA4150		Ruby Red (RAL 1503) (Additional 15 day lead time)	(NC 4066)	€39
DHA4170		Dynamite Orange (Additional 15 day lead time)	(NC 4066)	€39
DHA4180		Forest Green (RAL 6029) (Additional 15 day lead time)	(NC 4066)	€39
UHA4101		Diamond Silver/Mirror Polish with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066)		€513
UHA4102		Pastel Green with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066)		€413
UHA4103		Sky Blue with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066)		€413
UHA4104		Plum Purple with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066)		€413
UHA4105		Salmon with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066)		€413
UHA0008		Your personal frame colour. Add RAL colour of your choice in special comments box, for a personalised colour frame		€384
CERAMIC GLAZE FRAME FINISH				
UHA4066		Ceramic glaze finish on frame, for added durability		€432
FRAME DECAL TYPE*				
DHA4190		Bold Decal (Full coloured transfer on frame)	(NC 4300, 4310, 1111)	€39
DHA4195		Tag (Small Kuschall brand transfer on front of frame)	(NC 4200, 4210, 4220, 4230, 4240, 4250, 2222)	◆
DHA4197		Without Tag (Kuschall branding only)	(NC 4200, 4210, 4220, 4230, 4240, 4250, 4300, 4310)	○
BOLD DECAL COLOUR*				
DHA4200		Black		○
DHA4210		Dark Grey		○
DHA4220		Light Grey		○
DHA4230		Red		○
DHA4240		Blue		○
DHA4250		Gold		○
DHA1111		No bold colour, select only with Option DHA4195 Tag or DHA4197 Without Tag	(NC 4190)	○
TAG COLOUR*				
DHA4300		Black		○
DHA4310		Red		○
DHA2222		No tag colour, select only with Option DHA4190 Bold Decal or DHA4197 Without Tag	(NC 4195)	○

FRONT CASTORS*			
DHA5000		75mm x 25mm (3" x 1") Supersport	(NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) ○
DHA5050		100mm x 30mm (4" x 1.2") Everyday Black hub and tyre	(NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) ○
DHA5100		100mm x 30mm (4" x 1.2") Comfort Black hub and tyre	(NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) ○
DHA5150		100mm x 25mm (4" x 1") Starec Black hub and tyre	(NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) €102
DHA5200		100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frogs Legs) Use this option to access 4" anodised colour options	€114
DHA5250		100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs)	(NC 5502, 5504, 5506, 5507, 5508, 5509) €114
DHA5300		125mm x 30mm (5" x 1.2") Everyday Black hub and tyre	(NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) ○
DHA5350		125mm x 30mm (5" x 1.2") Comfort Black hub and tyre	(NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) ◆
DHA5400		125mm x 25mm (5" x 1") Starec Black hub and tyre	(NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) €102
DHA5450		125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) Use this option to access 5" anodised colour options	€114
DHA5500		125mm x 36mm (5" x 1.4") (Frog Legs) Aluminium Softroll Silver hub, Black tyre (Frogs Legs)	(NC 5502, 5504, 5506, 5507, 5508, 5509) €114
ANODISED COLOUR SOFTROLL CASTORS			
UHA5502		Pink	€64
UHA5504		Orange	€64
UHA5506		Blue	€64
UHA5507		Purple	€64
UHA5508		Red	€64
UHA5509		Gold	€64
CERAMIC CASTOR BEARINGS			
UHA5510		Ceramic castor bearings for Aluminium soft roll castors	€75
REAR WHEEL POSITION (CENTRE OF GRAVITY)*			
DHA6100		1st from front Most Active Position	○
DHA6110		2nd from front	○
DHA6120		3rd from front	○
DHA6130		4th from front	○
DHA6140		5th from front	◆
DHA6150		6th from front	○
DHA6160		7th from front	○
DHA6170		8th from front	○
DHA6180		9th from front	○
QUAD QUICK RELEASE AXLE			
DHA6300		Quad Quick Release/Tetra Clip	(NC 6550, 6551, 6552, 6553, 6554, 0907) €96
REAR WHEEL SIZE*			
DHA6410		24"	◆
DHA6420		25"	(NC 6510, 6600, 6630, 6660, 6700, 6710, 6720, 6760, 6770, 6775, 6900) ○
DHA6430		26"	(NC 6520, 6550, 6551, 6552, 6553, 6554, 6600, 6630, 6660, 6720, 6740, 6750, 6760, 6770, 6775) ○

REAR WHEEL TYPE*			
DHA6510		24", 26" Kuschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub.	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6579) ◆
DHA6520		24", 25" High Performance Wheel Black hub and rim	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6579, 6660) €311
DHA6550		24", 25" Spinergy LXX, Black hub, rim and 18 spokes in Black	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6660) €746
UHA6551		24", 25" Spinergy Lite Extreme LX, Weight 687 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)	(NC 6660) €1,009
UHA6552		24", 25" Spinergy Lite Extreme X-laced XLX, Weight 821 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)	(NC 6660) €1,084
UHA6553		24", 25" Spinergy X-Laced Sport Lite Extreme XSLX, Weight 757 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)	(NC 6660) €1,147
UHA6554		24", 25" Spinergy Carbon Fibre CLX, Weight 530 grams per wheel (Unique coloured hub and spokes available, supplied in black as standard)	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6660) €1,393
SWO0907		24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. (Unique coloured spokes available)	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6660) €543
ANODISED COLOUR SPINERGY RIM (DESIGN LAB)			
UHA6557		Blue	€230
UHA6558		Gold	€230
UHA6559		Green	€230
UHA6560		Orange	€230
UHA6561		Red	€230
UHA6562		Purple	€230
ANODISED COLOUR SPINERGY HUB (DESIGN LAB)			
UHA6563		Blue	€230
UHA6564		Gold	€230
UHA6565		Silver	€230
UHA6566		Green	€230
UHA6567		Orange	€230
UHA6568		Red	€230
UHA6569		Purple	€230
SPINERGY COLOURED SPOKES			
UHA6570		Blue	€108
UHA6571		Green	€108
UHA6572		Orange	€108
UHA6573		Pink	€108
UHA6574		Red	€108
UHA6575		Purple	€108
UHA6576		White	€108
UHA6577		Yellow	€108
CERAMIC BEARING REAR WHEEL			
UHA6579		Ceramic Bearing for Spinergy rear wheels	€218

REAR WHEEL TYRE*			
DHA6600		24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection	◆
DHA6625		24", 25", 26" x 1" Maxxis Alpha Protect. Lightweight, Pneumatic and highly puncture resistant.	€113
DHA6630		24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance)	€96
DHA6660		24" x 1 3/8" Solid Profile Black Tyre (NC 7230)	€82
HANDRIM*			
DHA6700		24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance	€252
DHA6710		26" Steel, short tab with 1 hole, mounted at 3 cm distance	€76
DHA6720		24" Titanium, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance	€200
DHA6740		24", 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. For complete or limited hand function offering low-effort propulsion. Ergonomic shape with fully integrated, replaceable silicon strip for extra grip, profile of 32mm diameter.	€280
DHA6750		24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. For complete or limited hand function offering low-effort propulsion. Ergonomic shape with fully integrated, replaceable silicon strip for extra grip, smaller profile of 24mm diameter for smaller hands.	€280
DHA6760		24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape, all metal black coated, low-friction surface for complete hand function for low-effort propulsion	€250
DHA6770		24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics without tricep function. Features a unique angle shaped profile and anti-slip coating providing a larger contact surface designed to be propelled using the heel of the hand or thumb pad.	€290
DHA6775		24" Carbolife Curve L Tetra handrim, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics with tricep function, but limited hand strength. Ergonomic shape with a flat, non-slip coated pushing surface for the thumb and a curved side for the fingers to wrap around.	€200
DHA6780		24", 25", 26" Light Aluminium Grey, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance	◆
DHA6790		24", 25", 26" Light Aluminium Black, short tab with 1 hole for 3cm mounting	€51
REAR WHEEL CAMBER*			
DHA6820		0° Camber	◆
DHA6830		1° Camber	○
DHA6840		3° Camber (NC 7230)	○
SPOKEGUARD			
DHA6900		24", 26" Transparent Spokeguards (When ordered with Carbolife handrims, handrim distance is set at 4cm)	€155
FIXATION CLIP*			
DHA7000		Fixation Strap (Hook and loop strap under the seat holds the frame together when folded for easy transportation, designed for clients with limited hand function)	○
DHA7010		Fixation Clip (Plastic clip on either side of the seat holds the frame together when folded for easy transportation)	◆
DHA3333		No fixation clip, select with Tapered frame options 1170, 1180	○

UNFOLDING HELP*			
DHA7100	<input type="checkbox"/>	Unfolding Help (Cord to help open the folding mechanism from the front of the chair)	◆
UNIQUE CUSTOMISATION COLOUR PACKAGE*			
UHA6960	<input type="checkbox"/>	UNIQUE KIT 1 - Includes custom colour forks and bearing cover. Select colour from the 10 choices below	€167
UHA6961	<input type="checkbox"/>	UNIQUE KIT 2 - Includes custom colour forks and bearing cover, castor fork supporter, rear wheel supporter, brake clamps. Select colour from the 10 choices below	€366
UHA1234	<input type="checkbox"/>	Without colour package	○
UNIQUE CUSTOMISATION COLOURS			
UHA6963	<input type="checkbox"/>	Red	○
UHA6964	<input type="checkbox"/>	Purple	○
UHA6965	<input type="checkbox"/>	Copper	○
UHA6966	<input type="checkbox"/>	Green	○
UHA6967	<input type="checkbox"/>	Orange	○
UHA6968	<input type="checkbox"/>	Gold	○
UHA6969	<input type="checkbox"/>	Champagne	○
UHA6970	<input type="checkbox"/>	Pink	○
UHA6971	<input type="checkbox"/>	Blue	○
UHA6972	<input type="checkbox"/>	Poly Mirror Finish	○
BRAKES*			
DHA7200	<input type="checkbox"/>	Pull Brakes	€85
DHA7210	<input type="checkbox"/>	Standard Push to Lock Brakes, Lever Bent	◆
DHA7220	<input type="checkbox"/>	Standard Push to Lock Brakes, Lever Straight	○
DHA7230	<input type="checkbox"/>	Performance Brakes, Push to Lock. Brake barrel flips away from tyre	(NC 7300, 7810, 7820, 7830) €85
DHA7240	<input type="checkbox"/>	Active/Scissor Brakes	(NC 7300) €85
DHA7250	<input type="checkbox"/>	Light/Compact Active, Aluminium Brakes	(NC 7300) €248
BRAKE LEVER EXTENSION			
DHA7300	<input type="checkbox"/>	Brake Lever Extension	€48
SAFETY OPTIONS			
DHA7500	<input type="checkbox"/>	Anti-Tipper Right	€149
DHA7510	<input type="checkbox"/>	Anti-Tipper Left	€149
DHA7520	<input type="checkbox"/>	Anti-Tipper Both sides	€299
DHA7800	<input type="checkbox"/>	Standard Pelvic belt	€74
BODYPOINT BELTS			
DHA7810	<input type="checkbox"/>	Bodypoint Posture Belt Evoflex Small	€171
DHA7820	<input type="checkbox"/>	Bodypoint Posture Belt Evoflex Medium	€171
DHA7830	<input type="checkbox"/>	Bodypoint Posture Belt Evoflex Large	€171

OTHER OPTIONS		
DHA8010		Front Frame Protectors €36
DHA8100		Tipper-Aid, Mounted Right Side Max. 2 adaptions per side (Adaptations = Anti tipper, Tipper aid and Cane holder) €86
DHA8200		Cane Holder, Mounted Right Side Max. 2 adaptions per side (Adaptations = Anti tipper, Tipper aid and cane holder) €86
DHA8500		Air Pump €17
DHA8550		Tool Kit ◆
DHA8650		Kuschall branded Wheel bag (Fits all wheel sizes) €58
DHA8680		24" Mountain Bike Wheels with bag. Fitted with Specific Black Tyres, Grey Aluminium Handrims; Quick Release Axles Included €335
DHA8700		Travel bag for Champion, Wheels can remain on the wheelchair €294
BODYPOINT ACCESSORIES		
DHA8800		Bodypoint Calf Strap €80
DHA8850		Bodypoint Mobility Bag €80

TOTAL PRICE	
--------------------	--

Non Conforming Options:

DHA1320	360 mm
DHA1325	380 mm
DHA1350	480 mm
DHA1946	430 mm
DHA1948	440 mm
DHA1950	450 mm
DHA1952	460 mm
DHA1954	470 mm
DHA1956	480 mm
DHA1958	490 mm
DHA1960	500 mm
DHA2200	Black Standard Light Upholstery, without Tension Adjustment
DHA2240	Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required
UHA2241	Neon Yellow
UHA2242	Neon Pink
UHA2243	Blue
UHA2244	White
UHA2245	Red
UHA0002	Your personal colour stitching. Add RAL colour of your choice in special comments box, for a personalised colour stitching on back upholstery.
UHA0003	Your personal design embroidery. Supply your design with your order, subject to copyright. Design size up to 250mm x 250mm
UHA0004	Your personal flock design. Supply your design with your order, subject to copyright. Design size up to 280mm x 400mm
DHA2400	Without Push Handles
DHA2430	Push Handles Height Adjustable, Rear Set
DHA2710	Flo-Tech Lite Cushion 78mm high- Comfort/Low Risk
DHA2720	Matrx Libra Cushion 83mm high, Black Smooth Cover - Very High Risk
DHA2748	Contour Visco NG Cushion 100mm high - High Risk
DCA2750	Vicair Active O2 lightweight cushion 60mm high - High risk
DHA2810	Fluid supplement sac (Pair)
DHA2820	Tapered wedge (Pair)
DHA2830	Abductor
DHA2840	Lateral side wedge
DHA2850	Anterior front wedge
DHA2860	Seat rigidizer
DHA3480	Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included
DHA3490	Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included
UHA4066	Ceramic glaze finish on frame, for added durability
DHA4190	Bold Decal (Full coloured transfer on frame)
DHA4195	Tag (Small Kuschall brand transfer on front of frame)
DHA4200	Black
DHA4210	Dark Grey
DHA4220	Light Grey
DHA4230	Red
DHA4240	Blue

DHA4250	Gold
DHA1111	No bold colour, select only with Option DHA4195 Tag or DHA4197 Without Tag
DHA4300	Black
DHA4310	Red
DHA2222	No tag colour, select only with Option DHA4190 Bold Decal or DHA4197 Without Tag
DHA5000	75mm x 25mm (3" x 1") Supersport
DHA5050	100mm x 30mm (4" x 1.2") Everyday Black hub and tyre
DHA5100	100mm x 30mm (4" x 1.2") Comfort Black hub and tyre
DHA5150	100mm x 25mm (4" x 1") Starec Black hub and tyre
DHA5200	100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frogs Legs) Use this option to access 4" anodised colour options
DHA5250	100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs)
UHA5502	Pink
UHA5504	Orange
UHA5506	Blue
UHA5507	Purple
UHA5508	Red
UHA5509	Gold
UHA5510	Ceramic castor bearings for Aluminium soft roll castors
DHA6180	9th from front
DHA6410	24"
DHA6420	25"
DHA6430	26"
DHA6510	24", 26" Kuschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub.
DHA6520	24", 25" High Performance Wheel Black hub and rim
DHA6550	24", 25" Spinergy LXX, Black hub, rim and 18 spokes in Black
UHA6551	24", 25" Spinergy Lite Extreme LX, Weight 687 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)
UHA6552	24", 25" Spinergy Lite Extreme X-laced XLX, Weight 821 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)
UHA6553	24", 25" Spinergy X-Laced Sport Lite Extreme XSLX, Weight 757 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)
UHA6554	24", 25" Spinergy Carbon Fibre CLX, Weight 530 grams per wheel (Unique coloured hub and spokes available, supplied in black as standard)
SWO0907	24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. (Unique coloured spokes available)
UHA6557	Blue
UHA6558	Gold
UHA6559	Green
UHA6560	Orange
UHA6561	Red
UHA6562	Purple
UHA6563	Blue
UHA6564	Gold
UHA6565	Silver
UHA6566	Green
UHA6567	Orange
UHA6568	Red

UHA6569	Purple
UHA6570	Blue
UHA6571	Green
UHA6572	Orange
UHA6573	Pink
UHA6574	Red
UHA6575	Purple
UHA6576	White
UHA6577	Yellow
UHA6579	Ceramic Bearing for Spinergy rear wheels
DHA6600	24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection
DHA6630	24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance)
DHA6660	24" x 1 3/8" Solid Profile Black Tyre
DHA6700	24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance
DHA6710	26" Steel, short tab with 1 hole, mounted at 3 cm distance
DHA6720	24" Titanium, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance
DHA6740	24", 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. For complete or limited hand function offering low-effort propulsion. Ergonomic shape with fully integrated, replaceable silicon strip for extra grip, profile of 32mm diameter.
DHA6750	24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. For complete or limited hand function offering low-effort propulsion. Ergonomic shape with fully integrated, replaceable silicon strip for extra grip, smaller profile of 24mm diameter for smaller hands.
DHA6760	24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape, all metal black coated, low-friction surface for complete hand function for low-effort propulsion
DHA6770	24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics without tricep function. Features a unique angle shaped profile and anti-slip coating providing a larger contact surface designed to be propelled using the heel of the hand or thumb pad.
DHA6775	24" Carbolife Curve L Tetra handrim, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics with tricep function, but limited hand strength. Ergonomic shape with a flat, non-slip coated pushing surface for the thumb and a curved side for the fingers to wrap around.
DHA6900	24", 26" Transparent Spokeguards (When ordered with Carbolife handrims, handrim distance is set at 4cm)
DHA7000	Fixation Strap (Hook and loop strap under the seat holds the frame together when folded for easy transportation, designed for clients with limited hand function)
DHA7010	Fixation Clip (Plastic clip on either side of the seat holds the frame together when folded for easy transportation)
DHA7230	Performance Brakes, Push to Lock. Brake barrel flips away from tyre
DHA7300	Brake Lever Extension
DHA7810	Bodypoint Posture Belt Evoflex Small
DHA7820	Bodypoint Posture Belt Evoflex Medium
DHA7830	Bodypoint Posture Belt Evoflex Large



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):