



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Manual Retail Prescription Form

August 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Küschall Champion 2.0 One Solution

Standard Features

- Max User Weight 120kg (19 Stone)
- Lightweight Folding Fully Hydroformed Aluminium Wheelchair
- Multiple adjustments, Front, Rear Seat Heights and Centre of Gravity
- Aluminium Front Frame with 75° or 90° Front Frame Angles
- 24" Self Propel Q/R Wheels with Grey Aluminium Handrims
- and Pneumatic Tyres
- Choice of 7 Seat Widths 360mm to 480mm (14" - 19")
- Choice of 8 Seat Depths 340mm to 480mm (13" - 19")
- Choice of 5 Standard Frame Colours
- Range of Front Seat Heights 450mm to 540mm (17" - 21")
- Range of Rear Seat Heights 390mm to 480mm (15" - 19")
- Folding Backrest, Angle Adjustable and Soft Lock Feature with Matrix backrest
- Angle and Height Adjustable Rigid Footrest, with Calf Strap
- Without Push Handles, Flip up Carbon Clothesguards
- Unfolding Help Cord. Black Comfort Front Castor
- Standard Brake Lever Bent
- Tool Kit



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description
------	-------------

BASE PRODUCT*

DHA0031	☐	Küschall Champion 2.0	☐
---------	--------------------------	-----------------------	-----------------------

UNIQUE SPECIALS (For special modifications only) Additional lead time of 4 weeks

DHA0000	☐	UNIQUE Special Request - For use when requesting any special options or builds	☐
---------	--------------------------	--	-----------------------

UNIQUE CUSTOMISATION AND PERSONALISATION OPTIONS Additional 4 weeks lead time for Unique options

UHA0001	☐	UNIQUE Personalisation - For use when selecting UHA (NC 1350, 4185) customisation options	☐
---------	--------------------------	---	-----------------------

FRAME*				
DHA1150		Frame Angle 75° Aluminium	(NC 1350, 4185)	£2,695
DHA1160		Frame Angle 90° Aluminium	(NC 1320, 1325, 7010, 4185)	£2,695
DHA1170		Frame Angle 75° Aluminium - Tapered/Adducted Front Frame	(NC 1320, 1325, 7010, 4185)	£2,695
DHA1180		Frame Angle 90° Aluminium - Tapered/Adducted Front Frame		£2,695
DHA1185		Frame Angle 75° Carbon	(NC 1260, 1350, 4000, 4010, 4020, 4030, 4040, 4060, 4100, 4130, 4145, 4150, 4170, 4180, 4101, 4102, 4103, 4104, 4190, 4197, 4300)	£3,100
DHA1190		Frame Angle 90° Carbon	(NC 1260, 1350, 4000, 4010, 4020, 4030, 4040, 4060, 4100, 4130, 4145, 4150, 4170, 4180, 4101, 4102, 4103, 4104, 4190, 4197, 4300)	£3,100
CHAMPION SL HYBRID RIGID FRAME				
UHA1260		Complete frame with fully welded front seat height & centre of gravity rear wheel supports	(NC 6180)	£450
SEAT WIDTH*				
DHA1320		360 mm	(NC 2364, 2366, 2368, 2370, 2372, 2710, 2748, 7000)	○
DHA1325		380 mm	(NC 2366, 2368, 2370, 2372, 7000)	○
DHA1330		400 mm	(NC 2360, 2366, 2368, 2370, 2372)	○
DHA1335		420 mm	(NC 2360, 2362, 2368, 2370, 2372)	○
DHA1340		440 mm	(NC 2360, 2362, 2364, 2370, 2372)	○
DHA1345		460 mm	(NC 2360, 2362, 2364, 2370, 2372)	○
DHA1350		480 mm	(NC 2300, 2305, 2310, 2315, 2345, 2360, 2362, 2364, 2366)	○
SEAT DEPTH*				
DHA1405		340 mm	(NC 2710, 2720, 2748, 2750, 2860, 7810, 7820, 7830)	○
DHA1415		360 mm	(NC 2710, 2720, 2748, 2860, 7810, 7820, 7830)	○
DHA1425		380 mm	(NC 2720)	○
DHA1435		400 mm		○
DHA1445		420 mm		○
DHA1460		440 mm		○
DHA1470		460 mm		○
DHA1480		480 mm	(NC 2710, 2748)	○
BACK HEIGHT*				
DHA1615		300 mm	(NC 2325, 2330, 2331, 2345, 2350, 2355, 2357, 2400)	○
DHA1620		315 mm	(NC 2325, 2330, 2331, 2350, 2355, 2357, 2400)	○
DHA1625		330 mm	(NC 2325, 2330, 2331, 2350, 2355, 2357)	○
DHA1630		345 mm	(NC 2330, 2331, 2355, 2357)	○
DHA1635		360 mm		○
DHA1640		375 mm		○
DHA1645		390 mm		○
DHA1650		405 mm		○
DHA1655		420 mm		○
DHA1660		435 mm		○
DHA1665		450 mm		○
DHA1670		465 mm		○

FRONT SEAT HEIGHT*				
DHA1740		450 mm	(NC 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960)	○
DHA1745		460 mm	(NC 1948, 1950, 1952, 1954, 1956, 1958, 1960)	○
DHA1750		470 mm	(NC 1950, 1952, 1954, 1956, 1958, 1960)	○
DHA1755		480 mm	(NC 1952, 1954, 1956, 1958, 1960)	○
DHA1760		490 mm	(NC 1954, 1956, 1958, 1960)	○
DHA1765		500 mm	(NC 1956, 1958, 1960)	○
DHA1770		510 mm	(NC 1958, 1960)	○
DHA1775		520 mm	(NC 1960)	○
DHA1780		530 mm	(NC 5000)	○
DHA1785		540 mm	(NC 5000, 5050, 5100, 5200, 5250)	○
REAR SEAT HEIGHT. On 75° frame, Rear Seat Height can be 0mm - 70mm lower then front seat height and 90° can be 30mm to 100mm lower*				
DHA1815		390 mm	(NC 6420, 6430)	○
DHA1820		400 mm	(NC 6420, 6430)	○
DHA1825		410 mm	(NC 6430)	○
DHA1830		420 mm	(NC 6430)	○
DHA1835		430 mm		○
DHA1840		440 mm		○
DHA1845		450 mm		○
DHA1850		460 mm		○
DHA1855		470 mm	(NC 6410)	○
DHA1860		480 mm	(NC 6410)	○
LOWER LEG LENGTH Seat canvas to footplate measurement, cushion height must be deducted.*				
DHA1924		320 mm		○
DHA1926		330 mm		○
DHA1928		340 mm		○
DHA1930		350 mm		○
DHA1932		360 mm		○
DHA1934		370 mm		○
DHA1936		380 mm		○
DHA1938		390 mm		○
DHA1940		400 mm		○
DHA1942		410 mm		○
DHA1944		420 mm		○
DHA1946		430 mm		○
DHA1948		440 mm		○
DHA1950		450 mm		○
DHA1952		460 mm		○
DHA1954		470 mm		○
DHA1956		480 mm		○
DHA1958		490 mm		○
DHA1960		500 mm		○
BACKREST*				
DHA2140		Angle Adjustable, Forward Folding with Soft Lock. Backrest remains locked for transportation		◆
BACKREST UPHOLSTERY*				
DHA2240		Without Backrest Upholstery Select if fitting a Matrx/Fixed backrest, no upholstery required		○

MATRX BACKREST*			
DHA2300		Matrx MX2, Height 150 mm with Q/R Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting. (NC 2370, 2372)	£751
DHA2305		Matrx MX2, Height 230 mm with Q/R Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting. (NC 2370, 2372)	£751
DHA2310		Matrx MX2, Height 305 mm with Q/R Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting. (NC 2370, 2372)	£751
DHA2315		Matrx MX2, Height 405 mm with Q/R Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting. (NC 2370, 2372)	£751
DHA2316		Matrx Elite Deep E2, Height 160 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.	£546
DHA2317		Matrx Elite Deep E2, Height 210 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.	£546
DHA2318		Matrx Elite Deep E2, Height 260 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.	£546
DHA2320		Matrx Elite Deep E2, Height 310 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.	£546
DHA2325		Matrx Elite Deep E2, Height 410 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.	£546
DHA2330		Matrx Elite Deep E2, Height 460 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.	£546
DHA2331		Matrx Elite Deep E2, Height 510 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.	£546
DHA2332		Matrx Elite E2, Height 160 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.	£494
DHA2334		Matrx Elite E2, Height 210 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.	£494
DHA2335		Matrx Elite E2, Height 260 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.	£494
DHA2340		Matrx Elite E2, Height 310 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.	£494
DHA2345		Matrx Elite E2, Height 360 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.	£494
DHA2350		Matrx Elite E2, Height 410 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.	£494
DHA2355		Matrx Elite E2, Height 460 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.	£494
DHA2357		Matrx Elite E2, Height 510 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.	£494

MATRX BACKREST WIDTH CHOICE*				
DHA2360		360mm Backrest width, to fit chest width 265mm		○
DHA2362		380mm Backrest width, to fit chest width 290mm		○
DHA2364		410mm Backrest width, to fit chest width 315mm		○
DHA2366		430mm Backrest width, to fit chest width 340mm		○
DHA2368		460mm Backrest width, to fit chest width 370mm		○
DHA2370		480mm Backrest width, to fit chest width 390mm		○
DHA2372		510mm Backrest width, to fit chest width 420mm		○
PUSH HANDLES*				
DHA2400		Without Push Handles		◆
DHA2410		Push Handles Fixed Short/Mini Grip		£133
DHA2420		Push Handles Fixed Long Grip		○
DHA2440		Push Handles Height Adjustable, Integrated in backrest		£219
DHA2450		Push Handles Foldable		£198
BACKREST ANGLE*				
DHA2510		76° angle from the seat		○
DHA2530		80.5° angle from the seat		○
DHA2550		85° angle from the seat		◆
DHA2570		89.5° angle from the seat		○
DHA2590		94° angle from the seat		○
SEAT UPHOLSTERY*				
DHA2600		Seat Upholstery with 1 central loop		○
DHA2601		Standard Seat Upholstery with 1 central loop and 2 side hook		◆
DHA2610		Tension Adjustable Seat Upholstery with one 5cm central band (Loop)		£45
SEAT CUSHION				
DHA2700		Standard Seat Cushion, 60mm high Double Foam with Visco Top Layer - Comfort/Low Risk	(NC 2810, 2820, 2830, 2840, 2850, 2815, 2825, 2835, 2845, 2855)	£66
DHA2710		Flo-Tech Lite Cushion 78mm high- Comfort/Low Risk	(NC 2810, 2830, 2815, 2835)	£121
DHA2720		Matrx Libra Cushion 70mm high, Black Smooth Cover - Very High Risk		£410
DHA2748		Contour Visco NG Cushion 120mm high - High Risk	(NC 2810, 2820, 2830, 2840, 2850, 2815, 2825, 2835, 2845, 2855)	£169
DHA2750		Vicair Active O2 lightweight cushion 60mm high - High risk	(NC 2810, 2820, 2830, 2840, 2850, 2815, 2825, 2835, 2845, 2855, 2860)	£521
SEAT BAG				
DHA2800		Seat Bag, Click in and removable, includes the fixation clips	(NC 7000, 7010)	£74
MATRX LIBRA CUSHION ACCESSORIES				
DHA2810		Fluid supplement sac (Single unit)		£44
DHA2820		Tapered wedge (Pair)		£30
DHA2830		Abductor		£23
DHA2840		Lateral side wedge		£20
DHA2850		Anterior front wedge		£27

2ND MATRX LIBRA CUSHION POSIIONING ACCESSORIES			
DHA2815		Fluid supplement sac (Single unit)	£44
DHA2825		Tapered wedge (Pair)	£30
DHA2835		Abductor	£23
DHA2845		Lateral side wedge	£20
DHA2855		Anterior front wedge	£27
SEAT RIGIDIZER			
DHA2860		Seat rigidizer	£39
FOOTREST*			
DHA3240		Rigid, Plastic Footplate, Set to front, Height and Angle Adjustable. Note: This does not flip up for transfer, but folds easy for transport and storage	◆
DHA3250		Rigid, Plastic Footplate, Set to rear, Height and Angle Adjustable. Note: This does not flip up for transfer, but folds easy for transport and storage	○
CALF STRAP*			
DHA3320		Calf Strap Mounted on Supporter	◆
DHA3330		Calf Strap Mounted on Frame	○
CLOTHESGUARDS/MUDGUARDS*			
DHA3405		Composite Clothesguards, Flip Back - (only possible with 24" and 25" rear wheels) (NC 6430)	○
DHA3420		Carbon Clothesguards, Flip Up	◆
DHA3425		Composite Mudguards, Fixed (only possible with 24" and 25" rear wheels) (NC 3490, 6430)	£65
DHA3430		Carbon Mudguards, Fixed (NC 3490)	£247
DHA3440		Carbon Mudguards, Removable (NC 3490)	£307
DHA3450		Composite Mudguards, Removable (only possible with 24" & 25" rear wheels) (NC 3490, 6430)	£65
ARMREST			
DHA3490		Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included	£106

FRAME COLOUR*				
DHA4000		Space Black (RAL 9005)	(NC 4066)	○
DHA4010		Anthracite Grey (RAL 7016)	(NC 4066)	○
DHA4020		Pure White (RAL 9016)	(NC 4066, 4220)	○
DHA4030		Carmine Red (RAL 3002)	(NC 4066)	○
DHA4040		Azure Blue (RAL 5017)	(NC 4066)	○
DHA4060		Flashy Yellow (RAL 1516) (Additional 15 day lead time)	(NC 4066)	£31
DHA4100		Black Glossy (RAL 9005) (Additional 15 day lead time)	(NC 4066)	£31
DHA4130		Silver Grey (RAL 9006) (Additional 15 day lead time)	(NC 4066)	£31
DHA4145		Pearl Night Blue (RAL 5026) (Additional 15 day lead time)	(NC 4066)	£31
DHA4150		Ruby Red (RAL 1503) (Additional 15 day lead time)	(NC 4066)	£31
DHA4170		Dynamite Orange (Additional 15 day lead time)	(NC 4066)	£31
DHA4180		Forest Green (RAL 6029) (Additional 15 day lead time)	(NC 4066)	£31
UHA4101		Diamond Silver/Mirror Polish with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066)		£399
UHA4102		Pastel Green with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066)		£299
UHA4103		Sky Blue with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066)		£299
UHA4104		Plum Purple with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066)		£299
UHA0008		Your personal frame colour. Add RAL colour of your choice in special comments box, for a personalised colour frame		£299
DHA4185		Carbon Finish		○
CERAMIC GLAZE FRAME FINISH				
UHA4066		Ceramic glaze finish on frame, for added durability		£336
FRAME DECAL TYPE*				
DHA4190		Bold Decal (Full coloured transfer on frame)	(NC 4300, 4310, 1111)	£31
DHA4195		Tag (Small Küschall brand transfer on front of frame)	(NC 4200, 4210, 4220, 4230, 4240, 4250, 2222)	◆
DHA4197		Without Tag (Kuschall branding only)	(NC 4200, 4210, 4220, 4230, 4240, 4250, 4300, 4310)	○
BOLD DECAL COLOUR*				
DHA4200		Black		○
DHA4210		Dark Grey		○
DHA4220		Light Grey		○
DHA4230		Red		○
DHA4240		Blue		○
DHA4250		Gold		○
DHA1111		No bold colour, select only with Option DHA4195 Tag or DHA4197 Without Tag	(NC 4190)	○
TAG COLOUR*				
DHA4300		Black		○
DHA4310		Red		○
DHA2222		No tag colour, select only with Option DHA4190 Bold Decal or DHA4197 Without Tag	(NC 4195)	○

FRONT CASTORS*				
DHA5000		75mm x 25mm (3" x 1") Supersport	(NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)	○
DHA5050		100mm x 30mm (4" x 1.2") Everyday Black hub and tyre	(NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)	○
DHA5100		100mm x 30mm (4" x 1.2") Comfort Black hub and tyre	(NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)	○
DHA5200		100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frogs Legs) Use this option to access 4" anodised colour options		£93
DHA5250		100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs)	(NC 5502, 5504, 5506, 5507, 5508, 5509)	£93
DHA5300		125mm x 30mm (5" x 1.2") Everyday Black hub and tyre	(NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)	○
DHA5350		125mm x 30mm (5" x 1.2") Comfort Black hub and tyre	(NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)	◆
DHA5450		125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) Use this option to access 5" anodised colour options		£93
DHA5500		125mm x 36mm (5" x 1.4") (Frog Legs) Aluminium Softroll Silver hub, Black tyre (Frogs Legs)	(NC 5502, 5504, 5506, 5507, 5508, 5509)	£93
ANODISED COLOUR SOFTROLL CASTORS				
UHA5502		Pink		£50
UHA5504		Orange		£50
UHA5506		Blue		£50
UHA5507		Purple		£50
UHA5508		Red		£50
UHA5509		Gold		£50
CERAMIC CASTOR BEARINGS				
UHA5510		Ceramic castor bearings for Aluminium soft roll castors		£58
REAR WHEEL POSITION (CENTRE OF GRAVITY 12MM INCREMENTS)*				
DHA6100		1st from front Most Active Position 102mm		○
DHA6110		2nd from front 90.5mm		○
DHA6120		3rd from front 79mm		○
DHA6130		4th from front 67.5mm		○
DHA6140		5th from front 56mm		◆
DHA6150		6th from front 44.5mm		○
DHA6160		7th from front 33mm		○
DHA6170		8th from front 21.5mm		○
DHA6180		9th from front 10mm		○
QUAD QUICK RELEASE AXLE				
DHA6300		Quad Quick Release/Tetra Clip	(NC 6550, 6551, 6552, 6553, 6554, 0907)	£85
REAR WHEEL SIZE*				
DHA6410		24"		◆
DHA6420		25"	(NC 6510, 6600, 6630, 6660, 6700, 6710, 6720, 6760, 6770, 6775, 6900)	○
DHA6430		26"	(NC 6520, 6550, 6551, 6552, 6553, 6554, 6600, 6630, 6660, 6720, 6740, 6750, 6760, 6770, 6775)	○

REAR WHEEL TYPE*			
DHA6510		24", 26" KÜschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub.	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 8888, 6579) ◆
DHA6520		24", 25" High Performance Wheel Black hub and rim	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 8888, 6579, 6660) £182
DHA6550		24", 25" Spinergy LXX, Black hub, rim and 18 spokes in Black	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6660, 8888) £611
UHA6551		24", 25" Spinergy Lite Extreme LX, Weight 687 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)	(NC 6660) £785
UHA6552		24", 25" Spinergy Lite Extreme X-laced XLX, Weight 821 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)	(NC 6660) £844
UHA6553		24", 25" Spinergy X-Laced Sport Lite Extreme XSLX, Weight 757 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)	(NC 6660) £893
UHA6554		24", 25" Spinergy Carbon Fibre CLX, Weight 530 grams per wheel (Unique coloured hub and spokes available, supplied in black as standard)	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6660) £1,084
SWO0907		24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKÜschall	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6660) £486
ANODISED COLOUR SPINERGY RIM (DESIGN LAB)			
UHA6557		Blue	£179
UHA6558		Gold	£179
UHA6559		Green	£179
UHA6560		Orange	£179
UHA6561		Red	£179
UHA6562		Purple	£179
ANODISED COLOUR SPINERGY HUB (DESIGN LAB)			
UHA6563		Blue	£179
UHA6564		Gold	£179
UHA6565		Silver	£179
UHA6566		Green	£179
UHA6567		Orange	£179
UHA6568		Red	£179
UHA6569		Purple	£179
SPINERGY COLOURED SPOKES			
UHA6570		Blue	○
UHA6571		Green	○
UHA6572		Orange	○
UHA6573		Pink	○
UHA6574		Red	○
UHA6575		Purple	○
UHA6576		White	○
UHA6577		Yellow	○
SWO8888		Custom mixed colour spokes. Please list colours required in special comments box	£149
CERAMIC BEARING REAR WHEEL			
UHA6579		Ceramic Bearing for Spinergy rear wheels	£170

REAR WHEEL TYRE*			
DHA6600		24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection	◆
DHA6625		24", 25", 26" x 1" Maxxis Alpha Protect. Lightweight, Pneumatic and highly puncture resistant.	£100
DHA6630		24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance)	£97
DHA6660		24" x 1 3/8" Solid Profile Black Tyre (NC 7230)	£66
HANDRIM*			
DHA6700		24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance	£223
DHA6710		26" Steel, short tab with 1 hole, mounted at 3 cm distance	£62
DHA6720		24" Titanium, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance	£192
DHA6740		24", 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. For complete or limited hand function offering low-effort propulsion. Ergonomic shape with fully integrated, replaceable silicon strip for extra grip, profile of 32mm diameter.	£215
DHA6750		24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. For complete or limited hand function offering low-effort propulsion. Ergonomic shape with fully integrated, replaceable silicon strip for extra grip, smaller profile of 24mm diameter for smaller hands.	£215
DHA6760		24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape, all metal black coated, low-friction surface for complete hand function for low-effort propulsion	£190
DHA6770		24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics without tricep function. Features a unique angle shaped profile and anti-slip coating providing a larger contact surface designed to be propelled using the heel of the hand or thumb pad.	£225
DHA6775		24" Carbolife Curve L Tetra handrim, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics with tricep function, but limited hand strength. Ergonomonic shape with a flat, non-slip coated pushing surface for the thumb and a curved side for the fingers to wrap around.	£155
DHA6780		24", 25", 26" Light Aluminium Grey, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance	◆
DHA6790		24", 25", 26" Light Aluminium Black, short tab with 1 hole for 3cm mounting	£42
REAR WHEEL CAMBER*			
DHA6820		0° Camber	◆
DHA6830		1° Camber	○
DHA6840		3° Camber (NC 7230)	○
SPOKEGUARD			
DHA6900		24", 26" Transparent Spokeguards (When ordered with Carbolife handrims, handrim distance is set at 4cm)	£126
FIXATION CLIP*			
DHA7000		Fixation Strap (Hook and loop strap under the seat holds the frame together when folded for easy transportation, designed for clients with limited hand function)	○
DHA7010		Fixation Clip (Plastic clip on either side of the seat holds the frame together when folded for easy transportation)	◆
DHA3333		No fixation clip, select with Tapered frame options 1170, 1180 & 2800	○

UNFOLDING HELP*			
DHA7100	☐	Unfolding Help (Cord to help open the folding mechanism from the front of the chair)	◆
UNIQUE CUSTOMISATION COLOUR PACKAGE*			
UHA6960	☐	UNIQUE KIT 1 - Includes custom colour forks and bearing cover. Select colour from the 10 choices below	£111
UHA6961	☐	UNIQUE KIT 2 - Includes custom colour forks and bearing cover, castor fork supporter, rear wheel supporter, brake clamps. Select colour from the 10 choices below	£285
UHA1234	☐	Without colour package	○
UNIQUE CUSTOMISATION COLOURS			
UHA6963	☐	Red	○
UHA6964	☐	Purple	○
UHA6965	☐	Copper	○
UHA6966	☐	Green	○
UHA6967	☐	Orange	○
UHA6968	☐	Gold	○
UHA6969	☐	Champagne	○
UHA6970	☐	Pink	○
UHA6971	☐	Blue	○
UHA6972	☐	Poly Mirror Finish	○
BRAKES*			
DHA7200	☐	Pull Brakes	£78
DHA7210	☐	Standard Push to Lock Brakes, Lever Bent	◆
DHA7220	☐	Standard Push to Lock Brakes, Lever Straight	○
DHA7230	☐	Performance Brakes, Push to Lock. Brake barrel flips away from tyre (NC 7300, 7810, 7820, 7830)	£78
DHA7240	☐	Active/Scissor Brakes (NC 7300)	£78
DHA7250	☐	Light/Compact Active, Aluminium Brakes (NC 7300)	£136
BRAKE LEVER EXTENSION			
DHA7300	☐	Brake Lever Extension	£49
SAFETY OPTIONS			
DHA7500	☐	Anti-Tipper Right	£122
DHA7510	☐	Anti-Tipper Left	£122
DHA7520	☐	Anti-Tipper Both sides	£244
DHA7800	☐	Standard Pelvic belt	£74
BODYPOINT BELTS			
DHA7810	☐	Bodypoint Posture Belt Evoflex Small	£81
DHA7820	☐	Bodypoint Posture Belt Evoflex Medium	£81
DHA7830	☐	Bodypoint Posture Belt Evoflex Large	£81

OTHER OPTIONS			
DHA8010		Front Frame Protectors	£37
DHA8100		Tipper-Aid, Mounted Right Side Max. 2 adaption's per side (Adaptations = Anti tipper, Tipper aid and Cane holder)	£46
DHA8200		Cane Holder, Mounted Right Side Max. 2 adaption's per side (Adaptations = Anti tipper, Tipper aid and cane holder)	£75
DHA8500		Air Pump	£13
DHA8550		Tool Kit	◆
DHA8650		Kuschall branded Wheel bag (Fits all wheel sizes)	£45
DHA8680		24" Mountain Bike Wheels with bag. Fitted with Specific Black Tyres, Grey Aluminium Handrims; Quick Release Axles Included	£320
DHA8700		Travel bag for Champion, Wheels can remain on the wheelchair	£254
BODYPOINT ACCESSORIES			
DHA8800		Bodypoint Calf Strap	£57
DHA8850		Bodypoint Mobility Bag	£55

			TOTAL PRICE	
--	--	--	--------------------	--

Non Conforming Options:	
UHA1260	Complete frame with fully welded front seat height & centre of gravity rear wheel supports
DHA1320	360 mm
DHA1325	380 mm
DHA1350	480 mm
DHA1946	430 mm
DHA1948	440 mm
DHA1950	450 mm
DHA1952	460 mm
DHA1954	470 mm
DHA1956	480 mm
DHA1958	490 mm
DHA1960	500 mm
DHA2300	Matrx MX2, Height 150 mm with Q/R Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.
DHA2305	Matrx MX2, Height 230 mm with Q/R Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.
DHA2310	Matrx MX2, Height 305 mm with Q/R Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.
DHA2315	Matrx MX2, Height 405 mm with Q/R Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.
DHA2325	Matrx Elite Deep E2, Height 410 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.
DHA2330	Matrx Elite Deep E2, Height 460 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.
DHA2331	Matrx Elite Deep E2, Height 510 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.
DHA2345	Matrx Elite E2, Height 360 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.
DHA2350	Matrx Elite E2, Height 410 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.
DHA2355	Matrx Elite E2, Height 460 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.
DHA2357	Matrx Elite E2, Height 510 mm with Easy-Set Hardware. Ensure the top of the Matrx backrest is no higher than 465mm, the highest backrest setting.
DHA2360	360mm Backrest width, to fit chest width 265mm
DHA2362	380mm Backrest width, to fit chest width 290mm
DHA2364	410mm Backrest width, to fit chest width 315mm
DHA2366	430mm Backrest width, to fit chest width 340mm
DHA2368	460mm Backrest width, to fit chest width 370mm
DHA2370	480mm Backrest width, to fit chest width 390mm
DHA2372	510mm Backrest width, to fit chest width 420mm
DHA2400	Without Push Handles
DHA2710	Flo-Tech Lite Cushion 78mm high- Comfort/Low Risk
DHA2720	Matrx Libra Cushion 70mm high, Black Smooth Cover - Very High Risk
DHA2748	Contour Visco NG Cushion 120mm high - High Risk
DHA2750	Vicair Active O2 lightweight cushion 60mm high - High risk
DHA2810	Fluid supplement sac (Single unit)
DHA2820	Tapered wedge (Pair)
DHA2830	Abductor
DHA2840	Lateral side wedge

DHA2850	Anterior front wedge
DHA2815	Fluid supplement sac (Single unit)
DHA2825	Tapered wedge (Pair)
DHA2835	Abductor
DHA2845	Lateral side wedge
DHA2855	Anterior front wedge
DHA2860	Seat rigidizer
DHA3490	Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included
DHA4000	Space Black (RAL 9005)
DHA4010	Anthracite Grey (RAL 7016)
DHA4020	Pure White (RAL 9016)
DHA4030	Carmin Red (RAL 3002)
DHA4040	Azure Blue (RAL 5017)
DHA4060	Flashy Yellow (RAL 1516) (Additional 15 day lead time)
DHA4100	Black Glossy (RAL 9005) (Additional 15 day lead time)
DHA4130	Silver Grey (RAL 9006) (Additional 15 day lead time)
DHA4145	Pearl Night Blue (RAL 5026) (Additional 15 day lead time)
DHA4150	Ruby Red (RAL 1503) (Additional 15 day lead time)
DHA4170	Dynamite Orange (Additional 15 day lead time)
DHA4180	Forest Green (RAL 6029) (Additional 15 day lead time)
UHA4101	Diamond Silver/Mirror Polish with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066)
UHA4102	Pastel Green with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066)
UHA4103	Sky Blue with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066)
UHA4104	Plum Purple with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066)
DHA4185	Carbon Finish
UHA4066	Ceramic glaze finish on frame, for added durability
DHA4190	Bold Decal (Full coloured transfer on frame)
DHA4195	Tag (Small Kuschall brand transfer on front of frame)
DHA4197	Without Tag (Kuschall branding only)
DHA4200	Black
DHA4210	Dark Grey
DHA4220	Light Grey
DHA4230	Red
DHA4240	Blue
DHA4250	Gold
DHA1111	No bold colour, select only with Option DHA4195 Tag or DHA4197 Without Tag
DHA4300	Black
DHA4310	Red
DHA2222	No tag colour, select only with Option DHA4190 Bold Decal or DHA4197 Without Tag
DHA5000	75mm x 25mm (3" x 1") Supersport
DHA5050	100mm x 30mm (4" x 1.2") Everyday Black hub and tyre
DHA5100	100mm x 30mm (4" x 1.2") Comfort Black hub and tyre
DHA5200	100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frogs Legs) Use this option to access 4" anodised colour options

DHA5250	100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs)
UHA5502	Pink
UHA5504	Orange
UHA5506	Blue
UHA5507	Purple
UHA5508	Red
UHA5509	Gold
UHA5510	Ceramic castor bearings for Aluminium soft roll castors
DHA6180	9th from front 10mm
DHA6410	24"
DHA6420	25"
DHA6430	26"
DHA6510	24", 26" Kuschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub.
DHA6520	24", 25" High Performance Wheel Black hub and rim
DHA6550	24", 25" Spinergy LXX, Black hub, rim and 18 spokes in Black
UHA6551	24", 25" Spinergy Lite Extreme LX, Weight 687 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)
UHA6552	24", 25" Spinergy Lite Extreme X-laced XLX, Weight 821 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)
UHA6553	24", 25" Spinergy X-Laced Sport Lite Extreme XSLX, Weight 757 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)
UHA6554	24", 25" Spinergy Carbon Fibre CLX, Weight 530 grams per wheel (Unique coloured hub and spokes available, supplied in black as standard)
SWO0907	24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKuschall
UHA6557	Blue
UHA6558	Gold
UHA6559	Green
UHA6560	Orange
UHA6561	Red
UHA6562	Purple
UHA6563	Blue
UHA6564	Gold
UHA6565	Silver
UHA6566	Green
UHA6567	Orange
UHA6568	Red
UHA6569	Purple
UHA6570	Blue
UHA6571	Green
UHA6572	Orange
UHA6573	Pink
UHA6574	Red
UHA6575	Purple
UHA6576	White
UHA6577	Yellow
SWO8888	Custom mixed colour spokes. Please list colours required in special comments box

UHA6579	Ceramic Bearing for Spinergy rear wheels
DHA6600	24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection
DHA6630	24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance)
DHA6660	24" x 1 3/8" Solid Profile Black Tyre
DHA6700	24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance
DHA6710	26" Steel, short tab with 1 hole, mounted at 3 cm distance
DHA6720	24" Titanium, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance
DHA6740	24", 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. For complete or limited hand function offering low-effort propulsion. Ergonomic shape with fully integrated, replaceable silicon strip for extra grip, profile of 32mm diameter.
DHA6750	24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. For complete or limited hand function offering low-effort propulsion. Ergonomic shape with fully integrated, replaceable silicon strip for extra grip, smaller profile of 24mm diameter for smaller hands.
DHA6760	24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape, all metal black coated, low-friction surface for complete hand function for low-effort propulsion
DHA6770	24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics without tricep function. Features a unique angle shaped profile and anti-slip coating providing a larger contact surface designed to be propelled using the heel of the hand or thumb pad.
DHA6775	24" Carbolife Curve L Tetra handrim, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics with tricep function, but limited hand strength. Ergonomonic shape with a flat, non-slip coated pushing surface for the thumb and a curved side for the fingers to wrap around.
DHA6900	24", 26" Transparent Spokeguards (When ordered with Carbolife handrims, handrim distance is set at 4cm)
DHA7000	Fixation Strap (Hook and loop strap under the seat holds the frame together when folded for easy transportation, designed for clients with limited hand function)
DHA7010	Fixation Clip (Plastic clip on either side of the seat holds the frame together when folded for easy transportation)
DHA7230	Performance Brakes, Push to Lock. Brake barrel flips away from tyre
DHA7300	Brake Lever Extension
DHA7810	Bodypoint Posture Belt Evoflex Small
DHA7820	Bodypoint Posture Belt Evoflex Medium
DHA7830	Bodypoint Posture Belt Evoflex Large



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):