

Küschall®
CHAMPION



UNLIMIT YOUR WORLD

We at Kuschall believe anything is possible when you have the power to achieve it. The Kuschall Champion offers a superior solution for you. Unlimit your adventures, unlimit your performance, **Unlimit your world.**

KÜSCHALL CHAMPION

Kuschall Champion is the **unique** wheelchair that will give you the driving experience of a rigid frame with the convenience of a folding. It's the ultimate combination of modern design with proven **hybrid** performance. The SK (Swiss knife) option will make your Kuschall Champion **super compact** to the size of a backpack, the perfect travel companion for your adventures.

HYDROFORMING

Strength, lightness and rigidity. The hydroformed front frame on the Kuschall Champion enhances your driving performance, increasing the stiffness and enabling effortless transfers with the ergonomic grip.



HYDROFORMED FRONT FRAME

- Optimised driving performance due to lower flex (higher stiffness) in knee area.

LIGHT ACTIVE BRAKE

- Minimalistic design, precision even in the smallest details.

SUPER COMPACT

- The SK (Swiss knife) option enables the Champion to fold down to the size of a backpack



Küschall® CHAMPION

CARBOLIFE HANDRIMS

- The innovative handrim profile enables you to gain control, the low effort grip reduces the stresses during propelling.

FOLDING MECHANISM

- Fold the wheelchair and the backrest effortlessly in one movement.

UNIQUE FOLDING UNIT

- Quick and easy centre of gravity adjustments ensure you're sitting at the optimum position for your preferred activity level.

KEEP ON TRACK

- Intuitive castor housing adjustment, keeps the Champion rolling smooth and steady.

KEY FEATURES

UNIQUE

As well as being small, practical and configurable, the Champion is also adjustable.

9 position centre of gravity

Rear and Front seat height

Backrest height and angle

Multiple adjustments ensures the Champion remains the perfect fit around you.



HYBRID

The folding mechanism is the reason the Champion has a superior driving performance vs other foldable wheelchairs. Perfectly engineered, to deliver the rigidity, appearance and driving performance of a rigid wheelchair, with the practicality of a folding.

SUPER COMPACT

Simply fold the main frame with the unique designed folding mechanism, for compact and easy transportation or travel. Get it super compact with the SK (Swiss knife-folding) option, enables the Champion to be folded even smaller.




HYDRO-FORMED FRONT FRAME

The newly designed hydroformed front frame provides an ergonomic grip to support comfortably while repositioning and during transfers.

CUSTOMISATION

1 CHOOSE THE FRAME COLOUR



Pure White



Anthracite Grey



Space Black



Azure Blue



Carmine Red

2 CHOOSE THE DECAL STYLE



3 CHOOSE THE DECAL COLOUR



Black



Red



Light Grey



Gold



Red



Blue

















Dark Grey




Black

TECHNICAL SPECIFICATION

SEAT WIDTH 360 - 480 mm In increments of 20 mm 	SEAT DEPTH 340 - 480 mm In increments of 20 mm 	KNEE-TO-HEEL LENGTH 320 - 500 mm Stepless adjustable 	BACKREST HEIGHT 300 - 465 mm In increments of 15 mm 	BACKREST ANGLE 76° / 80.5° / 85° / 89.5° / 94° 
SEAT HEIGHT Front: 450 - 540 mm Rear: 390 - 480 mm 	FRAME ANGLE 75° / 90° Straight or tapered 	CENTRE OF GRAVITY 9 positions 	CAMBER 0° / 1° / 3° 	TOTAL LENGTH 75° approx. 930 mm 90° approx. 830 mm 
TOTAL WIDTH SW + 160 - 240 mm 	TOTAL WIDTH FOLDED Approx. 280 - 340 mm 	TOTAL WEIGHT UNLOADED Approx. 9.2 kg 	TRANSPORT WEIGHT Approx. 6.7 kg 	MAX. USER WEIGHT 120 kg 

All measurements and weights indicated are based on seat width 400 mm and seat depth 400 mm in the lightest and most compact configuration. Weight and measurements can change with the different configurations. For further information please check the user manual on our website.



/ UNLIMIT YOUR WORLD

Invacare Limited
 Pencoed Technology Park - Pencoed
 Bridgend CF35 5AQ - United Kingdom
 Tel: +44 1656 776 222
 Fax: +44 1656 776 220
 E-mail: uk@invacare.com
 Sales Order E-mail: ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
 Unit 5 - Seatown Business Campus
 Seatown Road - SWORDS - County Dublin - Ireland
 Tel. +353 1 8107084
 Fax +353 1 8107085
 Email: Ireland@invacare.com
www.invacare.ie

Küschall, Küschall Unlimit your world, K, Unlimit your world and Champion are trademarks of Invacare International GmbH. Every effort has been made to ensure that the contents of this publication are fully up-to-date at the time of going to print. As part of its ongoing improvement of products, Invacare International GmbH reserves the right to modify existing models at any time. © July 2023, Invacare International GmbH. All rights reserved. UK - 1632213