



Yes, you can:

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Manual Invacare Ireland Prescription Form

07 JANUARY 2026
 Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

Kuschall Champion SK 2.0

Standard Features

- Max User Weight 120kg (19 Stone)
- Active, Lightweight Folding Fully Hydroformed Aluminium Wheelchair
- Multiple adjustments, Front and Rear Seat Heights and Centre of Gravity
- Aluminium Folding Front Frame with 75° or 90° Front Frame Angles
- 24" Self Propel Q/R Wheels with Grey Aluminium Handrims and Pneumatic Tyres
- Choice of 6 Seat Widths 360mm to 460mm (14" - 18")
- Choice of 8 Seat Depths 340mm to 480mm (13" - 19")
- Choice of 5 Standard Frame Colours
- Range of Front Seat Heights 450mm to 540mm (17" - 21")
- Range of Rear Seat Heights 390mm to 480mm (15" - 19")
- Folding Backrest, Angle Adjustable and Soft Lock Feature with Nylon Light Upholstery
- Angle and Height Adjustable Rigid Footrest, with Calf Strap
- Without Push Handles
- Flip up Carbon Clothesguards
- Unfolding Help Cord
- Black Comfort Front Castor
- Standard Brake Lever Bent
- Tool Kit



Standard Base price: €4,312

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

| Code | Description | |
|---|--|-----|
| BASE PRODUCT* | | |
| DHA0031 | Kuschall Champion SK 2.0 | ○ |
| UNIQUE SPECIALS (For special modifications only) Additional lead time of 4 weeks | | |
| DHA0000 | UNIQUE Special Request - For use when requesting any special options or builds | POA |
| UNIQUE CUSTOMISATION AND PERSONALISATION OPTIONS Additional 4 weeks lead time for Unique options | | |
| UHA0001 | UNIQUE Personalisation - For use when selecting UHA customisation options | ○ |

| FRAME* | | | |
|---------------------------------------|--|--|--------|
| DHA1150 | | Frame Angle 75° Aluminium | ○ |
| DHA1160 | | Frame Angle 90° Aluminium | ○ |
| SK FOLDING FRONT* | | | |
| DHA1250 | | SK Folding Front Frame | €4,312 |
| CHAMPION SL HYBRID RIGID FRAME | | | |
| UHA1260 | | Complete frame with fully welded front seat height & centre of gravity rear wheel supports (NC 6180) | €578 |
| SEAT WIDTH* | | | |
| DHA1320 | | 360 mm (NC 2710, 2748, 7000) | ○ |
| DHA1325 | | 380 mm (NC 7000) | ○ |
| DHA1330 | | 400 mm | ○ |
| DHA1335 | | 420 mm | ○ |
| DHA1340 | | 440 mm | ○ |
| DHA1345 | | 460 mm | ○ |
| SEAT DEPTH* | | | |
| DHA1405 | | 340 mm (NC 2710, 2720, 2748, 2860, 7810, 7820, 7830) | ○ |
| DHA1415 | | 360 mm (NC 2710, 2720, 2748, 2860, 7810, 7820, 7830) | ○ |
| DHA1425 | | 380 mm (NC 2720) | ○ |
| DHA1435 | | 400 mm | ○ |
| DHA1445 | | 420 mm | ○ |
| DHA1460 | | 440 mm | ○ |
| DHA1470 | | 460 mm | ○ |
| DHA1480 | | 480 mm (NC 2710, 2748) | ○ |
| BACK HEIGHT* | | | |
| DHA1615 | | 300 mm (NC 2200, 2400) | ○ |
| DHA1620 | | 315 mm (NC 2200, 2400) | ○ |
| DHA1625 | | 330 mm | ○ |
| DHA1630 | | 345 mm | ○ |
| DHA1635 | | 360 mm | ○ |
| DHA1640 | | 375 mm | ○ |
| DHA1645 | | 390 mm | ○ |
| DHA1650 | | 405 mm | ○ |
| DHA1655 | | 420 mm | ○ |
| DHA1660 | | 435 mm | ○ |
| DHA1665 | | 450 mm | ○ |
| DHA1670 | | 465 mm | ○ |
| FRONT SEAT HEIGHT* | | | |
| DHA1740 | | 450 mm (NC 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960) | ○ |
| DHA1745 | | 460 mm (NC 1948, 1950, 1952, 1954, 1956, 1958, 1960) | ○ |
| DHA1750 | | 470 mm (NC 1950, 1952, 1954, 1956, 1958, 1960) | ○ |
| DHA1755 | | 480 mm (NC 1952, 1954, 1956, 1958, 1960) | ○ |
| DHA1760 | | 490 mm (NC 1954, 1956, 1958, 1960) | ○ |
| DHA1765 | | 500 mm (NC 1956, 1958, 1960) | ○ |
| DHA1770 | | 510 mm (NC 1958, 1960) | ○ |
| DHA1775 | | 520 mm (NC 1960) | ○ |
| DHA1780 | | 530 mm (NC 5000) | ○ |
| DHA1785 | | 540 mm (NC 5000, 5050, 5100, 5150, 5200, 5250) | ○ |

REAR SEAT HEIGHT. On 75° frame, Rear Seat Height can be 0mm - 70mm lower then front seat height and 90° can be 30mm to 100mm lower*

| | | | | |
|---------|--|--------|-----------------|---|
| DHA1815 | | 390 mm | (NC 6420, 6430) | ○ |
| DHA1820 | | 400 mm | (NC 6420, 6430) | ○ |
| DHA1825 | | 410 mm | (NC 6430) | ○ |
| DHA1830 | | 420 mm | (NC 6430) | ○ |
| DHA1835 | | 430 mm | | ○ |
| DHA1840 | | 440 mm | | ○ |
| DHA1845 | | 450 mm | | ○ |
| DHA1850 | | 460 mm | | ○ |
| DHA1855 | | 470 mm | (NC 6410) | ○ |
| DHA1860 | | 480 mm | (NC 6410) | ○ |

LOWER LEG LENGTH Seat canvas to footplate measurement, cushion height must be deducted.*

| | | | | |
|---------|--|--------|--|---|
| DHA1936 | | 380 mm | | ○ |
| DHA1938 | | 390 mm | | ○ |
| DHA1940 | | 400 mm | | ○ |
| DHA1942 | | 410 mm | | ○ |
| DHA1944 | | 420 mm | | ○ |
| DHA1946 | | 430 mm | | ○ |
| DHA1948 | | 440 mm | | ○ |
| DHA1950 | | 450 mm | | ○ |
| DHA1952 | | 460 mm | | ○ |
| DHA1954 | | 470 mm | | ○ |
| DHA1956 | | 480 mm | | ○ |
| DHA1958 | | 490 mm | | ○ |
| DHA1960 | | 500 mm | | ○ |

BACKREST*

| | | | | |
|---------|--|--|--|---|
| DHA2140 | | Angle Adjustable, Forward Folding with Soft Lock. Backrest remains locked for transportation | | ◆ |
|---------|--|--|--|---|

BACKREST UPHOLSTERY*

| | | | | |
|---------|--|--|-----------|------|
| DHA2200 | | Black Standard Light Upholstery, without Tension Adjustment | (NC 2430) | ◆ |
| DHA2220 | | Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener | | €160 |
| DHA2230 | | Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener | | €253 |
| DHA2240 | | Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required | (NC 2430) | ○ |

BACK UPHOLSTERY STITCHING COLOUR

| | | | | |
|---------|--|--|--|-----|
| UHA2241 | | Neon Yellow | | €96 |
| UHA2242 | | Neon Pink | | €96 |
| UHA2243 | | Blue | | €96 |
| UHA2244 | | White | | €96 |
| UHA2245 | | Red | | €96 |
| UHA0002 | | Your personal colour stitching. Add RAL colour of your choice in special comments box, for a personalised colour stitching on back upholstery. | | €96 |

| BACK UPHOLSTERY EMBROIDERY | | | |
|--|--|---|--|
| UHA0003 | | Your personal design embroidery. Supply your design with your order, subject to copyright. Design size up to 250mm x 250mm | (NC 0004) €256 |
| BACK UPHOLSTERY FLOCKING | | | |
| UHA0004 | | Your personal flock design. Supply your design with your order, subject to copyright. Design size up to 280mm x 400mm | €127 |
| PUSH HANDLES* | | | |
| DHA2400 | | Without Push Handles | ◆ |
| DHA2410 | | Push Handles Fixed Short/Mini Grip | €149 |
| DHA2420 | | Push Handles Fixed Long Grip | ○ |
| DHA2430 | | Push Handles Height Adjustable, Rear Set | (NC 3480) €294 |
| DHA2440 | | Push Handles Height Adjustable, Integrated in backrest | €225 |
| DHA2450 | | Push Handles Foldable | €287 |
| BACKREST ANGLE* | | | |
| DHA2510 | | 76° angle from the seat | ○ |
| DHA2530 | | 80.5° angle from the seat | ○ |
| DHA2550 | | 85° angle from the seat | ◆ |
| DHA2570 | | 89.5° angle from the seat | ○ |
| DHA2590 | | 94° angle from the seat | ○ |
| SEAT UPHOLSTERY* | | | |
| DHA2600 | | Seat Upholstery with 1 central loop | ○ |
| DHA2601 | | Standard Seat Upholstery with 1 central loop and 2 side hook | ◆ |
| DHA2610 | | Tension Adjustable Seat Upholstery with one 5cm central band (Loop) | €54 |
| SEAT CUSHION | | | |
| DHA2700 | | Standard Seat Cushion, 60mm High Double Foam with Visco Top Layer - Comfort/Low Risk | (NC 2810, 2820, 2830, 2840, 2850) €74 |
| DHA2710 | | Flo-Tech Lite Cushion 78mm High - Comfort/Low Risk | (NC 2810, 2820, 2830, 2840, 2850) €168 |
| DHA2720 | | Matrx Libra Cushion, Black Smooth Cover 70mm High - Very High Risk | €645 |
| DHA2748 | | Contour Visco NG Cushion 100mm High - High Risk | (NC 2810, 2820, 2830, 2840, 2850) €236 |
| DCA2750 | | Vicair Active O2 lightweight cushion 60mm High - High risk | (NC 2810, 2820, 2830, 2840, 2850, 2860) €731 |
| SEAT BAG | | | |
| DHA2800 | | Seat Bag, Click in and removable, includes the fixation clips | (NC 7010) €85 |
| MATRX LIBRA CUSHION ACCESSORIES | | | |
| DHA2810 | | Fluid supplement sac (Pair) | €123 |
| DHA2820 | | Tapered wedge (Pair) | €28 |
| DHA2830 | | Abductor | €58 |
| DHA2840 | | Lateral side wedge | €46 |
| DHA2850 | | Anterior front wedge | €43 |
| SEAT RIGIDIZER | | | |
| DHA2860 | | Seat rigidizer | €79 |
| FOOTREST* | | | |
| DHA3240 | | Rigid, Plastic Footplate, Set to front, Height and Angle Adjustable. Note: This does not flip up for transfer, but folds easy for transport and storage | ◆ |
| DHA3250 | | Rigid, Plastic Footplate, Set to rear, Height and Angle Adjustable. Note: This does not flip up for transfer, but folds easy for transport and storage | ○ |

| CALF STRAP* | | | |
|----------------------------|--|--|------|
| DHA3320 | | Calf Strap Mounted on Supporter | ◆ |
| DHA3330 | | Calf Strap Mounted on Frame | ○ |
| CLOTHESGUARDS/MUDGUARDS* | | | |
| DHA3405 | | Composite Clothesguards, Flip Back - (only possible with 24" and 25" rear wheels) (NC 6430) | ○ |
| DHA3420 | | Carbon Clothesguards, Flip Up | ◆ |
| DHA3425 | | Composite Mudguards, Fixed (only possible with 24" and 25" rear wheels) (NC 3490, 6430) | €91 |
| DHA3430 | | Carbon Mudguards, Fixed (NC 3490) | €192 |
| DHA3440 | | Carbon Mudguards, Removable (NC 3490) | €327 |
| DHA3450 | | Composite Mudguards, Removable (only possible with 24" & 25" rear wheels) (NC 3490, 6430) | €91 |
| ARMREST | | | |
| DHA3480 | | Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included | €281 |
| DHA3490 | | Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included | €109 |
| FRAME COLOUR* | | | |
| DHA4000 | | Space Black (RAL 9005) (NC 4066) | ○ |
| DHA4010 | | Anthracite Grey (RAL 7016) (NC 4066) | ○ |
| DHA4020 | | Pure White (RAL 9016) (NC 4066) | ○ |
| DHA4030 | | Carmine Red (RAL 3002) (NC 4066) | ○ |
| DHA4040 | | Azure Blue (RAL 5017) (NC 4066) | ○ |
| DHA4060 | | Flashy Yellow (RAL 1516) (Additional 15 day lead time) (NC 4066) | €39 |
| DHA4100 | | Black Glossy (RAL 9005) (Additional 15 day lead time) (NC 4066) | €39 |
| DHA4130 | | Silver Grey (RAL 9006) (Additional 15 day lead time) (NC 4066) | €39 |
| DHA4145 | | Pearl Night Blue (RAL 5026) (Additional 15 day lead time) (NC 4066) | €39 |
| DHA4150 | | Ruby Red (RAL 1503) (Additional 15 day lead time) (NC 4066) | €39 |
| DHA4170 | | Dynamite Orange (Additional 15 day lead time) (NC 4066) | €39 |
| DHA4180 | | Forest Green (RAL 6029) (Additional 15 day lead time) (NC 4066) | €39 |
| UHA4101 | | Diamond Silver/Mirror Polish with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066) | €513 |
| UHA4102 | | Pastel Green with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066) | €413 |
| UHA4103 | | Sky Blue with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066) | €413 |
| UHA4104 | | Plum Purple with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066) | €413 |
| UHA4105 | | Salmon with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066) | €413 |
| UHA0008 | | Your personal frame colour. Add RAL colour of your choice in special comments box, for a personalised colour frame | €384 |
| CERAMIC GLAZE FRAME FINISH | | | |
| UHA4066 | | Ceramic glaze finish on frame, for added durability | €432 |
| FRAME DECAL TYPE* | | | |
| DHA4195 | | Tag (Small Kuschall brand transfer on front of frame) (NC 2222) | ◆ |
| DHA4197 | | Without Tag (Kuschall branding only) (NC 4300, 4310) | ○ |

| TAG COLOUR* | | | |
|--|--|--|---|
| DHA4300 | | Black | ○ |
| DHA4310 | | Red | ○ |
| DHA2222 | | No tag colour, select only with Option DHA4197 Without Tag | (NC 4195) ○ |
| FRONT CASTORS* | | | |
| DHA5000 | | 75mm x 25mm (3" x 1") Supersport | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) ○ |
| DHA5050 | | 100mm x 30mm (4" x 1.2") Everyday Black hub and tyre | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) ○ |
| DHA5100 | | 100mm x 30mm (4" x 1.2") Comfort Black hub and tyre | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) ○ |
| DHA5150 | | 100mm x 25mm (4" x 1") Starec Black hub and tyre | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) €102 |
| DHA5200 | | 100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frogs Legs) Use this option to access 4" anodised colour options | €114 |
| DHA5250 | | 100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs) | (NC 5502, 5504, 5506, 5507, 5508, 5509) €114 |
| DHA5300 | | 125mm x 30mm (5" x 1.2") Everyday Black hub and tyre | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) ○ |
| DHA5350 | | 125mm x 30mm (5" x 1.2") Comfort Black hub and tyre | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) ◆ |
| DHA5400 | | 125mm x 25mm (5" x 1") Starec Black hub and tyre | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) €102 |
| DHA5450 | | 125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) Use this option to access 5" anodised colour options | €114 |
| DHA5500 | | 125mm x 36mm (5" x 1.4") (Frog Legs) Aluminium Softroll Silver hub, Black tyre (Frogs Legs) | (NC 5502, 5504, 5506, 5507, 5508, 5509) €114 |
| ANODISED COLOUR SOFTROLL CASTORS | | | |
| UHA5502 | | Pink | €64 |
| UHA5504 | | Orange | €64 |
| UHA5506 | | Blue | €64 |
| UHA5507 | | Purple | €64 |
| UHA5508 | | Red | €64 |
| UHA5509 | | Gold | €64 |
| CERAMIC CASTOR BEARINGS | | | |
| UHA5510 | | Ceramic castor bearings for Aluminium soft roll castors | €64 |
| REAR WHEEL POSITION (CENTRE OF GRAVITY)* | | | |
| DHA6100 | | 1st from front Most Active Position | ○ |
| DHA6110 | | 2nd from front | ○ |
| DHA6120 | | 3rd from front | ○ |
| DHA6130 | | 4th from front | ○ |
| DHA6140 | | 5th from front | ◆ |
| DHA6150 | | 6th from front | ○ |
| DHA6160 | | 7th from front | ○ |
| DHA6170 | | 8th from front | ○ |
| DHA6180 | | 9th from front | ○ |
| QUAD QUICK RELEASE AXLE | | | |
| DHA6300 | | Quad Quick Release/Tetra Clip | (NC 6550, 6551, 6552, 6553, 6554, 0907) €96 |
| REAR WHEEL SIZE* | | | |
| DHA6410 | | 24" | ◆ |
| DHA6420 | | 25" | (NC 6510, 6600, 6630, 6660, 6700, 6710, 6720, 6760, 6770, 6775, 6900) ○ |
| DHA6430 | | 26" | (NC 6520, 6550, 6551, 6552, 6553, 6554, 6600, 6630, 6660, 6720, 6740, 6750, 6760, 6770, 6775) ○ |

| REAR WHEEL TYPE* | | | |
|--|--|---|--|
| DHA6510 | | 24", 26" Kuschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub. | (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6579) ◆ |
| DHA6520 | | 24", 25" High Performance Wheel Black hub and rim | (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6579, 6660) €311 |
| DHA6550 | | 24", 25" Spinergy LXX, Black hub, rim and 18 spokes in Black | (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6660) €746 |
| UHA6551 | | 24", 25" Spinergy Lite Extreme LX, Weight 687 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard) | (NC 6660) €1,009 |
| UHA6552 | | 24", 25" Spinergy Lite Extreme X-laced XLX, Weight 821 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard) | (NC 6660) €1,084 |
| UHA6553 | | 24", 25" Spinergy X-Laced Sport Lite Extreme XSLX, Weight 757 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard) | (NC 6660) €1,147 |
| UHA6554 | | 24", 25" Spinergy Carbon Fibre CLX, Weight 530 grams per wheel (Unique coloured hub and spokes available, supplied in black as standard) | (NC 6557, 6558, 6559, 6560, 6561, 6562, 6660) €1,393 |
| SWO0907 | | 24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. (Unique coloured spokes available) | (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6660) €543 |
| ANODISED COLOUR SPINERGY RIM (DESIGN LAB) | | | |
| UHA6557 | | Blue | €230 |
| UHA6558 | | Gold | €230 |
| UHA6559 | | Green | €230 |
| UHA6560 | | Orange | €230 |
| UHA6561 | | Red | €230 |
| UHA6562 | | Purple | €230 |
| ANODISED COLOUR SPINERGY HUB (DESIGN LAB) | | | |
| UHA6563 | | Blue | €230 |
| UHA6564 | | Gold | €230 |
| UHA6565 | | Silver | €230 |
| UHA6566 | | Green | €230 |
| UHA6567 | | Orange | €230 |
| UHA6568 | | Red | €230 |
| UHA6569 | | Purple | €230 |
| SPINERGY COLOURED SPOKES | | | |
| UHA6570 | | Blue | €108 |
| UHA6571 | | Green | €108 |
| UHA6572 | | Orange | €108 |
| UHA6573 | | Pink | €108 |
| UHA6574 | | Red | €108 |
| UHA6575 | | Purple | €108 |
| UHA6576 | | White | €108 |
| UHA6577 | | Yellow | €108 |
| CERAMIC BEARING REAR WHEEL | | | |
| UHA6579 | | Ceramic Bearing for Spinergy rear wheels | €218 |

| REAR WHEEL TYRE* | | | |
|---------------------------|--|---|------|
| DHA6600 | | 24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection | ◆ |
| DHA6625 | | 24", 25", 26" x 1" Maxxis Alpha Protect. Lightweight, Pneumatic and highly puncture resistant. | €113 |
| DHA6630 | | 24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance) | €96 |
| DHA6660 | | 24" x 1 3/8" Solid Profile Black Tyre (NC 7230) | €82 |
| HANDRIM* | | | |
| DHA6700 | | 24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance | €252 |
| DHA6710 | | 24", 26" Steel, short tab with 1 hole, mounted at 3 cm distance | €76 |
| DHA6720 | | 24" Titanium, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance | €200 |
| DHA6740 | | 24", 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. For complete or limited hand function offering low-effort propulsion. Ergonomic shape with fully integrated, replaceable silicon strip for extra grip, profile of 32mm diameter. | €280 |
| DHA6750 | | 24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. For complete or limited hand function offering low-effort propulsion. Ergonomic shape with fully integrated, replaceable silicon strip for extra grip, smaller profile of 24mm diameter for smaller hands. | €280 |
| DHA6760 | | 24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape, all metal black coated, low-friction surface for complete hand function for low-effort propulsion | €250 |
| DHA6770 | | 24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics without tricep function. Features a unique angle shaped profile and anti-slip coating providing a larger contact surface designed to be propelled using the heel of the hand or thumb pad. | €290 |
| DHA6775 | | 24" Carbolife Curve L Tetra handrim, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics with tricep function, but limited hand strength. Ergonomic shape with a flat, non-slip coated pushing surface for the thumb and a curved side for the fingers to wrap around. | €200 |
| DHA6780 | | 24", 25", 26" Light Aluminium Grey, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance | ◆ |
| DHA6790 | | 24", 25", 26" Light Aluminium Black, short tab with 1 hole for 3cm mounting | €51 |
| REAR WHEEL CAMBER* | | | |
| DHA6820 | | 0° Camber | ◆ |
| DHA6830 | | 1° Camber | ○ |
| DHA6840 | | 3° Camber (NC 7230) | ○ |
| SPOKEGUARD | | | |
| DHA6900 | | 24", 26" Transparent Spokeguards (When ordered with Carbolife handrims, handrim distance is set at 4cm) | €155 |
| FIXATION CLIP* | | | |
| DHA7000 | | Fixation Strap (Hook and loop strap under the seat holds the frame together when folded for easy transportation, designed for clients with limited hand function) | ○ |
| DHA7010 | | Fixation Clip (Plastic clip on either side of the seat holds the frame together when folded for easy transportation) | ◆ |
| DHA3333 | | No fixation clip, select with Tapered frame options 1170, 1180 | ○ |

| UNFOLDING HELP* | | | |
|---|--------------------------|---|---------------------------------|
| DHA7100 | <input type="checkbox"/> | Unfolding Help (Cord to help open the folding mechanism from the front of the chair) | ◆ |
| UNIQUE CUSTOMISATION COLOUR PACKAGE* | | | |
| UHA6960 | <input type="checkbox"/> | UNIQUE KIT 1 - Includes custom colour forks and bearing cover. Select colour from the 10 choices below | €167 |
| UHA6961 | <input type="checkbox"/> | UNIQUE KIT 2 - Includes custom colour forks and bearing cover, castor fork supporter, rear wheel supporter, brake clamps & SK front frame folding brackets. Select colour from the 10 choices below | €366 |
| UHA1234 | <input type="checkbox"/> | Without colour package | ○ |
| UNIQUE CUSTOMISATION COLOURS | | | |
| UHA6963 | <input type="checkbox"/> | Red | ○ |
| UHA6964 | <input type="checkbox"/> | Purple | ○ |
| UHA6965 | <input type="checkbox"/> | Copper | ○ |
| UHA6966 | <input type="checkbox"/> | Green | ○ |
| UHA6967 | <input type="checkbox"/> | Orange | ○ |
| UHA6968 | <input type="checkbox"/> | Gold | ○ |
| UHA6969 | <input type="checkbox"/> | Champagne | ○ |
| UHA6970 | <input type="checkbox"/> | Pink | ○ |
| UHA6971 | <input type="checkbox"/> | Blue | ○ |
| UHA6972 | <input type="checkbox"/> | Mirror Finish | ○ |
| BRAKES* | | | |
| DHA7200 | <input type="checkbox"/> | Pull Brakes | €85 |
| DHA7210 | <input type="checkbox"/> | Standard Push to Lock Brakes, Lever Bent | ◆ |
| DHA7220 | <input type="checkbox"/> | Standard Push to Lock Brakes, Lever Straight | ○ |
| DHA7230 | <input type="checkbox"/> | Performance Brakes, Push to Lock. Brake barrel flips away from tyre | (NC 7300, 7810, 7820, 7830) €85 |
| DHA7240 | <input type="checkbox"/> | Active/Scissor Brakes | (NC 7300) €85 |
| DHA7250 | <input type="checkbox"/> | Light/Compact Active, Aluminium Brakes | (NC 7300) €248 |
| BRAKE LEVER EXTENSION | | | |
| DHA7300 | <input type="checkbox"/> | Brake Lever Extension | €48 |
| SAFETY OPTIONS | | | |
| DHA7500 | <input type="checkbox"/> | Anti-Tipper Right | €149 |
| DHA7510 | <input type="checkbox"/> | Anti-Tipper Left | €149 |
| DHA7520 | <input type="checkbox"/> | Anti-Tipper Both sides | €299 |
| DHA7800 | <input type="checkbox"/> | Standard Pelvic belt | €74 |
| BODYPOINT BELTS | | | |
| DHA7810 | <input type="checkbox"/> | Bodypoint Posture Belt Evoflex Small | €171 |
| DHA7820 | <input type="checkbox"/> | Bodypoint Posture Belt Evoflex Medium | €171 |
| DHA7830 | <input type="checkbox"/> | Bodypoint Posture Belt Evoflex Large | €171 |

| OTHER OPTIONS | | | |
|-----------------------|--|--|------|
| DHA8010 | | Front Frame Protectors; 2 Piece. Removable with hook and loop | €36 |
| DHA8100 | | Tipper-Aid, Mounted Right Side Max. 2 adaptations per side (Adaptations = Anti tipper, Tipper aid and Cane holder) | €86 |
| DHA8200 | | Cane Holder, Mounted Right Side Max. 2 adaptations per side (Adaptations = Anti tipper, Tipper aid and cane holder) | €86 |
| DHA8500 | | Air Pump | €17 |
| DHA8550 | | Tool Kit | ◆ |
| DHA8650 | | Kuschall branded Wheel bag (Fits all wheel sizes) | €58 |
| DHA8680 | | 24" Mountain Bike Wheels with bag. Fitted with Specific Black Tyres, Grey Aluminium Handrims; Quick Release Axles Included | €335 |
| DHA8750 | | Travel bag for Champion SK, | €294 |
| BODYPOINT ACCESSORIES | | | |
| DHA8800 | | Bodypoint Calf Strap | €80 |
| DHA8850 | | Bodypoint Mobility Bag | €80 |

| | |
|--------------------|--|
| TOTAL PRICE | |
|--------------------|--|

Non Conforming Options:

| | |
|---------|--|
| DHA1946 | 430 mm |
| DHA1948 | 440 mm |
| DHA1950 | 450 mm |
| DHA1952 | 460 mm |
| DHA1954 | 470 mm |
| DHA1956 | 480 mm |
| DHA1958 | 490 mm |
| DHA1960 | 500 mm |
| DHA2200 | Black Standard Light Upholstery, without Tension Adjustment |
| UHA0004 | Your personal flock design. Supply your design with your order, subject to copyright. Design size up to 280mm x 400mm |
| DHA2400 | Without Push Handles |
| DHA2430 | Push Handles Height Adjustable, Rear Set |
| DHA2710 | Flo-Tech Lite Cushion 78mm High - Comfort/Low Risk |
| DHA2720 | Matrx Libra Cushion, Black Smooth Cover 70mm High - Very High Risk |
| DHA2748 | Contour Visco NG Cushion 100mm High - High Risk |
| DHA2810 | Fluid supplement sac (Pair) |
| DHA2820 | Tapered wedge (Pair) |
| DHA2830 | Abductor |
| DHA2840 | Lateral side wedge |
| DHA2850 | Anterior front wedge |
| DHA2860 | Seat rigidizer |
| DHA3480 | Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included |
| DHA3490 | Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included |
| UHA4066 | Ceramic glaze finish on frame, for added durability |
| DHA4195 | Tag (Small Kuschall brand transfer on front of frame) |
| DHA4300 | Black |
| DHA4310 | Red |
| DHA2222 | No tag colour, select only with Option DHA4197 Without Tag |
| DHA5000 | 75mm x 25mm (3" x 1") Supersport |
| DHA5050 | 100mm x 30mm (4" x 1.2") Everyday Black hub and tyre |
| DHA5100 | 100mm x 30mm (4" x 1.2") Comfort Black hub and tyre |
| DHA5150 | 100mm x 25mm (4" x 1") Starec Black hub and tyre |
| DHA5200 | 100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frogs Legs) Use this option to access 4" anodised colour options |
| DHA5250 | 100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs) |
| UHA5502 | Pink |
| UHA5504 | Orange |
| UHA5506 | Blue |
| UHA5507 | Purple |
| UHA5508 | Red |
| UHA5509 | Gold |
| UHA5510 | Ceramic castor bearings for Aluminium soft roll castors |
| DHA6180 | 9th from front |
| DHA6410 | 24" |

| | |
|---------|---|
| DHA6420 | 25" |
| DHA6430 | 26" |
| DHA6510 | 24", 26" Kuschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub. |
| DHA6520 | 24", 25" High Performance Wheel Black hub and rim |
| DHA6550 | 24", 25" Spinergy LXK, Black hub, rim and 18 spokes in Black |
| UHA6551 | 24", 25" Spinergy Lite Extreme LX, Weight 687 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard) |
| UHA6552 | 24", 25" Spinergy Lite Extreme X-laced XLX, Weight 821 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard) |
| UHA6553 | 24", 25" Spinergy X-Laced Sport Lite Extreme XSLX, Weight 757 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard) |
| UHA6554 | 24", 25" Spinergy Carbon Fibre CLX, Weight 530 grams per wheel (Unique coloured hub and spokes available, supplied in black as standard) |
| SWO0907 | 24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. (Unique coloured spokes available) |
| UHA6557 | Blue |
| UHA6558 | Gold |
| UHA6559 | Green |
| UHA6560 | Orange |
| UHA6561 | Red |
| UHA6562 | Purple |
| UHA6563 | Blue |
| UHA6564 | Gold |
| UHA6565 | Silver |
| UHA6566 | Green |
| UHA6567 | Orange |
| UHA6568 | Red |
| UHA6569 | Purple |
| UHA6570 | Blue |
| UHA6571 | Green |
| UHA6572 | Orange |
| UHA6573 | Pink |
| UHA6574 | Red |
| UHA6575 | Purple |
| UHA6576 | White |
| UHA6577 | Yellow |
| UHA6579 | Ceramic Bearing for Spinergy rear wheels |
| DHA6600 | 24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection |
| DHA6630 | 24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance) |
| DHA6660 | 24" x 1 3/8" Solid Profile Black Tyre |
| DHA6700 | 24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance |
| DHA6710 | 24", 26" Steel, short tab with 1 hole, mounted at 3 cm distance |
| DHA6720 | 24" Titanium, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance |
| DHA6740 | 24", 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. For complete or limited hand function offering low-effort propulsion. Ergonomic shape with fully integrated, replaceable silicon strip for extra grip, profile of 32mm diameter. |
| DHA6750 | 24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. For complete or limited hand function offering low-effort propulsion. Ergonomic shape with fully integrated, replaceable silicon strip for extra grip, smaller profile of 24mm diameter for smaller hands. |

| | |
|---------|---|
| DHA6760 | 24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape, all metal black coated, low-friction surface for complete hand function for low-effort propulsion |
| DHA6770 | 24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics without tricep function. Features a unique angle shaped profile and anti-slip coating providing a larger contact surface designed to be propelled using the heel of the hand or thumb pad. |
| DHA6775 | 24" Carbolife Curve L Tetra handrim, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics with tricep function, but limited hand strength. Ergonomic shape with a flat, non-slip coated pushing surface for the thumb and a curved side for the fingers to wrap around. |
| DHA6900 | 24", 26" Transparent Spokeguards (When ordered with Carbolife handrims, handrim distance is set at 4cm) |
| DHA7000 | Fixation Strap (Hook and loop strap under the seat holds the frame together when folded for easy transportation, designed for clients with limited hand function) |
| DHA7010 | Fixation Clip (Plastic clip on either side of the seat holds the frame together when folded for easy transportation) |
| DHA7230 | Performance Brakes, Push to Lock. Brake barrel flips away from tyre |
| DHA7300 | Brake Lever Extension |
| DHA7810 | Bodypoint Posture Belt Evoflex Small |
| DHA7820 | Bodypoint Posture Belt Evoflex Medium |
| DHA7830 | Bodypoint Posture Belt Evoflex Large |



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):