



# Küschall® CHAMPION SL

MADE TO  
TRAVEL

—  
BUILT TO  
PERFORM



Kuschall.com



*FOLDABLE  
AND COMPACT*

**CHAMPION**

# SL

**Meet the all new Kuschall Champion SL**, an exciting and stylish new Kuschall design that retains all the travel-conscious characteristics of the original Champion series but welcomes a lighter and more agile performance. Engineered with the increased stiffness and rigidity of 'The KSL', the new Champion SL gives optimum driver control and high levels of responsiveness for an even greater drive.

*SUPER LIGHT  
AND RIGID*





The Champion SL promises to give the most discerning of active wheelchair users an impeccable drive.



We poured our passion and expertise into every detail of the Champion SL.

Our goal was to design a compact, foldable chair that gives an exceptional drive. Not only did we commit to the highest standards of functionality and durability in our design, but we reached for a chair that captures the spirit and aesthetic of Kuschall.



# *MADE TO TRAVEL* *BUILT TO PERFORM*

Designed for an active lifestyle, the Champion SL helps take you on a journey to explore new adventures, pursue passions and help unlimit your world.





## *TRANSPORTS WITH EASE*

We would never change what you love about the Kuschall Champion. A travel-conscious design, the Champion SL is a winner when it comes to convenience. Simply fold the frame for easy transportation and compact storage. The new Champion SL takes you places you want to go – with confidence.





# *EXCEPTIONAL DRIVE PERFORMANCE*

The Champion SL is a highly versatile active chair built especially to perform. Super light and robust, we have focused on a design precisely engineered for energy efficiency. With increased resistance and responsiveness, the Champion SL pays attention to efficient propulsion and optimal biomechanics.





# TAKE A CLOSER LOOK

## **SUPER LIGHT**

Now 5.9kg transport weight for easier handling.

## **RESPONSIVE**

Gives easy and quick changes in direction, speed and navigation of tight corners.

## **TRANSPORTS WITH EASE**

Unique folding mechanism makes transporting and storage simple.



## **ERGONOMIC GRIP**

For easier repositioning and transfers.

## **HYDROFORMED TECHNOLOGY**

Strengthens the frame for increased rigidity and improved energy efficiency when propelling.



**WATCH  
THE  
VIDEO**



# CUSTOMISE YOUR CHAMPION SL

## 1. FRAME COLOUR

Standard frame colours



Pure White    Anthracite Grey    Space Black



Carmine Red    Azure Blue

Extra frame colours\*



Flashy Yellow    Black Glossy    Silver Grey    Pearl Night Blue



Ruby Red    Dynamite Orange    Forest Green

Please refer to local prescription forms for colour options.

\*Extra lead time


## 2. DECAL STYLE



## 3. DECAL COLOUR



Red    Black



Light Grey    Gold    Red    Blue    Dark Grey    Black

## 4. UNIQUE SERVICE

Create your new Küschall wheelchair from the wheels up with UNIQUE. This dedicated service will create your personalized chair, designed to fit you and your lifestyle.

**UNIQUE**  
Designed by you. Built by us.



Pictured: Spineray wheel with UNIQUE features. Scan the QR code on this page for more information.







# TECHNICAL SPECIFICATION

## SEAT WIDTH

360 - 480 mm  
In increments of 20 mm



## SEAT DEPTH

340 - 480 mm  
In increments of 20 mm



## KNEE-TO-HEEL LENGTH

320 - 500 mm  
Stepless adjustable



## BACKREST HEIGHT

300 - 465 mm  
In increments of 15 mm



## BACKREST ANGLE

76° / 80.5° / 85° / 89.5° / 94°



## SEAT HEIGHT

Front: 450 - 540 mm  
Rear: 390 - 490 mm



## FRAME ANGLE

75° / 90°  
Straight or tapered



## CENTER OF GRAVITY

8 positions



## CAMBER

0° / 1° / 3°



## TOTAL LENGTH

75° approx. 930 mm  
90° approx. 830 mm



## TOTAL WIDTH

SW + 160 - 240 mm



## TOTAL WIDTH FOLDED

Approx. 280 - 340 mm



## TOTAL WEIGHT UNLOADED

Approx. 8.4 kg



## TRANSPORT WEIGHT

Approx. 5.9 kg



## MAX. USER WEIGHT

100 kg



All measurements and weights indicated are based on seat width 400 mm and seat depth 400 mm in the lightest and most compact configuration. Weight and measurements can change with the different configurations. For further information please check the user manual on our website.



**UNLIMIT YOUR WORLD**

**T: +44 1656 776222**  
uk@invacare.com  
www.invacare.co.uk  
www.invacare.ie

Invacare UK & Ireland  
Unit 4 Pencoed  
Technology Park  
Bridgend CF35 5AQ

