



Küschall®

COMPACT

VERSATILITY
& COMFORT
IN ONE CHAIR



Küschall.com



FOLDABLE
& LIGHT

COMPACT

**Küschall Compact is the absolute standard
for reliability and quality in one foldable chair.**

It can easily be adapted to suit varying needs to enhance overall comfort. Available in a fixed front frame with impressive rigidity for excellent driving performance or a swing-away frame, making transfers, storage and transportation even easier.



ENHANCED
COMFORT





OUTSTANDING FUNCTIONALITY

A foldable backrest option makes the chair small and compact for easy transportation. Available in 70° or 80° legrest angles, in addition to the adjustable precision locking mechanism with the swing-away legrests for increased rigidity and optimal drive.



INTRODUCING *COMPACT CARBON*

Small detail. Noticeable difference. Optional carbon legrests bring a lighter, more refined finish to the Compact. A subtle upgrade that enhances both the feel and the overall aesthetic, without changing the reliability you depend on.



Watch
the video

—
Elevate your Compact and be
confident with a touch of carbon.

C

-200g*

*Lighter in total weight than the
standard aluminium Compact.





Küschall®
COMPACT



TAKE A CLOSER LOOK

ANGLE ADJUSTABLE BACKREST

Available as an option, the angle adjustable and foldable backrest folds neatly over the seat rails for compact storage and easy transportation.

EASY FOLDING

The unique design of the cross brace makes the Compact simple to fold. A key feature if you have reduced strength or limited hand function.

KEEP ON TRACK

Intuitive and unique castor housing adjustment, to ensure your wheelchair always performs at its most energy efficient, making propelling easier for longer.





FIXED FRAME

Increases rigidity and gives a more active aesthetic. 70° or 80° legrest angle.

SWING-AWAY LEGRESTS

Making transport, storage or transfers easier. 70° or 80° legrest angle.

Küschall® | e

Optional carbon legrests for a lighter more refined finish.

FOOTPLATE


The footplate range fits all seat widths and features an anti-slip surface for increased comfort and support.



CUSTOMISE YOUR **COMPACT**

1. FRAME COLOUR

Standard frame colours



Pure White



Anthracite Grey



Space Black



Carmine Red



Azure Blue

Extra frame colours*



Flashy Yellow



Black Glossy



Silver Grey



Pearl Night Blue



Ruby Red



Dynamite Orange



Forest Green

*Extra delivery time will apply.

2. DECAL STYLE



3. DECAL COLOUR



Red



Black



Light Grey



Gold



Red



Blue



Dark Grey



Black



4. UNIQUE SERVICE

Create your new Küschall wheelchair from the wheels up with UNIQUE. This dedicated service will create your personalized chair, designed to fit you and your lifestyle.

UNIQUE
Designed by you. Built by us.





TECHNICAL SPECIFICATION

**SEAT
WIDTH**
280 – 500 mm
in increments of 20 mm



**SEAT
DEPTH**
280 – 500 mm
in increments of 20 mm



**KNEE-TO-HEEL
LENGTH**
200 - 510 mm
in increments of 10 mm



**BACKREST
HEIGHT**
300 – 510 mm
in increments of 15 mm



**BACKREST
ANGLE**
82° / 86° / 90° /
94° / 98° / 102°



**SEAT
HEIGHT**
Front: 380 – 530 mm
Rear: 370 – 500 mm
in increments of 10 mm



**FRAME
ANGLE**
70° / 80°
Straight or tapered



**TOTAL
WIDTH**
SW + 160 mm



**TOTAL
WIDTH FOLDED**
Approx. 290 mm



**TOTAL
LENGTH**
70° approx. 870 – 1220 mm
80° approx. 850 – 1200 mm



CAMBER
0° / 3°



**TOTAL WEIGHT
UNLOADED**
Approx. 10.3 kg (swing away)
Approx. 9.6 kg (fixed frame)



**TRANSPORT
WEIGHT**
Approx. 7.8 kg (swing away)
Approx. 6.7 kg (fixed frame)



**MAX. USER
WEIGHT**
130 kg



CRASH TEST
Approved – Fulfills the
requirements of
ISO 7176-19



All measurements and weights indicated are based on seat width 400 mm and seat depth 400 mm in the lightest and most compact configuration.
Weight and measurements can change with the different configurations. For further information please check the user manual on our website.

/ UNLIMIT YOUR WORLD

T: +44 1656 776222
uk@invacare.com
www.invacare.co.uk
www.invacare.ie

Invacare UK & Ireland
Unit 4 Pencoed
Technology Park
Bridgend CF35 5AQ

