



**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

FEBRUARY 2024  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Küschall Compact One Solution

### Standard Features

- Max user weight 130kg/20.5 stone
- Crash tested ISO 7176-part 19
- Hydroformed aluminium, Active, Folding, with Swingaway or Fixed frame
- Multiple adjustments, front, rear seat heights, centre of gravity
- Hydroformed swingaway legrests, choice of 70° or 80° angles
- Hydroformed fixed frame, choice of 70° or 80° angles
- Height adjustable flip up plastic footrests
- Choice of seat widths 360mm - 500mm (14" - 20")
- Choice of 10 seat depths 320mm - 500mm (12" - 20")
- Fixed height backrest for Matrx backrests set at 90°
- Standard push handles
- Fixed composite clothesguards
- 5" comfort solid front castors with black tyres
- Kuschall Starec Q/R rear wheels with pneumatic tyres
- Grey aluminium handrims
- Pull to lock brakes
- Transportation kit
- Choice of frame colours

**Standard Base price: £2,025**



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
DCA0012	KÜSCHALL COMPACT 2.0 One Solution	£2,025
DCA0000	MYKÜSCHALL Special Request - Please include quotation reference in special description at the end of the form where applicable	POA

FRAME*				
DCA1100		Swing away Front frame	(NC 1900, 1902, 3040, 3140)	○
DCA1110		Fixed front frame 70°	(NC 1700, 1705, 1710, 1715, 1720, 1725, 1730, 1805, 3000, 3010, 3020, 3030, 3100, 3110, 3120, 3130)	○
DCA1120		Fixed front frame 80° (Minimum seat angle, front and rear seat heights is 40 mm with 5" front castors)	(NC 1700, 1705, 1710, 1715, 1720, 1725, 1730, 1805, 1962, 3000, 3010, 3020, 3030, 3100, 3110, 3120, 3130, 5550, 5650)	○
DCA1130		Fixed front frame tapered 70° Additional 10 day leadtime	(NC 1700, 1705, 1710, 1715, 1720, 1725, 1730, 1805, 3000, 3010, 3020, 3030, 3100, 3110, 3120, 3130)	○
DCA1140		Fixed front frame tapered 80° (minimum seat angle, front and rear seat heights is 40 mm with 5" front castors) Additional 10 day leadtime	(NC 1700, 1705, 1710, 1715, 1720, 1725, 1730, 1805, 1962, 3000, 3010, 3020, 3030, 3100, 3110, 3120, 3130, 5550, 5650)	○
SEAT WIDTH*				
DCA1320		360 mm Single cross brace	(NC 2364, 2366, 2368, 2370, 2372, 2710, 2748, 3020, 3030, 3120, 3130)	○
DCA1325		380 mm Single cross brace	(NC 2366, 2368, 2370, 2372, 3020, 3030, 3120, 3130)	○
DCA1330		400 mm Single cross brace	(NC 2360, 2366, 2368, 2370, 2372)	◆
DCA1335		420 mm Single cross brace	(NC 2360, 2362, 2368, 2370, 2372)	○
DCA1340		440 mm Double cross brace	(NC 2360, 2362, 2364, 2370, 2372)	○
DCA1345		460 mm Double cross brace	(NC 2360, 2362, 2364, 2372)	○
DCA1350		480 mm Double cross brace	(NC 2300, 2305, 2310, 2315, 2345, 2360, 2362, 2364, 2366)	○
DCA1355		500 mm Double cross brace	(NC 2300, 2305, 2310, 2315, 2345, 2360, 2362, 2364, 2366, 2368)	○
SEAT DEPTH*				
DCA1400		320 mm	(NC 2155, 2160, 2710, 2720, 2748, 2750, 2860, 5650)	○
DCA1405		340 mm	(NC 2155, 2160, 2710, 2720, 2748, 2860, 5650)	○
DCA1415		360 mm	(NC 2710, 2720, 2748, 2860, 5650)	○
DCA1425		380 mm	(NC 2720, 5650)	○
DCA1435		400 mm		◆
DCA1445		420 mm		○
DCA1460		440 mm		○
DCA1470		460 mm		○
DCA1480		480 mm	(NC 2748)	○
DCA1490		500 mm	(NC 2710, 2748)	○
BACK HEIGHT Transportation kit is only possible with backrest heights of 405mm and above*				
DCA1615		300 mm	(NC 2140, 7600)	○
DCA1620		315 mm	(NC 2140, 7600)	○
DCA1625		330 mm	(NC 2140, 7600)	○
DCA1630		345 mm	(NC 7600)	○
DCA1635		360 mm	(NC 7600)	○
DCA1640		375 mm	(NC 7600)	○
DCA1645		390 mm	(NC 7600)	○
DCA1650		405 mm		○
DCA1655		420 mm		◆
DCA1660		435 mm		○
DCA1665		450 mm		○
DCA1670		465 mm		○
DCA1675		480 mm		○
DCA1680		495 mm		○
DCA1685		510 mm		○

**FRONT SEAT HEIGHT. Seat heights 380mm TO 430mm offers a unique front frame design to achieve a lower front seat height, only available for the Compact Swing option\***

DCA1700		370 mm	(NC 1810, 1815, 1820, 1825, 1830, 1835, 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1932, 1934, 1936, 1938, 1940, 1942, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962, 3030, 3130, 3430, 6000, 6010, 6410, 6420, 7600)	○
DCA1705		380 mm	(NC 1815, 1820, 1825, 1830, 1835, 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1934, 1936, 1938, 1940, 1942, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962, 3030, 3130, 3430, 6000, 6010, 6410, 6420, 7600)	○
DCA1710		390 mm	(NC 1820, 1825, 1830, 1835, 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1936, 1938, 1940, 1942, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962, 3030, 3130, 3430, 6000, 6010, 6410, 6420, 7600)	○
DCA1715		400 mm	(NC 1825, 1830, 1835, 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1938, 1940, 1942, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962, 3030, 3130, 3430, 6000, 6010, 6410, 6420, 7600)	○
DCA1720		410 mm	(NC 1830, 1835, 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1940, 1942, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962, 3030, 3130)	○
DCA1725		420 mm	(NC 1835, 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1942, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962)	○
DCA1730		430 mm	(NC 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962)	○
DCA1735		440 mm	(NC 1845, 1850, 1855, 1860, 1865, 1870, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962)	○
DCA1740		450 mm	(NC 1850, 1855, 1860, 1865, 1870, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962)	○
DCA1745		460 mm	(NC 1855, 1860, 1865, 1870, 1950, 1952, 1954, 1956, 1958, 1960, 1962)	○
DCA1750		470 mm	(NC 1860, 1865, 1870, 1952, 1954, 1956, 1958, 1960, 1962)	○
DCA1755		480 mm	(NC 1805, 1865, 1870, 1954, 1956, 1958, 1960, 1962)	○
DCA1760		490 mm	(NC 1805, 1810, 1956, 1958, 1960, 1962)	○
DCA1765		500 mm	(NC 1805, 1810, 1815, 1958, 1960, 1962)	◆
DCA1770		510 mm	(NC 1805, 1810, 1815, 1820, 1960, 1962)	○
DCA1775		520 mm	(NC 1805, 1810, 1815, 1820, 1825, 1962)	○
DCA1780		530 mm	(NC 1805, 1810, 1815, 1820, 1825, 1830)	○

**REAR SEAT HEIGHT\***

DCA1805		370 mm	(NC 6000, 6010, 6410, 6420, 6430)	○
DCA1810		380 mm	(NC 6000, 6010, 6410, 6420, 6430)	○
DCA1815		390 mm	(NC 6420, 6430)	○
DCA1820		400 mm	(NC 6420, 6430)	○
DCA1825		410 mm	(NC 6420, 6430)	○
DCA1830		420 mm		○
DCA1835		430 mm		○
DCA1840		440 mm		○
DCA1845		450 mm		◆
DCA1850		460 mm		○
DCA1855		470 mm		○
DCA1860		480 mm	(NC 6400)	○
DCA1865		490 mm	(NC 6400)	○
DCA1870		500 mm	(NC 6000, 6010, 6400, 6410)	○

<b>LOWER LEG LENGTH</b> Seat canvas to footplate measurement, cushion height must be deducted.*				
DCA1900		200 mm	(NC 3030, 3130)	○
DCA1902		210 mm	(NC 3030, 3130)	○
DCA1904		220 mm	(NC 3030, 3130)	○
DCA1906		230 mm	(NC 3030, 3130)	○
DCA1908		240 mm	(NC 3030, 3130)	○
DCA1910		250 mm	(NC 3030, 3130)	○
DCA1912		260 mm	(NC 3030, 3130)	○
DCA1914		270 mm	(NC 3030, 3130)	○
DCA1916		280 mm	(NC 3030, 3130)	○
DCA1918		290 mm	(NC 3030, 3130)	○
DCA1920		300 mm	(NC 3030, 3130)	○
DCA1922		310 mm	(NC 3030, 3130)	○
DCA1924		320 mm	(NC 3030, 3130)	○
DCA1926		330 mm	(NC 3030, 3130)	○
DCA1928		340 mm	(NC 3030, 3130)	○
DCA1930		350 mm	(NC 3030, 3130)	○
DCA1932		360 mm	(NC 3030, 3130)	○
DCA1934		370 mm	(NC 3030, 3130)	○
DCA1936		380 mm		○
DCA1938		390 mm		○
DCA1940		400 mm		○
DCA1942		410 mm		○
DCA1944		420 mm		○
DCA1946		430 mm		◆
DCA1948		440 mm		○
DCA1950		450 mm		○
DCA1952		460 mm		○
DCA1954		470 mm		○
DCA1956		480 mm		○
DCA1958		490 mm		○
DCA1960		500 mm		○
DCA1962		510 mm	(NC 3000, 3100, 3010, 3110)	○
<b>BACKREST TYPE*</b>				
DCA2110		Fixed straight backrest tubes set at 90°	(NC 2155, 2160, 2165, 2170, 2175, 2180)	◆
DCA2130		Angle adjustable back (angle must be selected from the options below)	(NC 0000)	£202
DCA2140		Fold down, angle adjustable back, folds onto seat (angle must be selected from the options below)	(NC 0000, 2180, 3520, 3620)	£206
<b>BACKREST ANGLE*</b>				
DCA2155		Angle adjustable backrest adjusted 82° from the seat		○
DCA2160		Angle adjustable backrest adjusted 86° from the seat		○
DCA2165		Angle adjustable backrest adjusted 90° from the seat		○
DCA2170		Angle adjustable backrest adjusted 94° from the seat		○
DCA2175		Angle adjustable backrest adjusted 98° from the seat		○
DCA2180		Angle adjustable backrest adjusted 102° from the seat		○
DCA0000		Without angle adjustment. select with option 2110 fixed backrest		○

<b>BACKREST UPHOLSTERY*</b>			
DCA2240		Without Backrest Upholstery	◆
<b>MATRX BACKREST TYPE*</b>			
DCA2300		Matrx MX2, height 150 mm with quick release mini set hardware (NC 2370, 2372, 7600)	£751
DCA2305		Matrx MX2, height 230 mm with quick release mini set hardware (NC 2370, 2372, 7600)	£751
DCA2310		Matrx MX2, height 305 mm with quick release mini set hardware (NC 2370, 2372, 7600)	£751
DCA2315		Matrx MX2, height 405 mm with quick release mini set hardware (NC 2370, 2372, 7600)	£751
DCA2316		Matrx Elite Deep E2, height 160 mm with quick release mini set hardware (NC 7600)	£546
DCA2317		Matrx Elite Deep E2, height 210 mm with quick release mini set hardware (NC 7600)	£546
DCA2318		Matrx Elite Deep E2, height 260 mm with quick release mini set hardware (NC 7600)	£546
DCA2320		Matrx Elite Deep E2, height 310 mm with quick release mini set hardware (NC 7600)	£546
DCA2322		Matrx Elite Deep E2, height 360 mm with quick release mini set hardware (NC 7600)	£546
DCA2325		Matrx Elite Deep E2, height 410 mm with quick release mini set hardware	£546
DCA2330		Matrx Elite Deep E2, height 460 mm with easy-set mounting hardware	£546
DCA2331		Matrx Elite Deep E2, height 510 mm with easy-set mounting hardware	£546
DCA2332		Matrx Elite E2, height 160 mm with quick release mini set hardware (NC 7600)	£494
DCA2334		Matrx Elite E2, height 210 mm with quick release mini set hardware (NC 7600)	£494
DCA2335		Matrx Elite E2, height 260 mm with quick release mini set hardware (NC 7600)	£494
DCA2340		Matrx Elite E2, height 310 mm with quick release mini set hardware (NC 7600)	£494
DCA2345		Matrx Elite E2, height 360 mm with quick release mini set hardware (NC 7600)	£494
DCA2350		Matrx Elite E2, height 410 mm with quick release mini set hardware	£494
DCA2355		Matrx Elite E2, height 460 mm with easy-set mounting hardware	£494
DCA2357		Matrx Elite E2, height 510 mm with easy-set mounting hardware	£494
<b>MATRX BACKREST WIDTH CHOICE*</b>			
DCA2360		360mm Backrest width, to fit chest width 265mm	○
DCA2362		380mm Backrest width, to fit chest width 290mm	○
DCA2364		410mm Backrest width, to fit chest width 315mm	○
DCA2366		430mm Backrest width, to fit chest width 340mm	○
DCA2368		460mm Backrest width, to fit chest width 370mm	○
DCA2370		480mm Backrest width, to fit chest width 390mm	○
DCA2372		510mm Backrest width, to fit chest width 420mm	○
<b>PUSH HANDLES*</b>			
DCA2400		Without Push Handles	○
DCA2410		Push Handles Fixed Short/Mini Grip (NC 6580)	○
DCA2420		Push Handles Fixed Long Grip	◆
DCA2440		Push Handles Height Adjustable, Integrated in backrest	£219
DCA2450		Push Handles Foldable (NC 6580)	£198

<b>SEAT UPHOLSTERY*</b>			
DCA2610		Standard seat upholstery with 1 central band (loop)	○
DCA2620		Seat upholstery with 1 central (loop) and 2 side bands (hook)	◆
<b>SEAT CUSHION</b>			
DCA2700		Standard Seat Cushion, 60mm high Double Foam with Visco Top Layer - Comfort/Low Risk (NC 2810, 2820, 2830, 2840, 2850)	£66
DCA2710		Flo-Tech Lite Cushion 78mm high - Comfort/Low Risk (NC 2810, 2820, 2830, 2840, 2850)	£121
DCA2720		Matrix Libra Cushion, Black Smooth Cover 70mm High - Very High Risk	£410
DCA2748		Contour Visco NG Cushion 100mm High - High Risk (NC 2810, 2820, 2830, 2840, 2850)	£169
DCA2750		Vicair Active O2 lightweight cushion 60mm High - High risk (NC 2810, 2820, 2830, 2840, 2850, 2860)	£521
<b>SEAT BAG</b>			
DCA2800		Seat Bag, Click in and removable	£74
<b>MATRIX LIBRA CUSHION ACCESSORIES</b>			
DCA2810		Fluid supplement sac (Pair)	£44
DCA2820		Tapered Wedge Pair	£30
DCA2830		Abductor	£23
DCA2840		Lateral Side Wedge	£20
DCA2850		Anterior Front Wedge	£27
<b>SEAT RIGIDISER</b>			
DCA2860		Seat Rigidiser	£39
<b>LEGREST TYPE RH*</b>			
DCA3000		Legrest 70° RH (NC 3110, 9999)	◆
DCA3010		Legrest 80° RH (NC 3100, 5650, 9999)	○
DCA3020		Amputee legrest RH (NC 3280, 3300, 3305)	£138
DCA3030		Angle adjustable legrest - RH, including angle, width and depth adjustable plastic footplate (NC 3280, 3300, 3305)	£119
DCA3040		Without Legrests RH, Select with fixed front frames	○
<b>LEGREST TYPE LH*</b>			
DCA3100		Legrest 70° LH (NC 9999)	◆
DCA3110		Legrest 80° LH (NC 5650, 9999)	○
DCA3120		Amputee legrest LH (NC 3280, 3300, 3305)	£138
DCA3130		Angle adjustable legrest - LH, including angle, width and depth adjustable plastic footplate (NC 3280, 3300, 3305)	£119
DCA3140		Without Legrests LH, Select with fixed front frames	○
<b>FOOTREST*</b>			
DCA3260		2 pc. Plastic flip up	◆
DCA3270		2 pc. Carbon angle and depth adjustable (NC 3310)	£333
DCA3280		1 pc. Carbon flip-up to the left, angle and depth adjustable (NC 3310)	£167
DCA9999		No footrest (Select when ordering Amputee brackets both sides)	○
<b>CALF STRAP</b>			
DCA3300		Calf strap mounted on the frame	○
DCA3305		Calf strap mounted on bracket	○
<b>HEEL LOOPS</b>			
DCA3310		Heel loops	○

<b>CLOTHESGUARDS/MUDGUARDS We are currently experiencing an extended leadtime of 12 weeks on carbon mudguards*</b>			
DCA3405		Composite clothes guards, removable	(NC 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 7600, 8450) ◆
DCA3410		Carbon clothes guards, fixed	(NC 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 7600, 8450) ○
DCA3425		Composite mudguards, fixed	(NC 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 6000, 6010, 6150, 6580, 7600, 8450) £65
DCA3430		Carbon Mudguards, Fixed (Currently we have an extended leadtime of 12 weeks on carbon mudguards)	(NC 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 6000, 6010, 6150, 6580, 7600, 8450) £247
DCA3440		Carbon Mudguards, Removable (Currently we have an extended leadtime of 12 weeks on carbon mudguards)	(NC 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 6150, 8450) £307
DCA3450		Composite mudguards, removable	(NC 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 6150, 8450) £65
DCA6666		Without clothesguards or mudguards, select only with options 3500, 3510, 3530, 3540, 3600, 3610, 3630 or 3640 Armrests supplied with clothesguards	(NC 3520, 3620) ○
<b>ARMREST RH*</b>			
DCA3500		Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, Carbon Clothesguard included RH	(NC 6150) ○
DCA3510		Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, Carbon Clothesguard included RH	(NC 6150) ○
DCA3520		Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included	(NC 8450) £119
DCA3530		Küschall armrest, height adjustable, removable, swing-away, desk armpad with black neopreen sideguard upholstery, RH	(NC 6130, 6140, 6150) £178
DCA3540		Hemi armrest, height adjustable, removable, swing-away, with hemi armpad RH	(NC 6130, 6140, 6150, 8450) £103
DCA7777		Without armrest, select with Clothesguard/Mudguard options	(NC 3500, 3510, 3530, 3540) ○
<b>ARMREST LH*</b>			
DCA3600		Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, Carbon Clothesguard included LH	○
DCA3610		Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, Carbon Clothesguard included LH	○
DCA3620		Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included LH	(NC 8450) £119
DCA3630		Küschall armrest, height adjustable, removable, swing-away, desk armpad with black neopreen sideguard upholstery, LH	(NC 6130, 6140, 6150) £178
DCA3640		Küschall hemi armrest, height adjustable, removable, swing-away, with hemi armpad LH	(NC 6130, 6140, 6150, 8450) £103
DCA8888		Without armrest, select with Clothesguard/Mudguard options	(NC 3600, 3610, 3630, 3640) ○
<b>FRAME COLOUR: We are removing the Atomic Violet and Royal Blue frame colours, but will be replacing them shortly with an alternative selection of similar colours*</b>			
DCA4000		Space Black (RAL 9005)	○
DCA4010		Anthracite Grey (RAL 7016)	○
DCA4020		Pure White (RAL 9016)	(NC 4220) ○
DCA4030		Carmine Red (RAL 3002)	○
DCA4040		Azure Blue (RAL 5017)	○
DCA4060		Flashy Yellow (RAL 1516) (Additional 15 day lead time)	£31
DCA4100		Black Glossy (RAL 9005) (Additional 15 day lead time)	£31
DCA4130		Silver Grey (RAL 9006) (Additional 15 day lead time)	£31
DCA4150		Ruby Red (RAL 1503) (Additional 15 day lead time)	£31
DCA4170		Dynamite Orange (Additional 15 day lead time)	£31
DCA4180		Forest Green (RAL 6029) (Additional 15 day lead time)	£31
DCA4050		Personal frame colour. Provide standard RAL colour code of your choice	£311

FRAME DECAL TYPE*			
DCA4190		Bold Decal (Full coloured transfer on frame)	(NC 4300, 4310, 1111) £31
DCA4195		Tag (Small Kuschall brand transfer on front of frame)	(NC 4200, 4210, 4220, 4230, 4240, 4250, 2222) ◆
DCA4197		Without Tag (Kuschall branding only)	(NC 4200, 4210, 4220, 4230, 4240, 4250, 4300, 4310) ○
BOLD DECAL COLOUR*			
DCA4200		Black	○
DCA4210		Dark Grey	○
DCA4220		Light Grey	○
DCA4230		Red	○
DCA4240		Blue	○
DCA4250		Gold	○
DCA1111		No bold colour, select only with Option DCA4195 Tag or DCA4197 Without Tag	(NC 4190) ○
TAG COLOUR*			
DCA4300		Black	○
DCA4310		Red	○
DCA2222		No tag colour, select only with Option DCA4190 Bold Decal or DCA4197 Without Tag	(NC 4195) ○
FRONT CASTORS*			
DCA5050		100mm x 30mm (4" x 1.2") Everyday Black hub and tyre	○
DCA5100		100mm x 30mm (4" x 1.2") Comfort Black hub and tyre	○
DCA5150		100mm x 25mm (4" x 1") Starec Black hub and tyre	£83
DCA5200		100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs)	£93
DCA5250		100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs)	£93
DCA5300		125mm x 30mm (5" x 1.2") Everyday Black hub and tyre	○
DCA5350		125mm x 30mm (5" x 1.2") Comfort Black hub and tyre	◆
DCA5400		125mm x 25mm (5" x 1") Starec Black hub and tyre	£83
DCA5450		125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs)	£93
DCA5500		125mm x 36mm (5" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs)	£93
DCA5550		150mm x 30mm (6" x 1.2") Comfort Black hub and tyre	○
DCA5650		175mm x 30mm (7" x 1.2") Comfort Black hub and tyre	○
ONE ARM DRIVE			
DCA6000		One Arm Drive system Right side (with MAXXIS Lite guard pneumatic tyres, alu handrims)	(NC 6040, 6130, 6140, 6150, 6300, 6400, 6420, 6430, 6510, 6520, 6550, 6580, 6625, 6630, 6660, 6840, 6900, 6700, 6710, 6720, 6740, 6750, 6760, 6770, 6780, 6790, 8680) £482
DCA6010		One Arm Drive system Left side (with MAXXIS Lite guard pneumatic tyres, alu handrims)	(NC 6040, 6130, 6140, 6150, 6300, 6400, 6420, 6430, 6510, 6520, 6550, 6580, 6625, 6630, 6660, 6840, 6900, 6700, 6710, 6720, 6740, 6750, 6760, 6770, 6780, 6790, 8680) £482
REAR WHEEL FIXATION*			
DCA6020		Standard fixation, 4 rear wheel positions	(NC 6140, 6150) ◆
DCA6040		Rear wheel extension, 6 rear wheel positions	(NC 6100, 6580) £63
REAR WHEEL POSITION. Position 1 is most active and position 6 is most stable*			
DCA6100		1st Position	(NC 6580, 6840) ○
DCA6110		2nd Position	○
DCA6120		3rd Position	◆
DCA6130		4th Position	(NC 6580) ○
DCA6140		5th Position	○
DCA6150		6th Position	○



QUAD QUICK RELEASE AXLE			
DCA6300		Quad quick release/tetra clip	(NC 6550, 6580) £85
REAR WHEEL SIZE*			
DCA6400		22"	(NC 6520, 6550, 6600, 6630, 6660, 6760, 6710, 6720, 6740, 6750, 6760, 6770) ○
DCA6410		24"	◆
DCA6420		25"	(NC 6510, 6580, 6600, 6630, 6660, 6700, 6710, 6720, 6760, 6770, 6900) ○
DCA6430		26"	(NC 6520, 6550, 6600, 6630, 6660, 6720, 6740, 6750, 6760, 6770, 6840) ○
REAR WHEEL TYPE*			
DCA6510		22", 24", 26" Kuschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub. 22" Silver Cross-Spokes, Colourless hub and Black rim	◆
DCA6520		24", 25" High Performance Wheel Black hub and rim	(NC 6660) £182
DCA6550		24", 25" Spinerger LXX, Black hub, rim and 18 spokes in Black	(NC 6660) £611
SWO0907		24", 25" Spinerger SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKuschall	(NC 6000, 6010, 6660) £486
SWO0909		24", 25" Spinerger LX, Chrome hub and 12 spokes in Black. Coloured spokes are available via MyKuschall	(NC 6000, 6010, 6660) £735
DCA6580		22", 24", 26" Wheels with Attendant Drum Brake	(NC 8680) £366
DCA3333		select with options 6000 and 6010 One Arm Drive wheels, specific wheel supplied	○
REAR WHEEL TYRE*			
DCA6600		24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection	◆
DCA6625		24", 25", 26" x 1" Maxxis Alpha Protect. Lightweight, Pneumatic and highly puncture resistant.	£100
DCA6630		24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance)	£97
DCA6660		24" x 1 3/8" Solid Profile Black Tyre	(NC 7230) £66
HANDRIM*			
DCA6700		22", 24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance	£223
DCA6710		24", 26" Steel, short tab with 1 hole, mounted at 3cm distance	£62
DCA6720		24" Titanium, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance	£192
DCA6740		24", 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip	£320
DCA6750		24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands	£320
DCA6760		24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion	£284
DCA6770		24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim	£356
DCA6780		22", 24", 25", 26" Light Aluminium Grey, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance	◆
DCA6790		22", 24", 25", 26" Light Aluminium Black, short tab with 1 hole for 3cm mounting	£42
DCA5555		select with options 6000 and 6010 one arm drive wheels, supplied with specific handrims	○

<b>CAMBER ANGLE*</b>			
DCA6820		Camber 0°	◆
DCA6840		Camber 3° (NC 7230)	○
<b>SPOKEGUARD</b>			
DCA6900		22", 24", 26" Transparent Spokeguards	£126
<b>BRAKES*</b>			
DCA7200		Pull Brakes	◆
DCA7210		Standard Push to Lock Brakes, Lever Bent	○
DCA7220		Standard Push to Lock Brakes, Lever Straight	○
DCA7230		Performance Brakes, Push to Lock. Brake barrel flips away from tyre (NC 7300, 7810, 7820, 7830)	£78
DCA7240		Active/Scissor Brakes (NC 7300)	£78
DCA7250		Light/Compact Active, Aluminium Brakes (NC 7300)	£136
<b>BRAKE LEVER EXTENSION</b>			
DCA7300		Brake Lever Extension	£49
<b>SAFETY OPTIONS 2 adaptions can be mounted on the same side, this includes Anti tippers, Tipper aid and Cane holder</b>			
DCA7500		Anti-tipper Mounted on Right side. Height adjustable and swingaway	£122
DCA7510		Anti-tipper Mounted on Left side. Height adjustable and swingaway	£122
DCA7520		Anti-tipper Mounted on both sides. Height adjustable and swingaway	£244
DCA7600		Transportation Kit Includes Tiedown Labels. Seat belt 7800 must also be selected	○
DCA7800		Standard Pelvic belt. Must be selected if Transportation kit is selected	£74
<b>BODYPOINT BELTS</b>			
DCA7810		Bodypoint Evoflex posture belt, 2 point Small, flat-mount with push button	£81
DCA7820		Bodypoint Evoflex posture belt, 2 point Medium, flat-mount with push button	£81
DCA7830		Bodypoint Evoflex posture belt, 2 point Large, flat-mount with push button	£81
<b>OTHER OPTIONS</b>			
DCA8000		Rear frame protection, plastic guard attached to underside of frame	£19
DCA8010		Front frame protection, 2 piece removable with hook and loop	£37
DCA8100		Tipper-aid, Right (NC 8200)	£46
DCA8110		Tipper-aid, Left (NC 8210)	£46
DCA8120		Tipper-aid, Both sides (NC 8200, 8210)	£92
DCA8200		Cane holder, Right	£65
DCA8210		Cane holder, Left	£65
DCA8300		TRANSIT WHEELS 3"	£230
DCA8450		TABLE TRAY	£336
DCA8500		Air pump	£13
DCA8550		Tool kit	◆
DCA8680		24" Mountain bike wheels with bag. fitted with specific black tyres and grey aluminium handrims; quick release axles included	£320
<b>BODYPOINT ACCESSORIES</b>			
DCA8800		Bodypoint calf strap	£57
DCA8850		Bodypoint mobility bag	£55

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

DCA0000	MYKÜSCHALL Special Request - Please include quotation reference in special description at the end of the form where applicable
DCA1700	370 mm
DCA1705	380 mm
DCA1710	390 mm
DCA1715	400 mm
DCA1720	410 mm
DCA1725	420 mm
DCA1730	430 mm
DCA1805	370 mm
DCA1810	380 mm
DCA1815	390 mm
DCA1820	400 mm
DCA1825	410 mm
DCA1830	420 mm
DCA1835	430 mm
DCA1840	440 mm
DCA1845	450 mm
DCA1850	460 mm
DCA1855	470 mm
DCA1860	480 mm
DCA1865	490 mm
DCA1870	500 mm
DCA1900	200 mm
DCA1902	210 mm
DCA1932	360 mm
DCA1934	370 mm
DCA1936	380 mm
DCA1938	390 mm
DCA1940	400 mm
DCA1942	410 mm
DCA1944	420 mm
DCA1946	430 mm
DCA1948	440 mm
DCA1950	450 mm
DCA1952	460 mm
DCA1954	470 mm
DCA1956	480 mm
DCA1958	490 mm
DCA1960	500 mm
DCA1962	510 mm
DCA2140	Fold down, angle adjustable back, folds onto seat (angle must be selected from the options below)
DCA2155	Angle adjustable backrest adjusted 82° from the seat
DCA2160	Angle adjustable backrest adjusted 86° from the seat

DCA2165	Angle adjustable backrest adjusted 90° from the seat
DCA2170	Angle adjustable backrest adjusted 94° from the seat
DCA2175	Angle adjustable backrest adjusted 98° from the seat
DCA2180	Angle adjustable backrest adjusted 102° from the seat
DCA0000	Without angle adjustment. select with option 2110 fixed backrest
DCA2300	Matrx MX2, height 150 mm with quick release mini set hardware
DCA2305	Matrx MX2, height 230 mm with quick release mini set hardware
DCA2310	Matrx MX2, height 305 mm with quick release mini set hardware
DCA2315	Matrx MX2, height 405 mm with quick release mini set hardware
DCA2345	Matrx Elite E2, height 360 mm with quick release mini set hardware
DCA2360	360mm Backrest width, to fit chest width 265mm
DCA2362	380mm Backrest width, to fit chest width 290mm
DCA2364	410mm Backrest width, to fit chest width 315mm
DCA2366	430mm Backrest width, to fit chest width 340mm
DCA2368	460mm Backrest width, to fit chest width 370mm
DCA2370	480mm Backrest width, to fit chest width 390mm
DCA2372	510mm Backrest width, to fit chest width 420mm
DCA2710	Flo-Tech Lite Cushion 78mm high - Comfort/Low Risk
DCA2720	Matrx Libra Cushion, Black Smooth Cover 70mm High - Very High Risk
DCA2748	Contour Visco NG Cushion 100mm High - High Risk
DCA2750	Vicair Active O2 lightweight cushion 60mm High - High risk
DCA2810	Fluid supplement sac (Pair)
DCA2820	Tapered Wedge Pair
DCA2830	Abductor
DCA2840	Lateral Side Wedge
DCA2850	Anterior Front Wedge
DCA2860	Seat Rigidiser
DCA3000	Legrest 70° RH
DCA3010	Legrest 80° RH
DCA3020	Amputee legrest RH
DCA3030	Angle adjustable legrest - RH, including angle, width and depth adjustable plastic footplate
DCA3040	Without Legrests RH, Select with fixed front frames
DCA3100	Legrest 70° LH
DCA3110	Legrest 80° LH
DCA3120	Amputee legrest LH
DCA3130	Angle adjustable legrest - LH, including angle, width and depth adjustable plastic footplate
DCA3140	Without Legrests LH, Select with fixed front frames
DCA3280	1 pc. Carbon flip-up to the left, angle and depth adjustable
DCA9999	No footrest (Select when ordering Amputee brackets both sides)
DCA3300	Calf strap mounted on the frame
DCA3305	Calf strap mounted on bracket
DCA3310	Heel loops
DCA3430	Carbon Mudguards, Fixed (Currently we have an extended leadtime of 12 weeks on carbon mudguards)
DCA3500	Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, Carbon Clothesguard included RH

DCA3510	Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, Carbon Clothesguard included RH
DCA3520	Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included
DCA3530	Küschall armrest, height adjustable, removable, swing-away, desk armpad with black neopreen sideguard upholstery, RH
DCA3540	Hemi armrest, height adjustable, removable, swing-away, with hemi armpad RH
DCA3600	Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, Carbon Clothesguard included LH
DCA3610	Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, Carbon Clothesguard included LH
DCA3620	Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included LH
DCA3630	Küschall armrest, height adjustable, removable, swing-away, desk armpad with black neopreen sideguard upholstery, LH
DCA3640	Küschall hemi armrest, height adjustable, removable, swing-away, with hemi armpad LH
DCA4190	Bold Decal (Full coloured transfer on frame)
DCA4195	Tag (Small Küschall brand transfer on front of frame)
DCA4200	Black
DCA4210	Dark Grey
DCA4220	Light Grey
DCA4230	Red
DCA4240	Blue
DCA4250	Gold
DCA1111	No bold colour, select only with Option DCA4195 Tag or DCA4197 Without Tag
DCA4300	Black
DCA4310	Red
DCA2222	No tag colour, select only with Option DCA4190 Bold Decal or DCA4197 Without Tag
DCA5550	150mm x 30mm (6" x 1.2") Comfort Black hub and tyre
DCA5650	175mm x 30mm (7" x 1.2") Comfort Black hub and tyre
DCA6000	One Arm Drive system Right side (with MAXXIS Lite guard pneumatic tyres, alu handrims)
DCA6010	One Arm Drive system Left side (with MAXXIS Lite guard pneumatic tyres, alu handrims)
DCA6040	Rear wheel extension, 6 rear wheel positions
DCA6100	1st Position
DCA6130	4th Position
DCA6140	5th Position
DCA6150	6th Position
DCA6300	Quad quick release/tetra clip
DCA6400	22"
DCA6410	24"
DCA6420	25"
DCA6430	26"
DCA6510	22", 24", 26" Küschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub. 22" Silver Cross-Spokes, Colourless hub and Black rim
DCA6520	24", 25" High Performance Wheel Black hub and rim
DCA6550	24", 25" Spinergy L XK, Black hub, rim and 18 spokes in Black
DCA6580	22", 24", 26" Wheels with Attendant Drum Brake
DCA6600	24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection
DCA6625	24", 25", 26" x 1" Maxxis Alpha Protect. Lightweight, Pneumatic and highly puncture resistant.
DCA6630	24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance)
DCA6660	24" x 1 3/8" Solid Profile Black Tyre
DCA6700	22", 24", 26" Supergrripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance

DCA6710	24", 26" Steel, short tab with 1 hole, mounted at 3cm distance
DCA6720	24" Titanium, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance
DCA6740	24" , 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip
DCA6750	24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands
DCA6760	24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion
DCA6770	24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim
DCA6780	22", 24", 25", 26" Light Aluminium Grey, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance
DCA6790	22", 24", 25", 26" Light Aluminium Black, short tab with 1 hole for 3cm mounting
DCA6840	Camber 3°
DCA6900	22", 24", 26" Transparent Spokeguards
DCA7230	Performance Brakes, Push to Lock. Brake barrel flips away from tyre
DCA7300	Brake Lever Extension
DCA7600	Transportation Kit Includes Tiedown Labels. Seat belt 7800 must also be selected
DCA7810	Bodypoint Evoflex posture belt, 2 point Small, flat-mount with push button
DCA7820	Bodypoint Evoflex posture belt, 2 point Medium, flat-mount with push button
DCA7830	Bodypoint Evoflex posture belt, 2 point Large, flat-mount with push button
DCA8200	Cane holder, Right
DCA8210	Cane holder, Left
DCA8450	TABLE TRAY
DCA8680	24" Mountain bike wheels with bag. fitted with specific black tyres and grey aluminium handrims; quick release axles included



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):