



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

FEBRUARY 2024
 Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

Kuschall Compact 2.0

Standard Features

- Max user weight 130kg/20.5 stone
- Crash tested ISO 7176-part 19
- Hydroformed aluminium, active, folding wheelchair, with swingaway or fixed frame
- Multiple adjustments including front and rear seat heights and centre of gravity
- Hydroformed swingaway legrests, choice of 70° and 80° angles
- Hydroformed fixed frame, choice of 70° or 80° angles
- Height adjustable flip up plastic footrests
- Choice of seat widths 280mm - 500mm (11" - 20")
- Choice of 10 seat depths 320mm - 500mm (12" - 20")
- Height adjustable fixed angle backrest set at 90°
- Standard light black nylon upholstery, wipeable
- Standard push handles
- Fixed carbon clothesguards
- 5" comfort solid front castors with black tyres
- Kuschall starec q/r rear wheels with light pneumatic tyres and aluminium handrims
- Pull to lock brakes
- Choice of frame colours

Standard Base price: £2,025



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

| Code | Description | |
|----------------------|--|--------|
| BASE PRODUCT* | | |
| DCA0012 | KÜSCHALL COMPACT 2.0 | £2,025 |
| DCA0000 | MYKÜSCHALL Special Request - Please include quotation reference in special description at the end of the form where applicable | POA |

| FRAME* | | | |
|-------------|--|---|---|
| DCA1100 | | Swing away Front frame | (NC 1900, 1902, 3040, 3140) ○ |
| DCA1110 | | Fixed front frame 70° | (NC 1700, 1705, 1710, 1715, 1720, 1725, 1730, 1805, 3000, 3010, 3020, 3030, 3100, 3110, 3120, 3130) ○ |
| DCA1120 | | Fixed front frame 80° (Minimum seat angle, front and rear seat heights is 40 mm with 5" front castors) | (NC 1700, 1705, 1710, 1715, 1720, 1725, 1730, 1805, 1962, 3000, 3010, 3020, 3030, 3100, 3110, 3120, 3130, 5550, 5650) ○ |
| DCA1130 | | Fixed front frame tapered 70° Additional 10 day leadtime | (NC 1300, 1305, 1700, 1705, 1710, 1715, 1720, 1725, 1730, 1805, 3000, 3010, 3020, 3030, 3100, 3110, 3120, 3130) ○ |
| DCA1140 | | Fixed front frame tapered 80° (minimum seat angle, front and rear seat heights is 40 mm with 5" front castors) Additional 10 day leadtime | (NC 1300, 1305, 1700, 1705, 1710, 1715, 1720, 1725, 1730, 1805, 1962, 2748, 3000, 3010, 3020, 3030, 3100, 3110, 3120, 3130, 5550, 5650) ○ |
| SEAT WIDTH* | | | |
| DCA1300 | | 280 mm Single cross brace | (NC 2210, 2490, 2710, 2720, 2748, 2750, 2748, 2750, 2860, 3020, 3030, 3120, 3130, 3260, 6000, 6010, 8450) ○ |
| DCA1305 | | 300 mm Single cross brace | (NC 2210, 2490, 2710, 2720, 2748, 2750, 2748, 2750, 3020, 3030, 3120, 3130, 3260, 6000, 6010, 8450) ○ |
| DCA1310 | | 320 mm Single cross brace | (NC 2490, 2710, 2720, 2748, 2750, 2748, 3020, 3030, 3120, 3130, 3260, 6000, 6010, 8450) ○ |
| DCA1315 | | 340 mm Single cross brace | (NC 2490, 2710, 2720, 2748, 2750, 2748, 3020, 3030, 3120, 3130, 3260, 6000, 6010, 8450) ○ |
| DCA1320 | | 360 mm Single cross brace | (NC 2490, 2710, 2748, 3020, 3030, 3120, 3130) ○ |
| DCA1325 | | 380 mm Single cross brace | (NC 2490, 3020, 3030, 3120, 3130) ○ |
| DCA1330 | | 400 mm Single cross brace | ◆ |
| DCA1335 | | 420 mm Single cross brace | ○ |
| DCA1340 | | 440 mm Double cross brace | ○ |
| DCA1345 | | 460 mm Double cross brace | ○ |
| DCA1350 | | 480 mm Double cross brace | ○ |
| DCA1355 | | 500 mm Double cross brace | (NC 2210) ○ |
| SEAT DEPTH* | | | |
| DCA1400 | | 320 mm | (NC 2155, 2160, 2710, 2720, 2748, 2750, 2860, 5650) ○ |
| DCA1405 | | 340 mm | (NC 2155, 2160, 2710, 2720, 2748, 2860, 5650) ○ |
| DCA1415 | | 360 mm | (NC 2710, 2720, 2748, 2860, 5650) ○ |
| DCA1425 | | 380 mm | (NC 2720, 5650) ○ |
| DCA1435 | | 400 mm | ◆ |
| DCA1445 | | 420 mm | ○ |
| DCA1460 | | 440 mm | ○ |
| DCA1470 | | 460 mm | ○ |
| DCA1480 | | 480 mm | (NC 2748) ○ |
| DCA1490 | | 500 mm | (NC 2710, 2748) ○ |

BACK HEIGHT Transportation kit is only possible with backrest heights of 405mm and above*

| | | | | |
|---------|--|--------|-----------------------------|---|
| DCA1615 | | 300 mm | (NC 2140, 2210, 2595, 7600) | ○ |
| DCA1620 | | 315 mm | (NC 2140, 2210, 2595, 7600) | ○ |
| DCA1625 | | 330 mm | (NC 2140, 2210, 2595, 7600) | ○ |
| DCA1630 | | 345 mm | (NC 2210, 2595, 7600) | ○ |
| DCA1635 | | 360 mm | (NC 2595, 7600) | ○ |
| DCA1640 | | 375 mm | (NC 2595, 7600) | ○ |
| DCA1645 | | 390 mm | (NC 2595, 7600) | ○ |
| DCA1650 | | 405 mm | | ○ |
| DCA1655 | | 420 mm | | ◆ |
| DCA1660 | | 435 mm | | ○ |
| DCA1665 | | 450 mm | | ○ |
| DCA1670 | | 465 mm | | ○ |
| DCA1675 | | 480 mm | | ○ |
| DCA1680 | | 495 mm | | ○ |
| DCA1685 | | 510 mm | | ○ |

FRONT SEAT HEIGHT. Seat heights 380mm TO 430mm offers a unique front frame design to achieve a lower front seat height, only available for the Compact Swing option*

| | | | | |
|---------|--|--------|---|---|
| DCA1700 | | 370 mm | (NC 1810, 1815, 1820, 1825, 1830, 1835, 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1932, 1934, 1936, 1938, 1940, 1942, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962, 3030, 3130, 3430, 6000, 6010, 6410, 6420, 7600) | ○ |
| DCA1705 | | 380 mm | (NC 1815, 1820, 1825, 1830, 1835, 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1934, 1936, 1938, 1940, 1942, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962, 3030, 3130, 3430, 6000, 6010, 6410, 6420, 7600) | ○ |
| DCA1710 | | 390 mm | (NC 1820, 1825, 1830, 1835, 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1936, 1938, 1940, 1942, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962, 3030, 3130, 3430, 6000, 6010, 6410, 6420, 7600) | ○ |
| DCA1715 | | 400 mm | (NC 1825, 1830, 1835, 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1938, 1940, 1942, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962, 3030, 3130, 3430, 6000, 6010, 6410, 6420, 7600) | ○ |
| DCA1720 | | 410 mm | (NC 1830, 1835, 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1940, 1942, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962, 3030, 3130) | ○ |
| DCA1725 | | 420 mm | (NC 1835, 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1942, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962) | ○ |
| DCA1730 | | 430 mm | (NC 1840, 1845, 1850, 1855, 1860, 1865, 1870, 1944, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962) | ○ |
| DCA1735 | | 440 mm | (NC 1845, 1850, 1855, 1860, 1865, 1870, 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962) | ○ |
| DCA1740 | | 450 mm | (NC 1850, 1855, 1860, 1865, 1870, 1948, 1950, 1952, 1954, 1956, 1958, 1960, 1962) | ○ |
| DCA1745 | | 460 mm | (NC 1855, 1860, 1865, 1870, 1950, 1952, 1954, 1956, 1958, 1960, 1962) | ○ |
| DCA1750 | | 470 mm | (NC 1860, 1865, 1870, 1952, 1954, 1956, 1958, 1960, 1962) | ○ |
| DCA1755 | | 480 mm | (NC 1805, 1865, 1870, 1954, 1956, 1958, 1960, 1962) | ○ |
| DCA1760 | | 490 mm | (NC 1805, 1810, 1956, 1958, 1960, 1962) | ○ |
| DCA1765 | | 500 mm | (NC 1805, 1810, 1815, 1958, 1960, 1962) | ◆ |
| DCA1770 | | 510 mm | (NC 1805, 1810, 1815, 1820, 1960, 1962) | ○ |
| DCA1775 | | 520 mm | (NC 1805, 1810, 1815, 1820, 1825, 1962) | ○ |
| DCA1780 | | 530 mm | (NC 1805, 1810, 1815, 1820, 1825, 1830) | ○ |

| REAR SEAT HEIGHT* | | | | |
|-------------------|--|--------|-----------------------------------|---|
| DCA1805 | | 370 mm | (NC 6000, 6010, 6410, 6420, 6430) | ○ |
| DCA1810 | | 380 mm | (NC 6000, 6010, 6410, 6420, 6430) | ○ |
| DCA1815 | | 390 mm | (NC 6420, 6430) | ○ |
| DCA1820 | | 400 mm | (NC 6420, 6430) | ○ |
| DCA1825 | | 410 mm | (NC 6420, 6430) | ○ |
| DCA1830 | | 420 mm | | ○ |
| DCA1835 | | 430 mm | | ○ |
| DCA1840 | | 440 mm | | ○ |
| DCA1845 | | 450 mm | | ◆ |
| DCA1850 | | 460 mm | | ○ |
| DCA1855 | | 470 mm | | ○ |
| DCA1860 | | 480 mm | (NC 6400) | ○ |
| DCA1865 | | 490 mm | (NC 6400) | ○ |
| DCA1870 | | 500 mm | (NC 6000, 6010, 6400, 6410) | ○ |

LOWER LEG LENGTH Seat canvas to footplate measurement, cushion height must be deducted.*

| | | | | |
|---------|--|--------|-----------------------------|---|
| DCA1900 | | 200 mm | (NC 3030, 3130) | ○ |
| DCA1902 | | 210 mm | (NC 3030, 3130) | ○ |
| DCA1904 | | 220 mm | (NC 3030, 3130) | ○ |
| DCA1906 | | 230 mm | (NC 3030, 3130) | ○ |
| DCA1908 | | 240 mm | (NC 3030, 3130) | ○ |
| DCA1910 | | 250 mm | (NC 3030, 3130) | ○ |
| DCA1912 | | 260 mm | (NC 3030, 3130) | ○ |
| DCA1914 | | 270 mm | (NC 3030, 3130) | ○ |
| DCA1916 | | 280 mm | (NC 3030, 3130) | ○ |
| DCA1918 | | 290 mm | (NC 3030, 3130) | ○ |
| DCA1920 | | 300 mm | (NC 3030, 3130) | ○ |
| DCA1922 | | 310 mm | (NC 3030, 3130) | ○ |
| DCA1924 | | 320 mm | (NC 3030, 3130) | ○ |
| DCA1926 | | 330 mm | (NC 3030, 3130) | ○ |
| DCA1928 | | 340 mm | (NC 3030, 3130) | ○ |
| DCA1930 | | 350 mm | (NC 3030, 3130) | ○ |
| DCA1932 | | 360 mm | (NC 3030, 3130) | ○ |
| DCA1934 | | 370 mm | (NC 3030, 3130) | ○ |
| DCA1936 | | 380 mm | | ○ |
| DCA1938 | | 390 mm | | ○ |
| DCA1940 | | 400 mm | | ○ |
| DCA1942 | | 410 mm | | ○ |
| DCA1944 | | 420 mm | | ○ |
| DCA1946 | | 430 mm | | ◆ |
| DCA1948 | | 440 mm | | ○ |
| DCA1950 | | 450 mm | | ○ |
| DCA1952 | | 460 mm | | ○ |
| DCA1954 | | 470 mm | | ○ |
| DCA1956 | | 480 mm | | ○ |
| DCA1958 | | 490 mm | | ○ |
| DCA1960 | | 500 mm | | ○ |
| DCA1962 | | 510 mm | (NC 3000, 3100, 3010, 3110) | ○ |

| BACKREST TYPE* | | | | |
|-----------------------------------|--|--|---|------|
| DCA2110 | | Fixed straight backrest tubes set at 90° | (NC 2155, 2160, 2165, 2170, 2175, 2180) | ◆ |
| DCA2120 | | Fixed backrest with lumbal bend angled at 7° | (NC 2155, 2160, 2165, 2170, 2175, 2180) | ○ |
| DCA2130 | | Angle adjustable back (angle must be selected from the options below) | (NC 0000) | £202 |
| DCA2140 | | Fold down, angle adjustable back, folds onto seat (angle must be selected from the options below) | (NC 0000, 2180, 2430, 3520, 3620) | £206 |
| BACKREST ANGLE* | | | | |
| DCA2155 | | Angle adjustable backrest adjusted 82° from the seat | | ○ |
| DCA2160 | | Angle adjustable backrest adjusted 86° from the seat | | ○ |
| DCA2165 | | Angle adjustable backrest adjusted 90° from the seat | | ○ |
| DCA2170 | | Angle adjustable backrest adjusted 94° from the seat | | ○ |
| DCA2175 | | Angle adjustable backrest adjusted 98° from the seat | | ○ |
| DCA2180 | | Angle adjustable backrest adjusted 102° from the seat | | ○ |
| DCA0000 | | Without angle adjustment. select with options 2110 fixed and 2120 lumbal backrests, supplied with fixed angle | | ○ |
| BACKREST UPHOLSTERY* | | | | |
| DCA2200 | | Black Standard Light Upholstery, without Tension Adjustment | (NC 2400, 2430, 2595) | ◆ |
| DCA2210 | | Lightweight back upholstery, no tension adjustment | (NC 2430, 2595) | £69 |
| DCA2220 | | Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener | | £137 |
| DCA2230 | | Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener | | £223 |
| DCA2240 | | Without Backrest Upholstery Select if fitting a Matrx/Fixed backrest, no upholstery required | (NC 2430, 2595) | ○ |
| PUSH HANDLES* | | | | |
| DCA2400 | | Without Push Handles | (NC 2490) | ○ |
| DCA2410 | | Push Handles Fixed Short/Mini Grip | (NC 2490, 6580) | ○ |
| DCA2420 | | Push Handles Fixed Long Grip | | ◆ |
| DCA2430 | | Push handles height adjustable, rear set (maximum height of push handles is rear seat height plus 600mm) | (NC 2490, 3520, 3620) | £254 |
| DCA2440 | | Push Handles Height Adjustable, Integrated in backrest | (NC 2490, 2595) | £219 |
| DCA2450 | | Push Handles Foldable | (NC 2490, 2595, 6580) | £198 |
| HEADREST | | | | |
| DCA2490 | | Headrest | (NC 2595, 6580) | £168 |
| BACKREST STABILIZATION BAR | | | | |
| DCA2595 | | Backrest Stabilization Bar, Swingaway and Removable | | £164 |
| SEAT UPHOLSTERY* | | | | |
| DCA2610 | | Standard seat upholstery with 1 central band (loop) | | ○ |
| DCA2620 | | Seat upholstery with 1 central (loop) and 2 side bands (hook) | | ◆ |
| SEAT CUSHION | | | | |
| DCA2700 | | Standard Seat Cushion, 60mm high Double Foam with Visco Top Layer - Comfort/Low Risk | (NC 2810, 2820, 2830, 2840, 2850) | £63 |
| DCA2710 | | Flo-Tech Lite Cushion 78mm high - Comfort/Low Risk | (NC 2810, 2820, 2830, 2840, 2850) | £121 |
| DCA2720 | | Matrx Libra Cushion, Black Smooth Cover 83mm High - Very High Risk | | £410 |
| DCA2748 | | Contour Visco NG Cushion 100mm High - High Risk | (NC 2810, 2820, 2830, 2840, 2850) | £169 |
| DCA2750 | | Vicair Active O2 lightweight cushion 60mm High - High risk | (NC 2810, 2820, 2830, 2840, 2850, 2860) | £521 |

| SEAT BAG | | | |
|--|--|--|------|
| DCA2800 | | Seat Bag, Click in and removable | £74 |
| MATRX LIBRA CUSHION ACCESSORIES | | | |
| DCA2810 | | Fluid supplement sac (Pair) | £44 |
| DCA2820 | | Tapered Wedge Pair | £30 |
| DCA2830 | | Abductor | £23 |
| DCA2840 | | Lateral Side Wedge | £20 |
| DCA2850 | | Anterior Front Wedge | £27 |
| SEAT RIGIDISER | | | |
| DCA2860 | | Seat Rigidiser | £39 |
| LEGREST TYPE RH* | | | |
| DCA3000 | | Legrest 70° RH (NC 3110, 9999) | ◆ |
| DCA3010 | | Legrest 80° RH (NC 3100, 5650, 9999) | ○ |
| DCA3020 | | Amputee legrest RH (NC 3280, 3300, 3305) | £138 |
| DCA3030 | | Angle adjustable/Elevating legrest - RH, including angle, width and depth adjustable plastic footplate (NC 3280, 3300, 3305) | £119 |
| DCA3040 | | Without Legrests RH, Select with fixed front frames | ○ |
| LEGREST TYPE LH* | | | |
| DCA3100 | | Legrest 70° LH (NC 9999) | ◆ |
| DCA3110 | | Legrest 80° LH (NC 5650, 9999) | ○ |
| DCA3120 | | Amputee legrest LH (NC 3280, 3300, 3305) | £138 |
| DCA3130 | | Angle adjustable/Elevating legrest - LH, including angle, width and depth adjustable plastic footplate (NC 3280, 3300, 3305) | £119 |
| DCA3140 | | Without Legrests LH, Select with fixed front frames | ○ |
| FOOTREST* | | | |
| DCA3260 | | 2 pc. Plastic flip up | ◆ |
| DCA3270 | | 2 pc. Carbon angle and depth adjustable (NC 3310) | £333 |
| DCA3280 | | 1 pc. Carbon flip-up to the left, angle and depth adjustable (NC 3310) | £167 |
| DCA9999 | | No footrest (Select when ordering Amputee brackets both sides) | ○ |
| CALF STRAP | | | |
| DCA3300 | | Calf strap mounted on the frame | ○ |
| DCA3305 | | Calf strap mounted on bracket | ○ |
| HEEL LOOPS | | | |
| DCA3310 | | Heel loops | ○ |
| CLOTHESGUARDS/MUDGUARDS We are currently experiencing an extended leadtime of 12 weeks on carbon mudguards* | | | |
| DCA3405 | | Composite clothes guards, removable (NC 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 7600, 8450) | ◆ |
| DCA3410 | | Carbon clothes guards, fixed (NC 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 7600, 8450) | ○ |
| DCA3425 | | Composite mudguards, fixed (NC 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 6000, 6010, 6140, 6150, 6580, 7600, 8450) | £65 |
| DCA3430 | | Carbon Mudguards, Fixed (Currently we have an extended leadtime of 12 weeks on carbon mudguards) (NC 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 6000, 6010, 6140, 6150, 6580, 7600, 8450) | £247 |
| DCA3440 | | Carbon Mudguards, Removable (Currently we have an extended leadtime of 12 weeks on carbon mudguards) (NC 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 6140, 6150, 8450) | £307 |
| DCA3450 | | Composite mudguards, removable (NC 3500, 3510, 3530, 3540, 3600, 3610, 3630, 3640, 6000, 6010, 6140, 6150, 6580, 8450) | £65 |
| DCA6666 | | Without clothesguards or mudguards, select only with options 3500, 3510, 3530, 3540, 3600, 3610, 3630 or 3640 Armrests supplied with clothesguards (NC 3520, 3620) | ○ |

| ARMREST RH* | | | |
|-------------|--|--|----------------------------------|
| DCA3500 | | Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, Carbon Clothesguard included RH | (NC 6150) ○ |
| DCA3510 | | Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, Carbon Clothesguard included RH | (NC 6150) ○ |
| DCA3520 | | Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included | (NC 8450) £119 |
| DCA3530 | | Küschall armrest, height adjustable, removable, swing-away, desk armpad with black neopreen sideguard upholstery, RH | (NC 6130, 6140, 6150) £212 |
| DCA3540 | | Hemi armrest, height adjustable, removable, swing-away, with hemi armpad RH | (NC 6130, 6140, 6150, 8450) £103 |
| DCA7777 | | Without armrest, select with Clothesguard/Mudguard options | (NC 3500, 3510, 3530, 3540) ○ |

| ARMREST LH* | | | |
|-------------|--|--|----------------------------------|
| DCA3600 | | Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, Carbon Clothesguard included LH | ○ |
| DCA3610 | | Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, Carbon Clothesguard included LH | ○ |
| DCA3620 | | Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included LH | (NC 8450) £119 |
| DCA3630 | | Küschall armrest, height adjustable, removable, swing-away, desk armpad with black neopreen sideguard upholstery, LH | (NC 6130, 6140, 6150) £212 |
| DCA3640 | | Küschall hemi armrest, height adjustable, removable, swing-away, with hemi armpad LH | (NC 6130, 6140, 6150, 8450) £103 |
| DCA8888 | | Without armrest, select with Clothesguard/Mudguard options | (NC 3600, 3610, 3630, 3640) ○ |

FRAME COLOUR: We are removing the Atomic Violet and Royal Blue frame colours, but will be replacing them shortly with an alternative selection of similar colours*

| | | | |
|---------|--|--|-------------|
| DCA4000 | | Space Black (RAL 9005) | ○ |
| DCA4010 | | Anthracite Grey (RAL 7016) | ○ |
| DCA4020 | | Pure White (RAL 9016) | (NC 4220) ○ |
| DCA4030 | | Carmine Red (RAL 3002) | ○ |
| DCA4040 | | Azure Blue (RAL 5017) | ○ |
| DCA4060 | | Flashy Yellow (RAL 1516) (Additional 15 day lead time) | £31 |
| DCA4100 | | Black Glossy (RAL 9005) (Additional 15 day lead time) | £31 |
| DCA4130 | | Silver Grey (RAL 9006) (Additional 15 day lead time) | £31 |
| DCA4150 | | Ruby Red (RAL 1503) (Additional 15 day lead time) | £31 |
| DCA4170 | | Dynamite Orange (Additional 15 day lead time) | £31 |
| DCA4180 | | Forest Green (RAL 6029) (Additional 15 day lead time) | £31 |

| FRAME DECAL TYPE* | | | |
|-------------------|--|---|---|
| DCA4190 | | Bold Decal (Full coloured transfer on frame) | (NC 4300, 4310, 1111) £31 |
| DCA4195 | | Tag (Small Küschall brand transfer on front of frame) | (NC 4200, 4210, 4220, 4230, 4240, 4250, 2222) ◆ |
| DCA4197 | | Without Tag (Kuschall branding only) | (NC 4200, 4210, 4220, 4230, 4240, 4250, 4300, 4310) ○ |

| BOLD DECAL COLOUR* | | | |
|--------------------|--|--|-------------|
| DCA4200 | | Black | ○ |
| DCA4210 | | Dark Grey | ○ |
| DCA4220 | | Light Grey | ○ |
| DCA4230 | | Red | ○ |
| DCA4240 | | Blue | ○ |
| DCA4250 | | Gold | ○ |
| DCA1111 | | No bold colour, select only with Option DCA4195 Tag or DCA4197 Without Tag | (NC 4190) ○ |

| TAG COLOUR* | | | |
|---|--|--|------|
| DCA4300 | | Black | ○ |
| DCA4310 | | Red | ○ |
| DCA2222 | | No tag colour, select only with Option DCA4190 Bold Decal or DCA4197 Without Tag (NC 4195) | ○ |
| FRONT CASTORS* | | | |
| DCA5050 | | 100mm x 30mm (4" x 1.2") Everyday Black hub and tyre | ○ |
| DCA5100 | | 100mm x 30mm (4" x 1.2") Comfort Black hub and tyre | ○ |
| DCA5150 | | 100mm x 25mm (4" x 1") Starec Black hub and tyre | £83 |
| DCA5200 | | 100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) | £93 |
| DCA5250 | | 100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs) | £93 |
| DCA5300 | | 125mm x 30mm (5" x 1.2") Everyday Black hub and tyre | ○ |
| DCA5350 | | 125mm x 30mm (5" x 1.2") Comfort Black hub and tyre | ◆ |
| DCA5400 | | 125mm x 25mm (5" x 1") Starec Black hub and tyre | £83 |
| DCA5450 | | 125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) | £93 |
| DCA5500 | | 125mm x 36mm (5" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs) | £93 |
| DCA5550 | | 150mm x 30mm (6" x 1.2") Comfort Black hub and tyre | ○ |
| DCA5650 | | 175mm x 30mm (7" x 1.2") Comfort Black hub and tyre | ○ |
| ONE ARM DRIVE | | | |
| DCA6000 | | One Arm Drive system Right side (with MAXXIS Lite guard pneumatic tyres, alu handrims) (NC 6040, 6130, 6140, 6150, 6300, 6400, 6420, 6430, 6510, 6520, 6550, 6580, 6625, 6630, 6660, 6840, 6900, 6700, 6710, 6720, 6740, 6750, 6760, 6770, 6780, 6790, 8680) | £482 |
| DCA6010 | | One Arm Drive system Left side (with MAXXIS Lite guard pneumatic tyres, alu handrims) (NC 6040, 6130, 6140, 6150, 6300, 6400, 6420, 6430, 6510, 6520, 6550, 6580, 6625, 6630, 6660, 6840, 6900, 6700, 6710, 6720, 6740, 6750, 6760, 6770, 6780, 6790, 8680) | £482 |
| REAR WHEEL FIXATION* | | | |
| DCA6020 | | Standard fixation, 4 rear wheel positions (NC 6140, 6150) | ◆ |
| DCA6040 | | Rear wheel extension, 6 rear wheel positions (NC 6100, 6580) | £66 |
| REAR WHEEL POSITION. Position 1 is most active and position 6 is most stable* | | | |
| DCA6100 | | 1st Position (NC 6580, 6840) | ○ |
| DCA6110 | | 2nd Position | ○ |
| DCA6120 | | 3rd Position | ◆ |
| DCA6130 | | 4th Position (NC 6580) | ○ |
| DCA6140 | | 5th Position | ○ |
| DCA6150 | | 6th Position | ○ |
| QUAD QUICK RELEASE AXLE | | | |
| DCA6300 | | Quad quick release/tetra clip (NC 6550, 6580) | £85 |
| REAR WHEEL SIZE* | | | |
| DCA6400 | | 22" (NC 6520, 6550, 6600, 6630, 6660, 6760, 6710, 6720, 6740, 6750, 6760, 6770) | ○ |
| DCA6410 | | 24" | ◆ |
| DCA6420 | | 25" (NC 6510, 6580, 6600, 6630, 6660, 6700, 6710, 6720, 6760, 6770, 6900) | ○ |
| DCA6430 | | 26" (NC 6520, 6550, 6600, 6630, 6660, 6720, 6740, 6750, 6760, 6770, 6840) | ○ |

| REAR WHEEL TYPE* | | | |
|------------------|--|---|------|
| DCA6510 | | 22", 24", 26" Kuschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub. 22" Silver Cross-Spokes, Colourless hub and Black rim | ◆ |
| DCA6520 | | 24", 25" High Performance Wheel Black hub and rim (NC 6660) | £182 |
| DCA6550 | | 24", 25" Spinergy L XK, Black hub, rim and 18 spokes in Black (NC 6660) | £611 |
| SWO0907 | | 24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKuschall (NC 6000, 6010, 6660) | £486 |
| SWO0909 | | 24", 25" Spinergy LX, Chrome hub and 12 spokes in Black. Coloured spokes are available via MyKuschall (NC 6000, 6010, 6660) | £735 |
| DCA6580 | | 22", 24", 26" Wheels with Attendant Drum Brake (NC 8680) | £366 |
| DCA3333 | | select with options 6000 and 6010 One Arm Drive wheels, specific wheel supplied | ○ |
| REAR WHEEL TYRE* | | | |
| DCA6600 | | 24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection | ◆ |
| DCA6625 | | 22", 24", 25", 26" x 1" Maxxis Alpha Protect. Lightweight, Pneumatic and highly puncture resistant. | £100 |
| DCA6630 | | 24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance) | £97 |
| DCA6660 | | 24" x 1 3/8" Solid Profile Black Tyre (NC 7230) | £66 |
| HANDRIM* | | | |
| DCA6700 | | 22", 24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance | £223 |
| DCA6710 | | 24", 26" Steel, short tab with 1 hole, mounted at 3cm distance | £62 |
| DCA6720 | | 24" Titanium, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance | £192 |
| DCA6740 | | 24" , 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip | £320 |
| DCA6750 | | 24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands | £320 |
| DCA6760 | | 24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion | £284 |
| DCA6770 | | 24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim | £356 |
| DCA6780 | | 22", 24", 25", 26" Light Aluminium Grey, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance | ◆ |
| DCA6790 | | 22", 24", 25", 26" Light Aluminium Black, short tab with 1 hole for 3cm mounting | £42 |
| DCA5555 | | select with options 6000 and 6010 one arm drive wheels, supplied with specific handrims | ○ |
| CAMBER ANGLE* | | | |
| DCA6820 | | Camber 0° | ◆ |
| DCA6840 | | Camber 3° (NC 7230) | ○ |
| SPOKEGUARD | | | |
| DCA6900 | | 22", 24", 26" Transparent Spokeguards | £126 |

| BRAKES* | | | |
|---|--|---|------|
| DCA7200 | | Pull Brakes | ◆ |
| DCA7210 | | Standard Push to Lock Brakes, Lever Bent | ○ |
| DCA7220 | | Standard Push to Lock Brakes, Lever Straight | ○ |
| DCA7230 | | Performance Brakes, Push to Lock. Brake barrel flips away from tyre (NC 7300, 7810, 7820, 7830) | £78 |
| DCA7240 | | Active/Scissor Brakes (NC 7300) | £78 |
| DCA7250 | | Light/Compact Active, Aluminium Brakes (NC 7300) | £136 |
| BRAKE LEVER EXTENSION | | | |
| DCA7300 | | Brake Lever Extension | £49 |
| SAFETY OPTIONS 2 adaptions can be mounted on the same side, this includes Anti tippers, Tipper aid and Cane holder | | | |
| DCA7500 | | Anti-tipper Mounted on Right side. Height adjustable and swingaway | £122 |
| DCA7510 | | Anti-tipper Mounted on Left side. Height adjustable and swingaway | £122 |
| DCA7520 | | Anti-tipper Mounted on both sides. Height adjustable and swingaway | £244 |
| DCA7600 | | Transportation Kit Includes Tiedown Labels. Seat belt 7800 must also be selected with this option | ○ |
| DCA7800 | | Standard Pelvic belt. Must be selected if Transportation kit is selected | £74 |
| BODYPOINT BELTS | | | |
| DCA7810 | | BODYPOINT POSTURE BELT EVOFLEX SMALL | £81 |
| DCA7820 | | BODYPOINT POSTURE BELT EVOFLEX MEDIUM | £81 |
| DCA7830 | | BODYPOINT POSTURE BELT EVOFLEX LARGE | £81 |
| OTHER OPTIONS | | | |
| DCA8000 | | Rear frame protection, plastic guard attached to underside of frame | £19 |
| DCA8010 | | Front frame protection, 2 piece removable with hook and loop | £37 |
| DCA8100 | | Tipper-aid, Right (NC 8200) | £46 |
| DCA8110 | | Tipper-aid, Left (NC 8210) | £46 |
| DCA8120 | | Tipper-aid, Both sides (NC 8200, 8210) | £92 |
| DCA8200 | | Cane holder, Right | £65 |
| DCA8210 | | Cane holder, Left | £65 |
| DCA8300 | | TRANSIT WHEELS 3" | £230 |
| DCA8450 | | TABLE TRAY | £336 |
| DCA8500 | | Air pump | £13 |
| DCA8550 | | Tool kit | ◆ |
| DCA8680 | | 24" Mountain bike wheels with bag. fitted with specific black tyres and grey aluminium handrims; quick release axles included | £320 |
| BODYPOINT ACCESSORIES | | | |
| DCA8800 | | BODYPOINT CALF STRAP | £57 |
| DCA8850 | | BODYPOINT MOBILITY BAG | £55 |

| | |
|--------------------|--|
| TOTAL PRICE | |
|--------------------|--|

Non Conforming Options:

| | |
|---------|--|
| DCA0000 | MYKÜSCHALL Special Request - Please include quotation reference in special description at the end of the form where applicable |
| DCA1300 | 280 mm Single cross brace |
| DCA1305 | 300 mm Single cross brace |
| DCA1700 | 370 mm |
| DCA1705 | 380 mm |
| DCA1710 | 390 mm |
| DCA1715 | 400 mm |
| DCA1720 | 410 mm |
| DCA1725 | 420 mm |
| DCA1730 | 430 mm |
| DCA1805 | 370 mm |
| DCA1810 | 380 mm |
| DCA1815 | 390 mm |
| DCA1820 | 400 mm |
| DCA1825 | 410 mm |
| DCA1830 | 420 mm |
| DCA1835 | 430 mm |
| DCA1840 | 440 mm |
| DCA1845 | 450 mm |
| DCA1850 | 460 mm |
| DCA1855 | 470 mm |
| DCA1860 | 480 mm |
| DCA1865 | 490 mm |
| DCA1870 | 500 mm |
| DCA1900 | 200 mm |
| DCA1902 | 210 mm |
| DCA1932 | 360 mm |
| DCA1934 | 370 mm |
| DCA1936 | 380 mm |
| DCA1938 | 390 mm |
| DCA1940 | 400 mm |
| DCA1942 | 410 mm |
| DCA1944 | 420 mm |
| DCA1946 | 430 mm |
| DCA1948 | 440 mm |
| DCA1950 | 450 mm |
| DCA1952 | 460 mm |
| DCA1954 | 470 mm |
| DCA1956 | 480 mm |
| DCA1958 | 490 mm |
| DCA1960 | 500 mm |
| DCA1962 | 510 mm |
| DCA2140 | Fold down, angle adjustable back, folds onto seat (angle must be selected from the options below) |

| | |
|---------|---|
| DCA2155 | Angle adjustable backrest adjusted 82° from the seat |
| DCA2160 | Angle adjustable backrest adjusted 86° from the seat |
| DCA2165 | Angle adjustable backrest adjusted 90° from the seat |
| DCA2170 | Angle adjustable backrest adjusted 94° from the seat |
| DCA2175 | Angle adjustable backrest adjusted 98° from the seat |
| DCA2180 | Angle adjustable backrest adjusted 102° from the seat |
| DCA0000 | Without angle adjustment. select with options 2110 fixed and 2120 lumbal backrests, supplied with fixed angle |
| DCA2210 | Lightweight back upholstery, no tension adjustment |
| DCA2400 | Without Push Handles |
| DCA2430 | Push handles height adjustable, rear set (maximum height of push handles is rear seat height plus 600mm) |
| DCA2490 | Headrest |
| DCA2595 | Backrest Stabilization Bar, Swingaway and Removable |
| DCA2710 | Flo-Tech Lite Cushion 78mm high - Comfort/Low Risk |
| DCA2720 | Matrx Libra Cushion, Black Smooth Cover 83mm High - Very High Risk |
| DCA2748 | Contour Visco NG Cushion 100mm High - High Risk |
| DCA2750 | Vicair Active O2 lightweight cushion 60mm High - High risk |
| DCA2810 | Fluid supplement sac (Pair) |
| DCA2820 | Tapered Wedge Pair |
| DCA2830 | Abductor |
| DCA2840 | Lateral Side Wedge |
| DCA2850 | Anterior Front Wedge |
| DCA2860 | Seat Rigidiser |
| DCA3000 | Legrest 70° RH |
| DCA3010 | Legrest 80° RH |
| DCA3020 | Amputee legrest RH |
| DCA3030 | Angle adjustable/Elevating legrest - RH, including angle, width and depth adjustable plastic footplate |
| DCA3040 | Without Legrests RH, Select with fixed front frames |
| DCA3100 | Legrest 70° LH |
| DCA3110 | Legrest 80° LH |
| DCA3120 | Amputee legrest LH |
| DCA3130 | Angle adjustable/Elevating legrest - LH, including angle, width and depth adjustable plastic footplate |
| DCA3140 | Without Legrests LH, Select with fixed front frames |
| DCA3260 | 2 pc. Plastic flip up |
| DCA3280 | 1 pc. Carbon flip-up to the left, angle and depth adjustable |
| DCA9999 | No footrest (Select when ordering Amputee brackets both sides) |
| DCA3300 | Calf strap mounted on the frame |
| DCA3305 | Calf strap mounted on bracket |
| DCA3310 | Heel loops |
| DCA3430 | Carbon Mudguards, Fixed (Currently we have an extended leadtime of 12 weeks on carbon mudguards) |
| DCA3500 | Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, Carbon Clothesguard included RH |
| DCA3510 | Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, Carbon Clothesguard included RH |
| DCA3520 | Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included |
| DCA3530 | Küschar armrest, height adjustable, removable, swing-away, desk armpad with black neopreen sideguard upholstery, RH |
| DCA3540 | Hemi armrest, height adjustable, removable, swing-away, with hemi armpad RH |

| | |
|---------|--|
| DCA3600 | Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, Carbon Clothesguard included LH |
| DCA3610 | Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, Carbon Clothesguard included LH |
| DCA3620 | Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included LH |
| DCA3630 | Küschall armrest, height adjustable, removable, swing-away, desk armpad with black neopreen sideguard upholstery, LH |
| DCA3640 | Küschall hemi armrest, height adjustable, removable, swing-away, with hemi armpad LH |
| DCA4190 | Bold Decal (Full coloured transfer on frame) |
| DCA4195 | Tag (Small Küschall brand transfer on front of frame) |
| DCA4200 | Black |
| DCA4210 | Dark Grey |
| DCA4220 | Light Grey |
| DCA4230 | Red |
| DCA4240 | Blue |
| DCA4250 | Gold |
| DCA1111 | No bold colour, select only with Option DCA4195 Tag or DCA4197 Without Tag |
| DCA4300 | Black |
| DCA4310 | Red |
| DCA2222 | No tag colour, select only with Option DCA4190 Bold Decal or DCA4197 Without Tag |
| DCA5550 | 150mm x 30mm (6" x 1.2") Comfort Black hub and tyre |
| DCA5650 | 175mm x 30mm (7" x 1.2") Comfort Black hub and tyre |
| DCA6000 | One Arm Drive system Right side (with MAXXIS Lite guard pneumatic tyres, alu handrims) |
| DCA6010 | One Arm Drive system Left side (with MAXXIS Lite guard pneumatic tyres, alu handrims) |
| DCA6040 | Rear wheel extension, 6 rear wheel positions |
| DCA6100 | 1st Position |
| DCA6130 | 4th Position |
| DCA6140 | 5th Position |
| DCA6150 | 6th Position |
| DCA6300 | Quad quick release/tetra clip |
| DCA6400 | 22" |
| DCA6410 | 24" |
| DCA6420 | 25" |
| DCA6430 | 26" |
| DCA6510 | 22", 24", 26" Küschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub. 22" Silver Cross-Spokes, Colourless hub and Black rim |
| DCA6520 | 24", 25" High Performance Wheel Black hub and rim |
| DCA6550 | 24", 25" Spinergy L XK, Black hub, rim and 18 spokes in Black |
| DCA6580 | 22", 24", 26" Wheels with Attendant Drum Brake |
| DCA6600 | 24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection |
| DCA6625 | 22", 24", 25", 26" x 1" Maxxis Alpha Protect. Lightweight, Pneumatic and highly puncture resistant. |
| DCA6630 | 24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance) |
| DCA6660 | 24" x 1 3/8" Solid Profile Black Tyre |
| DCA6700 | 22", 24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance |
| DCA6710 | 24", 26" Steel, short tab with 1 hole, mounted at 3cm distance |
| DCA6720 | 24" Titanium, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance |
| DCA6740 | 24" , 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip |

| | |
|---------|---|
| DCA6750 | 24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands |
| DCA6760 | 24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion |
| DCA6770 | 24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim |
| DCA6780 | 22", 24", 25", 26" Light Aluminium Grey, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance |
| DCA6790 | 22", 24", 25", 26" Light Aluminium Black, short tab with 1 hole for 3cm mounting |
| DCA6840 | Camber 3° |
| DCA6900 | 22", 24", 26" Transparent Spokeguards |
| DCA7230 | Performance Brakes, Push to Lock. Brake barrel flips away from tyre |
| DCA7300 | Brake Lever Extension |
| DCA7600 | Transportation Kit Includes Tiedown Labels. Seat belt 7800 must also be selected with this option |
| DCA7810 | BODYPOINT POSTURE BELT EVOFLEX SMALL |
| DCA7820 | BODYPOINT POSTURE BELT EVOFLEX MEDIUM |
| DCA7830 | BODYPOINT POSTURE BELT EVOFLEX LARGE |
| DCA8200 | Cane holder, Right |
| DCA8210 | Cane holder, Left |
| DCA8450 | TABLE TRAY |
| DCA8680 | 24" Mountain bike wheels with bag. fitted with specific black tyres and grey aluminium handrims; quick release axles included |



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):