



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Manual Retail Prescription Form

August 2026

Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Kuschall K-Series 2.0 One Solution

Standard Features

- Max User Weight 130kg (20.5 stone)
- Crash Tested ISO 7176-Part 19
- Active Rigid Wheelchair with Height Adjustable Rigid Footrest
- Multiple adjustments, Front and Rear Seat Heights and Centre of Gravity
- Full Hydroformed Aluminium Front Frame 75° Front Angle
- 24" Self Propel Q/R Wheels with Grey Aluminium Handrims and Pneumatic Tyres
- Choice of 10 Seat Widths 360mm to 500mm (14" -20")
- Choice of 8 Seat Depths 350mm to 525mm (13.8" -20.5")
- Range of Front Seat Heights from 450mm to 520mm
- Range of Rear Seat Heights from 390mm to 490mm
- Lockdown Folding, Angle Adjustable Fixed Height Backrest
- Without Push Handles, Carbon Clothesguards
- Black Aluminium Axle Tube with 3° Camber
- Black Comfort Front Castor
- Standard Brake Lever Bent
- Tool Kit



Standard Base price: £2,540

Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
DKA0067	KÜSCHALL K-SERIES	£2,540
DKA0074	KÜSCHALL K-SERIES CARBON (NC 1040, 1050, 1060, 1070, 1080, 1090)	£2,995
UNIQUE SPECIALS (For special modifications only) Additional lead time of 4 weeks		
DKA0000	UNIQUE Special Request - For use when requesting any special options or builds	POA
UNIQUE CUSTOMISATION AND PERSONALISATION OPTIONS Additional 4 weeks lead time for Unique options		
UKA0001	UNIQUE Personalisation - For use when selecting UKA customisation and personalisation options	☐

FRAME*				
DKA1040		V Style Standard Frame Angle 75° Aluminium (Footrest width = seat width minus 120mm)	(NC 1475, 1490, 1495, 4185, 7240)	◆
DKA1050		Tapered/Adducted Style Standard Frame Angle 75° Aluminium (Footrest width = seat width minus 100mm)	(NC 1475, 1490, 1495, 1740, 1745, 1750, 3230, 4185, 7240)	○
DKA1060		V Style +5cm Frame Angle 75° Aluminium. Frame to accommodate a longer seat depth and seat height for taller clients (Footrest width = seat width minus 120mm)	(NC 1410, 1420, 1435, 1740, 1745, 1750, 1755, 1760, 4185)	○
DKA1070		V Style Standard Frame Angle 90° Aluminium (Footrest width = seat width minus 120mm)	(NC 1475, 1490, 1495, 4185, 7240)	○
DKA1080		Tapered/Adducted Style Standard Frame Angle 90° Aluminium. Taper before the frame bend creating a parallel style front frame with a wider footrest space (Footrest width = seat width minus 100mm)	(NC 1475, 1490, 1495, 1740, 1745, 1750, 3230, 4185, 5300, 5350, 5450, 5500, 7240)	○
DKA1090		V Style +5cm Frame Angle 90° Aluminium. Frame to accommodate a longer seat depth and seat height for taller clients (Footrest width = seat width minus 120mm)	(NC 1410, 1420, 1435, 1740, 1745, 1750, 1755, 1760, 4185)	○
DKA1095		V Style Standard Frame Angle 75° Carbon (Footrest width = seat width minus 120mm)	(NC 1475, 1490, 1495, 4000, 4010, 4020, 4030, 4040, 4060, 4100, 4130, 4150, 4170, 4180, 4101, 4102, 4103, 4104, 4190, 4197, 4300, 7240)	○
DKA1100		Tapered/Adducted Style Standard Frame Angle 75° Carbon (Footrest width = seat width minus 100mm)	(NC 1475, 1490, 1495, 1740, 1745, 1750, 3230, 4000, 4010, 4020, 4030, 4040, 4060, 4100, 4130, 4150, 4170, 4180, 4101, 4102, 4103, 4104, 4190, 4197, 4300, 7240)	○
DKA1105		V Style +5cm Frame Angle 75° Carbon. Frame to accommodate a longer seat depth and seat height for taller clients (Footrest width = seat width minus 120mm)	(NC 1410, 1420, 1435, 1740, 1745, 1750, 1755, 1760, 4000, 4010, 4020, 4030, 4040, 4060, 4100, 4130, 4150, 4170, 4180, 4101, 4102, 4103, 4104, 4190, 4197, 4300)	○
DKA1110		V Style Standard Frame Angle 90° Carbon (Footrest width = seat width minus 120mm)	(NC 1475, 1490, 1495, 4000, 4010, 4020, 4030, 4040, 4060, 4100, 4130, 4150, 4170, 4180, 4101, 4102, 4103, 4104, 4190, 4197, 4300, 7240)	○
DKA1115		Tapered/Adducted Style Standard Frame Angle 90° Carbon (Footrest width = seat width minus 100mm)	(NC 1475, 1490, 1495, 1740, 1745, 1750, 3230, 4000, 4010, 4020, 4030, 4040, 4060, 4100, 4130, 4150, 4170, 4180, 4101, 4102, 4103, 4104, 4190, 4197, 4300, 5300, 5350, 5450, 5500, 7240)	○
DKA1120		V Style +5cm Frame Angle 90° Carbon. Frame to accommodate a longer seat depth and seat height for taller clients (Footrest width = seat width minus 120mm)	(NC 1410, 1420, 1435, 1740, 1745, 1750, 1755, 1760, 4000, 4010, 4020, 4030, 4040, 4060, 4100, 4130, 4150, 4170, 4180, 4101, 4102, 4103, 4104, 4190, 4197, 4300)	○
SUSPENSION				
DKA1200		Seat Suspension, Weight Range 40-65 KG	(NC 1350, 1355, 1410, 1520, 1815, 1820, 1825, 1830, 1835, 6030, 6040, 6420, 6430)	£434
DKA1210		Seat Suspension, Weight Range 65-85 KG	(NC 1350, 1355, 1410, 1520, 1815, 1820, 1825, 1830, 1835, 6030, 6040, 6420, 6430)	£434
DKA1220		Seat Suspension, Weight Range 85-100 KG	(NC 1350, 1355, 1410, 1520, 1815, 1820, 1825, 1830, 1835, 6030, 6040, 6420, 6430)	£434
SEAT WIDTH*				
DKA1320		360 mm	(NC 2364, 2366, 2368, 2370, 2372, 2710, 2748, 2750)	○
DKA1325		380 mm	(NC 2366, 2368, 2370, 2372, 2750)	○
DKA1330		400 mm	(NC 2360, 2366, 2368, 2370, 2372)	◆
DKA1335		420 mm	(NC 2360, 2362, 2368, 2370, 2372)	○
DKA1340		440 mm	(NC 2360, 2362, 2364, 2370, 2372, 2630)	○
DKA1345		460 mm	(NC 2360, 2362, 2364, 2372, 2630)	○
DKA1350		480 mm	(NC 2300, 2305, 2310, 2315, 2345, 2360, 2362, 2364, 2366, 2630, 3230)	○
DKA1355		500 mm	(NC 2300, 2305, 2310, 2315, 2345, 2360, 2362, 2364, 2366, 2368, 2630, 3210, 3230)	○
SEAT DEPTH*				
DKA1410		350 mm	(NC 1500, 1510, 2630, 2710, 2720, 2748, 2860, 7810, 7820, 7830)	○
DKA1420		375 mm	(NC 2710, 2720, 2748, 7810, 7820, 7830)	○
DKA1435		400 mm		◆
DKA1450		425 mm		○
DKA1465		450 mm	(NC 2630)	○
DKA1475		475 mm	(NC 2630, 2748)	○
DKA1490		500 mm	(NC 2630, 2710, 2748)	○
DKA1495		525 mm	(NC 2630, 2610, 2710, 2748)	○

SEAT POSITION. Seat module can be moved forward and backward along the frame to achieve a more active or stable rear wheel position*			
DKA1500		Position 1 = 140 mm MOST ACTIVE POSITION	○
DKA1510		Position 2 = 125 mm	○
DKA1520		Position 3 = 110 mm	◆
DKA1530		Position 4 = 95 mm	○
DKA1540		Position 5 = 80 mm LEAST ACTIVE POSITION	○
BACK HEIGHT. Transportation kit is only possible with backrest heights of 405mm and above*			
DKA1605		270 mm (NC 2325, 2330, 2331, 2350, 2355, 2357, 7600)	○
DKA1610		285 mm (NC 2325, 2330, 2331, 2350, 2355, 2357, 7600)	○
DKA1615		300 mm (NC 2325, 2330, 2331, 2350, 2355, 2357, 7600)	○
DKA1620		315 mm (NC 7600)	○
DKA1625		330 mm (NC 7600)	○
DKA1630		345 mm (NC 7600)	○
DKA1635		360 mm (NC 7600)	○
DKA1640		375 mm (NC 7600)	○
DKA1645		390 mm (NC 7600)	○
DKA1650		405 mm	○
DKA1655		420 mm	○
FRONT SEAT HEIGHT. For 75° frames, the maximum lower leg must be front seat height minus 30mm and minus 40mm for 90° frames*			
DKA1740		450 mm (NC 5300, 5350, 5450, 5500)	○
DKA1745		460 mm (NC 5300, 5350, 5450, 5500)	○
DKA1750		470 mm	○
DKA1755		480 mm	○
DKA1760		490 mm	○
DKA1765		500 mm	◆
DKA1770		510 mm	○
DKA1775		520 mm	○
REAR SEAT HEIGHT*			
DKA1815		390 mm (NC 3490, 3495, 6030, 6040, 6430)	○
DKA1820		400 mm (NC 3490, 3495, 6030, 6040, 6430)	○
DKA1825		410 mm (NC 3490, 3495)	○
DKA1830		420 mm	○
DKA1835		430 mm	○
DKA1840		440 mm	◆
DKA1845		450 mm	○
DKA1850		460 mm	○
DKA1855		470 mm	○
DKA1860		480 mm	○
DKA1865		490 mm (NC 6410)	○

LOWER LEG LENGTH Seat canvas to footplate measurement, cushion height must be deducted.*			
DKA1918		290 mm	○
DKA1920		300 mm	○
DKA1922		310 mm	○
DKA1924		320 mm	○
DKA1926		330 mm	○
DKA1928		340 mm	○
DKA1930		350 mm	○
DKA1932		360 mm	○
DKA1934		370 mm	○
DKA1936		380 mm	○
DKA1938		390 mm	○
DKA1940		400 mm	○
DKA1942		410 mm	○
DKA1944		420 mm	○
DKA1946		430 mm	◆
DKA1948		440 mm	○
DKA1950		450 mm	○
DKA1952		460 mm	○
DKA1954		470 mm	○
DKA1956		480 mm	○
BACKREST MATERIAL*			
DKA2000		Aluminium	◆
BACKREST TYPE*			
DKA2150		Backrest with Lower Backbrace for use with Matrix fixed back. Lockdown folding for easy transportation	◆
BACKREST UPHOLSTERY*			
DKA2240		Without Backrest Upholstery	◆

MATRIX BACKREST: Please ensure the top of the Matrix backrest is not positioned higher than 480mm, which is the highest possible setting on the K-Series.*

DKA2300		Matrx MX2, height 150 mm with quick release mini set hardware	(NC 2370, 2372, 7600)	£751
DKA2305		Matrx MX2, height 230 mm with quick release mini set hardware	(NC 2370, 2372, 7600)	£751
DKA2310		Matrx MX2, height 305 mm with quick release mini set hardware	(NC 2370, 2372)	£751
DKA2315		Matrx MX2, height 405 mm with quick release mini set hardware	(NC 2370, 2372)	£751
DKA2316		Matrx Elite Deep E2, height 160 mm with Mini-Set Mounting Hardware	(NC 7600)	£546
DKA2317		Matrx Elite Deep E2, Height 210 mm with Mini-Set Mounting Hardware	(NC 7600)	£546
DKA2318		Matrx Elite Deep E2, Height 260 mm with Mini-Set Mounting Hardware	(NC 7600)	£546
DKA2320		Matrx Elite Deep E2, Height 310 mm with Mini-Set Mounting Hardware		£546
DKA2322		Matrx Elite Deep E2, Height 360 mm with Mini-Set Mounting Hardware		£546
DKA2325		Matrx Elite Deep E2, Height 410 mm with Mini-Set Mounting Hardware		£546
DKA2330		Matrx Elite Deep E2, Height 460 mm with Easy-Set Mounting Hardware		£546
DKA2331		Matrx Elite Deep E2, Height 510 mm with Easy-Set Mounting Hardware		£546
DKA2332		Matrx Elite E2, Height 160 mm with Mini-Set Mounting Hardware	(NC 7600)	£494
DKA2334		Matrx Elite E2, Height 210 mm with Mini-Set Mounting Hardware	(NC 7600)	£494
DKA2335		Matrx Elite E2, Height 260 mm with Mini-Set Mounting Hardware	(NC 7600)	£494
DKA2340		Matrx Elite E2, Height 310 mm with Mini-Set Mounting Hardware		£494
DKA2345		Matrx Elite E2, Height 360 mm with Mini-Set Mounting Hardware		£494
DKA2350		Matrx Elite E2, Height 410 mm with Mini-Set Mounting Hardware		£494
DKA2355		Matrx Elite E2, Height 460 mm with Easy-Set Mounting Hardware		£494
DKA2357		Matrx Elite E2, Height 510 mm with Easy-Set Mounting Hardware		£494

MATRIX BACKREST WIDTH CHOICE*

DCA2360		360mm Backrest width, to fit chest width 265mm	○
DCA2362		380mm Backrest width, to fit chest width 290mm	○
DCA2364		410mm Backrest width, to fit chest width 315mm	○
DCA2366		430mm Backrest width, to fit chest width 340mm	○
DCA2368		460mm Backrest width, to fit chest width 370mm	○
DCA2370		480mm Backrest width, to fit chest width 390mm	○
DCA2372		510mm Backrest width, to fit chest width 420mm	○

PUSH HANDLES*

DKA2400		Without Push Handles	◆
DKA2410		Push Handles Fixed Short/Mini Hand Grip	£133
DKA2420		Push Handles Fixed Long Hand Grip	£133
DKA2440		Push Handles Height Adjustable, Integrated in backrest (Height range with Matrx backs 340mm to 460mm from the seat)	£219
DKA2450		Push Handles Foldable	£198

BACKREST ANGLE*				
DKA2520		78° angle from the seat		○
DKA2540		82° angle from the seat		○
DKA2560		86° angle from the seat		○
DKA2580		90° angle from the seat		◆
DKA2590		94° angle from the seat		○
SEAT UPHOLSTERY*				
DKA2600		Nylon Seat Upholstery with one 5cm central band (Loop)		○
DKA2601		Standard Seat Upholstery with 1 central loop and 2 side hook		◆
DKA2610		Tension Adjustable Seat Upholstery with one 5cm central band (Loop)		£45
DKA2630		Light Seat Upholstery, mesh material to reduce weight		£69
SEAT CUSHION				
DKA2700		Standard Seat Cushion, 60mm high Double Foam with Visco Top Layer - Comfort/Low Risk	(NC 2810, 2820, 2830, 2840, 2850, 2815, 2825, 2835, 2845, 2855)	£66
DKA2710		Flo-Tech Lite Cushion 78mm high - Comfort/Low Risk	(NC 2810, 2830, 2815, 2835)	£121
DKA2720		Matrx Libra Cushion, Black Smooth Cover 70mm High - Very High Risk		£410
DKA2748		Contour Visco NG Cushion 100mm High - High Risk	(NC 2810, 2820, 2830, 2840, 2850, 2815, 2825, 2835, 2845, 2855)	£169
DKA2750		Vicair Active O2 lightweight cushion 60mm High - High risk	(NC 2810, 2820, 2830, 2840, 2850, 2810, 2820, 2830, 2840, 2850, 2860)	£521
SEAT BAG				
DKA2800		Seat Bag, Click in and removable		£74
1ST MATRX LIBRA CUSHION POSIIONING ACCESSORIES				
DKA2810		Fluid supplement sac (Single unit)		£44
DKA2820		Tapered Wedge Pair		£30
DKA2830		Abductor		£23
DKA2840		Lateral Side Wedge		£20
DKA2850		Anterior Front Wedge		£27
2ND MATRX LIBRA CUSHION POSIIONING ACCESSORIES				
DKA2815		Fluid supplement sac (Single unit)		£44
DKA2825		Tapered Wedge Pair		£30
DKA2835		Abductor		£23
DKA2845		Lateral Side Wedge		£20
DKA2855		Anterior Front Wedge		£27
SEAT RIGIDISER				
DKA2860		Seat Rigidiser		£39
FOOTREST*				
DKA3200		Footrest Tube, Rigid - Aluminium Hydroformed profile. Weight approx. 260g		◆
DKA3210		Footrest Tube, Rigid - Titanium Rounded profile tube. Weight approx. 210g	(NC 6962)	£185
DKA3220		Carbon Footplate, Angle Adjustable	(NC 3290, 6962)	£225
DKA3230		Footplate Flip-to Side, Angle Adjustable	(NC 3290, 6962, 6963, 6964, 6965, 6966, 6967, 6968, 6969, 6970, 6971, 6972)	£143
FOOTREST COVER				
DKA3290		Footrest Cover		£37
CALF STRAP*				
DKA3300		Calf strap mounted on the frame with hook and loop		◆
DKA3305		Calf strap mounted with bracket support		○

CLOTHESGUARDS/MUDGUARDS*			
DKA3405		Composite Clothesguards, Flip up, Fold in	○
DKA3420		Carbon clothes guards, Flip up, Fold in.	◆
DKA3440		Carbon Mudguards, Removable (NC 3490, 3495)	£307
DKA3450		Composite mudguards, removable (NC 3490, 3495)	£65
ARMREST			
DKA3480		Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included	£241
DKA3490		Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included	£164
DKA3495		Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, no Clothesguard included	£164
FRAME COLOUR*			
DKA4000		Space Black (RAL 9005) (NC 4066)	○
DKA4010		Anthracite Grey (RAL 7016) (NC 4066)	○
DKA4020		Pure White (RAL 9016) (NC 4066, 4220)	○
DKA4030		Carmine Red (RAL 3002) (NC 4066)	○
DKA4040		Azure Blue (RAL 5017) (NC 4066)	○
DKA4060		Flashy Yellow (RAL 1516) (Additional 15 day lead time) (NC 4066)	£31
DKA4100		Black Glossy (RAL 9005) (Additional 15 day lead time) (NC 4066)	£31
DKA4130		Silver Grey (RAL 9006) (Additional 15 day lead time) (NC 4066)	£31
DKA4150		Ruby Red (RAL 1503) (Additional 15 day lead time) (NC 4066)	£31
DKA4170		Dynamite Orange (Additional 15 day lead time) (NC 4066)	£31
DKA4180		Forest Green (RAL 6029) (Additional 15 day lead time) (NC 4066)	£31
DKA4185		Carbon Finish	○
UKA4101		Diamond Silver/Mirror Polish with standard varnish finish (Premium ceramic glaze finish available, for this option select UKA4066)	£399
UKA4102		Pastel Green with standard varnish finish (Premium ceramic glaze finish available, for this option select UKA4066)	£299
UKA4103		Sky Blue with standard varnish finish (Premium ceramic glaze finish available, for this option select UKA4066)	£299
UKA4104		Plum Purple with standard varnish finish (Premium ceramic glaze finish available, for this option select UKA4066)	£299
UKA0008		Your personal frame colour. Add RAL colour of your choice in special comments box, for a personalised colour frame	£299
CERAMIC GLAZE FRAME FINISH			
UKA4066		Ceramic glaze finish on frame, for added durability	£336
FRAME DECAL TYPE*			
DKA4190		Bold Decal (Full coloured transfer on frame) (NC 4300, 4310, 1111)	£31
DKA4195		Tag (Small Küschall brand transfer on front of frame) (NC 4200, 4210, 4220, 4230, 4240, 4250, 2222)	◆
DKA4197		Without Tag (Kuschall branding only) (NC 4200, 4210, 4220, 4230, 4240, 4250, 4300, 4310)	○
BOLD DECAL COLOUR*			
DKA4200		Black	○
DKA4210		Dark Grey	○
DKA4220		Light Grey	○
DKA4230		Red	○
DKA4240		Blue	○
DKA4250		Gold	○
DKA1111		No bold colour, select only with Option DKA4195 Tag or DKA4197 Without Tag (NC 4190)	○

TAG COLOUR*			
DKA4300		Black	○
DKA4310		Red	○
DKA2222		No tag colour, select only with Option DKA4190 Bold Decal or DKA4197 Without Tag (NC 4195)	○
FRONT CASTORS*			
DKA5000		75mm x 25mm (3" x 1") Supersport (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)	○
DKA5050		100mm x 30mm (4" x 1.2") Everyday Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)	○
DKA5100		100mm x 30mm (4" x 1.2") Comfort Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)	○
DKA5200		100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frogs Legs) Use this option to access 4" anodised colour options	£93
DKA5250		100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs) (NC 5502, 5504, 5506, 5507, 5508, 5509)	£93
DKA5300		125mm x 30mm (5" x 1.2") Everyday Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)	○
DKA5350		125mm x 30mm (5" x 1.2") Comfort Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)	○
DKA5450		125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) Use this option to access 5" anodised colour options	£93
DKA5500		125mm x 36mm (5" x 1.4") (Frog Legs) Aluminium Softroll Silver hub, Black tyre (Frogs Legs) (NC 5502, 5504, 5506, 5507, 5508, 5509)	£93
ANODISED COLOUR SOFTROLL CASTORS			
UKA5502		Pink	£50
UKA5504		Orange	£50
UKA5506		Blue	£50
UKA5507		Purple	£50
UKA5508		Red	£50
UKA5509		Gold	£50
CERAMIC CASTOR BEARINGS			
UKA5510		Ceramic castor bearings for Aluminium soft roll castors	£58
REAR WHEEL POSITION*			
DKA6020		Standard, Fixed Rear Wheel Position (NC 6100, 6110, 6120, 6130, 6140, 6230, 6240)	◆
DKA6030		Light Rear Wheel Extension (3 Positions) (NC 6070, 6130, 6140, 6230, 6240)	£113
DKA6040		Rear Wheel Extension (5 Positions) (NC 6070)	£131
CAMBER AXLE TYPE*			
DKA6050		Vario-Ax, Adjustable Camber can achieve 1°, 3°, 7° and 10°. Only active/scissor brake will lock 10° camber position (NC 6200, 6210, 6220, 6270, 7230)	£287
DKA6060		Standard Aluminium Axle (NC 4444)	◆
DKA6070		Carbon Axle (NC 4444)	£212
POSITION OF AXLE			
DKA6100		1st Position from the front. If 2nd axle is selected, 1st axle must be placed in 1st position and 2nd in the 5th position	○
DKA6110		2nd Position from the front (NC 6230, 6240)	○
DKA6120		3rd Position from the front (NC 6230, 6240)	○
DKA6130		4th Position from the front (NC 6230, 6240)	○
DKA6140		5th Position from the front (NC 6230, 6240)	○

CAMBER ANGLE*				
DKA6200		0° Camber	(NC 6270)	○
DKA6210		3° Camber		◆
DKA6220		6° Camber	(NC 6240, 6250)	○
DKA4444		Without camber angle, select only with DKA6050 Vario-Axle		○
2ND CAMBER AXLE TYPE. If a 2nd axle is selected, 1st axle must be placed in 1st position and the 2nd axle in the 5th position				
DKA6230		Standard Aluminium Axle		£232
DKA6240		Vario-Ax, Adjustable Camber can achieve 1°, 3°, 7° and 10°	(NC 6250, 6260, 6270, 7230)	£287
CAMBER ANGLE 2ND AXLE. Max of 3° difference in camber from 1st axle. If only 1 of the axles is a Vario-Axle, 2nd axle is only possible in 0° or 3° camber				
DKA6250		0° Camber		○
DKA6260		3° Camber		○
DKA6270		6° Camber		○
DKA5555		Without camber angle, select only with DKA6240 Vario-Axle		○
QUAD QUICK RELEASE AXLE				
DKA6300		Quad Quick Release/Tetra Clip	(NC 6550, 6551, 6552, 6553, 6554, 0907)	£85
REAR WHEEL SIZE*				
DKA6410		24"		◆
DKA6420		25"	(NC 6510, 6600, 6630, 6660, 6700, 6710, 6720, 6760, 6770, 6775, 6900, 6910)	○
DKA6430		26"	(NC 6520, 6550, 6551, 6552, 6553, 6554, 6600, 6630, 6660, 6720, 6740, 6750, 6760, 6770, 6775)	○
REAR WHEEL TYPE*				
DKA6510		24", 26" Kuschall Starec Wheel, 24" Black Cross Spokes, Rim & Hub. 26" Silver Cross-Spokes, Colourless hub and Black rim	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6579)	◆
DKA6520		24", 25" High Performance Wheel Black hub and rim	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6579, 6660)	£182
DKA6550		24", 25" Spinergy LXX, Black hub, rim and 18 spokes in Black	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6660, 6910)	£611
UKA6551		24", 25" Spinergy Lite Extreme LX, Weight 687 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)	(NC 6660, 6910)	£785
UKA6552		24", 25" Spinergy Lite Extreme X-laced XLX, Weight 821 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)	(NC 6660, 6910)	£844
UKA6553		24", 25" Spinergy X-Laced Sport Lite Extreme XSLX, Weight 757 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)	(NC 6660, 6910)	£893
UKA6554		24", 25" Spinergy Carbon Fibre CLX, Weight 530 grams per wheel (Unique coloured hub and spokes available, supplied in black as standard)	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6660, 6910)	£1,084
SWO0907		24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. (Unique coloured spokes available)	(NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6660, 6910)	£486
ANODISED COLOUR SPINERGY RIM (DESIGN LAB)				
UKA6557		Blue		£179
UKA6558		Gold		£179
UKA6559		Green		£179
UKA6560		Orange		£179
UKA6561		Red		£179
UKA6562		Purple		£179

ANODISED COLOUR SPINERGY HUB (DESIGN LAB)			
UKA6563		Blue	£179
UKA6564		Gold	£179
UKA6565		Silver	£179
UKA6566		Green	£179
UKA6567		Orange	£179
UKA6568		Red	£179
UKA6569		Purple	£179
SPINERGY COLOURED SPOKES			
UKA6570		Blue	£84
UKA6571		Green	£84
UKA6572		Orange	£84
UKA6573		Pink	£84
UKA6574		Red	£84
UKA6575		Purple	£84
UKA6576		White	£84
UKA6577		Yellow	£84
CERAMIC BEARING REAR WHEEL			
UKA6579		Ceramic Bearing for Spinergy rear wheels	£170
REAR WHEEL TYRE*			
DKA6600		24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection	◆
DKA6625		24", 25", 26" x 1" Maxxis Alpha Protect. Lightweight, Pneumatic and highly puncture resistant.	£100
DKA6630		24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance)	£97
DKA6660		24" x 1 3/8" Solid Profile Black Tyre (NC 7230)	£66
UNIQUE CUSTOMISATION COLOUR PACKAGE*			
UKA6960		UNIQUE KIT 1 - Includes custom colour forks and bearing cover. Select colour from 10 choices below	£111
UKA6961		UNIQUE KIT 2 - Includes custom colour forks, bearing cover, backrest bracket, seat end caps and brake clamps. Select colour from the 10 choices below	£285
UKA6962		UNIQUE KIT 3 - Includes custom colour forks, bearing cover, backrest bracket, seat end caps, brake clamps, axle clamps, calf strap bracket, rear wheel extension. Select colour from the 10 choices below	£329
UKA1234		Without colour package	○
UNIQUE CUSTOMISATION COLOURS			
UKA6963		Red	○
UKA6964		Purple	○
UKA6965		Copper	○
UKA6966		Green	○
UKA6967		Orange	○
UKA6968		Gold	○
UKA6969		Champagne	○
UKA6970		Pink	○
UKA6971		Blue	○
UKA6972		Mirror Finish	○

HANDRIM*			
DKA6700		24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance	£223
DKA6710		24", 26" Steel, short tab with 1 hole, mounted at 3cm distance	£62
DKA6720		24" Titanium, long tab with 2 holes for 3 cm or 4 cm mountings, mounted at 3cm distance	£192
DKA6740		24", 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. For complete or limited hand function offering low-effort propulsion. Ergonomic shape with fully integrated, replaceable silicon strip for extra grip, profile of 32mm diameter.	£215
DKA6750		24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. For complete or limited hand function offering low-effort propulsion. Ergonomic shape with fully integrated, replaceable silicon strip for extra grip, smaller profile of 24mm diameter for smaller hands.	£215
DKA6760		24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape, all metal black coated, low-friction surface for complete hand function for low-effort propulsion	£190
DKA6770		24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics without tricep function. Features a unique angle shaped profile and anti-slip coating providing a larger contact surface designed to be propelled using the heel of the hand or thumb pad.	£225
DKA6775		24" Carbolife Curve L Tetra handrim, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics with tricep function, but limited hand strength. Ergonomic shape with a flat, non-slip coated pushing surface for the thumb and a curved side for the fingers to wrap around.	£155
DKA6780		24", 25", 26" Light Aluminium Grey, long tab with 2 holes for 3 cm or 4cm mountings, mounted at 3cm distance	◆
DKA6790		24", 25", 26" Light Aluminium Black, short tab with 1 hole for 3cm mounting	£42
SPOKEGUARD			
DKA6900		24", 26" Transparent Spokeguards (When ordered with Carbolife handrims, handrim distance is set at 4cm)	£126
DKA6910		24", 26" Black Spokeguard (When ordered with Carbolife handrims, handrim distance is set at 4cm)	£126
BRAKES*			
DKA7210		Standard Push to Lock Brakes, Lever Bent	◆
DKA7230		Performance Brakes, Push to Lock. Brake barrel flips away from tyre (NC 7300)	£78
DKA7240		Active/Scissor Brakes (NC 7300)	£78
DKA7250		Light/Compact Active, Aluminium Brakes (NC 7300)	£136
BRAKE LEVER EXTENSION			
DKA7300		Brake Lever Extension	£49
SAFETY OPTIONS			
DKA7500		Anti-Tipper Right (NC 7520)	£122
DKA7510		Anti-Tipper Left (NC 7520)	£122
DKA7520		Anti-Tipper Both sides (NC 8300)	£244
DKA7600		Transportation Kit Includes Tiedown Labels	○
DKA7800		Standard lap belt. We recommend the fitment of a seat belt if the Transportation kit is selected	£5

BODYPOINT BELTS		
DKA7810		Bodypoint Posture Belt Evoflex Small £81
DKA7820		Bodypoint Posture Belt Evoflex Medium £81
DKA7830		Bodypoint Posture Belt Evoflex Large £81
OTHER OPTIONS		
DKA8010		Front Frame Protectors; 2 Piece. Removable with hook and loop £37
DKA8100		Tipper-aid, Right £46
DKA8110		Tipper-aid, Left £46
DKA8120		Tipper-aid, Both sides £92
DKA8500		Air Pump £13
DKA8650		Kuschall branded Wheel bag (Fits all wheel sizes) £45
DKA8680		24" Mountain Bike Wheels with bag. Fitted with Specific Black Tyres, Grey Aluminium Handrims; Quick Release Axles Included £320
CANE HOLDER		
DKA8200		Cane Holder, Mounted Right Side £75
DKA8210		Cane Holder, Mounted Left Side £75
TRANSIT WHEELS		
DKA8300		Transit Wheels 3" £230
BODYPOINT ACCESSORIES		
DKA8800		Bodypoint Calf Strap £57
DKA8850		Bodypoint Mobility Bag £55
TOTAL PRICE		

Non Conforming Options:	
DKA1040	V Style Standard Frame Angle 75° Aluminium (Footrest width = seat width minus 120mm)
DKA1050	Tapered/Adducted Style Standard Frame Angle 75° Aluminium (Footrest width = seat width minus 100mm)
DKA1060	V Style +5cm Frame Angle 75° Aluminium. Frame to accommodate a longer seat depth and seat height for taller clients (Footrest width = seat width minus 120mm)
DKA1070	V Style Standard Frame Angle 90° Aluminium (Footrest width = seat width minus 120mm)
DKA1080	Tapered/Adducted Style Standard Frame Angle 90° Aluminium. Taper before the frame bend creating a parallel style front frame with a wider footrest space (Footrest width = seat width minus 100mm)
DKA1090	V Style +5cm Frame Angle 90° Aluminium. Frame to accommodate a longer seat depth and seat height for taller clients (Footrest width = seat width minus 120mm)
DKA1350	480 mm
DKA1355	500 mm
DKA1410	350 mm
DKA1420	375 mm
DKA1435	400 mm
DKA1475	475 mm
DKA1490	500 mm
DKA1495	525 mm
DKA1500	Position 1 = 140 mm MOST ACTIVE POSITION
DKA1510	Position 2 = 125 mm
DKA1520	Position 3 = 110 mm
DKA1740	450 mm
DKA1745	460 mm
DKA1750	470 mm
DKA1755	480 mm
DKA1760	490 mm
DKA1815	390 mm
DKA1820	400 mm
DKA1825	410 mm
DKA1830	420 mm
DKA1835	430 mm
DKA2300	Matrx MX2, height 150 mm with quick release mini set hardware
DKA2305	Matrx MX2, height 230 mm with quick release mini set hardware
DKA2310	Matrx MX2, height 305 mm with quick release mini set hardware
DKA2315	Matrx MX2, height 405 mm with quick release mini set hardware
DKA2325	Matrx Elite Deep E2, Height 410 mm with Mini-Set Mounting Hardware
DKA2330	Matrx Elite Deep E2, Height 460 mm with Easy-Set Mounting Hardware
DKA2331	Matrx Elite Deep E2, Height 510 mm with Easy-Set Mounting Hardware
DKA2345	Matrx Elite E2, Height 360 mm with Mini-Set Mounting Hardware
DKA2350	Matrx Elite E2, Height 410 mm with Mini-Set Mounting Hardware
DKA2355	Matrx Elite E2, Height 460 mm with Easy-Set Mounting Hardware
DKA2357	Matrx Elite E2, Height 510 mm with Easy-Set Mounting Hardware
DCA2360	360mm Backrest width, to fit chest width 265mm
DCA2362	380mm Backrest width, to fit chest width 290mm
DCA2364	410mm Backrest width, to fit chest width 315mm
DCA2366	430mm Backrest width, to fit chest width 340mm

DCA2368	460mm Backrest width, to fit chest width 370mm
DCA2370	480mm Backrest width, to fit chest width 390mm
DCA2372	510mm Backrest width, to fit chest width 420mm
DKA2610	Tension Adjustable Seat Upholstery with one 5cm central band (Loop)
DKA2630	Light Seat Upholstery, mesh material to reduce weight
DKA2710	Flo-Tech Lite Cushion 78mm high - Comfort/Low Risk
DKA2720	Matrx Libra Cushion, Black Smooth Cover 70mm High - Very High Risk
DKA2748	Contour Visco NG Cushion 100mm High - High Risk
DKA2750	Vicair Active O2 lightweight cushion 60mm High - High risk
DKA2810	Fluid supplement sac (Single unit)
DKA2820	Tapered Wedge Pair
DKA2830	Abductor
DKA2840	Lateral Side Wedge
DKA2850	Anterior Front Wedge
DKA2815	Fluid supplement sac (Single unit)
DKA2825	Tapered Wedge Pair
DKA2835	Abductor
DKA2845	Lateral Side Wedge
DKA2855	Anterior Front Wedge
DKA2860	Seat Rigidiser
DKA3210	Footrest Tube, Rigid - Titanium Rounded profile tube. Weight approx. 210g
DKA3230	Footplate Flip-to Side, Angle Adjustable
DKA3290	Footrest Cover
DKA3490	Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included
DKA3495	Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, no Clothesguard included
DKA4000	Space Black (RAL 9005)
DKA4010	Anthracite Grey (RAL 7016)
DKA4020	Pure White (RAL 9016)
DKA4030	Carmine Red (RAL 3002)
DKA4040	Azure Blue (RAL 5017)
DKA4060	Flashy Yellow (RAL 1516) (Additional 15 day lead time)
DKA4100	Black Glossy (RAL 9005) (Additional 15 day lead time)
DKA4130	Silver Grey (RAL 9006) (Additional 15 day lead time)
DKA4150	Ruby Red (RAL 1503) (Additional 15 day lead time)
DKA4170	Dynamite Orange (Additional 15 day lead time)
DKA4180	Forest Green (RAL 6029) (Additional 15 day lead time)
DKA4185	Carbon Finish
UKA4101	Diamond Silver/Mirror Polish with standard varnish finish (Premium ceramic glaze finish available, for this option select UKA4066)
UKA4102	Pastel Green with standard varnish finish (Premium ceramic glaze finish available, for this option select UKA4066)
UKA4103	Sky Blue with standard varnish finish (Premium ceramic glaze finish available, for this option select UKA4066)
UKA4104	Plum Purple with standard varnish finish (Premium ceramic glaze finish available, for this option select UKA4066)
UKA4066	Ceramic glaze finish on frame, for added durability
DKA4190	Bold Decal (Full coloured transfer on frame)
DKA4195	Tag (Small Kuschall brand transfer on front of frame)

DKA4197	Without Tag (Kuschall branding only)
DKA4200	Black
DKA4210	Dark Grey
DKA4220	Light Grey
DKA4230	Red
DKA4240	Blue
DKA4250	Gold
DKA1111	No bold colour, select only with Option DKA4195 Tag or DKA4197 Without Tag
DKA4300	Black
DKA4310	Red
DKA2222	No tag colour, select only with Option DKA4190 Bold Decal or DKA4197 Without Tag
DKA5300	125mm x 30mm (5" x 1.2") Everyday Black hub and tyre
DKA5350	125mm x 30mm (5" x 1.2") Comfort Black hub and tyre
DKA5450	125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) Use this option to access 5" anodised colour options
DKA5500	125mm x 36mm (5" x 1.4") (Frog Legs) Aluminium Softroll Silver hub, Black tyre (Frogs Legs)
UKA5502	Pink
UKA5504	Orange
UKA5506	Blue
UKA5507	Purple
UKA5508	Red
UKA5509	Gold
UKA5510	Ceramic castor bearings for Aluminium soft roll castors
DKA6030	Light Rear Wheel Extension (3 Positions)
DKA6040	Rear Wheel Extension (5 Positions)
DKA6070	Carbon Axle
DKA6100	1st Position from the front. If 2nd axle is selected, 1st axle must be placed in 1st position and 2nd in the 5th position
DKA6110	2nd Position from the front
DKA6120	3rd Position from the front
DKA6130	4th Position from the front
DKA6140	5th Position from the front
DKA6200	0° Camber
DKA6210	3° Camber
DKA6220	6° Camber
DKA4444	Without camber angle, select only with DKA6050 Vario-Axle
DKA6230	Standard Aluminium Axle
DKA6240	Vario-Ax, Adjustable Camber can achieve 1°, 3°, 7° and 10°
DKA6250	0° Camber
DKA6260	3° Camber
DKA6270	6° Camber
DKA6410	24"
DKA6420	25"
DKA6430	26"
DKA6510	24", 26" Kuschall Starec Wheel, 24" Black Cross Spokes, Rim & Hub. 26" Silver Cross-Spokes, Colourless hub and Black rim
DKA6520	24", 25" High Performance Wheel Black hub and rim

DKA6550	24", 25" Spinergy LXX, Black hub, rim and 18 spokes in Black
UKA6551	24", 25" Spinergy Lite Extreme LX, Weight 687 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)
UKA6552	24", 25" Spinergy Lite Extreme X-laced XLX, Weight 821 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)
UKA6553	24", 25" Spinergy X-Laced Sport Lite Extreme XSLX, Weight 757 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard)
UKA6554	24", 25" Spinergy Carbon Fibre CLX, Weight 530 grams per wheel (Unique coloured hub and spokes available, supplied in black as standard)
SWO0907	24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. (Unique coloured spokes available)
UKA6557	Blue
UKA6558	Gold
UKA6559	Green
UKA6560	Orange
UKA6561	Red
UKA6562	Purple
UKA6563	Blue
UKA6564	Gold
UKA6565	Silver
UKA6566	Green
UKA6567	Orange
UKA6568	Red
UKA6569	Purple
UKA6570	Blue
UKA6571	Green
UKA6572	Orange
UKA6573	Pink
UKA6574	Red
UKA6575	Purple
UKA6576	White
UKA6577	Yellow
UKA6579	Ceramic Bearing for Spinergy rear wheels
DKA6600	24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection
DKA6630	24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance)
DKA6660	24" x 1 3/8" Solid Profile Black Tyre
UKA6962	UNIQUE KIT 3 - Includes custom colour forks, bearing cover, backrest bracket, seat end caps, brake clamps, axle clamps, calf strap bracket, rear wheel extension. Select colour from the 10 choices below
UKA6963	Red
UKA6964	Purple
UKA6965	Copper
UKA6966	Green
UKA6967	Orange
UKA6968	Gold
UKA6969	Champagne
UKA6970	Pink
UKA6971	Blue
UKA6972	Mirror Finish

DKA6700	24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance
DKA6710	24", 26" Steel, short tab with 1 hole, mounted at 3cm distance
DKA6720	24" Titanium, long tab with 2 holes for 3 cm or 4 cm mountings, mounted at 3cm distance
DKA6740	24", 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. For complete or limited hand function offering low-effort propulsion. Ergonomic shape with fully integrated, replaceable silicon strip for extra grip, profile of 32mm diameter.
DKA6750	24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. For complete or limited hand function offering low-effort propulsion. Ergonomic shape with fully integrated, replaceable silicon strip for extra grip, smaller profile of 24mm diameter for smaller hands.
DKA6760	24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape, all metal black coated, low-friction surface for complete hand function for low-effort propulsion
DKA6770	24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics without tricep function. Features a unique angle shaped profile and anti-slip coating providing a larger contact surface designed to be propelled using the heel of the hand or thumb pad.
DKA6775	24" Carbolife Curve L Tetra handrim, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Suitable for Tetraplegics with tricep function, but limited hand strength. Ergonomonic shape with a flat, non-slip coated pushing surface for the thumb and a curved side for the fingers to wrap around.
DKA6900	24", 26" Transparent Spokeguards (When ordered with Carbolife handrims, handrim distance is set at 4cm)
DKA6910	24", 26" Black Spokeguard (When ordered with Carbolife handrims, handrim distance is set at 4cm)
DKA7230	Performance Brakes, Push to Lock. Brake barrel flips away from tyre
DKA7240	Active/Scissor Brakes
DKA7300	Brake Lever Extension
DKA7520	Anti-Tipper Both sides
DKA7600	Transportation Kit Includes Tiedown Labels
DKA7810	Bodypoint Posture Belt Evoflex Small
DKA7820	Bodypoint Posture Belt Evoflex Medium
DKA7830	Bodypoint Posture Belt Evoflex Large
DKA8300	Transit Wheels 3"



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):