



Yes, you can:

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

JANUARY 2025
 Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

Kuschall The KSL 2.0

Standard Features

- Max User Weight 100kg (16 stone)
- Active, Fully Welded, Lightweight, Rigid Wheelchair with Height Adjustable Rigid Footrest
- Fully Hydroformed Aluminium Front Frame 75° or 90° Front Angle
- 24" Self Propel Q/R Wheels with Grey Aluminium Handrims & Pneumatic Tyres
- Choice of 6 Seat Widths 340mm to 440mm (13" - 17")
- Choice of 5 Seat Depths 385mm to 485mm (15" - 19")
- Range of Front Seat Heights from 450mm to 530mm (17" - 21")
- Range of Rear Seat Heights from 400mm to 480mm (15" - 19")
- Lockdown Folding, Angle & Height Adjustable Backrest
- Black Nylon Upholstery
- Without Push Handles, Carbon Clothesguards
- Carbon Fibre Axle Tube with 3° Camber
- Black Comfort Front Castor
- Standard Brake Lever Bent
- Tool Kit



Standard Base price: £3,150

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

| Code | Description | |
|---|--|-----------------------|
| BASE PRODUCT* | | |
| DSA0085 | <input type="checkbox"/> KÜSCHALL THE KSL 2.0 | £3,150 |
| UNIQUE SPECIALS (For special modifications only) Additional lead time of 4 weeks | | |
| DSA0000 | <input type="checkbox"/> UNIQUE Special Request - For use when requesting any special options or builds | POA |
| UNIQUE CUSTOMISATION AND PERSONALISATION OPTIONS Additional 4 weeks lead time for Unique options | | |
| USA0001 | <input type="checkbox"/> UNIQUE Personalisation - For use when selecting USA customisation and personalisation options | <input type="radio"/> |

| FRAME* | | | |
|--|--|---|------|
| DSA1000 | | Aluminium Hydroformed V Style Frame Angle 75° to the ground | ◆ |
| DSA1010 | | Aluminium Hydroformed Tapered Frame Angle 75° to the ground | ○ |
| DSA1020 | | Aluminium Hydroformed V Style Frame Angle 90° to the ground | ○ |
| DSA1030 | | Aluminium Hydroformed Tapered Frame Angle 90° to the ground (NC 5300, 5350, 5400, 5450, 5500) | ○ |
| WELDED UNDERSEAT RIGID BAR* | | | |
| SWO3333 | | Under seat rigidiser bar (Welded) | £100 |
| USA5555 | | Without under seat rigidiser bar (NC 4444) | ◆ |
| SEAT WIDTH* | | | |
| DSA1315 | | 340 mm (NC 2710, 2720, 2750, 2748) | ○ |
| DSA1320 | | 360 mm (NC 2748) | ○ |
| DSA1325 | | 380 mm | ○ |
| DSA1330 | | 400 mm | ◆ |
| DSA1335 | | 420 mm | ○ |
| DSA1340 | | 440 mm | ○ |
| SEAT DEPTH* | | | |
| DSA1430 | | 385 mm (NC 2720, 7810, 7820, 7830) | ○ |
| DSA1440 | | 410 mm | ◆ |
| DSA1455 | | 435 mm | ○ |
| DSA1470 | | 460 mm | ○ |
| DSA1485 | | 485 mm (NC 2710, 2748) | ○ |
| BACK HEIGHT* | | | |
| DSA1605 | | 270 mm (NC 2240) | ○ |
| DSA1610 | | 285 mm (NC 2240) | ○ |
| DSA1615 | | 300 mm (NC 2240) | ○ |
| DSA1620 | | 315 mm (NC 2240) | ○ |
| DSA1625 | | 330 mm (NC 2240) | ○ |
| DSA1630 | | 345 mm (NC 2240) | ○ |
| DSA1635 | | 360 mm (NC 2240) | ○ |
| DSA1640 | | 375 mm (NC 2240) | ○ |
| DSA1645 | | 390 mm (NC 2240) | ○ |
| DSA1650 | | 405 mm (NC 2240) | ○ |
| DSA1655 | | 420 mm (NC 2240) | ○ |
| FRONT SEAT HEIGHT. Lower leg must be front height -100mm To -30mm on 75° frames and -80mm To -40mm on 90° frames* | | | |
| DSA1740 | | 450 mm (NC 1840, 1845, 1850, 1855, 1860, 6420) | ○ |
| DSA1745 | | 460 mm (NC 1845, 1850, 1855, 1860) | ○ |
| DSA1750 | | 470 mm (NC 1850, 1855, 1860) | ○ |
| DSA1755 | | 480 mm (NC 1855, 1860) | ○ |
| DSA1760 | | 490 mm (NC 1860) | ○ |
| DSA1765 | | 500 mm | ◆ |
| DSA1770 | | 510 mm (NC 1820) | ○ |
| DSA1775 | | 520 mm (NC 1820, 1825) | ○ |
| DSA1780 | | 530 mm (NC 1820, 1825, 1830, 6410) | ○ |

| REAR SEAT HEIGHT. Rear height = "Front height -100 MM" TO "Front height -20 MM"* | | | |
|---|--|---|--|
| DSA1820 | | 400 mm | (NC 6420) ○ |
| DSA1825 | | 410 mm | ○ |
| DSA1830 | | 420 mm | ○ |
| DSA1835 | | 430 mm | ◆ |
| DSA1840 | | 440 mm | ○ |
| DSA1845 | | 450 mm | ○ |
| DSA1850 | | 460 mm | ○ |
| DSA1855 | | 470 mm | ○ |
| DSA1860 | | 480 mm | ○ |
| LOWER LEG LENGTH Seat canvas to footplate measurement, cushion height must be deducted.* | | | |
| DSA1930 | | 350 mm | ○ |
| DSA1932 | | 360 mm | ○ |
| DSA1934 | | 370 mm | ○ |
| DSA1936 | | 380 mm | ○ |
| DSA1938 | | 390 mm | ○ |
| DSA1940 | | 400 mm | ○ |
| DSA1942 | | 410 mm | ○ |
| DSA1944 | | 420 mm | ○ |
| DSA1946 | | 430 mm | ◆ |
| DSA1948 | | 440 mm | ○ |
| DSA1950 | | 450 mm | ○ |
| DSA1952 | | 460 mm | ○ |
| DSA1954 | | 470 mm | ○ |
| DSA1956 | | 480 mm | ○ |
| DSA1958 | | 490 mm | ○ |
| DSA1960 | | 500 mm | ○ |
| BACKREST TYPE* | | | |
| DSA2110 | | Fixed Aluminium Hydroformed backrest in frame colour | (NC 6666, 7777) ◆ |
| SWO4444 | | K-Series style lockdown folding backrest, supplied with carbon fibre flip up, fold in clothesguards and tension adjustable upholstery | (NC 2200, 2220, 2240, 3410, 3430) £80 |
| SWO9999 | | K-Series style lockdown folding backrest, supplied with carbon fibre flip up, removable mudguards and tension adjustable upholstery | (NC 2200, 2220, 2240, 3410, 3430) £307 |
| BACKREST UPHOLSTERY* | | | |
| DSA2200 | | Standard Light Black Nylon upholstery (No tension adjustment) | (NC 0002, 0003, 0004, 2241, 2242, 2243, 2244, 2245) ◆ |
| DSA2220 | | Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener | (NC 0002, 0003, 0004, 2241, 2242, 2243, 2244, 2245) £137 |
| DSA2230 | | Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener | £223 |
| DSA2240 | | Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required | (NC 0002, 0003, 0004, 2241, 2242, 2243, 2244, 2245) ○ |
| SWO7777 | | Nylon tension adjustable upholstery for options 4444 & 9999 | (NC 0002, 0003, 0004, 2241, 2242, 2243, 2244, 2245) ○ |

| BACK UPHOLSTERY STITCHING COLOUR | | | |
|----------------------------------|--|--|------|
| USA2241 | | Neon Yellow | £75 |
| USA2242 | | Neon Pink | £75 |
| USA2243 | | Blue | £75 |
| USA2244 | | White | £75 |
| USA2245 | | Red | £75 |
| USA0002 | | Your personal colour stitching. Add RAL colour of your choice in special comments box, for a personalised colour stitching on back upholstery. | £75 |
| BACK UPHOLSTERY EMBROIDERY | | | |
| USA0003 | | Your personal design embroidery. Supply your design with your order, subject to copyright. Design size up to 250mm x 250mm (NC 0004) | £199 |
| BACK UPHOLSTERY FLOCKING | | | |
| USA0004 | | Your personal flock design. Supply your design with your order, subject to copyright. Design size up to 280mm x 400mm | £99 |
| PUSH HANDLES* | | | |
| DSA2400 | | Without Push Handles | ◆ |
| DSA2410 | | Push Handles Fixed Short/Mini Hand Grip | £133 |
| DSA2420 | | Push Handles Fixed Long Hand Grip | £133 |
| DSA2440 | | Push Handles Height Adjustable, Integrated in backrest | £219 |
| DSA2450 | | Push Handles Foldable | £198 |
| BACKREST ANGLE* | | | |
| DSA2500 | | 74° angle from the seat | ○ |
| DSA2520 | | 78° angle from the seat | ○ |
| DSA2540 | | 82° angle from the seat | ◆ |
| DSA2560 | | 86° angle from the seat | ○ |
| DSA2580 | | 90° angle from the seat | ○ |
| SEAT UPHOLSTERY* | | | |
| DSA2600 | | Standard Nylon Seat Upholstery with one 5cm central band (Loop) | ◆ |
| DSA2630 | | Light Seat Upholstery, mesh material to reduce weight | £69 |
| SEAT CUSHION | | | |
| DSA2700 | | Standard Seat Cushion, 60mm high Double Foam with Visco Top Layer - Comfort/Low Risk (NC 2810, 2820, 2830, 2840, 2850) | £66 |
| DSA2710 | | Flo-Tech Lite Cushion 78mm high - Comfort/Low Risk (NC 2810, 2820, 2830, 2840, 2850) | £121 |
| DSA2720 | | Matrx Libra Cushion, Black Smooth Cover 70mm High - Very High Risk | £410 |
| DSA2748 | | Contour Visco NG Cushion 100mm High - High Risk (NC 2810, 2820, 2830, 2840, 2850) | £169 |
| DSA2750 | | Vicair Active O2 lightweight cushion 60mm High - High risk (NC 2810, 2820, 2830, 2840, 2850, 2860) | £521 |
| MATRX LIBRA CUSHION ACCESSORIES | | | |
| DSA2810 | | Fluid supplement sac (Pair) | £44 |
| DSA2820 | | Tapered Wedge Pair | £30 |
| DSA2830 | | Abductor | £23 |
| DSA2840 | | Lateral Side Wedge | £20 |
| DSA2850 | | Anterior Front Wedge | £27 |
| SEAT RIGIDISER | | | |
| DSA2860 | | Seat Rigidiser | £39 |

| FOOTREST* | | | |
|-----------------------------------|--|---|------|
| DSA3200 | | Footrest Tube, Rigid - Aluminium Hydroformed profile. Weight approx. 260g | ◆ |
| DSA3210 | | Footrest Tube, Rigid - Titanium Rounded profile tube. (NC 6962) Weight approx. 210g | £185 |
| DSA3220 | | Carbon Footplate, Angle Adjustable (NC 3290, 6962) | £225 |
| FOOTREST COVER | | | |
| DSA3290 | | Footrest Cover | £37 |
| CALF STRAP* | | | |
| DSA3300 | | Calf strap Mounted on the frame with hook and loop | ◆ |
| DSA3305 | | Calf strap Mounted on the supporter with bracket | ○ |
| CLOTHESGUARDS/MUDGUARDS* | | | |
| DSA3410 | | Carbon Clothesguards, fixed between the back rest and seat rails | ◆ |
| DSA3430 | | Carbon Mudguards, fixed between the backrest and seat rails (NC 7230) | £255 |
| SWO6666 | | Carbon fibre flip up clothesguards, for use with option SWO5555 K-Series style folding backrest | ○ |
| FRAME COLOUR* | | | |
| DSA4000 | | Space Black (RAL 9005) (NC 4066) | ○ |
| DSA4010 | | Anthracite Grey (RAL 7016) (NC 4066) | ○ |
| DSA4020 | | Pure White (RAL 9016) (NC 4066, 4220) | ○ |
| DSA4030 | | Carmine Red (RAL 3002) (NC 4066) | ○ |
| DSA4040 | | Azure Blue (RAL 5017) (NC 4066) | ○ |
| DSA4060 | | Flashy Yellow (RAL 1516) (Additional 15 day lead time) (NC 4066) | £31 |
| DSA4100 | | Black Glossy (RAL 9005) (Additional 15 day lead time) (NC 4066) | £31 |
| DSA4130 | | Silver Grey (RAL 9006) (Additional 15 day lead time) (NC 4066) | £31 |
| DSA4145 | | Pearl Night Blue (Additional 15 day lead time) (NC 4066) | £31 |
| DSA4150 | | Ruby Red (RAL 1503) (Additional 15 day lead time) (NC 4066) | £31 |
| DSA4170 | | Dynamite Orange (Additional 15 day lead time) (NC 4066) | £31 |
| DSA4180 | | Forest Green (RAL 6029) (Additional 15 day lead time) (NC 4066) | £31 |
| USA4101 | | Mirror Polish with standard varnish finish (Premium ceramic glaze finish available, for this option select USA4066) | £399 |
| USA0008 | | Your personal frame colour. Add RAL colour of your choice in special comments box, for a personalised colour frame | £299 |
| CERAMIC GLAZE FRAME FINISH | | | |
| USA4066 | | Ceramic glaze finish on frame, for added durability | £336 |
| FRAME DECAL TYPE* | | | |
| DSA4190 | | Bold Decal (Full coloured transfer on frame) (NC 4300, 4310, 1111) | £31 |
| DSA4195 | | Tag (Small Kuschall brand transfer on front of frame) (NC 4200, 4210, 4220, 4230, 4240, 4250, 2222) | ◆ |
| DSA4197 | | Without Tag (Kuschall branding only) (NC 4200, 4210, 4220, 4230, 4240, 4250, 4300, 4310) | ○ |
| BOLD DECAL COLOUR* | | | |
| DSA4200 | | Black | ○ |
| DSA4210 | | Dark Grey | ○ |
| DSA4220 | | Light Grey | ○ |
| DSA4230 | | Red | ○ |
| DSA4240 | | Blue | ○ |
| DSA4250 | | Gold | ○ |
| DSA1111 | | No bold colour, select only with Option DSA4195 Tag or DSA4197 Without Tag (NC 4190) | ○ |

| TAG COLOUR* | | | |
|---|--------------------------|---|------|
| DSA4300 | <input type="checkbox"/> | Black | ○ |
| DSA4310 | <input type="checkbox"/> | Red | ○ |
| DSA2222 | <input type="checkbox"/> | No tag colour, select only with Option DSA4190 Bold Decal or DSA4197 Without Tag (NC 4195) | ○ |
| FRONT CASTORS* | | | |
| DSA5000 | <input type="checkbox"/> | 75mm x 25mm (3" x 1") Supersport (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) | ○ |
| DSA5050 | <input type="checkbox"/> | 100mm x 30mm (4" x 1.2") Everyday Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) | ○ |
| DSA5100 | <input type="checkbox"/> | 100mm x 30mm (4" x 1.2") Comfort Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) | ○ |
| DSA5150 | <input type="checkbox"/> | 100mm x 25mm (4" x 1") Starec Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) | £83 |
| DSA5200 | <input type="checkbox"/> | 100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frogs Legs) Use this option to access 4" anodised colour options | £93 |
| DSA5250 | <input type="checkbox"/> | 100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs) (NC 5502, 5504, 5506, 5507, 5508, 5509) | £93 |
| DSA5300 | <input type="checkbox"/> | 125mm x 30mm (5" x 1.2") Everyday Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) | ○ |
| DSA5350 | <input type="checkbox"/> | 125mm x 30mm (5" x 1.2") Comfort Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) | ○ |
| DSA5400 | <input type="checkbox"/> | 125mm x 25mm (5" x 1") Starec Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) | £83 |
| DSA5450 | <input type="checkbox"/> | 125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) Use this option to access 5" anodised colour options | £93 |
| DSA5500 | <input type="checkbox"/> | 125mm x 36mm (5" x 1.4") (Frog Legs) Aluminium Softroll Silver hub, Black tyre (Frogs Legs) (NC 5502, 5504, 5506, 5507, 5508, 5509) | £93 |
| ANODISED COLOUR SOFTROLL CASTORS | | | |
| USA5502 | <input type="checkbox"/> | Pink | £50 |
| USA5504 | <input type="checkbox"/> | Orange | £50 |
| USA5506 | <input type="checkbox"/> | Blue | £50 |
| USA5507 | <input type="checkbox"/> | Purple | £50 |
| USA5508 | <input type="checkbox"/> | Red | £50 |
| USA5509 | <input type="checkbox"/> | Gold | £50 |
| CERAMIC CASTOR BEARINGS | | | |
| USA5510 | <input type="checkbox"/> | Ceramic castor bearings | £58 |
| ADJUSTABLE REAR WHEEL POSITION (SELECT POSITION FROM THE 5 OPTIONS BELOW) | | | |
| USA6050 | <input type="checkbox"/> | Adjustable centre of gravity (adjustable +/- 1 position) Select from positions 1-5 below (15mm increments adjustment) | £445 |
| SEAT POSITION (Measured from centre of backpost to centre of axle)* | | | |
| DSA6100 | <input type="checkbox"/> | Position 1 = 140 mm MOST ACTIVE POSITION | ○ |
| DSA6110 | <input type="checkbox"/> | Position 2 = 125 mm | ○ |
| DSA6120 | <input type="checkbox"/> | Position 3 = 110 mm | ○ |
| DSA6130 | <input type="checkbox"/> | Position 4 = 95 mm | ◆ |
| DSA6140 | <input type="checkbox"/> | Position 5 = 80 mm LEAST ACTIVE POSITION | ○ |
| CAMBER OF CARBON AXLE* | | | |
| DSA6200 | <input type="checkbox"/> | 0° Camber | ○ |
| DSA6210 | <input type="checkbox"/> | 3° Camber | ◆ |
| DSA6220 | <input type="checkbox"/> | 6° Camber | ○ |
| QUAD QUICK RELEASE AXLE | | | |
| DSA6300 | <input type="checkbox"/> | Quad Quick Release/Tetra Clip (NC 6550, 6551, 6552, 6553, 6554, 0907) | £86 |
| REAR WHEEL SIZE* | | | |
| DSA6410 | <input type="checkbox"/> | 24" | ◆ |
| DSA6420 | <input type="checkbox"/> | 25" (NC 6510, 6600, 6630, 6700, 6710, 6720, 6760, 6770, 6900, 6910) | ○ |

| REAR WHEEL TYPE* | | | |
|--|--|--|--|
| DSA6510 | | 24" KÜschall Starec Wheel, Black Cross Spokes, Rim & Hub. | (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6579) ◆ |
| DSA6520 | | 24", 25" High Performance Wheel Black hub and rim | (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6579) £182 |
| DSA6550 | | 24", 25" Spinergy LXX, Black hub, rim and 18 spokes in Black | (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6910) £611 |
| USA6551 | | 24", 25" Spinergy Lite Extreme LX, Weight 687 grams per wheel (Unique coloured rims, hub and spokes available) | (NC 6910) £785 |
| USA6552 | | 24", 25" Spinergy Lite Extreme X-laced XLX, Weight 821 grams per wheel (Unique coloured rims, hub and spokes available) | (NC 6910) £844 |
| USA6553 | | 24", 25" Spinergy X-Laced Sport Lite Extreme XSLX, Weight 757 grams per wheel (Unique coloured rims, hub and spokes available) | (NC 6910) £893 |
| USA6554 | | 24", 25" Spinergy Carbon Fibre CLX, Weight 530 grams per wheel (Unique coloured hubs & spokes available) | (NC 6557, 6558, 6559, 6560, 6561, 6562, 6910) £1,084 |
| SWO0907 | | 24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKÜschall | (NC 6557, 6558, 6559, 6560, 6561, 6562, 6571, 6572, 6573, 6575) £486 |
| ANODISED COLOUR SPINERGY RIM (DESIGN LAB) | | | |
| USA6557 | | Blue | £179 |
| USA6558 | | Gold | £179 |
| USA6559 | | Green | £179 |
| USA6560 | | Orange | £179 |
| USA6561 | | Red | £179 |
| USA6562 | | Purple | £179 |
| ANODISED COLOUR SPINERGY HUB (DESIGN LAB) | | | |
| USA6563 | | Blue | £179 |
| USA6564 | | Gold | £179 |
| USA6565 | | Silver | £179 |
| USA6566 | | Green | £179 |
| USA6567 | | Orange | £179 |
| USA6568 | | Red | £179 |
| USA6569 | | Purple | £179 |
| SPINERGY COLOURED SPOKES | | | |
| USA6570 | | Blue | £84 |
| USA6571 | | Green | £84 |
| USA6572 | | Orange | £84 |
| USA6573 | | Pink | £84 |
| USA6574 | | Red | £84 |
| USA6575 | | Purple | £84 |
| USA6576 | | White | £84 |
| USA6577 | | Yellow | £84 |
| CERAMIC BEARING REAR WHEEL | | | |
| USA6579 | | Ceramic Bearing for Spinergy rear wheels | £170 |
| REAR WHEEL TYRE* | | | |
| DSA6600 | | 24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection | ◆ |
| DSA6625 | | 24", 25" x 1" Maxxis Alpha Protect. Pneumatic and puncture resistant. | £100 |
| DSA6630 | | 24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance) | £97 |

| UNIQUE CUSTOMISATION COLOUR PACKAGE | | | |
|-------------------------------------|--|--|------|
| USA6960 | | UNIQUE KIT 1 - Includes custom colour forks and bearing cover. Select colour from 10 choices below | £111 |
| USA6961 | | UNIQUE KIT 2 - Includes custom colour forks, bearing cover, backrest bracket with K-Series backrest, brake clamps. Select colour from the 10 choices below | £285 |
| USA6962 | | UNIQUE KIT 3 - Includes custom colour forks, bearing cover, backrest bracket with K-Series backrest, brake clamps, rear wheel axle clamps, calf strap bracket. Select colour from the 10 choices below | £329 |
| UNIQUE CUSTOMISATION COLOURS | | | |
| USA6963 | | Red | ○ |
| USA6964 | | Purple | ○ |
| USA6965 | | Copper | ○ |
| USA6966 | | Green | ○ |
| USA6967 | | Orange | ○ |
| USA6968 | | Gold | ○ |
| USA6969 | | Champagne | ○ |
| USA6970 | | Pink | ○ |
| USA6971 | | Blue | ○ |
| USA6972 | | Mirror Finish | ○ |
| HANDRIM* | | | |
| DSA6700 | | 24" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance | £223 |
| DSA6710 | | 24" Steel, mounted at 3 cm distance | £62 |
| DSA6720 | | 24" Titanium, long tab with 2 holes for 3 cm or 4 cm mountings, mounted at 3cm distance | £192 |
| DSA6740 | | 24" , 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip | £320 |
| DSA6750 | | 24" , 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands | £320 |
| DSA6760 | | 24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion | £284 |
| DSA6770 | | 24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim | £356 |
| DSA6780 | | 24", 25" Light Aluminium Grey, long tab with 2 holes for 3 cm or 4cm mountings, mounted at 3cm distance | ◆ |
| DSA6790 | | 25", 26" Light Aluminium Black, short tab with 1 hole for 3cm mounting | £42 |
| SPOKEGUARD | | | |
| DSA6900 | | 24" Transparent Spokeguard | £126 |
| DSA6910 | | 24" Black Spokeguard | £126 |
| BRAKES* | | | |
| DSA7210 | | Standard Push to Lock Brakes, Lever Bent | ◆ |
| DSA7230 | | Performance Brakes, Push to Lock. Brake barrel flips away from tyre | £78 |
| DSA7240 | | Active/Scissor Brakes | £78 |
| DSA7250 | | Light/Compact Active, Aluminium Brakes | £136 |

| SAFETY OPTIONS | | | |
|-----------------------|--|--|------|
| DSA7500 | | Anti-Tipper Right | £122 |
| DSA7510 | | Anti-Tipper Left | £122 |
| DSA7520 | | Anti-Tipper Both sides | £244 |
| DSA7810 | | Bodypoint Posture Belt Evoflex 2 points, Small, Flat-Mount with Push Button. User Weight: 25-50KG | £81 |
| DSA7820 | | Bodypoint Posture Belt Evoflex 2 points, Medium, Flat-Mount with Rehab Latch. User Weight: 50-75KG | £81 |
| DSA7830 | | Bodypoint Posture Belt Evoflex, 2 points, Large, Flat-Mount with Rehab Latch. User Weight: 75-100KG | £81 |
| OTHER OPTIONS | | | |
| DSA8010 | | Front Frame Protectors; 2 Piece. Removable with hook and loop | £37 |
| DSA8500 | | Air Pump | £13 |
| DSA8550 | | Tool Kit | ◆ |
| DSA8680 | | 24" Mountain Bike Wheels with bag. Fitted with Specific Black Tyres, Grey Aluminium Handrims; Quick Release Axles Included | £320 |
| BODYPOINT ACCESSORIES | | | |
| DSA8800 | | Bodypoint Calf Strap | £57 |
| DSA8850 | | Bodypoint Mobility Bag | £55 |

| | |
|--------------------|--|
| TOTAL PRICE | |
|--------------------|--|

Non Conforming Options:

| | |
|---------|--|
| DSA1820 | 400 mm |
| DSA1825 | 410 mm |
| DSA1830 | 420 mm |
| DSA1840 | 440 mm |
| DSA1845 | 450 mm |
| DSA1850 | 460 mm |
| DSA1855 | 470 mm |
| DSA1860 | 480 mm |
| SWO4444 | K-Series style lockdown folding backrest, supplied with carbon fibre flip up, fold in clothesguards and tension adjustable upholstery |
| DSA2200 | Standard Light Black Nylon upholstery (No tension adjustment) |
| DSA2220 | Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener |
| DSA2240 | Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required |
| SWO7777 | Nylon tension adjustable upholstery for options 4444 & 9999 |
| USA2241 | Neon Yellow |
| USA2242 | Neon Pink |
| USA2243 | Blue |
| USA2244 | White |
| USA2245 | Red |
| USA0002 | Your personal colour stitching. Add RAL colour of your choice in special comments box, for a personalised colour stitching on back upholstery. |
| USA0003 | Your personal design embroidery. Supply your design with your order, subject to copyright. Design size up to 250mm x 250mm |
| USA0004 | Your personal flock design. Supply your design with your order, subject to copyright. Design size up to 280mm x 400mm |
| DSA2710 | Flo-Tech Lite Cushion 78mm high - Comfort/Low Risk |
| DSA2720 | Matrix Libra Cushion, Black Smooth Cover 70mm High - Very High Risk |
| DSA2748 | Contour Visco NG Cushion 100mm High - High Risk |
| DSA2750 | Vicair Active O2 lightweight cushion 60mm High - High risk |
| DSA2810 | Fluid supplement sac (Pair) |
| DSA2820 | Tapered Wedge Pair |
| DSA2830 | Abductor |
| DSA2840 | Lateral Side Wedge |
| DSA2850 | Anterior Front Wedge |
| DSA2860 | Seat Rigidiser |
| DSA3290 | Footrest Cover |
| DSA3410 | Carbon Clothesguards, fixed between the back rest and seat rails |
| DSA3430 | Carbon Mudguards, fixed between the backrest and seat rails |
| SWO6666 | Carbon fibre flip up clothesguards, for use with option SWO5555 K-Series style folding backrest |
| USA4066 | Ceramic glaze finish on frame, for added durability |
| DSA4190 | Bold Decal (Full coloured transfer on frame) |
| DSA4195 | Tag (Small Kuschall brand transfer on front of frame) |
| DSA4200 | Black |
| DSA4210 | Dark Grey |
| DSA4220 | Light Grey |
| DSA4230 | Red |

| | |
|---------|---|
| DSA4240 | Blue |
| DSA4250 | Gold |
| DSA1111 | No bold colour, select only with Option DSA4195 Tag or DSA4197 Without Tag |
| DSA4300 | Black |
| DSA4310 | Red |
| DSA2222 | No tag colour, select only with Option DSA4190 Bold Decal or DSA4197 Without Tag |
| DSA5300 | 125mm x 30mm (5" x 1.2") Everyday Black hub and tyre |
| DSA5350 | 125mm x 30mm (5" x 1.2") Comfort Black hub and tyre |
| DSA5400 | 125mm x 25mm (5" x 1") Starec Black hub and tyre |
| DSA5450 | 125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) Use this option to access 5" anodised colour options |
| DSA5500 | 125mm x 36mm (5" x 1.4") (Frog Legs) Aluminium Softroll Silver hub, Black tyre (Frogs Legs) |
| USA5502 | Pink |
| USA5504 | Orange |
| USA5506 | Blue |
| USA5507 | Purple |
| USA5508 | Red |
| USA5509 | Gold |
| USA5510 | Ceramic castor bearings |
| DSA6410 | 24" |
| DSA6420 | 25" |
| DSA6510 | 24" Kuschall Starec Wheel, Black Cross Spokes, Rim & Hub. |
| DSA6550 | 24", 25" Spinergy L XK, Black hub, rim and 18 spokes in Black |
| USA6551 | 24", 25" Spinergy Lite Extreme LX, Weight 687 grams per wheel (Unique coloured rims, hub and spokes available) |
| USA6552 | 24", 25" Spinergy Lite Extreme X-laced XLX, Weight 821 grams per wheel (Unique coloured rims, hub and spokes available) |
| USA6553 | 24", 25" Spinergy X-Laced Sport Lite Extreme XSLX, Weight 757 grams per wheel (Unique coloured rims, hub and spokes available) |
| USA6554 | 24", 25" Spinergy Carbon Fibre CLX, Weight 530 grams per wheel (Unique coloured hubs & spokes available) |
| SWO0907 | 24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKuschall |
| USA6557 | Blue |
| USA6558 | Gold |
| USA6559 | Green |
| USA6560 | Orange |
| USA6561 | Red |
| USA6562 | Purple |
| USA6563 | Blue |
| USA6564 | Gold |
| USA6565 | Silver |
| USA6566 | Green |
| USA6567 | Orange |
| USA6568 | Red |
| USA6569 | Purple |
| USA6570 | Blue |
| USA6571 | Green |
| USA6572 | Orange |

| | |
|---------|--|
| USA6573 | Pink |
| USA6574 | Red |
| USA6575 | Purple |
| USA6576 | White |
| USA6577 | Yellow |
| USA6579 | Ceramic Bearing for Spinergy rear wheels |
| DSA6600 | 24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection |
| DSA6630 | 24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance) |
| USA6962 | UNIQUE KIT 3 - Includes custom colour forks, bearing cover, backrest bracket with K-Series backrest, brake clamps, rear wheel axle clamps, calf strap bracket. Select colour from the 10 choices below |
| DSA6700 | 24" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance |
| DSA6710 | 24" Steel, mounted at 3 cm distance |
| DSA6720 | 24" Titanium, long tab with 2 holes for 3 cm or 4 cm mountings, mounted at 3cm distance |
| DSA6760 | 24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion |
| DSA6770 | 24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim |
| DSA6900 | 24" Transparent Spokeguard |
| DSA6910 | 24" Black Spokeguard |
| DSA7230 | Performance Brakes, Push to Lock. Brake barrel flips away from tyre |
| DSA7810 | Bodypoint Posture Belt Evoflex 2 points, Small, Flat-Mount with Push Button. User Weight: 25-50KG |
| DSA7820 | Bodypoint Posture Belt Evoflex 2 points, Medium, Flat-Mount with Rehab Latch. User Weight: 50-75KG |
| DSA7830 | Bodypoint Posture Belt Evoflex, 2 points, Large, Flat-Mount with Rehab Latch. User Weight: 75-100KG |



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):