



**Küschall®**  
**THE KSL®**

ICONIC DESIGN  
FOR LIFE'S  
ADVENTURES



Küschall.com

PRECISE  
FIT



SUPER LIGHT  
AND RIGID



# THE KSL

Introducing 'The KSL' our lightest and most iconic award-winning Küschall chair designed for experienced users. Discover how this chair has it all, from a super light, rigid design, to a precise fit with even more options. It's time to Unlimit your World and achieve your ambitions.







## ***A MASTER IN REFINEMENT***

'The KSL' is a master in refinement. Elegant lines, discreet styling and iconic minimalism makes it an unrivalled choice. Designed to understand your individuality, 'The KSL' puts you first no matter where you go.





## *PRECISE FIT*

'The KSL' provides a tailored and personalised experience. It is built with a powerful geometry that flows seamlessly through its design. The precise fit ensures high levels of quality and functionality every time.





## ***SUPER LIGHT & EXTRA RESPONSIVE***

'The KSL' is our lightest and most agile Küschall designed with intuitive and responsive technologies. It confidently does the hard work for you and delivers a best-in-class performance.



# LATEST INNOVATIONS

## FW

The KSL with FW (Fully Welded) gives even better agility and performance in the chair. This option is a complete box frame from the front to the back and features a built-in grab bar for even easier handling.



## COG

The KSL with COG (Centre of Gravity) is a new dynamic option giving you more adjustment in 'The KSL'. Simply choose the right wheel position for you at anytime.







# TAKE A CLOSER LOOK

## ***SUPER LIGHT***

Fitted with carbon guards and axle for optimum weight reduction and drive performance. **Just 4.3 kg to transport.**

## ***RESPONSIVE***

The super-lightweight rigid frame provides effortless motion with great responsiveness indoors and outdoors.

## ***LATEST INNOVATIONS***

New COG (Centre of Gravity) and FW (Fully Welded) options.



## ***ERGONOMIC GRIP***

For easier repositioning and transfers.

## ***PRECISE FIT***

Shaped around you and your needs with no adjustments required.

## ***EASY TO HANDLE***

With the lightweight frame and optional welded grab bar.

# CUSTOMISE YOUR KSL

## 1. FRAME COLOUR

Standard frame colours



Pure White   Anthracite Grey   Space Black



Carmine Red   Azure Blue

Extra frame colours\*




Flashy Yellow   Black Glossy   Silver Grey   Pearl Night Blue



Ruby Red   Dynamite Orange   Forest Green

\*Extra delivery time will apply.


## 2. DECAL STYLE



## 3. DECAL COLOUR



Red   Black



Light Grey   Gold   Red   Blue   Dark Grey   Black

## 4. UNIQUE SERVICE

Create your new Küschall wheelchair from the wheels up with UNIQUE. This dedicated service will create your personalized chair, designed to fit you and your lifestyle.

**UNIQUE**  
Designed by you. Built by us.



Pictured: Spinergy wheel with personalised features from our UNIQUE service.







# TECHNICAL SPECIFICATION

## SEAT WIDTH

340 - 440 mm  
In increments of 20 mm



## SEAT DEPTH

385 - 485 mm  
In increments of 25 mm



## KNEE-TO-HEEL LENGTH

350 - 500 mm  
In increments of 10 mm



## BACKREST HEIGHT

270 - 420 mm  
In increments of 15 mm



## BACKREST ANGLE

74° / 78° / 82° / 86° / 90°



## SEAT HEIGHT

Front: 450 - 530 mm  
Rear: 400 - 480 mm  
In increments of 10 mm



## FRAME ANGLE

75° / 90°  
V-frame or tapered



## TOTAL WIDTH

SW + 170 mm



## TOTAL LENGTH

75° approx. 800 mm  
90° approx. 730 mm



## CENTER OF GRAVITY

80mm / 95mm / 110mm /  
125mm / 140mm



## CAMBER

0° / 3° / 6°



## TOTAL WEIGHT UNLOADED

Approx. 6.8 kg



## TRANSPORT WEIGHT

Approx. 4.3 kg



## MAX. USER WEIGHT

100 kg



All measurements and weights indicated are based on seat width 400 mm and seat depth 400 mm in the lightest and most compact configuration.  
Weight and measurements can change with the different configurations. For further information please check the user manual on our website.



**UNLIMIT YOUR WORLD**

**Tel: +41 61 487 70 70**  
info@kuschall.com  
www.kuschall.com

Invacare International GmbH  
Neuhofweg 51  
4147 Aesch BL  
Switzerland



Küschall, Küschall Unlimit your world, K, Unlimit your world and The KSL are trademarks of Invacare International GmbH. Every effort has been made to ensure that the contents of this publication are fully up-to-date at the time of going to print. As part of its ongoing improvement of products, Invacare International GmbH reserves the right to modify existing models at any time.

©April 2024, Invacare International GmbH. All rights reserved. EU – 01/2026