

# Invacare Leo<sup>®</sup>

Feel the Freedom



**Compact, safe with an exceptional drive**

Bringing independence to those moments that matter



**Yes, you can.**

# Feel the difference

Our scooters have been designed together with healthcare experts, to provide unrivalled standards of quality, safety and comfort.

Through continuous dedicated research and an unrivalled experience in designing for mobility, at Invacare we pride ourselves in scooters that improve quality of life. Being able to regain confidence, trust and independence from your scooter is important.

From transportable designs, cushioned seats and adjustable seat risers, every detail of our scooters is considered and refined to make getting about feel great.



## Quality



Market-leading build quality with exceptional attention to detail.

## Safety



Compliant with the latest European Medical Device Standards.

## Comfort



Adaptive and adjustable design to provide lasting comfort.

# Feel in control with Invacare Leo®

The Invacare Leo® is a compact scooter that gives a smooth, safe drive. Simply dismantle it for easy transportation and take it with you for lots of exciting days out and those moments that matter.

Leo® is packed with many safety and performance features. Its intelligent lighting package includes indicators that automatically switch off and a rear brake light for safe driving. A visual warning for the battery level display, alerts

with an audible warning (beep) when the battery needs recharging. Incredibly easy to use the Leo® puts you in control with confidence.

up to  
8km/h  
max speed

up to  
38km  
max distance  
on a single  
charge

136 kg  
max user  
weight



### Easy to use and transport

- Dismantles into five parts for easy transportation
- Ergonomic steering for added control
- Simple to steer and drive



### Safety in mind

- Two-step disengaging lever to prevent free-wheeling
- Brake lights, headlights and indicators
- Additional optional extras available



### Built around comfort

- High quality cushioned seat
- Adjustable seat height
- Flip up armrests and swivel seat makes getting on and off easy



### Adjustable steering for greater comfort



### Easy battery charge

with off-board  
charger kit



### Anti-splash guards

protects the scooter from  
dirt and water



Swivel seat for easy transfer



Easy to transport in your car



Two-step disengaging lever to prevent free-wheeling



### Add to your scooter

Tailor the Invacare Leo® to your taste and needs. Turn over to explore the optional extras available.

## Add to your scooter

We have a wide range of optional extras available, so you can customize your Invacare Leo® scooter. Ask your retailer to discover the full range.



Large mirror



Crutch holder



Seatbelt



Seat suspension module



Lockable front box



Windscreen

## Technical data

For more information about this product, including the user manual, please visit our website [www.invacare.eu.com](http://www.invacare.eu.com)



Backrest width ▼

470 mm



Seat depth ▼

410 mm



Seat to floor height ▼

460 - 560 mm



Backrest height ▼

475 mm



Driving unit width ▼

590 mm



Total length ▼

1220 mm



Empty weight ▼

59 kg (Excl. batteries)  
83 kg (Excl. batteries)



Max. user weight ▼

136 kg



Battery capacity ▼

2 x 36 Ah



Turning diameter ▼

2620 mm



Max. climbable obstacle height ▼

60 mm



Max. safe slope ▼

10°



Motor capacity ▼

240 W



Max. driving distance ▼

38 km



Max. driving speed ▼

8 km/h

## Colours\*



Moonstone Silver



Onyx Blue



Ruby Red



Black Vinyl

\* Scooter colours may vary slightly. Please check with your retailer.

## Upholstery

**Tel: +41 61 487 70 70**  
hqeurope@invacare.com  
[www.invacare.eu.com](http://www.invacare.eu.com)

Invacare International GmbH  
Neuhofweg 51  
4147 Aesch BL  
Switzerland

