

Maxxcare Heel Pro Evolution replacement SmartTubes

Air-filled inserts for advanced heel offloading and pressure redistribution

The **Maxxcare SmartTubes** are designed to maintain the clinical performance of the Maxxcare Heel Pro Evolution heel boot. These air-filled tubes play a crucial role in achieving effective heel offloading by lifting the heel and redistributing pressure away from vulnerable areas.



Based on proven **Vicair Technology**, the SmartTubes provide a lightweight and comfortable solution that supports both prevention and treatment of heel pressure injuries.

Engineered for stability. Designed for comfort.

Each replacement set consists of four air-filled tubes that fit precisely into the Maxxcare Heel Pro Evolution:

- 1 **Two side SmartTubes positioned along the side of the foot**
These tubes are size specific and vary depending on the boot size
- 2 **Two rear SmartTubes located at the back of the lower leg**
These tubes are universal in size and remain across all boot sizes



Effective heel offloading



Lifts the heel to reduce pressure on vulnerable areas.

Pressure redistribution



Supports skin integrity and helps prevent and treat pressure injuries.

Lightweight design



Enhances patient comfort and supports compliance.

Proven technology



Developed with the same proven technology trusted across the Vicair portfolio.

Seamless compatibility



Designed specifically for the Maxxcare Heel Pro Evolution heel boot.