


Please refer to www.invacare.ie for Local Distributor

Customer name:	Account No :
Contact no:	Order No :
Delivery address:	Order Qty :
	Order Date :

Matrx Libra
Ireland order Form
Maximum User Weight 21 stone (136kg) / HD 36 stone (227kg)


		WIDTH										
cm		36	38	41	43	46	48	51	53	56		
inch		14	15	16	17	18	19	20	21	22		
DEPTH	30	12										
	36	14										
	41	16										
	43	17										
	46	18										
	51	20										
56	22											



Ordering the correct size

Please note the complete product dimensions.
The example shown is of a 16"x16" cushion part number LC1616.
Please allow an additional 1 1/2" to the required seat depth for the rear cushion radius.

The rear cushion radius is a curved extension designed to maximise surface area and minimise the gap between seat and backrest; while stabilising fluid in tilted positions, with coccyx cut out for added relief.



Features:
Very High Risk
3" height
3 Year Warranty
Anatomically contoured foam base Innovative dual layer fluid sac
Choice of 3 covers
Positioning accessories available




Code	X	Description	Price
Libra Cushion			
LC ____ (Width x Depth)		Matrx Libra standard sizes with black smooth cover	€645
LCW ____ (Width x Depth)		Matrx Libra standard sizes with black wipedown cover	€645
LCS ____ (Width x Depth)		Matrx Libra standard sizes with black spacer cover	€645
Heavy Duty sizes			
LC2117		Heavy Duty Libra 21x17" (53x43cm)	Black Smooth cover €846
LC2118		Heavy Duty Libra 21x18" (53x46cm)	Black Smooth cover €846
LC2120		Heavy Duty Libra 21x20" (53x51cm)	Black Smooth cover €846
LC2122		Heavy Duty Libra 21x22" (53x56cm)	Black Smooth cover €846
LC2217		Heavy Duty Libra 22x17" (56x43cm)	Black Smooth cover €846
LC2218		Heavy Duty Libra 22x18" (56x46cm)	Black Smooth cover €846
LC2220		Heavy Duty Libra 22x20" (56x51cm)	Black Smooth cover €846
LC2222		Heavy Duty Libra 22x22" (56x56cm)	Black Smooth cover €846

LCW2117	Heavy Duty Libra 21x17" (53x43cm)	Black Wipedown Cover	€846
LCW2118	Heavy Duty Libra 21x18" (53x46cm)	Black Wipedown Cover	€846
LCW2120	Heavy Duty Libra 21x20" (53x51cm)	Black Wipedown Cover	€846
LCW2122	Heavy Duty Libra 21x22" (53x56cm)	Black Wipedown Cover	€846
LCW2217	Heavy Duty Libra 22x18" (56x43cm)	Black Wipedown Cover	€846
LCW2218	Heavy Duty Libra 22x18" (56x46cm)	Black Wipedown Cover	€846
LCW2220	Heavy Duty Libra 22x20" (56x51cm)	Black Wipedown Cover	€846
LCW2222	Heavy Duty Libra 22x22" (56x56cm)	Black Wipedown Cover	€846
LCS2117	Heavy Duty Libra 21x17" (53x43cm)	Spacer Fabric Cover	€846
LCS2118	Heavy Duty Libra 21x18" (53x46cm)	Spacer Fabric Cover	€846
LCS2120	Heavy Duty Libra 21x20" (53x51cm)	Spacer Fabric Cover	€846
LCS2122	Heavy Duty Libra 21x22" (53x56cm)	Spacer Fabric Cover	€846
LCS2217	Heavy Duty Libra 22x17" (56x43cm)	Spacer Fabric Cover	€846
LCS2218	Heavy Duty Libra 22x18" (56x46cm)	Spacer Fabric Cover	€846
LCS2220	Heavy Duty Libra 22x20" (56x51cm)	Spacer Fabric Cover	€846
LCS2222	Heavy Duty Libra 22x22" (56x56cm)	Spacer Fabric Cover	€846

Accessories



LABD	Matrix Abductor For comfortable mid-line positioning of thighs and hips (Provides 1 ^{1/8"} (29mm) increased height abductor contour)		€51
1487157	Anterior Seat Wedge Increases slope, prevents sliding and reduces shear (1" Anterior wedge increases pre-ischial contour - Underside mounted)		€40
1487158	Tapered Wedges (pair) For obliquity and pelvic rotation when used in rear (3/4" Tapered wedge - Underside mounted)		€51
1487155	Lateral Wedge Reinforces lateral stability for enhanced midline positioning (1" Lateral wedge increases hip/thigh loading, - Underside mounted)		€44
1578357	Libra Fluid Supplement For targeted protection/obliquity accommodation (May be used to accommodate mild pelvic obliquity/additional pressure reduction)		€70
1589366	Libra Accessory Kit Comprises of all the above items in an assessment presentation box with instruction guide. Parts included: LABD, 1487157, 1487158, 1487155, 1578357		€156
	Spare Platilon Inner Liner Provides added protection to the foam cushion from infection control issues.		€15
	Additional Cover (Widths 14" to 20")	<input type="checkbox"/> Black Smooth <input type="checkbox"/> Black Spacer <input type="checkbox"/> Black Wipedown	€108
	Additional Cover (Widths 21" and over)	<input type="checkbox"/> Black Smooth <input type="checkbox"/> Black Spacer <input type="checkbox"/> Black Wipedown	€128
1494494	14" x 16" Rigidiser		€79
1487165	16" x 16" Rigidiser		€79
1494495	16" x 18" Rigidiser		€79
1494496	16" x 20" Rigidiser		€79
1494497	18" x 16" Rigidiser		€79
1494498	18" x 18" Rigidiser		€79
1494499	18" x 20" Rigidiser		€79
1494500	20" x 16" Rigidiser		€79
1494501	20" x 18" Rigidiser		€79
1494502	20" x 20" Rigidiser		€79
1494503	22" x 18" Rigidiser		€79
1494504	22" x 20" Rigidiser		€79

TOTAL PRICE:

Special Instructions - Extended lead times will apply

Name: _____ **Email:** _____ **Phone:** _____



Description of special request (please include relevant option code):