

Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
www.invacare.co.uk

Matrx MX2 Back
January 2025

Account Name :	Account No :
Address :	Order No :
Tel No:	Order Qty :
Delivery address :	Order Date :

Matrx MX2
Retail order form
Max.user weight: 113kg / 17 stone





		WIDTH				
HEIGHT	mm	355	380	405	430	455
	inch	14	15	16	17	18
	150	6				
	230	9				
	310	12				
	410	16				



Standard Features:

- Lightweight carbon-fibre back shell
- Anatomical contouring with 76 mm depth for enhanced lateral stability
- Angle and depth adjustable
- Cutouts for breathability to help manage microclimate
- Slotted design to allow height adjustments even when in the chair
- Breathable Spacer cover
- Backrests 230mm or higher come with leatherette side panels
- Oversized foam for enhanced comfort and protection

Code	X	Description	Price
MX2 Back shell			
1629002		MX2 Back (35.5w x 15h)	14x6" £519.00
1629003		MX2 Back (35.5w x 23h)	14x9" £519.00
1629004		MX2 Back (35.5w x 31h)	14x12" £519.00
1629005		MX2 Back (35.5w x 41h)	14x16" £519.00
1629006		MX2 Back (38w x 15h)	15x6" £519.00
1629007		MX2 Back (38w x 23h)	15x9" £519.00
1629008		MX2 Back (38w x 31h)	15x12" £519.00
1629009		MX2 Back (38w x 41h)	15x16" £519.00
1629010		MX2 Back (40.5w x 15h)	16x6" £519.00
1629011		MX2 Back (40.5w x 23h)	16x9" £519.00
1629012		MX2 Back (40.5w x 31h)	16x12" £519.00
1629013		MX2 Back (40.5w x 41h)	16x16" £519.00
1629014		MX2 Back (43w x 15h)	17x6" £519.00
1629015		MX2 Back (43w x 23h)	17x9" £519.00
1629016		MX2 Back (43w x 31h)	17x12" £519.00
1629017		MX2 Back (43w x 41h)	17x16" £519.00

1629018		MX2 Back (45.5w x 15h)	18x6"	£519.00
1629019		MX2 Back (45.5w x 23h)	18x9"	£519.00
1629020		MX2 Back (45.5w x 31h)	18x12"	£519.00
1629021		MX2 Back (45.5w x 41h)	18x16"	£519.00
Hardware				
		<p>Hardware is width-adjustable allowing the backrest to be fitted to wheelchairs that are up to 25mm narrower or wider than the actual backrest width for a precise fit and positioning in rigid or foldable wheelchairs with backrest tube diameter of 19 mm, 22 mm and 25 mm.</p> <div>  <p>Fixed hardware</p>  <p>Mini set hardware</p> </div>		
1629027		Fixed carbon fibre hardware		£153.00
1629028		Mini set (quick release) hardware		£206.00
Privacy flap				
1551799		Extra Small	Fits 35.5cm width back	£49.00
1551800		Small	Fits 40.5cm width back	£49.00
1551903		Medium	Fits 43 and 45.5cm width back	£49.00