

Please refer to www.invacare.ie for Local Distributor

Account Name :	Account No :
Address :	Order No :
Tel No:	Order Qty :
Delivery address :	Order Date :

Matrx MX2
Ireland order form
Max. user weight: 113kg (17 stone)



		WIDTH					
		mm	355	380	405	430	455
		inch	14	15	16	17	18
HEIGHT	150	6					
	230	9					
	310	12					
	410	16					



Standard Features:
 Lightweight carbon-fibre back shell
 Anatomical contouring with 76 mm depth for enhanced lateral stability
 Angle and depth adjustable
 Cutouts for breathability to help manage microclimate
 Slotted design to allow height adjustments even when in the chair
 Breathable Spacer cover
 Backrests 230mm or higher come with leatherette side panels
 Oversized foam for enhanced comfort and protection
 Hardware not included

Code	X	Description	Price
MX2 Back shell			
1629002		MX2 Back (35.5w x 15h)	14x6" €1,116.00
1629003		MX2 Back (35.5w x 23h)	14x9" €1,116.00
1629004		MX2 Back (35.5w x 31h)	14x12" €1,116.00
1629005		MX2 Back (35.5w x 41h)	14x16" €1,116.00
1629006		MX2 Back (38w x 15h)	15x6" €1,116.00
1629007		MX2 Back (38w x 23h)	15x9" €1,116.00
1629008		MX2 Back (38w x 31h)	15x12" €1,116.00
1629009		MX2 Back (38w x 41h)	15x16" €1,116.00
1629010		MX2 Back (40.5w x 15h)	16x6" €1,116.00
1629011		MX2 Back (40.5w x 23h)	16x9" €1,116.00
1629012		MX2 Back (40.5w x 31h)	16x12" €1,116.00
1629013		MX2 Back (40.5w x 41h)	16x16" €1,116.00
1629014		MX2 Back (43w x 15h)	17x7" €1,116.00

1629015	MX2 Back (43w x 23h)	17x9"	€1,116.00
1629016	MX2 Back (43w x 31h)	17x12"	€1,116.00
1629017	MX2 Back (43w x 41h)	17x16"	€1,116.00
1629018	MX2 Back (45.5w x 15h)	18x6"	€1,116.00
1629019	MX2 Back (45.5w x 23h)	18x9"	€1,116.00
1629020	MX2 Back (45.5w x 31h)	18x12"	€1,116.00
1629021	MX2 Back (45.5w x 41h)	18x16"	€1,116.00

Hardware

Hardware is width-adjustable allowing the backrest to be fitted to wheelchairs that are up to 25mm narrower or wider than the actual backrest width for a precise fit and positioning in rigid or foldable wheelchairs with backrest tube diameter of 19 mm, 22 mm and 25 mm.



Fixed hardware



Mini set hardware

1629027	Fixed carbon fibre hardware		€221.00
1629028	Mini set (quick release) hardware		€177.00
1551799	Extra Small	Fits 35.5cm width back	€80.00
1551800	Small	Fits 40.5cm width back	€80.00
1551903	Medium	Fits 43 and 45.5cm width back	€80.00

Special Instructions - Extended lead times will apply



Name: _____ **Email:** _____ **Phone:** _____

Description of special request (please include relevant option code):